



EPOCH is a level 1, phase III randomized controlled trial using transarterial radiation therapy for mCRC liver metastases that demonstrated statistically significant improvements in both Progression-Free Survival (PFS) and Hepatic Progression-Free Survival (hPFS) in patients who progressed on first-line chemotherapy.

Mulcahy, M. et al, Radioembolization With Chemotherapy for Colorectal Liver Metastases: A Randomized, Open-Label, International, Multicenter, Phase III Trial (EPOCH). Journal of Clinical Oncology, 20 Sept 2021.

TRIAL OBJECTIVE

To evaluate the safety and efficacy of TheraSphere Y-90 Glass Microspheres combined with second-line therapy (oxaliplatin- or irinotecan-based chemotherapy) in patients with mCRC of the liver.


TRIAL DESIGN

An open-label, prospective, multicenter, phase III trial of 428 patients randomized 1:1 to treatment arm (TheraSphere + second-line chemotherapy) vs. control arm (second-line chemotherapy alone) across 95 centers in 12 countries, including North America, Europe and Asia.

PRIMARY ENDPOINTS

Progression-free survival (PFS) and hepatic PFS (hPFS)

- Time from randomization to progression by RECIST 1.1 or death
- Blinded independent central review (BICR)



1. TARE with Y-90 glass microspheres (TheraSphere™, Boston Scientific Corporation). Cycle 1= chemotherapy, Y-90 TARE replaces Cycle 2, Cycle 3 resume chemotherapy ± targeted therapy. ClinicalTrials.gov Identifier: NCT01483027. Chauhan N, Mulcahy MF, Salem R, et al. JMIR Res Protoc. 2019;8(1):e11545. doi: 10.2196/11545.

EPOCH TRIAL

Trial Objective & Design

Trial Conclusions

PRIMARY ENDPOINTS

Progression-Free Survival

Hepatic Progression-Free Survival

Subgroup Analyses for PFS & hPFS

SECONDARY ENDPOINTS

Overall Survival

Tumor Response

Time to Deterioration of Quality of Life

ADDITIONAL ANALYSES

Time to Progression

Hepatic Time to Progression

Time to Start of Subsequent Therapy

Post-Hoc Subgroup Analyses

KEY PATIENT & DISEASE CHARACTERISTICS

TREATMENT CHARACTERISTICS, DOSIMETRY & SAFETY



EPOCH TRIAL CONCLUSIONS

- Both primary endpoints were successfully met. Patients receiving TheraSphere Y-90 with second-line chemo were:
 - 31% less likely to experience disease progression or death (due to any cause) vs. chemo alone
 - 41% less likely to experience hepatic disease progression or death (due to any cause) vs. chemo alone
- The addition of TheraSphere Y-90 to second-line chemotherapy increased median Time to Progression (TTP) by 2.1 months* and increased median Hepatic Time to Progression (hTTP) by 4.9 months*
- Patients receiving TheraSphere Y-90 with second-line chemotherapy showed an Objective Response Rate (ORR) of 34.0% vs. 21.1% for the control arm*
- The addition of TheraSphere Y-90 to second-line chemotherapy:
 - Extended median Time to Subsequent Therapy by 10.9 months*
 - Did not compromise patients' ability to receive chemotherapy ± biologics
 - Did not compromise Quality of Life
 - Did not increase chemotherapy-related adverse events and no new safety signals were identified
- Subgroups receiving TheraSphere Y-90 with second-line chemo showed improved benefit in PFS, hPFS, and additional time to deterioration of QoL vs. chemo alone, and also showed greater magnitude in benefit compared to the overall ITT population (Subgroup A: excludes ECOG 1 and CEA ≥35 ng/mL. Subgroup B: excludes ECOG 1, CEA ≥35 ng/mL, and KRAS-m)
- Post-hoc safety analyses of overall ITT population showed patients with <10% tumor volume replacement and/or >10 lesions treated with TheraSphere Y-90 + second-line chemo experienced more liver-related TEAEs. Sequential lobar treatment, as opposed to same day whole liver treatment, may mitigate liver-related TEAEs.

EPOCH TRIAL

Trial Objective & Design

Trial Conclusions

PRIMARY ENDPOINTS

Progression-Free Survival

Hepatic Progression-Free Survival

Subgroup Analyses for PFS & hPFS

SECONDARY ENDPOINTS

Overall Survival

Tumor Response

Time to Deterioration of Quality of Life

ADDITIONAL ANALYSES

Time to Progression

Hepatic Time to Progression

Time to Start of Subsequent Therapy

Post-Hoc Subgroup Analyses

KEY PATIENT & DISEASE CHARACTERISTICS

TREATMENT CHARACTERISTICS, DOSIMETRY & SAFETY

*Indicates improvement in the treatment arm (TheraSphere + chemotherapy) compared to the control arm (chemotherapy alone) corresponding to 1-sided p <0.025






PRIMARY ENDPOINTS¹

EPOCH demonstrated statistically significant improvements in both primary endpoints of PFS and hPFS in patients with colorectal liver metastases.

Progression Free Survival

31%

Patients receiving TheraSphere with second-line chemo were **31% less likely** to experience disease progression or death (due to any cause) vs. chemo alone.



EPOCH TRIAL

Trial Objective & Design

Trial Conclusions

PRIMARY ENDPOINTS

Progression-Free Survival

Hepatic Progression-Free Survival

Subgroup Analyses for PFS & hPFS

SECONDARY ENDPOINTS

Overall Survival

Tumor Response

Time to Deterioration of Quality of Life

ADDITIONAL ANALYSES

Time to Progression

Hepatic Time to Progression

Time to Start of Subsequent Therapy

Post-Hoc Subgroup Analyses

KEY PATIENT & DISEASE CHARACTERISTICS

TREATMENT CHARACTERISTICS, DOSIMETRY & SAFETY

1. Time from randomization to progression according to RECIST 1.1 by Blinded Independent Central Review (BICR) or death (due to any cause), whichever occurred first.






PRIMARY ENDPOINTS¹

EPOCH demonstrated statistically significant improvements in both primary endpoints of PFS and hPFS in patients with colorectal liver metastases.

Hepatic Progression Free Survival

41%

Patients receiving TheraSphere with second-line chemo were **41% less likely** to experience hepatic disease progression or death (due to any cause) vs. chemo alone.



EPOCH TRIAL

Trial Objective & Design

Trial Conclusions

PRIMARY ENDPOINTS

Progression-Free Survival

Hepatic Progression-Free Survival

Subgroup Analyses for PFS & hPFS

SECONDARY ENDPOINTS

Overall Survival

Tumor Response

Time to Deterioration of Quality of Life

ADDITIONAL ANALYSES

Time to Progression

Hepatic Time to Progression

Time to Start of Subsequent Therapy

Post-Hoc Subgroup Analyses


KEY PATIENT & DISEASE CHARACTERISTICS

TREATMENT CHARACTERISTICS, DOSIMETRY & SAFETY

1. Time from randomization to hepatic progression according to RECIST 1.1 by Blinded Independent Central Review (BICR) or death (due to any cause), whichever occurred first.



SUBGROUP ANALYSES FOR PFS & hPFS¹



Continued

EPOCH TRIAL

Trial Objective & Design

Trial Conclusions

PRIMARY ENDPOINTS

Progression-Free Survival

Hepatic Progression-Free Survival

Subgroup Analyses for PFS & hPFS

SECONDARY ENDPOINTS

Overall Survival

Tumor Response

Time to Deterioration of Quality of Life

ADDITIONAL ANALYSES

Time to Progression

Hepatic Time to Progression

Time to Start of Subsequent Therapy

Post-Hoc Subgroup Analyses

KEY PATIENT & DISEASE CHARACTERISTICS


TREATMENT CHARACTERISTICS, DOSIMETRY & SAFETY

1. According to RECIST 1.1 by Blinded Independent Central Review (BICR).





SUBGROUP ANALYSES FOR PFS & hPFS¹



[Previous](#) [Continued](#)

EPOCH TRIAL

Trial Objective & Design

Trial Conclusions

PRIMARY ENDPOINTS

Progression-Free Survival

Hepatic Progression-Free Survival

Subgroup Analyses for PFS & hPFS

SECONDARY ENDPOINTS

Overall Survival

Tumor Response

Time to Deterioration of Quality of Life

ADDITIONAL ANALYSES

Time to Progression

Hepatic Time to Progression

Time to Start of Subsequent Therapy

Post-Hoc Subgroup Analyses

KEY PATIENT & DISEASE CHARACTERISTICS


TREATMENT CHARACTERISTICS, DOSIMETRY & SAFETY

1. According to RECIST 1.1 by Blinded Independent Central Review (BICR).
2. LLN = lower limit of normal.
3. CEA = carcinoembryonic antigen.





SUBGROUP ANALYSES FOR PFS & hPFS¹



[Previous](#)

EPOCH TRIAL

Trial Objective & Design

Trial Conclusions

PRIMARY ENDPOINTS

Progression-Free Survival

Hepatic Progression-Free Survival

Subgroup Analyses for PFS & hPFS

SECONDARY ENDPOINTS

Overall Survival

Tumor Response

Time to Deterioration of Quality of Life

ADDITIONAL ANALYSES

Time to Progression

Hepatic Time to Progression

Time to Start of Subsequent Therapy

Post-Hoc Subgroup Analyses

KEY PATIENT & DISEASE CHARACTERISTICS

TREATMENT CHARACTERISTICS, DOSIMETRY & SAFETY


1. According to RECIST 1.1 by Blinded Independent Central Review (BICR).





OVERALL SURVIVAL (INTENT TO TREAT)¹

There was no statistically significant difference in OS between the treatment and control arms in the intent to treat (ITT) population.



TheraSphere + Chemo:
Median 14.0 months (95% CI: 11.8, 15.5)

Chemo:
Median 14.4 months (95% CI: 12.8, 16.4)

Hazard Ratio 1.07 (95% CI: 0.86, 1.32)
1-sided *p* = 0.7229

N at Risk

	0	6	12	18	24	30	36	42	48	54	60	66
TheraSphere + Chemo	215	183	112	71	43	17	8	7	3	3	2	0
Chemo	213	164	115	62	38	25	12	10	3	0	0	0

OVERALL SURVIVAL (PER PROTOCOL)^{1,2}

TheraSphere + Chemo: Median 15.2 months (95% CI: 12.7, 17.7) **Chemo:** Median 14.3 months (95% CI: 12.6, 16.4) **Hazard Ratio 0.96** (95% CI: 0.74, 1.24)
1-sided *p* = 0.3841

EPOCH TRIAL

Trial Objective & Design

Trial Conclusions

PRIMARY ENDPOINTS

Progression-Free Survival

Hepatic Progression-Free Survival

Subgroup Analyses for PFS & hPFS

SECONDARY ENDPOINTS

Overall Survival

Tumor Response

Time to Deterioration of Quality of Life

ADDITIONAL ANALYSES

Time to Progression

Hepatic Time to Progression

Time to Start of Subsequent Therapy

Post-Hoc Subgroup Analyses

KEY PATIENT & DISEASE CHARACTERISTICS

TREATMENT CHARACTERISTICS, DOSIMETRY & SAFETY

1. Time from randomization to death or last date known alive in absence of death.

2. OS Per Protocol: TheraSphere + Chemo (N=100) and Chemo (N=40) patients excluded from Per Protocol analysis due to major deviations.






TUMOR RESPONSE¹


Patients receiving TheraSphere Y-90 with second-line chemotherapy showed an Objective Response Rate (ORR) of 34.0% vs. 21.1% for the control arm; a difference of 12.8%.

BEST OVERALL RESPONSE




OBJECTIVE RESPONSE RATE (CR+PR)

1-sided $p = 0.0019$




DISEASE CONTROL RATE (CR+PR+SD)

1-sided $p = 0.0626$



MEDIAN DURATION OF OBJECTIVE RESPONSE OR DISEASE CONTROL²



Duration of Disease Control was longer in the TheraSphere + Chemo group; however, Duration of Response in responders was not different between the two groups.

1. According to RECIST 1.1 by Blinded Independent Central Review (BICR).
2. By Kaplan-Meier analysis.
3. Time from first date of overall response of CR or PR by BICR until date of PD by BICR or death, whichever occurred first.
4. Time from first date of overall response of CR, PR, or SD by BICR until date of PD by BICR or death, whichever occurred first.

EPOCH TRIAL

Trial Objective & Design

Trial Conclusions

PRIMARY ENDPOINTS

Progression-Free Survival

Hepatic Progression-Free Survival

Subgroup Analyses for PFS & hPFS

SECONDARY ENDPOINTS

Overall Survival

Tumor Response

Time to Deterioration of Quality of Life

ADDITIONAL ANALYSES

Time to Progression

Hepatic Time to Progression

Time to Start of Subsequent Therapy

Post-Hoc Subgroup Analyses


KEY PATIENT & DISEASE CHARACTERISTICS

TREATMENT CHARACTERISTICS, DOSIMETRY & SAFETY



TIME TO DETERIORATION OF QUALITY OF LIFE (TTDQoL)¹

The addition of TheraSphere Y-90 to second-line chemotherapy did not compromise quality of life.



TheraSphere + Chemo
Chemo

N at Risk

TheraSphere + Chemo	215	43	11	7	2
Chemo	213	22	1	1	0

EPOCH TRIAL

Trial Objective & Design

Trial Conclusions

PRIMARY ENDPOINTS

Progression-Free Survival

Hepatic Progression-Free Survival

Subgroup Analyses for PFS & hPFS

SECONDARY ENDPOINTS

Overall Survival

Tumor Response

Time to Deterioration of Quality of Life

ADDITIONAL ANALYSES

Time to Progression

Hepatic Time to Progression

Time to Start of Subsequent Therapy

Post-Hoc Subgroup Analyses

KEY PATIENT & DISEASE CHARACTERISTICS

TREATMENT CHARACTERISTICS, DOSIMETRY & SAFETY


1. Time from randomization to the change from baseline in FACT-c total score \leq -7 points or death, whichever occurred first.





TIME TO PROGRESSION (TTP)¹

The addition of TheraSphere Y-90 to second-line chemotherapy increased median TTP by 2.1 months.



	0	6	12	18	24	30	36	42	48
N at Risk									
TheraSphere + Chemo	215	126	39	19	8	1	1	1	0
Chemo	213	93	26	8	5	3	1	1	0

EPOCH TRIAL

Trial Objective & Design

Trial Conclusions

PRIMARY ENDPOINTS

Progression-Free Survival

Hepatic Progression-Free Survival

Subgroup Analyses for PFS & hPFS

SECONDARY ENDPOINTS

Overall Survival

Tumor Response

Time to Deterioration of Quality of Life

ADDITIONAL ANALYSES

Time to Progression

Hepatic Time to Progression

Time to Start of Subsequent Therapy

Post-Hoc Subgroup Analyses

KEY PATIENT & DISEASE CHARACTERISTICS

TREATMENT CHARACTERISTICS, DOSIMETRY & SAFETY


1. Time from randomization to progression according to RECIST 1.1 by Blinded Independent Central Review (BICR) or death (due to any cause), whichever occurred first.





HEPATIC TIME TO PROGRESSION (hTTP)¹

The addition of TheraSphere Y-90 to second-line chemotherapy increased median hTTP by 4.9 months.



TheraSphere + Chemo
Chemo

N at Risk

TheraSphere + Chemo	215	134	43	20	10	1	1	1	1	0
Chemo	213	93	26	8	5	3	1	1	0	0

EPOCH TRIAL

Trial Objective & Design

Trial Conclusions

PRIMARY ENDPOINTS

Progression-Free Survival

Hepatic Progression-Free Survival

Subgroup Analyses for PFS & hPFS

SECONDARY ENDPOINTS

Overall Survival

Tumor Response

Time to Deterioration of Quality of Life

ADDITIONAL ANALYSES

Time to Progression

Hepatic Time to Progression

Time to Start of Subsequent Therapy

Post-Hoc Subgroup Analyses

KEY PATIENT & DISEASE CHARACTERISTICS


TREATMENT CHARACTERISTICS, DOSIMETRY & SAFETY

1. Time from randomization to hepatic progression according to RECIST 1.1 by Blinded Independent Central Review (BICR) or death (due to any cause), whichever occurred first.



TIME TO START OF SUBSEQUENT THERAPY¹

The addition of TheraSphere Y-90 to second-line chemotherapy extended median time to subsequent therapy by 10.9 months.



	N at Risk										
	0	6	12	18	24	30	36	42	48	54	60
TheraSphere + Chemo	215	170	84	45	23	8	4	3	2	2	1
Chemo	213	138	51	22	15	9	2	2	1	0	0

EPOCH TRIAL

Trial Objective & Design

Trial Conclusions

PRIMARY ENDPOINTS

Progression-Free Survival

Hepatic Progression-Free Survival

Subgroup Analyses for PFS & hPFS

SECONDARY ENDPOINTS

Overall Survival

Tumor Response

Time to Deterioration of Quality of Life

ADDITIONAL ANALYSES

Time to Progression

Hepatic Time to Progression

Time to Start of Subsequent Therapy

Post-Hoc Subgroup Analyses

KEY PATIENT & DISEASE CHARACTERISTICS

TREATMENT CHARACTERISTICS, DOSIMETRY & SAFETY


1. Time from randomization to start of the subsequent mCRC therapy (i.e. a complete change in the treatment regimen or addition of another locoregional therapy, including ablation or resection).



POST-HOC SUBGROUP ANALYSES: PFS, hPFS, TTDQoL, & OS

Patients receiving TheraSphere Y-90 with second-line chemotherapy showed improved PFS and hPFS benefit, and additional time to deterioration of quality of life (TTDQoL) in subgroups vs. chemo alone, and a greater magnitude in benefit compared to the overall intent to treat (ITT) population.

Two subgroup populations were identified based on three prognostic factors that impact TTDQoL



Outcome (median, months)	Overall ITT Population ¹		Subgroup A ²		Subgroup B ²	
	TheraSphere + Chemo N=215	Chemo N=213	TheraSphere + Chemo N=143	Chemo N=160	TheraSphere + Chemo N=77	Chemo N=91
PFS	8.0	7.2	9.4	7.6	11.6	8.5
	Difference: +0.8 months HR: 0.69 (95% CI: 0.54, 0.88) 1-sided p = 0.0013		Difference: +1.8 months HR: 0.64 (95% CI: 0.47, 0.87) 1-sided p = 0.0020		Difference: +3.1 months HR: 0.60 (95% CI: 0.39, 0.92) 1-sided p = 0.0089	
hPFS	9.1	7.2	10.8	7.6	12.5	8.5
	Difference: +1.9 months HR: 0.59 (95% CI: 0.46, 0.77) 1-sided p < 0.0001		Difference: +3.2 months HR: 0.53 (95% CI: 0.39, 0.73) 1-sided p < 0.0001		Difference: +4.0 months HR: 0.51 (95% CI: 0.33, 0.79) 1-sided p = 0.0011	
TTDQoL	3.8	3.8	5.7	3.9	7.8	3.9
	Difference: +0.0 months HR: 0.86 (95% CI: 0.65, 1.14) 1-sided p = 0.1513		Difference: +1.8 months HR: 0.65 (95% CI: 0.46, 0.91) 1-sided p = 0.0063		Difference: +3.9 months HR: 0.48 (95% CI: 0.30, 0.76) 1-sided p = 0.0008	

Overall Survival: No statistically significant difference in OS across Overall ITT Population or either Subgroups.

*35 ng/mL baseline CEA cutoff represents the median for the Overall ITT population in the study

1. Mulcahy, M. et al, Radioembolization With Chemotherapy for Colorectal Liver Metastases: A Randomized, Open-Label, International, Multicenter, Phase III Trial (EPOCH). Journal of Clinical Oncology, 20 Sept 2021.
2. Harris, W. The EPOCH Trial: Identifying Key Patient Subgroups to Optimize Treatment Planning. Poster presented at: ASCO-GI; January 21, 2023; San Francisco, CA.

EPOCH TRIAL

Trial Objective & Design

Trial Conclusions

PRIMARY ENDPOINTS

Progression-Free Survival

Hepatic Progression-Free Survival

Subgroup Analyses for PFS & hPFS

SECONDARY ENDPOINTS

Overall Survival

Tumor Response

Time to Deterioration of Quality of Life

ADDITIONAL ANALYSES

Time to Progression

Hepatic Time to Progression

Time to Start of Subsequent Therapy

Post-Hoc Subgroup Analyses

KEY PATIENT & DISEASE CHARACTERISTICS

TREATMENT CHARACTERISTICS, DOSIMETRY & SAFETY



KEY PATIENT & DISEASE CHARACTERISTICS

Patient and disease characteristics were well-balanced between the treatment and control arms.

	TheraSphere + Chemo (N = 215)	Chemo (N = 213)
Median Age (y)	63.0	60.0
Male	135 (62.8%)	138 (64.8%)
Region		
North America	63 (29.3%)	56 (26.3%)
Europe	131 (60.9%)	145 (68.1%)
Asia	21 (9.8%)	12 (5.6%)
ECOG 0	119 (55.3%)	133 (62.4%)
Albumin ≥ Site LLN ¹	182 (84.7%)	177 (83.1%)
CEA ² ≥ 35 ng/mL	116 (54.0%)	105 (49.3%)
KRAS Status		
Mutant	100 (46.5%)	101 (47.4%)
Wild Type	115 (53.5%)	112 (52.6%)
Bilobar disease	176 (81.9%)	173 (81.2%)
Liver tumor burden at baseline by BICR		
< 10%	124 (57.7%)	121 (56.8%)
≥ 10% to < 25%	54 (25.1%)	47 (22.1%)
≥ 25%	29 (13.5%)	28 (13.1%)
Missing	8 (3.7%)	17 (8.0%)
Maximum liver lesion size ≥ 4 cm	162 (75.3%)	142 (66.7%)
Primary tumor in situ	83 (38.6%)	69 (32.4%)
Left side primary tumor location	150 (69.8%)	136 (63.8%)
Extrahepatic metastases at baseline	113 (52.6%)	95 (44.6%)

1. LLN = lower limit of normal.
2. CEA = carcinoembryonic antigen.

EPOCH TRIAL

Trial Objective & Design

Trial Conclusions

PRIMARY ENDPOINTS

Progression-Free Survival

Hepatic Progression-Free Survival

Subgroup Analyses for PFS & hPFS

SECONDARY ENDPOINTS

Overall Survival

Tumor Response

Time to Deterioration of Quality of Life

ADDITIONAL ANALYSES

Time to Progression

Hepatic Time to Progression

Time to Start of Subsequent Therapy

Post-Hoc Subgroup Analyses

KEY PATIENT & DISEASE CHARACTERISTICS

TREATMENT CHARACTERISTICS, DOSIMETRY & SAFETY






TREATMENT CHARACTERISTICS

Treatment characteristics were well-balanced between the treatment and control arms.

	TheraSphere + Chemo (N = 215)	Chemo (N = 213)
Received Assigned Therapy	187 (87.0%)	191 (89.7%)
2nd Line Chemo Administered	203 (94.4%)	191 (89.7%)
Irinotecan-based / Mean Number of Cycles	130 (60.5%) / 9.0	123 (57.7%) / 9.5
Oxaliplatin-based / Mean Number of Cycles	73 (34.0%) / 8.5	68 (31.9%) / 8.8
Biological Agent	88 (40.9%)	93 (43.7%)
TheraSphere Y-90 Glass Microspheres Treatment		
Median time to TheraSphere Y-90 treatment, days (range)	25 (12, 90)	NA

DOSIMETRY APPROACH

In the 185 patients treated with TheraSphere Y-90 prior to progression¹:



Treatment Median: **Day 4** (1st week Thursday)

Median Specific Activity: **1,400 Bq** (single sphere)

Median dose absorbed by perfused volume: **117.0 Gy** (range: 61.7-156)

SAFETY

The addition of TheraSphere Y-90 to second-line chemotherapy did not increase chemo-related adverse events and no new safety signals were identified.²

Post-Hoc Analyses:

Liver-related treatment emergent adverse events (TEAEs) occurred more frequently in patients with <10% liver volume replaced by tumor and/or in patients with >10 lesions, likely due to increased proportion of irradiated normal liver tissue. Sequential lobar treatment, as opposed to same day whole liver treatment (as required by EPOCH protocol), may mitigate liver-related TEAEs.³

EPOCH TRIAL

Trial Objective & Design

Trial Conclusions

PRIMARY ENDPOINTS

Progression-Free Survival

Hepatic Progression-Free Survival

Subgroup Analyses for PFS & hPFS

SECONDARY ENDPOINTS

Overall Survival

Tumor Response

Time to Deterioration of Quality of Life

ADDITIONAL ANALYSES

Time to Progression

Hepatic Time to Progression

Time to Start of Subsequent Therapy

Post-Hoc Subgroup Analyses

KEY PATIENT & DISEASE CHARACTERISTICS

TREATMENT CHARACTERISTICS, DOSIMETRY & SAFETY

1. As assessed by investigator. 2. Mulcahy, M. et al, Radioembolization With Chemotherapy for Colorectal Liver Metastases: A Randomized, Open-Label, International, Multicenter, Phase III Trial (EPOCH). Journal of Clinical Oncology, 20 Sept 2021. 3. Salem, R. Optimizing patient selection for treating colorectal liver metastases with glass radioembolization plus chemotherapy: The EPOCH study. Lecture presented at: Society of Interventional Oncology; January 22, 2023; Washington, DC.



TheraSphere™ Y-90 Glass Microspheres | EPOCH Trial

EPOCH TRIAL

Trial Objective & Design

Trial Conclusions

PRIMARY ENDPOINTS

Progression-Free Survival

Hepatic Progression-Free Survival

Subgroup Analyses for PFS & hPFS

SECONDARY ENDPOINTS

Overall Survival

Tumor Response

Time to Deterioration of Quality of Life

ADDITIONAL ANALYSES

Time to Progression

Hepatic Time to Progression

Time to Start of Subsequent Therapy

Post-Hoc Subgroup Analyses

KEY PATIENT & DISEASE CHARACTERISTICS

TREATMENT CHARACTERISTICS, DOSIMETRY & SAFETY

TheraSphere™ Yttrium-90 Glass Microspheres

INDICATION FOR USE: TheraSphere is indicated for use as selective internal radiation therapy (SIRT) for local tumor control of solitary tumors (1-8 cm in diameter), in patients with unresectable hepatocellular carcinoma (HCC), Child-Pugh Score A cirrhosis, well-compensated liver function, no macrovascular invasion, and good performance status. **CONTRAINDICATIONS:** TheraSphere is contraindicated in patients: whose Tc-99m macroaggregated albumin (MAA) hepatic arterial perfusion scintigraphy shows any deposition to the gastrointestinal tract that may not be corrected by angiographic techniques • who show shunting of blood to the lungs that could result in delivery of greater than 16.5 mCi (0.61 GBq) of Y-90 to the lungs. Radiation pneumonitis has been seen rarely in patients receiving doses to the lungs greater than 30 Gy in a single treatment. • in whom hepatic artery catheterization is contraindicated, such as patients with vascular abnormalities or bleeding diathesis • who have pulmonary insufficiency (conventionally defined by an arterial oxygen pressure (Pa,O₂) of < 60 mmHg, or oxygen saturation (Sa,O₂) of < 90%) or severe liver dysfunction, including hepatic encephalopathy, clinically evident ascites or treatment with diuretics for ascites • with portal vein thrombosis (PVT) Type 4 involvement and lack of Tc-99m MAA deposition on the PVT seen on the Tc-99m MAA imaging with >70% tumor replacement in the liver • with comorbidities or poor overall health (e.g., ECOG performance status rating > 2) which may make the patient a poor candidate for locoregional radiation treatment. • who are pregnant. **WARNINGS:** The following pre-treatment, high-risk factors (disease characteristics) have been associated with serious adverse events deemed possibly related to use of the device: infiltrative tumor type • tumor nodules too numerous to count • AST or ALT > 5 times ULN • bilirubin > 2 mg/dL • tumor volume > 50% combined with albumin < 3 g/dL. Keep the TheraSphere dose vial upright and stored in its lead pot before and during patient treatment, except as required for radiation measurement. Do not open the dose vial acrylic shield prior to patient treatment. Post-treatment, waste materials require caution to prevent contamination and beta shielding due to residual glass microspheres. **PRECAUTIONS: GENERAL PRECAUTIONS:** As in any intra-arterial procedure, aseptic technique should be practiced, and care should be taken to ensure minimum patient anesthesia exposure extraneous to therapeutic objective. • Consideration of patient comorbidities should be used when determining the type and volume of fluid to infuse via catheter to avoid electrolyte imbalance, fluid shift, and hyperglycemia. • It is important to avoid any aggressive arterial procedure that may lead to arterial spasm that impairs TheraSphere distribution into the perfused liver target volume which may lead to underdosing or non-target deposition of TheraSphere. **PRECAUTION IN PATIENTS WITH IMPAIRED LIVER FUNCTION:** No efficacy or safety data from the LEGACY study are available to support the use of the device in patients with Child-Pugh score B or C cirrhosis. **PRECAUTION IN VULNERABLE PATIENTS:** No effectiveness or safety data are available to support the use of the device in children or breast-feeding women. **ENDOCRINE DISRUPTION, CARCINOGENICITY, MUTAGENICITY, TOXICITY TO REPRODUCTION:** Ideally the use of this radioactive device in women of childbearing capability should be performed during the first few (approximately 10) days following the onset of menses. **RADIATION SAFETY:** Radioactive products should be used only by healthcare professionals who are qualified by specific training in the safe use and handling of radionuclides and whose experience and training have been approved by the appropriate government agency authorized to license the use of radionuclides. • As in the use of any radioactive material, ensure minimum radiation exposure to the patient extraneous to the therapeutic objective, and to minimize radiation exposure to workers and others in contact with the patient. **RELEASE AND POST-TREATMENT PRECAUTIONS:** Post treatment patient care: use universal precautions for body fluid contact. Trace Y-90 may be detectable in blood and urine; handle with gloves and dispose as normal body fluids. The radiation field is expected to be less than 1 mrem/h (10 μSv/h) at 3 ft (1 m) from the patient's abdomen. Supplemental shielding and segregation of the patient are not required to maintain exposure to others below regulated limits. • Release instructions: The patient should follow good hygiene (e.g., proper hand washing). Caregivers, family, and others do not require restrictions on patient contact; however, they can minimize their radiation exposure by avoiding prolonged time (>12 hours per day) within 1 ft (0.3 m) of the patient's abdomen for the first week post therapy. Patients should be advised that radiation emitted from the patient may be detectable at security screening (e.g., international travel). • Special precautions post-administration: If the patient requires hospitalization, surgery, medical assessment or treatment regarding any part of their thorax or abdomen within first 2 weeks of treatment, the patient should advise the hospital and treating physician of the Y-90 TheraSphere implant. The physician should consult their radiation safety staff for handling and disposal of liver tissue. • Special liver tissue handling: Special liver tissue handling may be required for post-treatment surgery, explant, or transplant since the glass microspheres remain permanently implanted in the liver tissue. Disclosure of the treatment will be required if cremation is considered. **POTENTIAL ADVERSE EVENTS:** The use of this product leads to irradiation of both tumorous and normal liver tissue. As a result, patients with compromised liver function may be at greater risk of liver function impairment and hence could experience complications. Clinical side effects usually occur within the first 4 to 6 weeks after treatment. Based on clinical trial data, literature reviews and post market surveillance, adverse events potentially associated with treatment using Y-90 microspheres, including TheraSphere, may include the following: Allergic reaction • Altered liver function, acute or chronic • Anorexia • Anxiety • Ascites • Bile Duct injury • Bleeding/hemorrhage • Chills / rigors • Cholecystitis (inflammatory or infectious) • Colitis • Death • Dehydration • Diarrhea • Dizziness • Dyspnea • Edema (any location) • Electrolyte abnormalities • Elevated BUN/creatinine • Fall • Fatigue • Fever • Gastrointestinal bleeding / hemorrhage • Gastrointestinal ulcer or ulceration • Hepatic encephalopathy • Hepatorenal failure • Hiccups • Hypertension • Hypotension • Infection (any location) • Liver failure, acute or chronic • Lymphopenia • Malaise • Mood alteration • Muscle weakness • Nausea • Neutropenia • Pain (any location) • Pancreatitis • Platelet count abnormalities • Pleural effusion • Portal hypertension • Pre-existing chronic liver disease decompensation • Pulmonary edema • Pulmonary fibrosis • Radiation hepatitis • Radiation induced disease, acute • Radio Embolization Induced Liver Disease (REILD) • Sepsis • Supraventricular arrhythmia • Thrombosis (arterial or venous) • Tumor inflammation (including tumor edema) • Tumor-lysis syndrome • Vomiting • Weight loss. Complications related to the administration procedure itself may include: Allergic reaction: Arterial injury including vessel dissection • Aspiration pneumonia • Bruising/bleeding/hematoma at site • Constipation/abdominal distension • Fatigue • Flushing • Infection • Nausea • Nerve damage. **CAUTION:** Federal (USA) law restricts this device to sale by or on order of a physician. PI-992004-AA. **Note:** Dose to the liver does not exceed 150 Gy.

TheraSphere is a registered trademark of Theragenics Corporation used under license by Boston Scientific Medical Device Limited, a wholly owned indirect subsidiary of Boston Scientific Corporation. **CAUTION:** TheraSphere is under an investigational device exemption for treatment of patients with metastatic colorectal cancer. The safety and effectiveness for this treatment has not been established. For full adverse event rates and complete data set, reference Journal of Clinical Oncology manuscript.



Peripheral Interventions
300 Boston Scientific Way
Marlborough, MA 01752-1234
www.bostonscientific.com

To order product or for more information contact customer service at 1.888.272.1001.

© 2024 Boston Scientific Corporation or its affiliates. All rights reserved.

PI-1226604-AB

