- For single patient use only. Do not reuse, reprocess, or re-sterilize as these may compromise the structural integrity of the device and/or lead to device failure, and may result in patient injury, illness, or death. Reuse, reprocessing, or re-sterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness, or death of the patient.
- After use, dispose of product and packaging in accordance with hospital, administrative, and/or local government policy.

PRODUCT DESCRIPTION

The Claret Medical Sentinel Cerebral Protection System is a percutaneously delivered embolic protection device, designed to capture and remove debris dislodged during endovascular procedures. The Sentinel System utilizes an embolic filter delivered to the brachiocephalic artery (Proximal Filter), and a second embolic filter delivered to the left common carotid artery (Distal Filter). At the completion of the procedure, the filters and debris are recaptured into the catheter and removed from the patient.

The Sentinel System consists of a 6 French catheter with deployable Proximal and Distal Filters, an Articulating Sheath and an integral handle assembly. The Articulating Sheath tip, Proximal Sheath tip, Proximal Filter hoop, Distal Filter hoop and Distal Filter tip are radiopaque to enable visualization during use. See Figures 1 and 2.

Table 1: Filter-Vessel Sizing Guide

REF (Model) Number for Ordering	Proximal Filter Size (mm)	Target Proximal Vessel Size (mm)	Distal Filter Size (mm)	Target Distal Vessel Size (mm)	Articulating Sheath Curve Shape
CMS15-7A	15	9.0 – 15.0	7	5.0 – 7.0	A
CMS15-7B					B
CMS15-7C					C
CMS15-10A					A
CMS15-10B	10			6.5 – 10.0	B
CMS15-10C					C

Table 2: Sentinel System Specifications

Delivery Profile	6F
Working Length	95 cm
Articulating Sheath Length	4 cm
Guidewire Compatibility	0.014" (0.36 mm) diameter floppy tip coronary guidewire, 175 cm minimum length

Package contains one (1) Sentinel System*

*Reference package label for specific model. Review patient anatomy by CT or angiogram prior to selecting appropriate curve shape per Figure 3.

INDICATIONS FOR USE

The Sentinel System is indicated for use as an embolic protection device to capture and remove embolic material (thrombus/debris) that may enter the cerebral vascular system during endovascular procedures. The diameters of the arteries at the sites of filter placement should be measured and the filters sized to the Proximal and Distal Target Vessels as listed in Table 1.

CONTRAINDICATIONS FOR USE

- Do not use with patients in whom anticoagulant and antiplatelet therapy is contraindicated.
- Do not use in vessels with excessive tortuosity.
- Do not use in patients with uncorrected bleeding disorders.
- Do not use this product on patients with compromised blood flow to the right upper extremity.
- Do not use in a patient who has arterial stenosis >70% in either the left common carotid artery or the brachiocephalic artery.
- Do not use in a patient whose brachiocephalic or left carotid artery reveals significant stenosis, ectasia, dissection, or aneurysm at the aortic ostium or within 3 cm of the aortic ostium.

WARNINGS

The Sentinel System should only be used by physicians who have received appropriate training and are familiar with the principles, clinical applications, complications, side effects, and hazards commonly associated with endovascular procedures.

- Carefully read all instructions and labeling prior to use. Observe all warnings, cautions, and precautions noted throughout these instructions. Failure to do so may result in complications.
- Refer to the instructions for use supplied with any interventional device to be used in conjunction with the Sentinel System for their intended uses, sizing, warnings, and precautions.
- The appropriate antiplatelet/anticoagulation therapy should be administered pre- and post-procedure in accordance with standard medical practice.
- Prior to use, the packaging and product should be inspected for signs of damage. Never use damaged product or product from a damaged package.
- Never advance or withdraw the Sentinel System without proper fluoroscopic guidance or against resistance until the cause is determined. Advancing with such resistance may lead to embolization of debris, and vessel and/or device damage.
- It is recommended that the patency of the right radial or brachial artery be assessed prior to the introduction of the Sentinel System.
- It is recommended that the patient be administered an Allen Test should the radial artery be used for device introduction.
- Do not use the device in left radial or left brachial access.
- Do not use the Sentinel System to deliver any type of fluid to the patient e.g. contrast media, heparinized saline, etc. due to risk of air embolization or compromise to device performance.
- Minimize movement of the Sentinel System after initial placement. Excessive movement of filters may lead to embolization of debris, vessel and/or device damage.
- Do not deploy the filters within a previously repaired artery.
- Observe the Sentinel System under fluoroscopy and monitor the patient to verify the filters have not become occluded with debris resulting in slow or no flow. The filters should be recovered if they become occluded or if flow is compromised (See Procedural Use – Retrieval).
- Indwell time of the Sentinel System is not to exceed 90 minutes. Occlusion could occur, resulting in slow or no flow.
- Failure to adequately close off the Flush Ports (Front Handle, Rear Handle) may result in air embolism.
- Do not undersize or oversize the filters in relation to the selected vessel diameter. This may result in inadequate vessel wall apposition or incomplete deployment of the filters. (Refer to Sizing Guide, Table 1).
- Do not apply excessive force to the Sentinel System. This may lead to distal embolization of debris, and vessel and or device damage.
- Contents supplied STERILE using a radiation process. Do not use if sterile barrier is damaged. If damage is found, call your Claret Medical representative.

PRECAUTIONS

- Do not forcefully bend or reshape the Articulating Sheath of the Sentinel System. This may cause device damage.
- Do not use the product if the packaging sterile barrier has been damaged or compromised.
- Improper bending of the Sentinel System may damage the catheter.
- Do not re-sterilize or reuse on another vessel or patient.

POTENTIAL COMPLICATIONS

Access Site Complication/Injury, Angina, Bleeding, Death, Dissection, Embolism, Emergent Surgery, Hematoma, Ischemia, Infection, Myocardial Infarction, Renal Insufficiency, Stroke, and Vessel Injury.

INSTRUCTIONS FOR USE

Preparing the Sentinel System for Use

- Administer anticoagulation medications and monitor activated clotting time per standard institutional guidelines. Anticoagulant therapy sufficient to maintain an Activated Clotting Time of at least 250 seconds for the duration of the procedure is recommended.
- Perform angiogram of the aortic arch.
- Identify the location within the vessels where the filters will be deployed. Measure the vessel diameter at these locations in millimeters and select the appropriate Sentinel System Proximal and Distal Filter diameter.

WARNING: Do not undersize or oversize the filters in relation to the selected vessel diameter. This may result in inadequate vessel wall apposition, incomplete deployment of the filters, or vessel damage.

- Ensure the introducer sheath size will accommodate the Sentinel System.
- Using sterile techniques remove the Sentinel System from the packaging and place in a sterile work area.

CAUTION: Do not use the product if the packaging sterile barriers have been damaged or compromised.

WARNING: Inspect the device for any damage. Never use damaged product or product from a damaged package.

Flushing the Sentinel System

CAUTION: Do not prepare the Sentinel System or sheath the Proximal and Distal Filters until immediately prior to use.

Note: The device handle has two handle locks, the Rear Handle Lock, and the Front Handle Lock. Refer to Figure 1. Closing these locks facilitates flushing, prevents back-bleeding, and prevents motion of the device handle components and Distal Filter. The locks should be temporarily opened to facilitate movement of the handle components as required.

Note: The primary controls used to deploy the device, the Proximal Filter Slider (#1), the Articulation Knob (#2), and the Distal Filter Slider (#3), are all marked with the number "1", "2" and "3" indicating the order in which they are used. In this document, these names will be shown with the control number appended to the name.

- Remove the packaging stylet from the distal guidewire lumen and discard.
- Ensure that both the Front Handle Lock and the Rear Handle Lock are tightened.
- Flush through the Flush Port in the Distal Filter Slider (#3) with heparinized saline until all air is removed and fluid passes from Distal Filter Tip guidewire lumen. See Figure 1.
- Flush through the Rear Handle Flush Port with heparinized saline until all air is removed and fluid passes from the tip of the Articulating Sheath. See Figure 1. Ensure that the flush port stopcock is closed following flushing.
- Flush through the Front Handle Flush Port with heparinized saline until all air is removed and fluid passes from the tip of the Proximal Sheath. See Figure 1. Ensure that the flush port stopcock is closed following flushing.
- Submerge the distal end of the device in heparinized saline, and loosen the Rear Handle Lock. With the distal tip submerged, slowly retract the Distal Filter by pulling back on the Distal Filter Slider (#3) until the filter is fully collapsed into the Articulating Sheath. The submerged filter may be agitated during sheathing in order to facilitate removal of bubbles. Tighten the Rear Handle Lock.

Note: Flushing and sheathing of the Distal Filter may be repeated to ensure all air has been removed from the system.

CAUTION: Do not over-retract the Distal Filter as damage may occur.

- Ensure the Articulating Sheath is fully advanced until the Articulation Knob (#2) is in contact with the Front Handle Lock to ensure it does not interfere with sheathing the Proximal Filter. Tighten the Front Handle Lock. While submerged, sheath the Proximal Filter by slowly advancing the Proximal Filter Slider (#1) relative to the Front Handle until the Proximal Filter is fully sheathed. The submerged filter may be agitated during sheathing in order to facilitate removal of bubbles. See Figures 1 and 9.
- Note:** Flushing and sheathing of the Proximal Filter may be repeated to ensure all air has been removed from the system.
- While submerged, again flush through the Front Handle Flush Port with heparinized saline until all air is removed and fluid passes from the tip of the Proximal Sheath. See Figure 1. Ensure that the flush port stopcock is closed following flushing.

WARNING: Do not use a Sentinel System that has not been properly flushed. Failure to prepare and flush the device before use may introduce air and patient injury may result.

Note: Tighten both the Rear and Front Handle Locks prior to delivering device to prevent inadvertent movement.

Note: Use a minimum of 10 cc of heparinized saline to flush through the Front Handle Flush Port to ensure all air has been removed from the system.

Note: Refer to the instructions for use supplied with any interventional device to be used in conjunction with the Sentinel System for their intended uses, sizing, warnings, and precautions.

Procedural Use - Delivery and Deployment

WARNING: Do not use a Sentinel System that has not been properly flushed. Failure to prep and flush the device before use may introduce air and patient injury may result.

WARNING: Never advance, manipulate, or withdraw the Sentinel System without proper fluoroscopic guidance to prevent damage to the System and/or harm to the patient.

WARNING: The Sentinel System is not to be used to deliver any type of fluid to the patient e.g. contrast media, heparinized saline, etc.

- Using standard interventional technique, place a 6 French introducer sheath into the radial or brachial artery of the patient's right arm.
- Backload a floppy tip 0.014" coronary guidewire into the Distal Filter Tip located at the distal end of the Sentinel System until the guidewire tip is located just inside the distal tip of the Sentinel catheter.
Note: A stiff guidewire may impact the shape of the Articulating Sheath curve and make cannulation of the left common carotid difficult.
Note: Guidewires with intermediate coils may become deformed while cannulating the left common carotid resulting in compromised guidewire movement.
- Introduce the Sentinel System into the introducer sheath.
- In the patient's right arm, advance the guidewire relative to the Sentinel System until the distal tip of the guidewire is a minimum of 10 cm beyond the distal tip of the Sentinel System using fluoroscopic guidance.
- Advance the Sentinel System distally until it contacts the introducer sheath hemostasis valve. Gently advance the Sentinel System until it is fully inserted into the introducer hemostasis valve.
- Advance the Sentinel System and the guidewire together using standard interventional technique until the Proximal Filter is in the intended target location in the brachiocephalic artery with the Articulating Sheath section of the catheter extending down the ascending aorta. Should the catheter tip extend down the descending aorta, pull the system back and rotate to advance down the ascending aorta.

WARNING: Do not advance the Sentinel System without a guidewire extending distally past the tip of the catheter a minimum of 10 cm.

WARNING: Do not use excessive force on the Sentinel System while introducing or advancing through the introducer sheath or blood vessels. Excessive force may cause damage to the device and/or patient harm.

- Note:** The Articulating Sheath will protrude into the aorta during proximal filter deployment.
- Deploy the Proximal Filter by holding the Front Handle in a fixed position and slowly retracting the Proximal Filter Slider (#1) fully.
 - Confirm proper Proximal Filter position using fluoroscopy. The Proximal Filter should be positioned in the brachiocephalic artery to prevent any debris from reaching the right carotid artery. See Figures 4 and 5.
 - If the filter position is not optimal, the filter may be retrieved and repositioned up to two times. This may be done by holding the Front Handle in a stationary position and advancing the Proximal Filter Slider (#1) until the Proximal Filter is re-sheathed. The Proximal Filter may then be repositioned by advancing or retracting the catheter until optimal positioning is achieved. Finally the Proximal Filter is redeployed by retracting the Proximal Filter Slider (#1) while holding the Front Handle in a fixed position.

CAUTION: Repositioning, if required, should only occur during initial placement.

- Confirm filter-to-vessel wall apposition using fluoroscopy, and ensure that the Proximal Filter and Proximal Sheath do not move after placement.
- Withdraw the guidewire until the tip is located just within the distal tip of Sentinel catheter.
- Loosen the Front Handle Lock to facilitate positioning of the Articulating Sheath.
- Position the Articulating Sheath by manipulating the Rear Handle relative to the Front Handle in order to position the catheter tip. Rotate the Articulation Knob (#2) on the Rear Handle in the direction of the arrows in order to deflect the tip of the Articulating Sheath as necessary toward the left common carotid artery ostium.

CAUTION: Do not move the Front Handle, and thus the Proximal Filter, while manipulating the Rear Handle.

- Advance the 0.014" guidewire beyond the distal tip of the Articulating Sheath in order to place the guidewire in the left common carotid artery.

CAUTION: Do not to advance the guidewire more than 5 cm into the left common carotid artery.

- Position the Articulating Sheath so that the curvature matches the Brachiocephalic Artery – Aorta – Left Common Carotid Artery junction and is pulled up to the carina between the two vessels, see Figure 6.
Note: Ensure that the Articulating Sheath is well apposed to the carina, and does not protrude into the aortic space. See Figure 6 for correct positioning and Figure 8 for incorrect positioning.
- Secure the position of the Articulating Sheath by tightening the Front Handle Lock.
- Loosen the Rear Handle Lock and advance the Distal Filter under fluoroscopy by pushing the Distal Filter Slider (#3) forward until the Distal Filter frame is fully expanded and apposed to the vessel wall. The Distal Filter should be positioned just beyond the Articulating Sheath tip and movement should be minimized once it is fully expanded in the vessel. See Figure 2.

WARNING: Minimize movement of the Sentinel System after filter deployment. Excessive movement may lead to embolization of debris, and vessel and/or device damage.

- Confirm filter-to-vessel wall apposition of the distal filter using fluoroscopy. See Figure 7.
- Tighten the Rear Handle Lock. See Figure 1.

CAUTION: Verify that the Front Handle Lock and the Rear Handle Lock are tight and secure before any subsequent procedures.

WARNING: The Sentinel System is not to be used to deliver any type of fluid to the patient e.g. contrast media, heparinized saline, etc.

CAUTION: Repositioning, if required, should only occur during initial placement.

- Cover the exposed portion of the Sentinel System with a drape to prevent movement during subsequent endovascular procedures.

CAUTION: Care must be taken NOT to kink the exposed catheter.

WARNING: Minimize movement of the Sentinel System and its filters after filter deployment. Excessive movement may lead to embolization of debris, and vessel and/or device damage.

WARNING: If gross movement of either the Proximal or Distal Filter is noted, check to ensure filters remain apposed to the vessel walls by fluoroscopy.

WARNING: If the arterial flow is believed to be compromised (slow / no flow), the filters should be retrieved. See Retrieval below.

Procedural Use – Retrieval

WARNING: Do not pull excessively on the Sentinel System to avoid filter membrane tears, filter hoop detachment, system damage or patient harm during use.

WARNING: Never advance or withdraw the Sentinel System without proper fluoroscopic guidance.

WARNING: Never withdraw or move an intravascular device against any resistance until the resistance cause is determined. Advancing or retracting with resistance may lead to embolization of debris, and vessel and / or device damage.

There are two methods for Distal Filter recovery: Partial and Full Enclosure Recovery

- Loosen the Rear Handle Lock. See Figure 1.
- Recover the Distal Filter using one of the following two methods:
 - Full Enclosure Recovery:** Gently withdraw the Distal Filter Slider (#3) relative to the Rear Handle until the radiopaque Distal Filter Tip is flush with the Radiopaque Articulating Sheath Tip Marker as visualized on fluoroscopy. Tighten the Rear Handle Lock. If resistance is felt during Distal Filter recovery, or if it is believed that the Distal Filter is excessively full, follow the Partial Enclosure Recovery method detailed below.
 - Partial Enclosure Recovery:** Gently withdraw the Distal Filter Slider (#3) relative to the Rear Handle until the Distal Filter Radiopaque Hoop is collapsed inside the Articulating Sheath tip as visualized on fluoroscopy. Tighten the Rear Handle Lock.

WARNING: Exercise caution when using the partial enclosure recovery method. If resistance is felt during catheter withdrawal, advance the Distal Filter and Articulating Sheath together distally and withdraw the Distal Filter more fully into the Articulating Sheath before re-attempting withdrawal of the catheter.

- Loosen the Front Handle Lock and withdraw the Articulating Sheath tip from the left common carotid artery by manipulating, straightening, rotating, and advancing or withdrawing the Rear Handle and rotating the Articulation Knob (#2) until the Articulating Sheath tip is straight and is within the aorta.
- Advance the Articulating Sheath completely by advancing the Rear Handle until the Articulation Knob (#2) contacts the Front Handle Lock to prevent interference with the Proximal Sheath or Proximal Filter during Proximal Filter retrieval. Tighten the Front Handle Lock. See Figure 9.
- Re-sheath the Proximal Filter by holding the Front Handle in a stationary position and slowly advancing the Proximal Filter Slider (#1) until the Proximal Sheath Radiopaque Marker meets the Articulating Sheath as visualized on fluoroscopy. See Figure 9. Minimize retracting or advancing the Front Handle during this step. Vessel damage may occur or debris may be lost should the Proximal Filter be moved when in the deployed state.
- Advance the guidewire prior to withdrawal of the Sentinel System. Withdraw the catheter system while using fluoroscopy.
Note: If there is any resistance to removing the Sentinel System from the introducer, remove the introducer and the Sentinel System together.

CAUTION: Do not re-sterilize or reuse this device.

Storage:

- Do not store in direct sunlight
- Keep dry

Note: After use, dispose of product and packaging in accordance with hospital, administrative, and/or local government policy.

WARRANTY DISCLAIMER

Although the product has been manufactured under careful controlled conditions, Claret Medical, Inc. has no control over the conditions under which the product is used. Claret Medical, Inc. therefore disclaims all warranties, both expressed and implied, with respect to product including, but not limited to, any implied warranty of merchantability or fitness for a particular purpose. Claret Medical, Inc. shall not be liable to any person or entity for any medical expenses or any direct, incidental or consequential damages caused by any use, defect, failure or malfunction of the product, whether a claim for such damages is based upon warranty, contract, tort or otherwise. No person has any authority to bind Claret Medical, Inc. to any representation or warranty with respect to the product. The exclusions and limitations set out above, are not intended to, and should not be construed so as to contravene mandatory provisions of applicable law. If any part of this Disclaimer of Warranty is held to be illegal, unenforceable or in conflict with applicable law by a court of competent jurisdiction, the validity of the remaining portions of this Disclaimer of Warranty shall not be affected, and the rights and obligations shall be construed and enforced as if this Disclaimer of Warranty did not contain the particular parts or term held to be invalid.

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