

# CoreDx™

## Pulmonary Mini-Forceps

For intranodal histologic sample collection—a complementary tool for EBUS procedures



# About CoreDx

Introducing CoreDx, a mini-forceps enabling pulmonologists to collect intranodal histologic samples. Quantity with quality, potentially alleviating the need for patients to undergo more invasive biopsy procedures.<sup>1</sup>

## Histopathological Core Samples

For molecular analysis and diagnosis of sarcoidosis and lymphoma

## Flexibility & Column Strength

CoreDx is designed with flexibility to allow access to all regions within the transbronchial tree


High column strength allows for ease of entry through access hole when used in conjunction with EBUS-TBNA

## Scope Compatibility

CoreDx is compatible with all currently available EBUS scopes (Olympus, Pentax, and Fuji)

Compatible with super slim bronchoscope (1.2mm working channel)


A pushable yet flexible shaft allows scope articulation while maintaining column strength




# Why Histology?

## Quality Histology

Needle histology can be bloody and dilute, but CoreDx maintains the tissue structure necessary for accurate diagnosis.<sup>2</sup>



QUALITY  
HISTOLOGY  
SAMPLE



POOR  
HISTOLOGY  
SAMPLE

## When does histology matter?

### Lymphoma

Lymphoma can be a complex disease, with over 100 different varieties many with their own test requiring large amounts of quality tissue.<sup>2</sup>

### Sarcoidosis


Diagnosing Sarcoidosis can be challenging, histologic samples give pathologists a view of structure, which helps ensure accurate an diagnosis.<sup>2</sup>

### PD-L1


Genetic sequencing protocols can require histologic samples. Histology enables the testing that leads to targeted therapies for lung cancer treatment.<sup>3</sup>

The increase in diagnostic yield for mini-forceps is most pronounced in patients with Sarcoidosis 88% vs 36% for TBNA and Lymphoma 81% vs 35%.<sup>1</sup>

# CoreDx™ Procedure Overview




Following a standard EBUS procedure, CoreDx is used as an adjunctive technique to obtain histologic samples.





● ULTRASOUND REFERENCE POINTS



When using in conjunction with EBUS, take note of the needle alignment on the ultrasound image and choose reference points to help maintain the scope position during instrument exchange.



With CoreDx at the end of the working channel, align the ultrasound image to your reference points. Then using short deliberate strokes advance the forceps to the location of planned biopsy.



Open jaws carefully under ultrasound visualization. Move forceps forward against the tissue to be sampled and close jaws. Slowly withdraw forceps keeping jaws closed.

# Ordering Information

## CoreDx™ Pulmonary Mini-Forceps

| Endobronchial | GTIN           | Description           | Minimum Working Channel | Length | Packaging |
|---------------|----------------|-----------------------|-------------------------|--------|-----------|
| M00515220     | 08714729971054 | Pulmonary min-forceps | 1.2mm                   | 130cm  | Each      |

- 1 Endobronchial Ultrasound-Guided Miniforceps Biopsy in the Biopsy of Subcarinal Masses in Patients with Low Likelihood of Non-Small Cell Lung Cancer. Herth, Felix J.F. et al. The Annals of Thoracic Surgery, Volume 85, Issue 6, 1874 – 1878
- 2 Endobronchial-Ultrasound Guided Miniforceps Biopsy of Mediastinal and Hilar Lesions. Chrissian, Ara et al. The Annals of Thoracic Surgery, Volume 92, Issue 1, 284 – 288
- 3 Endobronchial Ultrasound: Can It Give Us What We Need For Next Generation Sequencing? M. Samant, A. C. Burks, A. Chen. *Am J Respir Crit Care Med* 2017;195:A1651

**CAUTION:** The law restricts these devices to sale by or on the order of a physician. Indications, contraindications, warnings and instructions for use can be found in the product labeling supplied with each device. Information for the use only in countries with applicable health authority registrations. Material not intended for use in France.

Rx Only.

All trademarks are the property of their respective owners.

All photographs taken by Boston Scientific.

**Boston  
Scientific**  
Advancing science for life™

Boston Scientific Corporation  
300 Boston Scientific Way  
Marlborough, MA 01750  
[www.BostonScientific.com/Pulmonary](http://www.BostonScientific.com/Pulmonary)  
[www.BronchSuite.com](http://www.BronchSuite.com)

© 2019 by Boston Scientific Corporation  
or its affiliates. All rights reserved.

ENDO-549610-AA