MANTIS™

Clip

en	Instructions for Use	2
----	----------------------	---

51691375-01

2023-10

TABLE OF CONTENTS

REUSE WARNING	3
DEVICE DESCRIPTION	3
Contents	3
Operating Principle	3
Materials	3
User Information	4
INTENDED USE/ INDICATIONS FOR USE	4
Clinical Benefit Statement	4
Summary of Safety and Clinical Performance	4
CONTRAINDICATIONS	4
WARNINGS	4
NOTES	5
PRECAUTIONS	5
ADVERSE EVENTS	5
HOW SUPPLIED	5
Device Details	5
Handling and Storage	5
OPERATIONAL INSTRUCTIONS	6
Figure 1. MANTIS Clip	6
Preparation	6
Device Insertion	6
Procedure	6
Figure 2. Open the MANTIS Clip Jaws	6
Figure 3. Rotating the Clip	7
Figure 4. MANTIS Clip Closed	7
Figure 5. MANTIS Clip Permanently Deployed	8
DEVICE REMOVAL	8
Disposal	8
Post-Procedure	8
MRI SAFETY INFORMATION	9
PATIENT COUNSELING INFORMATION	9
Implantable Device Patient Information	9
Implant Card Instructions	10
WARRANTY	10
SYMBOL DEFINITIONS	10

MANTIS™

Clip

Rx ONLY

Caution: Federal Law (USA) restricts this device to sale by or on the order of a physician.

REUSE WARNING

For single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

DEVICE DESCRIPTION

The MANTIS Clip consists of a radiopaque, single-use clip with an 11 mm clip opening, pre-loaded on a flexible, rotatable delivery system.

The MANTIS Clip is designed to be compatible with forward viewing endoscopes with working channels equal to or greater than 2.8 mm.

The radiopaque MANTIS Clip is engineered to enable opening and closing no more than five times prior to deployment, aiding in repositioning of the clip at the lesion site. Reopening, closing and rotation capability may be limited by clinical circumstances and patient anatomy, among other factors.

Contents

(1) MANTIS Clip, pre-loaded on a flexible, rotatable delivery system.

Operating Principle

MANTIS Clip is a single use rotatable mechanical clip with enhanced grasping capabilities which can be introduced in the GI tract through an endoscope to be deployed in a target site. The clip can then be positioned for deployment via manipulation of the scope and/or the device itself which may include rotation of the clip. Once the clip is positioned at the treatment site, the handle is actuated to deploy the clip and the delivery catheter is then withdrawn through the scope. The deployed clip remains at the treatment site and facilitates the treatment; the clip eventually passes naturally through the Gastrointestinal (GI) tract.

Materials

Materials and substances for which the patient can be exposed to, by the implantable portion of the medical device, are the following:

Implantable Material	% (wt) (Total assembly weight = 0.125g)
Stainless steel	75%
Cobalt-Chrome	22%
Styrene acrylonitrile copolymers	3%

Contains cobalt: CAS No. 7440-48-4; EN No. 231-158-0. Defined as a 1B carcinogen and reproductive toxicant according to the European Commission in a concentration above 0.1% weight by weight.

Note: Current scientific evidence supports that metal alloys containing cobalt used in medical devices do not cause an increased risk of cancer or adverse reproductive effects.

User Information

The MANTIS Clip should only be used by or under the supervision of a physician with a thorough understanding of the technical principles, clinical applications, and risks associated with MANTIS Clip.

INTENDED USE/ INDICATIONS FOR USE

The MANTIS Clip is indicated for clip placement within the Gastrointestinal (GI) tract for the purpose of:

1. Endoscopic marking
2. Hemostasis for: Mucosal/sub-mucosal defects < 3 cm, Bleeding ulcers, Arteries < 2 mm, Polyps < 1.5 cm in diameter, Diverticula in the colon, Prophylactic clipping to reduce the risk of delayed bleeding post lesion resection
3. Anchoring to affix jejunal feeding tubes to the wall of the small bowel
4. As a supplementary method, closure of GI tract luminal perforations < 20 mm that can be treated conservatively

Clinical Benefit Statement

The clinical benefit of the MANTIS Clip is to enable endoscopic marking, provide acute and prophylactic hemostasis, anchor jejunal feeding tubes, and aid in closure of defects in the GI tract.

Summary of Safety and Clinical Performance

For customers in the European Union, use the device name found in the labeling to search for the device's Summary of Safety and Clinical Performance, which is available on the European database on medical devices (EUDAMED) website: (<https://ec.europa.eu/tools/eudamed>).

CONTRAINDICATIONS

- Do not use this device when hemostasis cannot be verified visually with an endoscopic field of view.
- Arteries greater than 2 mm.
- Polyps greater than 1.5 cm in diameter.
- Mucosal/Submucosal defects greater than 3 cm.

WARNINGS

- DO NOT FORCIBLY PULL BACK ON A CLIP THAT IS DEPLOYED AND HAS NOT DETACHED FROM THE COIL. THIS WILL TEAR THE TISSUE AND LIKELY RESULT IN SEVERE BLEEDING OR PERFORATION.
A wire-cutter should be available on the endoscopy cart and used to cut the coil near the device handle if needed. The endoscope can then be removed leaving the clip and coil intact. The patient may require URGENT SURGERY to separate the clip from the coil without tearing tissue or to manage bleeding resulting from manipulating the imbedded clip.
- If the clip detaches prematurely, there is a risk of re-bleeding, unsuccessful anchoring to affix jejunal feeding tubes, unsuccessful marking, or insufficient perforation closure.
- Although rates of occurrence are low, recurrent bleeding, ineffective clipping or endoscopic complications could result in the need for surgery.
- The use of clips in the presence of bacterial contamination may potentiate or prolong infection.
- Clipping hard or severely fibrotic lesions may result in ineffective hemostasis, unsuccessful anchoring to affix jejunal feeding tubes, unsuccessful marking, or insufficient perforation closure. Additional interventions may be required to control bleeding, successfully anchor to affix jejunal feeding tube, successfully mark, or sufficiently close the perforation.
- Contains nickel, which may cause an allergic reaction in individuals with nickel sensitivity.

- The MANTIS Clip is designed to be compatible with forward viewing endoscopes with working channels equal to or greater than 2.8 mm. Use with a side viewing scope may result in difficulty or inability to deploy, leading to patient injury.

NOTES

Limited studies indicate that:

- Lesions located in the esophagus and the lesser curvature of the stomach may be difficult to treat with a forward viewing endoscope.
 - The number of clips required for hemostasis, anchoring, marking, or closure may vary depending upon the anatomical site, histology, lesion type, and patient condition and history.
 - There is no clinical evidence to support the use of this device for clipping the neck of the diverticulum to treat bleeds.
 - There is no clinical evidence to support the use of this device for clipping GI tract luminal perforations >20 mm.
 - There is no clinical evidence to support the use of the device for esophageal varices.
-

PRECAUTIONS

- It is recommended that healthcare providers distribute patient implant cards with the name of the clip and date it was placed.
- Passage of the MANTIS Clip through a retroflexed or tortuous path, may result in the clip separating from the catheter and potentially kinking or damaging the device.
- Applying tangential pressure to an opened or closed clip may result in the clip separating from the catheter and potentially kinking or damaging the device.
- In a difficult scope position, it may be necessary to straighten the endoscope to facilitate the device passage, then reposition scope for treatment.
- If the catheter kinks or becomes damaged during insertion or passage, do not use it. Call Boston Scientific Customer Service and return the product.

ADVERSE EVENTS

- Allergic reaction
- Burn
- Hemorrhage
- Infection
- Inflammation
- Pain
- Perforation
- Surgery
- Tissue Damage

HOW SUPPLIED

The device is supplied sterile by ethylene oxide gas sterilization and is intended for single use only.

Device Details

Carefully examine the unit to verify that neither the contents nor the sterile package has been damaged.

Do not use if package is damaged or unintentionally opened before use.

Do not use if labeling is incomplete or illegible.


Immediately return damaged product to Boston Scientific Corporation.

Handling and Storage

Rotate inventory so that products are used prior to the expiration date on the package label.

This product has no special handling or storage requirements.

OPERATIONAL INSTRUCTIONS

[1] Thumb Ring [2] Slider [3] Handle [4] Rotation Control Knob
[5] Clip [6] Clip Jaws

Figure 1. MANTIS Clip

Preparation

1. Open the pouch and remove the device.
2. Inspect the device for kink or damage.

Note: If the device shows any sign of damage, do not use it, call Boston Scientific Customer Service and return the product. Do not attempt to repair nonfunctional or damaged devices.

Device Insertion

1. Carefully insert the MANTIS Clip through the biopsy channel of the endoscope with short, deliberate 2 to 3 centimeter strokes.


Caution: Do not advance an open clip through the endoscope working channel, otherwise endoscope working channel damage may result.

Note: Passage of the MANTIS Clip through a retroflexed or tortuous path may result in the clip separating from the catheter and potentially kinking or damaging the device. In a difficult scope position, it may be necessary to straighten the endoscope to facilitate the device passage, then reposition scope for treatment. If the device kinks or becomes damaged during insertion or passage, do not use it. Call Boston Scientific Customer Service and return the product. Do not attempt to repair nonfunctional or damaged devices.

Procedure

Note: Applying tangential pressure to an opened or closed clip may result in the clip separating from the catheter and potentially kinking or damaging the device. Prior to permanently deploying the MANTIS Clip, visually confirm that the device has not kinked, separated from the catheter, or become damaged in any way. If the device shows any sign of damage, DO NOT USE IT. Call Boston Scientific Customer Service and return the product.

1. When the MANTIS Clip is at the desired location, gently move the slider distally (away from the thumb ring) to open the MANTIS Clip jaws, as shown in Figure 2.


[1] Thumb Ring [2] Slider [3] Handle [4] Rotation Control Knob
[5] Clip [6] Clip Jaws

Figure 2. Open the MANTIS Clip Jaws

2. If additional positioning is desired, the clip may be rotated in either of the following two methods (see Figure 3):
 - a. The clip can be rotated by turning the rotation control knob in either direction. Typically, this is done by the nurse, technician, or assistant. Because the rotation control knob is separate from the handle, the user can keep his/her hand on the handle while rotating the knob with the other hand.
 - b. The clip can be rotated by spinning the catheter between your fingers in either direction where it enters the scope channel. Typically, this would be done by the physician.

Note: It is not recommended to utilize both rotation methods at the same time, as the rotation performance may be compromised.

Note: If there is no rotation response, do not continue to rotate the device more than three full rotations. Failure to do so may result in difficulty removing and/or scope damage.


[1] Thumb Ring [2] Slider [3] Handle [4] Rotation Control Knob
[5] Catheter [6] Clip [7] Clip Jaws

Figure 3. Rotating the Clip

3. To close the MANTIS Clip on the desired location, move the slider proximally until tactile resistance is felt in the handle, as shown in Figure 4. Clip position may now be assessed prior to deployment.

Caution: Do not continue moving the slider proximally beyond the tactile resistance until you are ready to deploy the clip, otherwise you may not be able to reopen the clip. If you hear or feel a click, the clip cannot be reopened, go to step 4 Option 2 to complete clip deployment.

[1] Thumb Ring [2] Slider [3] Handle [4] Rotation Control Knob
[5] Clip [6] Clip Jaws


Figure 4. MANTIS Clip Closed

4. At this point, there are 2 options:
 - Option 1: The MANTIS Clip can be reopened, rotated, and repositioned to the desired location (See steps 1-3).

Note: The MANTIS Clip is engineered to enable opening and closing up to five times prior to deployment, aiding in repositioning of the clip at the lesion site. Reopening, closing, and rotation capability may be limited by clinical circumstances and patient anatomy, among other factors.

- Option 2: The MANTIS Clip can be permanently deployed. To permanently deploy the MANTIS Clip, continue moving the slider proximally beyond the tactile resistance point, at which point a first click may be heard or felt. Continue moving the slider proximally until a second tactile resistance point and/or click is heard or felt. Continue moving the slider proximally until it reaches the thumb ring, as shown in Figure 5.

Note: Do not attempt to reopen the clip once the initial tactile resistance point has been passed and/or the first click has been heard or felt. Reopening the clip may result in the clip separating from the catheter and potentially kinking or damaging the device. After the first tactile resistance point is passed and/or click is observed, do not attempt to move the slider distally until both clicks are heard or felt. Do not continue to move the slider towards the thumb ring after both clicks are heard.

[1] Thumb Ring [2] Slider [3] Handle [4] Rotation Control Knob
[5] Clip [6] Clip Jaws

Figure 5. MANTIS Clip Permanently Deployed

5. Once the MANTIS Clip has been deployed, gently move the slider distally to separate the clip from the delivery device. Once the clip separates from the delivery device, release the slider.

Warning: Failure to release the slider after separation could result in patient injury.

If the MANTIS Clip has not been deployed, close the jaws, and withdraw the device slowly through the endoscope.

DEVICE REMOVAL

Withdraw the device slowly through the endoscope.

For partially deployed clip, attempt to fully deploy the clip. If unable to deploy clip, withdraw the scope and clip together.

Disposal

To minimize the risk of infection or microbial hazards after use, dispose device and packaging as follows:

After use, device and packaging may contain biohazardous substances. Any device and packaging that came into contact with biohazardous substances should be treated and disposed of as biohazardous waste or be treated and disposed of in accordance with any applicable hospital, administrative, and/or local government regulations. Use of a biohazardous container with biological hazard symbol is recommended. Untreated biohazardous waste should not be disposed of in the municipal waste system.

Post-Procedure

Any serious incident that occurs in relation to this device should be reported to the manufacturer and to the relevant local regulatory authority.

MRI SAFETY INFORMATION

 MRI SAFETY INFORMATION A person with the MANTIS Clip may be safely scanned under the following conditions. Failure to follow these conditions may result in injury.	
Device Name	MANTIS Clip
Static Magnetic Field Strength (B₀)	1.5T or 3.0T
Maximum Spatial Field Gradient	30 T/m (3,000 gauss/cm)
RF Excitation	Circularly Polarized (CP)
RF Transmit Coil Type	Cylindrical Whole-body Coil Cylindrical Head Coil
Operating Mode	Normal Operating Mode
Maximum Whole-Body SAR	2 W/kg (Normal Operating Mode)
Maximum Head SAR	3.2 W/kg (Normal Operating Mode)
Scan Duration / Temperature Rise	Under scan conditions defined above, MANTIS Clips are expected to have a temperature rise of less than 4°C and can be used for 60 minutes of continuous RF
MR Image Artifact	Image artifact caused by device may extend approximately 90 mm from the implant

Warning: Failure to follow the recommended MR Conditional labeling may result in a deployed GI hemoclip dislodging from tissue or heating of tissue at the GI hemoclip location. A GI hemoclip dislodgement may result in rebleeding requiring additional intervention or surgery, serious injury, or death.

Warning: This clip contains ferromagnetic material. Follow your institutional protocols to determine whether or not an x-ray should be performed prior to an MRI exam. There may be a small potential risk of clip dislodgement and rebleeding if the clip is used in friable or healing tissues due to magnetic forces acting on the clip when in or near an MRI scanner.

PATIENT COUNSELING INFORMATION

Provide the patient with the completed implant card to carry and explain that the Boston Scientific website has additional patient information with the summary of clip safety and clinical performance.

Inform the patient to present the implant card to their Healthcare professionals (doctors, dentist, technicians) including at MRI scans so they can take the necessary precautions.

Inform the patient if any serious incident that occurs in relation to this device should be reported to the manufacturer, and to the relevant local regulatory authority.

Brief the patient on any relevant post-procedure instructions, contraindications, warnings, precautions and/or adverse events found within this Instructions for Use (IFU), pertaining to the patient.

Implantable Device Patient Information

Advise the patient that additional information may be available to them on the Boston Scientific website (www.bostonscientific.com/patientlabeling).

Implant Card Instructions

- Apply the peel-off label from the product onto the supplied patient implant card
- Fill in the implantation date, patient name, name and address of healthcare institution and/or physician information
- Give the completed card to the patient

WARRANTY

For device warranty information, visit www.bostonscientific.com/warranty.

The following are trademarks of Boston Scientific Corporation or its affiliates: MANTIS.

All other trademarks are the property of their respective owners.

SYMBOL DEFINITIONS

Commonly used medical device symbols that appear on the labeling are defined at www.bostonscientific.com/SymbolsGlossary.

Additional symbols are defined at the end of this document.

INSTRUCTIONS FOR USE

en

The Instructions for Use (IFU) for this product are supplied in electronic form over the Internet.

Visit www.IFU-BSCI.com to access the IFU using your internet browser. Information from the product labeling such as Product Code, UDI-DI or Product Name may be used to search for the specific IFU.

If you have difficulty accessing the IFU online, or would prefer to receive a paper copy, please contact Boston Scientific Customer Service or your local country contact. A copy will be sent to you at no charge and should arrive within seven days.

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or its affiliates. All rights reserved.
Boston Scientific (Australia) Pty Ltd
Boston Scientific Limited
© 2023 Boston Scientific Corporation