OverStitch™

English page 2
[Endoscopic Suturing System] (ESS)
INSTRUCTIONS FOR USE

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ESS-G02-160
PLY-G02-020-A
PLY-G02-020-APL
CNH-G01-000
THX-165-028

Caution:

Please read all instructions prior to use

SINGLE PATIENT USE
DISPOSABLE

GRF-00002-00R23

English

Single use only. Disposable. Do not resterilize.

Caution: Federal law (USA) restricts this device to sale by or on the order of a physician.

Patent Pending.

DISCLAIMER OF WARRANTY AND LIMITATION OF REMEDY

THERE IS NO EXPRESS OR IMPLIED WARRANTY, INCLUDING WITHOUT LIMITATION ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE, ON THE APOLLO ENDOSURGERY, INC. PRODUCT(S) DESCRIBED IN THIS PUBLICATION. TO THE FULLEST EXTENT PERMITTED BY APPLICABLE LAW, APOLLO ENDOSURGERY, INC. DISCLAIMS ALL LIABILITY FOR ANY INDIRECT, SPECIAL, INCIDENTAL, OR CONSEQUENTIAL DAMAGES, REGARDLESS OF WHETHER SUCH LIABILITY IS BASED ON CONTRACT, TORT, NEGLIGENCE, STRICT LIABILITY, PRODUCTS LIABILITY OR OTHERWISE. THE SOLE AND ENTIRE MAXIMUM LIABILITY OF APOLLO ENDOSURGERY, INC., FOR ANY REASON, AND BUYER'S SOLE AND EXCLUSIVE REMEDY FOR ANY CAUSE WHATSOEVER, SHALL BE LIMITED TO THE AMOUNT PAID BY THE CUSTOMER FOR THE PARTICULAR ITEMS PURCHASED. NO PERSON HAS THE AUTHORITY TO BIND APOLLO ENDOSURGERY, INC. TO ANY REPRESENTATION OR WARRANTY EXCEPT AS SPECIFICALLY SET FORTH HEREIN. DESCRIPTIONS OR SPECIFICATIONS IN APOLLO ENDOSURGERY, INC. PRINTED MATTER, INCLUDING THIS PUBLICATION, ARE MEANT SOLELY TO GENERALLY DESCRIBE THE PRODUCT AT THE TIME OF MANUFACTURE AND DO NOT CONSTITUTE ANY EXPRESS WARRANTIES OR RECOMMENDATIONS FOR USE OF THE PRODUCT IN SPECIFIC CIRCUMSTANCES. APOLLO ENDOSURGERY, INC. EXPRESSLY DISCLAIMS ANY AND ALL LIABILITY, INCLUDING ALL LIABILITY FOR ANY DIRECT, INDIRECT, SPECIAL, INCIDENTAL, OR CONSEQUENTIAL DAMAGES, RESULTING FROM REUSE OF THE PRODUCT.

Table of Symbols

Description	Symbol	Description	Symbol	Description	Symbol
Consult Instructions for Use		Manufacturer		Use By	
Do Not Re-use		Caution: Federal law (USA) restricts this device to sale by or on the order of a physician		Lot Number	
Sterilized Using Ethylene Oxide		Date of manufacture		Do not use if package is damaged	
Reference Number		Do Not Resterilize		Authorized Representative in the European Community	
Medical Device		Caution, contact accompanying documents		Contains hazardous substance	

1. Intended Use

The Apollo Endosurgery OverStitch™ Endoscopic Suture System (ESS) is intended for endoscopic placement of suture(s) and approximation of soft tissue.

1.1 Contraindications

Contraindications include those specific to use of an endoscopic suturing system, and any endoscopic procedure, which may include, but not limited to, the following:

- This system is not for use where endoscopic techniques are contraindicated.
- This system is not for use with malignant tissue.

1.2 Warnings

- Do not use a device where the integrity of the sterile packaging has been compromised or if the device appears damaged.
- Only physicians possessing sufficient skill and experience in similar or the same techniques should perform endoscopic procedures.
- Contact of electrosurgical components with other components may result in injury to the patient and/or operator as well as damage to the device and/or endoscope.
- Verify compatibility of endoscopic instruments and accessories and ensure performance is not compromised.
- NOTE: Refurbished scopes may no longer confirm to original specifications.
- Ensure that there is sufficient space for the Needle to open.
- Ensure that the Handle Grip of the Endoscopic Suturing System is closed and locked during intubation and extubation.
- Users should be familiar with surgical procedures and techniques involving absorbable sutures before employing Synthetic Absorbable Sutures for wound closure, as the risk of wound dehiscence may vary with the site of application and the suture material used.
- Reuse or reprocessing of the OverStitch system could result in device malfunction or patient consequences to include:
 - Infection or the transmission of disease
 - Failure of the handle mechanism causing the device to become locked on tissue that may require surgical intervention
 - Reduced retention on the Endoscope, causing the Endcap to detach during use that may require surgical intervention to retrieve
 - Reduced retention of the Anchor to the Needle Body, resulting in an inadvertent Anchor drop causing procedural delay or requiring subsequent intervention
 - Bending of the Needle Body, preventing the physician from driving the Needle correctly or performing the intended procedure
 - Failure of the Helix to extend fully, limiting the ability to acquire tissue and perform intended procedure
- If subject device is used to oversee foreign objects, such as staples, stents, clips or mesh, it is possible for the needle to become trapped in the foreign body, requiring surgical intervention.

- In situations where the operative site poses a risk of harm to adjacent anatomic structures, use of endoscopic accessories such as the OverStitch Tissue Helix is recommended to retract the tissue intended to be sutured away from these unseen structures.

1.3 Precautions

- The System may only be used if purchased from Apollo Endosurgery, Inc. or one of its authorized agents.
- With the Endoscopic Suturing System installed, the endoscope's primary channel effectively becomes a 3.2 mm channel.
- An overtube with an internal diameter of at least 16.7 mm may be used with the system to protect the esophagus.

1.4 System Compatibility

The OverStitch ESS is compatible with the sutures identified on Page 1.

The OverStitch is compatible with the following endoscopes:

- Olympus 2T160, 2TH180, or 2T240
- Fuji EI-740D/S

Availability of compatible devices may vary depending on geography.


1.5 Adverse Events

Possible complications that may result from using the Endoscopic Suturing System include, but may not be limited to:

- Pharyngitis / Sore throat
- Nausea and / or Vomiting
- Abdominal pain and / or Bloating
- Hemorrhage
- Hematoma
- Conversion to laparoscopic or open procedure
- Stricture
- Infection / Sepsis
- Pharyngeal, colonic and/or esophageal perforation
- Esophageal, colonic and/or pharyngeal laceration
- Intra-abdominal (hollow or solid) visceral injury
- Aspiration
- Wound dehiscence
- Acute inflammatory tissue reaction
- Death


NOTE: Any serious incident that has occurred in relation to the device should be reported to Apollo Endosurgery (see contact information at the end of this document) and any appropriate government entity.

2. Packaged System


Package 1


Needle Driver & Anchor Exchange

Package 2
Suture Assembly

Package 3
Cinch

Package 4
Helix
(AVAILABLE SEPARATELY)

Polypropylene (nonabsorbable) - BLUE

3. Nomenclature


Endoscope Handle


Endoscope


Needle Driver Handle

Needle Driver

Endcap

Anchor Release Button


Anchor Exchange


Suture Assembly

Anchor Body

Suture

Helix Handle


Helix

Helix Tip

Suture Loading Loop

Collar Plug

Cinch Handle

Cinch


English

4. Assembly


- 4.1. Place the compatible endoscope on a suitable surface for assembly.

NOTE: Ensure standard channel valves are fitted to working channels.

- 4.2. Remove Needle Driver from packaging.


- 4.3. Close Needle Body.

NOTE: Ensure Endcap is not dropped or otherwise damaged.

- 4.4. Push alignment tube into the primary working channel until Endcap is flush with Endoscope face.


NOTE: Resistance should be felt when attaching the Endcap as this indicates correct attachment.

NOTE: Ensure that the secondary working channel is not obstructed on the Endoscope.


WARNING: Ensure Endcap is fully seated on Endoscope.

- 4.5. Assemble Needle Driver handle onto Endoscope.

NOTE: Ensure external cable is not twisted around Endoscope.

- 4.6. Connect Endoscope to Tower.


Load Anchor

- 4.7. Select appropriate Suture type / size.

- 4.8. Remove Suture Assembly and Anchor Exchange from packaging.

- 4.9. Load the Anchor.


- 4.10. Remove Suture from Suture Cartridge by holding and pulling the Suture, not the Anchor or Anchor Exchange.

NOTE: Ensure Suture is not tangled after removal from the cartridge.

- 4.11. Open valve cover and insert Anchor Exchange into primary channel of Endoscope.

CAUTION: Do not use when valve covers are closed as suture drag will be increased.

- 4.12. Using a 'pencil-grip' on the catheter for optimum control, advance Anchor Exchange until Anchor is positioned near distal end of Endoscope.


CAUTION: If resistance is encountered when advancing the Anchor Exchange through the working channel of the Endoscope, reduce the Endoscope angulation until the device passes smoothly.

- 4.13. Insert Endoscope into patient.

WARNING: Do not introduce the device with the Needle Body in its open position.


5. Navigating to target anatomy

- 5.1. Advance until target anatomy is located.


CAUTION: When intubating or extubating with the Endoscopic Suturing System, ensure that the working length of the Endoscope and the actuation catheter are advanced and retracted together.


- 5.2. Open Needle Body.


- 5.3. Advance Anchor Exchange and/or manipulate Endoscope to create Suture slack.


- 5.4. Once sufficient slack has been created, retract Anchor Exchange into the Endoscope.

- 5.5. Close Needle Body.

- 5.6. Use a 'pencil-grip' on the **white** portion of the cable to advance Anchor Exchange until Anchor fully seats onto Needle Body.


- 5.7. A tactile 'click' or firm stop may be felt as the Anchor seats fully onto Needle Body.

- 5.8. Check the monitor image to ensure Anchor is properly installed on Needle Body.


- 5.9. Ensure that the proximal end of the Suture is visible beyond the Endoscope channel valves.

- 5.10. Fully depress Anchor Release Button to release Anchor.


- 5.11. With Anchor Release Button still fully depressed, slightly retract Anchor Exchange.

6. Tissue and Suture Management

- 6.1. Open Needle Body.

CAUTION: If the Needle Body does not open, ensure that the Anchor has been released from the Anchor Exchange.


WARNING: Ensure appropriate Suture slack has been created for desired Suture path and pattern. Advance Anchor Exchange and/or manipulate Endoscope to create Suture slack.

- 6.2. Position tissue in appropriate location for suturing, using Apollo Helix or Endoscope compatible accessory if required.


7. Use of Tissue Helix (Optional)

- 7.1. Advance Helix into secondary channel of Endoscope in the retracted position, until distal tip is visible on the monitor.


CAUTION: If resistance is encountered when advancing the Helix through the working channel of the Endoscope, reduce the Endoscope angulation until the device passes smoothly, and ensure that the secondary working channel is not obstructed on the Endoscope.

- 7.2. Fully depress Helix Handle button to expose Helix Tip.


CAUTION: Do not depress Helix Handle button while advancing Helix through Endoscope.

- 7.3. Acquire tissue by rotating Helix Handle clockwise until correct tissue depth is achieved.


- 7.4. Maintain gentle forward pressure during tissue acquisition.
7.5. Advance / retract Helix to position tissue into desired location.

8. Suturing


- 8.1. Drive Needle through tissue by closing the Needle Body.

- 8.2. Ensure that the Needle Arm does not inadvertently close on any foreign object or device.


Retrieve Anchor

- 8.3. Use a 'pencil-grip' to advance the Anchor Exchange until the Anchor is engaged and resistance is felt.

NOTE: Resistance may vary due to the position of the Endoscope.

- 8.4. Retract Anchor Exchange to acquire Anchor.


- 8.5. Use a 'pencil-grip' on the white section of the cable and place the remaining fingers of the same hand on the Endoscope housing to prevent damage to Suture or tissue once the Anchor 'pops' off the Needle Body.


CAUTION: Do not press the Anchor Release Button, as this may cause inadvertent drop of the Anchor.
WARNING: If sufficient slack has not been created prior to driving the Anchor through tissue, retraction of the Anchor Exchange may be difficult and the Anchor may not release correctly from the Needle Body.

Release Tissue


- 8.6. Rotate Helix handle counter-clockwise until the device is free from tissue.

- 8.7. Return the Helix handle button.


- 8.8. Check monitor to ensure tip is fully retracted into shaft before withdrawing into Endoscope.

- 8.9. Withdraw Helix a short distance into Endoscope.

- 8.10. Open Needle Body.

CAUTION: Do not tension the Suture with the Anchor in the Needle Body.


- 8.11. To continue placing stitches with this Anchor, repeat sections 5-8.
If stitching is complete for this Anchor, proceed to section 9 to approximate tissue, secure and cut the Suture.

Multiple Anchors can be utilized with each Needle Driver and Anchor Exchange.


9. Securing and Cutting Suture

- 9.1. Remove accessories (Apollo Helix or similar) from Endoscope.


- 9.2. Ensure Anchor is in Anchor Exchange and open Needle Body.

- 9.3. Advance Anchor distal of the alignment tube.


- 9.4. Fully depress Anchor Release Button and retract Anchor Exchange to release Anchor.

CAUTION: Do not release the Anchor inside the Endoscope's working channel.

- 9.5. Remove Anchor Exchange from Endoscope.

NOTE: The Anchor Exchange can be utilized for additional sutures.

- 9.6. Feed proximal end of Suture into removable Suture Loading Loop.


- 9.7. After threading, release the proximal end of the Suture to enable loading.


- 9.8. Pull Suture Loading Loop parallel to the device, to pull Suture into Cinch.

- 9.9. Holding the proximal end of the Suture, feed the Cinch down the working channel until 'Plug and Collar' can be seen on the monitor.


- 9.10. Pull the Suture and apply counter traction to the Cinch until the tissue is approximated and the desired Suture tension is achieved between the Anchor and the Cinch Collar.

NOTE: It is the Collar that sets the final position of the Cinch, not the Plug.


WARNING: Excessive tension may damage tissue.

- 9.11. Remove safety spacer from Cinch handle.

CAUTION: The safety spacer must only be removed immediately prior to deploying Cinch.


- 9.12. Firmly squeeze Cinch Handle to deploy Cinch and cut Suture.


CAUTION: Suture tension must be maintained.

NOTE: Significant force is required to pull and lock the Plug into the Collar, a 'pop' can often be felt once the Suture is cut.

- 9.13. Remove the Cinch.

- 9.14. Close Needle Body.

10. Multiple Sutures

- 10.1. Anchor Exchange and Needle Driver can be utilized with multiple sutures. Removal of the Endoscope after cinching is not required if further sutures are to be deployed.

WARNING: If the Endoscope is removed between stitching for cleaning, ensure the Endcap is fully seated before next intubation. Reseat if necessary.


NOTE: If removed and reinstalled, ensure that the secondary working channel is not obstructed on the Endoscope.

- 10.2. To use another Suture Assembly, return to step 4.7 for Anchor loading and follow all subsequent steps. If suturing is complete, proceed to section 11 for removal of the device.

English

11. Removal of Device

- 11.1 Remove any accessory devices from Endoscope.
- 11.2 Ensure the Needle Body is closed and retract Endoscope from patient, making sure the external cable is retracted together with the Endoscope.
- 11.3 Remove Endcap from Endoscope.

NOTE: Do not hold Needle Arm to remove the Endcap.

- 11.4 Remove Needle Driver from Endoscope by flexing the scope attachment bracket around the working channels.

12. MRI Safety Information

MR Conditional

Non-clinical testing has demonstrated that the Sutures, Cinches and Anchors (collectively termed Anchoring System) deployed by the OverStitch Endoscopic Suturing System are MR Conditional.

A patient with this Anchoring System can be safely scanned immediately after placement in an MR system meeting the following conditions:

Static Magnetic Field

- Static magnetic field of 1.5 T or 3.0 T.
- Maximum spatial field gradient of 2,000 gauss/cm (20 T/m).
- Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 2 W/kg.

Under the scan conditions defined above, the Anchoring System is expected to produce a maximum temperature rise of less than 2 °C after 15 minutes of continuous scanning.

In non-clinical testing, the image artifact caused by the Anchoring System extends approximately 10 mm from this device when imaged with a gradient echo pulse sequence and a 3.0 T MRI system.

13. Hazardous Substance Information

Contains Hazardous Substance

CMR Statement – The stainless steel and cobalt alloy components within this device contain the following substance(s) defined as a CMR (carcinogenic, mutagenic or toxic to reproduction) 1A/1B and/or endocrine disrupting in a concentration above 0.1% weight by weight:

Cobalt (CAS No. 7740-48-4; EC No. 231-158-0)

Current scientific evidence supports medical devices manufactured from these cobalt alloys or stainless steels containing cobalt do not cause increased risk of cancer or adverse reproductive effects.

14. Troubleshooting

14.1. Needle Body will not open:

Cause	Resolution
14.1.1. Needle obstructed:	i. Assess the space in which you are working and maneuver the catheter and Endoscope as a system, straighten the Endoscope to the non-Retroflexed position.
14.1.2. Suture movement restricted:	i. If the Anchor is on the Needle Body, ensure the suture is not held proximally near the handle during the opening operation. ii. Transfer the Anchor to the Anchor Exchange. Open the Needle Body. Slowly retract the Anchor Exchange proximally and then advance the Needle Body distally to free the Suture.
14.1.3. Needle Driver catheter looped or kinked:	i. Check to ensure that the catheter running down the outside of the Endoscope is not looped or kinked. Straighten the Endoscope to the non-Retroflexed position. Advance the Endoscope forward and slightly pull any slack from the actuation catheter proximally until minimal resistance is felt. Grasp both the actuation catheter and Endoscope and adjust by advancing and retracting as a system.
14.1.4. Foreign body obstruction:	i. Remove the Anchor Exchange only: a. Transfer the Anchor to the Needle Body and remove the Anchor Exchange from the Endoscope. b. Load a grasper through the primary channel and push the Needle Body open. ii. Remove the Anchor and Anchor Exchange, cutting the Suture if necessary: a. Transfer the Anchor to the Needle Body and remove the Anchor Exchange from the Endoscope. b. Through either channel use an appropriate accessory to cut the Suture. c. Use an accessory to push open the Needle Body d. Use appropriate means to remove the cut Suture. iii. Once standard endoscopic techniques have been exhausted, utilize laparoscopic techniques to remove the device.

14.2. Needle Body will not close:

Cause	Resolution
14.2.1. General obstruction:	i. Follow steps 14.1.1, 14.1.2, 14.1.3 above (Needle Body will not open). ii. Ensure the Needle Driver handle is locked closed and: a. Stretch the Needle Driver catheter to change the effective length of the Needle Body drive cable. b. Remove the Anchor Exchange and use graspers (through the primary channel) to grab the Needle Body. Cut the Suture if necessary.
14.2.2. Needle Driver cable broken:	i. Deploy the Anchor and Cinch. Advance grasper through endoscope and secure Needle Driver. Pull Needle Driver closed while removing the device. If using an overtube, advance the overtube as far distally as possible and withdraw the endoscope and device into the overtube, using the distal tip of the overtube to close the Needle Body.

14.3. Anchor Exchange will not exchange:

Cause	Resolution
14.3.1. Anchor Exchange will not install Anchor onto the Needle Body:	i. Ensure there is sufficient Suture slack and the Suture outside Endoscope is not entangled. ii. Ensure Anchor Exchange is properly positioned in alignment tube of Needle Driver. iii. If Anchor and Suture are through tissue, either drop Anchor and deploy Cinch per Section 9 of Instructions for Use or drop Anchor and use a suitable accessory to cut and remove the Suture. iv. If Anchor and Suture are not through tissue, close handle of Needle Driver. Remove Endoscope. Replace Anchor and/or Anchor Exchange.
14.3.2. Anchor Exchange will not release an Anchor:	i. Ensure there is sufficient Suture slack and the Suture outside Endoscope is not entangled. ii. Ensure the Anchor Exchange Release Button is FULLY depressed while retracting Anchor Exchange.

Cause	Resolution
	<ul style="list-style-type: none"> iii. Reduce the articulation / tortuosity of the Endoscope (if possible) and try to release the Anchor. iv. Use accessories compatible with Endoscope secondary working channel to cut and remove the Suture. v. Replace the Anchor Exchange.
14.3.3. Anchor Exchange will not retrieve Anchor from Needle Body:	<ul style="list-style-type: none"> i. Ensure there is sufficient Suture slack and Needle Driver handle is in closed position. ii. Use 'pencil-grip' to advance Anchor Exchange until the Anchor is engaged and resistance is felt. Retract Anchor Exchange to acquire Anchor. iii. If the Anchor cannot be retrieved, replace the Anchor Exchange. Alternatively, use a suitable accessory to cut and remove the Suture. iv. Replace the Anchor and resume suturing per <i>Section 8</i> of Instructions for Use.

14.4. Cinch does not cut the Suture when fired:

Cause	Resolution
14.4.1. Suture not cut:	<ul style="list-style-type: none"> i. Use a suitable accessory through the secondary working channel to cut the Suture and remove the Cinch. ii. Use standard endoscopic techniques to remove the cut Suture.

14.5. Inadvertent drop of Anchor:

Cause	Resolution
14.5.1. Anchor Exchange button depressed out of sequence:	<ul style="list-style-type: none"> i. Retrieve the Anchor as a foreign body, or follow the procedure to Cinch in place. If the Anchor is dropped inside the working channel, use the Anchor Exchange or 3.2 mm compatible grasper to push the Anchor out through the Endoscope. <p>Do not attempt to pull the Anchor back through the Endoscope as it may get stuck in the channel or y-junction at the Endoscope handle.</p>

14.6. Suture entanglement:

Cause	Resolution
14.6.1. Suture outside field of view:	<ul style="list-style-type: none"> i. Close the Needle Driver Handle and manipulate the Endoscope back to release.
14.6.2. Suture behind Tissue Guard:	<ul style="list-style-type: none"> i. Slightly close the Needle Body while retracting the Endoscope. ii. If necessary, transfer the Anchor to the Anchor Exchange. iii. Open the Needle Body and advance the Anchor Exchange beyond the Endcap to push the Suture free.
14.6.3. Suture twisted:	<ul style="list-style-type: none"> i. If the Suture is twisted, move the Endoscope and transfer the Anchor between the Needle Body and Anchor Exchange, on the opposite side of the suture thread, as needed to untwist. <p>NOTE: If the twist was noted immediately after articulation of the Endoscope, first try to articulate in reverse order to remove.</p> <ul style="list-style-type: none"> ii. If the Anchor has been deployed, use the Cinch to push and guide the Suture free.

14.7. Endcap detachment from Endoscope:

Cause	Resolution
14.7.1. Pushed off during use:	<ul style="list-style-type: none"> i. Close the Needle Body. Using a rat-toothed grasper, capture the alignment tube "zigzag" cutout and remove any slack in the Actuation Catheter and slowly remove the device from the patient. If using an overtube, using a rat-toothed grasper, capture the alignment tube "zigzag" cutout and remove any slack in the Actuation Catheter and slowly remove the device from the patient.

14.8. Helix does not unscrew:

Cause	Resolution
14.8.1. Helix stuck in tissue:	<ul style="list-style-type: none"> i. Use a suitable accessory through the primary channel to apply counter traction to the tissue around the Helix, and pull the Helix free. ii. Once endoscopic techniques have been exhausted, utilize laparoscopic techniques to remove the Helix.

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All rights reserved.