



ACURATE™ Aortic Valve System Publication Compendium

CAUTION: The ACURATE Prime Aortic Valve System is an investigational device. Not available for sale. In Europe, ACURATE neo2 Aortic Valve System is CE-marked. In the USA, ACURATE neo2 is an investigational device and restricted under federal law to investigational use only. Not available for sale. © 2023 Boston Scientific Corporation or its affiliates. All rights reserved. SH-1376306-AD

TABLE OF CONTENTS:

BOSTON SCIENTIFIC SPONSORED CLINICAL TRIALS	4
INVESTIGATOR SPONSORED RESEARCH	9
EXTERNALLY DRIVEN EVIDENCE (BY TOPIC)	10
HEAD-TO-HEAD COMPARISONS	10
LENGTH OF STAY	12
CONDUCTION DISORDERS	12
SMALL ANNULI	14
CORONARY OBSTRUCTION	14
CORONARY ACCESS / COMMISSURAL ALIGNMENT	15
PARAVALVULAR LEAKAGE AND AORTIC REGURGITATION	16
PRE-DILATATION AND POST-DILATATION	18
TAVI HEMODYNAMICS AND DURABILITY	18
HORIZONTAL AORTA	19
EX VIVO ANALYSIS	19
GENERAL TOPICS	22
OFF-LABEL USE (SAFETY AND EFFECTIVNESS HAVE NOT BEEN ESTABLISHED IN THESE PATIENT POPULATIONS)	24
VALVE-IN-VALVE	24
PURE AORTIC INSUFFICIENCY / AORTIC REGURGITATION	24
BICUSPID VALVE	25
ALTERNATIVE ACCESS ROUTE	25

INTRODUCTION

PUBLICATIONS LISTED IN THIS COMPENDIUM ARE BOSTON SCIENTIFIC SPONSORED, INVESTIGATOR SPONSORED OR INDEPENDENTLY LED ANALYSES WHICH INCLUDE THE ACURATE™ AORTIC VALVE SYSTEM (ACURATE Prime, ACURATE neo2 AND ACURATE neo).





THE SAFETY AND EFFECTIVENESS OF THE ACURATE™ AORTIC VALVE PLATFORM FOR TRANSCATHETER AORTIC VALVE REPLACEMENT (TAVR) IN SUBJECTS WITH SEVERE NATIVE AORTIC STENOSIS IS ASSESSED VIA THE ACURATE IDE TRIAL – SEE FULL ACURATE IDE DETAILS AT [CLINICALTRIALS.GOV](https://clinicaltrials.gov).

EACH PUBLICATION ENTRY INCLUDES (IF AVAILABLE) ITS FULL REFERENCE, AND A LINK TO THE FULL ARTICLE (THIS LINK MAY REQUIRE THE READER TO PAY FOR THE FULL CONTENT). BOSTON SCIENTIFIC SPONSORED TRIALS INCLUDE PUBLICATION ABSTRACT.

THIS DOCUMENT WILL BE UPDATED AS FUTURE PUBLICATIONS BECOME AVAILABLE.



ACURATE PUBLICATIONS ARE COLOR CODED BY PLATFORM THROUGHOUT THIS DOCUMENT

-  ACURATE Prime™
Aortic Valve System Publication
-  ACURATE neo2™
Aortic Valve System Publication
-  ACURATE neo2™ and ACURATE neo™
Aortic Valve System Publication
-  ACURATE neo™
Aortic Valve System Publication

BOSTON SCIENTIFIC SPONSORED CLINICAL TRIALS

Outcomes and performance of the ACURATE neo2 transcatheter heart valve in clinical practice: one-year results of the ACURATE neo2 PMCF Study.

Kim WK, Tamburino C, Möllmann H, Montorfano M et al., EuroIntervention. 2023 Nov 19;EIJ-D-23-00823. doi: 10.4244/EIJ-D-23-00823. Epub ahead of print. PMID: 37982152.

Abstract:

Background: Transcatheter aortic valve implantation is an effective treatment for patients with aortic stenosis; however, complications related to paravalvular leakage (PVL) persist, including increased risk of mortality, cardiovascular mortality, and rehospitalisation.

Aims: We sought to evaluate the clinical outcomes and valve performance at 1 year in patients with severe aortic stenosis treated with the ACURATE neo2 valve in a post-market clinical setting.

Methods: Valve Academic Research Consortium-2 safety events were assessed up to 1 year. Independent core laboratories evaluated echocardiographic measures of valve performance and hypoattenuated leaflet thickening (HALT; as measured by four-dimensional computed tomography).

Results: The study enrolled 250 patients (64% female; mean age: 81 years; baseline Society of Thoracic Surgeons risk score: 2.9±2.0%); 246 patients were implanted with ACURATE neo2. All-cause mortality was 0.8% at 30 days and 5.1% at 1 year. The 1-year rates for stroke and disabling stroke were 3.0% and 1.3%, respectively. Overall, HALT of >50% leaflet involvement of at least one leaflet was present in 9% of patients at 30 days and in 12% of patients at 1 year. No association was observed between the presence of HALT and 1-year clinical or haemodynamic outcomes. Early haemodynamic improvements were maintained up to 1 year (mean aortic valve gradient: 47.6±14.5 mmHg at baseline, 7.6±3.2 mmHg at 1 year; mean aortic valve area: 0.7±0.2 cm² at baseline, 1.7±0.4 cm² at 1 year). At 1 year, 99% of patients had mild or no/trace PVL (<1% had moderate PVL; no patient had severe PVL).

Conclusions: The study outcomes confirm favourable performance and safety up to 1 year in patients treated with ACURATE neo2 in routine clinical practice.

Clinical outcomes of the ACURATE neo2 transcatheter heart valve: a prospective, multicenter, observational, post-market surveillance study.

Kim WK, Tamburino C, Möllmann H, Montorfano M et al., EuroIntervention. 2022 Nov 27;19(1):83–92. doi: 10.4244/EIJ-D-22-00914. Epub ahead of print. PMID: 36440588; PMCID: PMC10173758.

Abstract:

Background: The next-generation ACURATE neo2 transcatheter aortic valve was designed for simplified implantation and to mitigate the risk of paravalvular leak (PVL) compared to the earlier device.

Aims: To collect clinical outcomes and device performance data, including echocardiography and 4-dimensional computed tomography (4D-CT) data, with the ACURATE neo2 transcatheter heart valve in patients with severe aortic stenosis (AS).



Methods: ACURATE neo2 PMCF is a single-arm, multicentre study of patients with severe AS treated in routine clinical practice. The primary safety endpoint was all-cause mortality at 30-days. The primary imaging endpoint was hypo-attenuated leaflet thickening (HALT), measured by core laboratory-adjudicated 4D-CT at 30 days. Secondary endpoints included VARC safety endpoints, procedural success, and evaluation of valve performance via core laboratory-adjudicated echocardiography.

Results: The study enrolled 250 patients at 18 European centres (mean age: 80.8 years; 63.6% female; mean STS score: $2.9 \pm 2.0\%$); 246 (98.4%) were successfully treated with ACURATE neo2. The 30-day rates for mortality and disabling stroke were 0.8% and 0%, respectively. The new permanent pacemaker implantation rate was 6.5%. HALT >50% was present in 9.3% of patients at 30 days. Valve haemodynamics improved from baseline to 30 days (mean aortic valve gradient: from 47.6 ± 14.5 mmHg to 8.6 ± 3.9 mmHg; mean aortic valve area: from 0.7 ± 0.2 cm² to 1.6 ± 0.4 cm²). At 30 days, PVL was evaluated as none/trace in 79.2% of patients, mild in 18.9%, moderate in 1.9%, and severe in 0%.

Conclusions: The study results support the safety and efficacy of TAVI with ACURATE neo2 in patients in routine clinical practice.

The ACURATE neo2 valve system for transcatheter aortic valve implantation: 30-day and 1-year outcomes.

Möllmann H, Holzhey DM, Hilker M et al., Clin Res Cardiol. 2021 Dec;110(12):1912-1920. doi: 10.1007/s00392-021-01882-3. Epub 2021 Jun 20. PMID: 34148125; PMCID: PMC8639565.

Abstract:

Background: Transcatheter aortic valve implantation (TAVI) has become standard treatment for elderly patients with symptomatic severe aortic valve stenosis. The ACURATE neo AS study evaluates 30-day and 1-year clinical and hemodynamic outcomes in patients treated with the ACURATE neo2 valve.

Methods: The primary endpoint of this single-arm multicenter study is 30-day all-cause mortality. Other key endpoints include device performance, echocardiographic measures assessed by an independent core laboratory, and VARC-2 clinical efficacy and safety endpoints through 12 months.

Results: The study enrolled 120 patients (mean age 82.1 ± 4.0 years; 67.5% female, mean baseline STS score $4.8 \pm 3.8\%$). The VARC-2 composite safety endpoint at 30 days occurred in 13.3% of patients. All-cause mortality was 3.3% at 30 days and 11.9% at 1 year. The 30-day stroke rate was 2.5% (disabling stroke 1.7%); there were no new strokes between 30 days and 12 months. The rate of permanent pacemaker implantation was 15.0% (18/120) at 30 days and 17.8% (21/120) at 1 year. No patients required re-intervention for valve-related dysfunction and there were no cases of valve thrombosis or endocarditis. Patients demonstrated significant improvement in mean aortic valve gradient (baseline 38.9 ± 13.1 mmHg, 1 year 7.8 ± 3.5 mmHg; $P < 0.001$ in a paired analysis). In the overall population, paravalvular leak was evaluated at 1 year as none/trace in 60.5%, mild in 37.0%, and moderate in 2.5%; no patients had severe PVL.

Conclusions: One-year outcomes from the ACURATE neo AS study support the safety and performance of TAVI with the ACURATE neo2 valve.

Keywords: Aortic valve stenosis; Paravalvular regurgitation; Transcatheter aortic valve replacement; Transfemoral aortic valve implantation



The SAVI-TF Registry: 1-Year Outcomes of the European Post-Market Registry Using the ACURATE neo Transcatheter Heart Valve Under Real-World Conditions in 1,000 Patients.

Kim WK, Hengstenberg C, Hilker M et al., JACC Cardiovasc Interv. 2018 Jul 23;11(14):1368-1374. doi: 10.1016/j.jcin.2018.03.023. PMID: 30025731

Abstract:

Objectives: The SAVI-TF (Symetis ACURATE neo Valve Implantation Using Transfemoral Access) registry was initiated to study the ACURATE neo transcatheter heart valve in a large patient population treated under real-world conditions.

Background: The self-expanding, supra-annular ACURATE neo prosthesis is a transcatheter heart valve that gained the Conformité Européenne mark in 2014, but only limited clinical data are available so far.

Methods: This prospective, multicenter registry enrolled 1,000 patients at 25 European centers who were followed for 1 year post-procedure.

Results: Mean patient age was 81.1 ± 5.2 years; mean logistic European System for Cardiac Operative Risk Evaluation I score, European System for Cardiac Operative Risk Evaluation II score, and Society of Thoracic Surgeons score were $18.1 \pm 12.5\%$, $6.6 \pm 7.5\%$, and $6.0 \pm 5.6\%$, respectively. At 1 year, 8.0% (95% confidence interval [CI]: 6.3% to 9.7%) of patients had died, 2.3% (95% CI: 1.3% to 3.2%) had disabling strokes, and 9.9% (95% CI: 8.1% to 11.8%) had permanent pacemaker implantations. Through 1 year, 5 reinterventions (0.5%; 95% CI: 0.1% to 1.0%) were performed: 3 valve-in-valve and 2 surgical aortic valve replacements. Mean effective orifice area was 1.84 ± 0.43 cm², mean gradient was 7.3 ± 3.7 mm Hg, and greater than mild paravalvular leakage was observed in 3.6% of patients.

Conclusions: Transfemoral implantation of the ACURATE neo prosthesis resulted in favorable 1-year clinical and echocardiographic outcomes with very low mortality and new pacemaker rates.

Real-world experience using the ACURATE neo prosthesis: 30-day outcomes of 1000 patients enrolled in the SAVI TF registry.

Möllmann H, Hengstenberg C, Hilker M et al., EuroIntervention. 2018 Feb 2;13(15): e1764-e1770. doi: 10.4244/EIJ-D-17-00628.

Abstract:

Aims: The aim of the SAVI TF registry was to assess the safety and performance of the self-expanding ACURATE neo transfemoral transcatheter heart valve in a large patient population with severe aortic stenosis and to investigate whether the outcomes obtained in the CE-mark cohort can be replicated in an unselected all-comers population.

Methods and results: From October 2014 until April 2016, 1,000 patients were enrolled in this prospective, European multicentre registry. Patients were 81.1 ± 5.2 years and had a logistic EuroSCORE II and STS score of $6.6 \pm 7.5\%$ and $6.0 \pm 5.6\%$, respectively. Pre-dilatation was performed in 96.1% of patients and post-dilatation in 44.8%. Procedural and device success were both obtained in 98.7%; failure comprised nine valve-in-valve procedures, three conversions to surgery, and one aborted procedure. The primary endpoint was 30-day mortality, which was observed in 14 patients (1.4% [95% CI: 0.7-2.1]).



Disabling stroke was seen in 1.2% (95% CI: 0.5-1.9) and new pacemaker implantation in 8.3% (95% CI: 6.6-10.0). At discharge, mean effective orifice area was 1.77 ± 0.46 cm² and mean gradient 8.4 ± 4.0 mmHg; 4.1% of patients had a more than mild paravalvular leak.

Conclusions: In this initial experience, treatment with the ACURATE neo prosthesis resulted in good clinical outcomes with very low complication rates.

Transcatheter Aortic Valve Implantation with ACURATE neo: Results from the PROGRESS PVL Registry.

Won-Keun Kim, Holger Thiele, Axel Linke et al., J Interv Cardiol. 2022 Jun 25;2022:9138403. doi: 10.1155/2022/9138403.

Abstract:

Objectives: The PROGRESS PVL registry evaluated transcatheter aortic valve implantation (TAVI) in patients treated with ACURATE neo, a supra-annular self-expanding bioprosthetic aortic valve.

Background: While clinical outcomes with TAVI are comparable with those achieved with surgery, residual aortic regurgitation (AR) and paravalvular leak (PVL) are common complications. The ACURATE neo valve has a pericardial sealing skirt designed to minimize PVL.

Methods: The primary endpoint was the rate of total AR over time, as assessed by a core echocardiographic laboratory. The study enrolled 500 patients (mean age: 81.8 ± 5.1 years; 61% female; mean baseline STS score: $6.0 \pm 4.5\%$) from 22 centers in Europe and Canada; 498 patients were treated with ACURATE neo.

Results: The rate of \geq moderate AR was 4.6% at discharge and 3.1% at 12 months; the rate of \geq moderate PVL was 4.6% at discharge and 2.6% at 12 months. Paired analyses showed significant improvement in overall PVL between discharge and 12 months ($P < 0.001$); 64.6% of patients had no change in PVL grade, 24.9% improved, and 10.5% worsened. Patients also exhibited significant improvement in transvalvular gradient ($P < 0.001$) and effective orifice area ($P=0.01$). The mortality rate was 2.2% at 30 days and 11.3% at 12 months. The permanent pacemaker implantation (PPI) rate was 10.2% at 30 days and 12.2% at 12 months.

Conclusions: Results from PROGRESS PVL support the sustained safety and performance of TAVI with the ACURATE neo valve, showing excellent valve hemodynamics, good clinical outcomes, and significant interindividual improvement in PVL from discharge to 12-month follow-up.

Transfemoral TAVI using the self-expanding ACURATE neo prosthesis: one-year outcomes of the multicentre "CE-approval cohort".

Möllmann H, Walther T, Siqueira D et al., EuroIntervention. 2017 Oct 13;13(9): e1040-e1046. doi: 10.4244/EIJ-D-17-00187. PMID: 28804056

Abstract:

Aims: The aim of this study was to assess the safety and performance of the ACURATE neo transcatheter heart valve and its transfemoral delivery system.



Methods and results: The prospective, multicentre "CE-approval cohort" consists of a prospective series of the first 89 patients implanted with the ACURATE neo prosthesis. The primary endpoint was all-cause mortality at 30 days. Mean patient age was 83.7±4.4 years and logistic EuroSCORE I was 26.5±7.7%. Procedural success was obtained in 84 patients (94.4%). At 30 days, three patients had died, and two major strokes and one reintervention for a ventricular septal defect occurred, leading to a major adverse cardiac and cerebrovascular event (MACCE) rate of 6.7%. Eight patients (10.3%) received a permanent pacemaker. At one year, 20 patients (22.5%) had died and the MACCE rate was 27%. Effective orifice area was 1.76±0.34 cm², and mean gradient 7.5±2.8 mmHg. Only three patients (4.5%) had moderate paravalvular regurgitation. NYHA Class III/IV was present in 94.4% of patients at baseline, in 9.9% at 30 days and in 4.5% at one year post procedure.

Conclusions: This first-in-human experience with a novel self-expanding heart valve showed low rates of procedural mortality, major stroke and pacemaker implantation, and good performance outcomes.

First results from the ACURATE Prime XL human feasibility study.

Gooley R, Murdoch D, Modolo R et al., *Cardiovasc Revasc Med.* 2023 Jun 30:S1553-8389(23)00673-5. doi: 10.1016/j.carrev.2023.06.028. Epub ahead of print. PMID: 37429793.

Abstract:

Background/purpose: This prospective, open-label, single-arm study evaluated transcatheter aortic valve replacement (TAVR) in patients with severe aortic stenosis with ACURATE Prime XL, an iteration of the ACURATE neo2 device designed with improved radial force and adaptations for compatibility with a larger annulus diameter (≥ 26.5 mm and ≤ 29 mm based on pre-procedure diagnostic imaging).

Methods: The composite primary device success endpoint was based on Valve Academic Research Consortium (VARC)-2 criteria. The primary safety endpoint was a composite of all-cause mortality and all stroke at 30 days. Aortic valve (AV) performance, including mean AV gradient, AV area, and grade of paravalvular leak (PVL), was assessed by an independent core laboratory.

Results: 13 male patients were enrolled at 3 Australian centers (mean age: 83.1 years; 10 of 13 were considered high/extreme operative risk). 61.5 % of patients met the primary device success endpoint. At 30 days, no patients experienced death or stroke; one patient received a permanent pacemaker. Mean AV gradient improved from baseline (42.7 ± 11.0 mmHg) to discharge (7.7 ± 2.5 mmHg) and through 30 days (7.2 ± 2.3 mmHg). Mean AV area was 0.8 ± 0.1 cm² at baseline, 1.9 ± 0.3 cm² at discharge, and 1.7 ± 0.3 cm² at 30 days. Per core-laboratory adjudication, no patient had moderate or severe PVL at 30 days; 91.7 % had no/trace PVL and 8.3 % had mild PVL.

Conclusions and relevance: In this first-in-human feasibility study of the ACURATE Prime XL valve, there were no safety concerns, and no deaths or strokes within 30 days. Valve hemodynamics were favorable, and no patient had >mild PVL.



INVESTIGATOR SPONSORED RESEARCH

The Early neo2 Registry: Transcatheter Aortic Valve Implantation With ACURATE neo2 in a European Population.

Ruck A, Kim WK, Abdel-Wahab M et al., J Am Heart Assoc. 2023 Aug;12(15):e029464. doi: 10.1161/JAHA.122.029464. Epub 2023 Jul 25. PMID: 37489732; PMCID: PMC10493001.

Paravalvular Leak After Implantation of the ACURATE neo and the ACURATE neo2 Transcatheter Heart Valve.

Stefan Toggweiler, Miriam Brinkert, Mathias Wolfrum, et al., Open Heart. 2020 Nov;7(2):e001391. doi: 10.1136/openhrt-2020-001391. PMID: 33243930; PMCID: PMC7692991.

Safety and efficacy of a self-expanding versus a balloon-expandable bioprosthesis for transcatheter aortic valve replacement in patients with symptomatic severe aortic stenosis: a randomised non-inferiority trial. (SCOPE I)

Lanz J, Kim W-K, Walther T et al., Lancet 2019 DOI:https://doi.org/10.1016/S0140-6736(19)32220-2.

Final 3-Year Outcomes of a Randomized Trial Comparing a Self-Expanding to a Balloon-Expandable Transcatheter Aortic Valve.

Lanz J, Mollmann H, Kim WK et al., Circ Cardiovasc Interv. 2023 Jul;16(7):e012873. doi: 10.1161/CIRCINTERVENTIONS.123.012873. Epub 2023 Jul 7. PMID: 37417229.

Comparison of Self-Expanding Bioprostheses for Transcatheter Aortic Valve Replacement in Patients With Symptomatic Severe Aortic Stenosis: SCOPE 2 Randomized Clinical Trial.

Tamburino C, Bleiziffer S, Thiele H et al., Circulation. 2020 Dec 22;142(25):2431-2442. doi: 10.1161/CIRCULATIONAHA.120.051547. Epub 2020 Oct 15. PMID: 33054367.

Permanent pacemaker implantation and left bundle branch block with self-expanding valves - a SCOPE 2 subanalysis.

Pellegrini C, Garot P, Morice MC et al., EuroIntervention. 2022 Sep 21:EIJ-D-22-00558. doi: 10.4244/EIJ-D-22-00558. Epub ahead of print. PMID: 36128956.

Transcatheter aortic valve replacement using the iSleeve expandable sheath in small femoral arteries.

Glaser N, O'Sullivan CJ, Saleh N et al., Open Heart. 2021 Oct;8(2):e001703. doi: 10.1136/openhrt-2021-001703. PMID: 34642241; PMCID: PMC8513271.



Accurate commissural alignment during ACURATE neo TAVI procedure. Proof of concept.

Redondo A, Valencia-Serrano F, Santos-Martínez S et al., Rev Esp Cardiol (Engl Ed). 2021 Mar 26:S1885-5857(21)00068-2. English, Spanish. doi: 10.1016/j.rec.2021.02.004. Epub ahead of print. PMID: 33781722.

EXTERNALLY DRIVEN EVIDENCE (BY TOPIC)

HEAD-TO-HEAD COMPARISONS

ACURATE neo2 versus SAPIEN 3 Ultra for transcatheter aortic valve implantation.

Pellegrini C, Rheude T, Renker M et al., EuroIntervention. 2023 Jan 23;18(12):987-995. doi: 10.4244/EIJ-D-22-00164. PMID: 36250307; PMCID: PMC9853033.

Comparison of transcatheter aortic valve replacement with the ACURATE neo2 versus Evolut PRO/PRO+ devices.

Baggio S, Pagnesi M, Kim WK et al., EuroIntervention. 2022 Sep 12:EIJ-D-22-00498. doi: 10.4244/EIJ-D-22-00498. Epub ahead of print. PMID: 36093795.

Multicenter comparison of transcatheter aortic valve implantation with the self-expanding ACURATE neo2 versus Evolut PRO transcatheter heart valves.

Rheude T, Rellegri C, Landt M et al., Clin Res Cardiol. 2023 Apr 28. doi: 10.1007/s00392-023-02194-4. Epub ahead of print. PMID: 37115228.

Comparison of contemporary transcatheter heart valve prostheses: data from the German Aortic Valve Registry (GARY).

Ruldolph TK, Herrmann E, Bon Dimitria et al., Clin Res Cardiol. 2023 Jul 18. doi: 10.1007/s00392-023-02242-z. Epub ahead of print. PMID: 37462856.

Comparing the Safety and Effectiveness of Five Leading New-Generation Devices for Transcatheter Aortic Valve Implantation: Twelve-Month Results From the RISPEVA Study.

Corcione N, Morello A, Ferraro P et al., Registro Italiano GISE sull'Impiantodi Valvola Aortica Percutanea (RISPEVA) Study Investigators. J Invasive Cardiol. 2021 May;33(5):E320-E329. Epub 2021 Mar 19. PMID: 33739300.



Comparative one-month safety and effectiveness of five leading new-generation devices for transcatheter aortic valve implantation.

Giordano A, Corcione N, Ferraro et al., Registro Italiano GISE sull'impianto di Valvola Aortica Percutanea (RISPEVA) Study Investigators. *Sci Rep.* 2019 Nov 19;9(1):17098. doi: 10.1038/s41598-019-53081-w. PMID: 31745198

Transcatheter Aortic Valve Replacement With Next-Generation Self-Expanding Devices: A Multicenter, Retrospective, Propensity-Matched Comparison of Evolut PRO Versus Acurate neo Transcatheter Heart Valves.

Pagnesi M, Kim WK, Conradi L et al., *JACC Cardiovasc Interv.* 2019 Mar 11;12(5):433-443. doi: 10.1016/j.jcin.2018.11.036

Comparison of Self-Expanding RDV Perceval S versus TAVI ACURATE neo/TF.

Gerfer S, Mauri V, Kuhn E et al., *Thorac Cardiovasc Surg.* 2021 Aug;69(5):420-427. doi: 10.1055/s-0040-1722692. Epub 2021 Mar 24. PMID: 33761569.

Comparison of latest generation transfemoral self-expandable and balloon-expandable transcatheter heart valves.

Schaefer A, Linder M, Seiffert M et al., *Interact Cardiovasc Thorac Surg.* 2017 Dec 1;25(6):905-911. doi: 10.1093/icvts/ivx194. PMID: 28655156

Multicenter Comparison of Novel Self-Expanding Versus Balloon-Expandable Transcatheter Heart Valves.

Husser O, Kim WK, Pellegrini C et al., *JACC Cardiovasc Interv.* 2017 Oct 23;10(20):2078-2087. doi: 10.1016/j.jcin.2017.06.026. PMID: 29050625

Multi-Center Comparison of Two Self-Expanding Transcatheter Heart Valves: A Propensity Matched Analysis.

Blumenstein J, Eckel C, Husser O et al., *J Clin Med.* 2022 Jul 21;11(14):4228. doi: 10.3390/jcm11144228. PMID: 35887990; PMCID: PMC9318122.

Transcatheter aortic valve implantation with different self-expanding devices-a propensity score-matched multicenter comparison.

Wienemann H, Hof A, Ludwig S et al., *J Cardiol.* 2023 Mar-Apr;70:1-9. doi: 10.1016/j.hjc.2022.12.006. Epub 2022 Dec 17. PMID: 36538975.



Real-World Multiple Comparison of Transcatheter Aortic Valves: Insights From the Multicenter OBSERVANT II Study.

Costa G, Barbanti M, Rosato S, et al., Cardiovascular Interventions. 2022;15 2022<https://doi.org/10.1161/CIRCINTERVENTIONS.122.012294>Circulation:

Comparison of latest generation supra-annular and intra-annular self-expanding transcatheter heart valves.

Pellegrini C, Rheude T, Michel J et al., J Thorac Dis. 2020 Nov;12(11):6769-6779. doi: 10.21037/jtd-20-1700.PMID: 33282378

LENGTH OF STAY

Safety of Next-Day Discharge After Transfemoral Transcatheter Aortic Valve Replacement with a Self-Expandable Versus Balloon-Expandable Valve Prosthesis.

Moriyama N, Vento A, Laine M et al., Circ Cardiovasc Interv. 2019 Jun;12(6):e007756. doi:10.1161/CIRCINTERVENTIONS.118.007756. Epub 2019 Jun 6. PMID: 31167602

CONDUCTION DISORDERS

Conduction Disturbance After Transcatheter Aortic Valve Implantation With Self- or Balloon-Expandable Valve According to the Implantation Depth.

Miyashita H, Moriyama N, Sugiyama Y et al., Am J Cardiol. 2023 Sep 15;203:17-22. doi: 10.1016/j.amjcard.2023.05.025. Epub 2023 Jul 19. PMID: 37478637.

Transcatheter Valve SELECTION in Patients With Right Bundle Branch Block and Impact on Pacemaker Implantations.

Husser O, Pellegrini C, Kim WK et al., JACC Cardiovasc Interv. 2019 Sep 23;12(18):1781-1793. doi: 10.1016/j.jcin.2019.05.055.

Relevance of New Conduction Disorders After Implantation of the ACURATE Neo Transcatheter Heart Valve in the Aortic Valve Position.

Brinkert M, Wolfrum M, Moccetti F et al., Am J Cardiol. 2020 Mar 1;125(5):783-787. doi: 10.1016/j.amjcard.2019.11.036. Epub 2019 Dec 14. PMID: 31898969



Predictors of permanent pacemaker implantation after ACURATE neo transcatheter heart valve implantation.

Kim WK, Möllmann H, Walther T et al., CW.Pacing Clin Electrophysiol. 2020 Dec 29. doi: 10.1111/pace.14155. Online ahead of print. PMID: 33373045

Paravalvular regurgitation and permanent pacemaker implantation after TAVR with a next-generation self-expanding device.

Mauri V, Deuschl F, Frohn T et al., Clin Res Cardiol. 2018 Aug;107(8):688-697. doi: 10.1007/s00392-018-1235-1. Epub 2018 Apr 17. PMID: 29667013

Predictors of Need for Permanent Pacemaker Implantation and Conduction Abnormalities with a Novel Self-expanding Transcatheter Heart Valve.

Pellegrini C, Husser O, Kim WK et al., Rev Esp Cardiol (Engl Ed). 2019 Feb;72(2):145-153. doi: 10.1016/j.rec.2018.01.011. Epub 2018 Mar 16. PMID: 29551701 English, Spanish.

Very low pacemaker rate following ACURATE neo transcatheter heart valve implantation.

Toggweiler S, Nissen H, Mogensen B et al., EuroIntervention. 2017 Dec 20;13(11):1273-1280. doi: 10.4244/EIJ-D-17-00252. PMID: 28870877

Impact of Interventricular membranous septum length on pacemaker need with different Transcatheter aortic valve implantation systems.

Thijmen W.H, Maarten P.vanW, J F. O et al., International Journal of Cardiology March 02, 2021DOI:https://doi.org/10.1016/j.ijcard.2021.02.080

Development of atrioventricular and intraventricular conduction disturbances in patients undergoing transcatheter aortic valve replacement with new generation self-expanding valves: A real world multicenter analysis.

Castro-Mejía AF, Amat-Santos I, Ortega-Armas ME et al., Int J Cardiol. 2022 Sep 1;362:128-136. doi: 10.1016/j.ijcard.2022.05.014. Epub 2022 May 10. PMID: 35550389.



SMALL ANNULI

Implantation of contemporary transcatheter aortic valves in small aortic annuli: the international multicentre TAVI-SMALL 2 registry.

Leone PP, Regazzoli D, Pagnesi M et al., EuroIntervention. 2023 Jun 19;19(3):256-266. doi: 10.4244/EIJ-D-22-00843. PMID: 36950893; PMCID: PMC10266383.

Short-Term Outcome and Hemodynamic Performance of Next-Generation Self-Expanding Versus Balloon-Expandable Transcatheter Aortic Valves in Patients With Small Aortic Annulus: A Multicenter Propensity-Matched Comparison.

Mauri V, Kim WK, Abumayyaleh M, et al., Circ Cardiovasc Interv. 2017;10(10)

Procedural Outcomes of a Self-Expanding Transcatheter Heart Valve in Small Annuli.

Eckel C, Sötemann D, Kim WK, et al., J Clin Med. 2022 Sep 9;11(18):5313. doi: 10.3390/jcm11185313. PMID: 36142960; PMCID: PMC9502952.

Transcatheter aortic valve implantation in patients with a small aortic annulus: performance of supra-, intra- and infra-annular transcatheter heart valves.

Voigtländer L, Kim WK, Mauri V, et al., Clin Res Cardiol. 2021 Dec;110(12):1957-1966. doi: 10.1007/s00392-021-01918-8. Epub 2021 Aug 13. PMID: 34387736; PMCID: PMC8639544.

Transcatheter Self-Expandable Valve Implantation for Aortic Stenosis in Small Aortic Annuli: The TAVI-SMALL Registry.

Regazzoli D, Chiarito M, Cannata et al., TAVI-SMALL Investigators. JACC Cardiovasc Interv. 2020 Jan 27;13(2):196-206. doi: 10.1016/j.jcin.2019.08.041. Epub 2019 Dec 25. PMID: 31883714

CORONARY OBSTRUCTION

Risk Assessment Of Coronary Obstruction For The ACURATE Transcatheter Heart Valve.

Kim WK, Eckel C, Charitos EI, et al., Am J Cardiol. 2023 Nov 30;213:45-46. doi: 10.1016/j.amjcard.2023.11.058. Epub ahead of print. PMID: 38042270.



CORONARY ACCESS / COMMISSURAL ALIGNMENT

Patient-Specific Implantation Technique to Obtain Neo-Commissural Alignment With Self-Expanding Transcatheter Aortic Valves.

Bieliauskas G, Wong I, Bajoras V, Wang X et al., JACC Cardiovasc Interv. 2021 Oct 11;14(19):2097-2108. doi: 10.1016/j.jcin.2021.06.033. Epub 2021 Sep 15. PMID: 34538602.

Commissural Alignment With ACURATE neo2 Valve in an Unselected Population.

Meduri C, Rück A, Linder R et al., J Am Coll Cardiol Intv. 2023 Mar, 16 (6) 670–677. <https://doi.org/10.1016/j.jcin.2023.01.018>

A patient-specific algorithm to achieve commissural alignment with Acurate Neo: The sextant technique.

De Marco F, Casenghi M, Spagnolo P et al., Catheter Cardiovasc Interv. 2021 Nov 15;98(6):E847-E854. doi: 10.1002/ccd.29737. Epub 2021 May 7. PMID: 33960624.

Coronary Access After Transcatheter Aortic Valve Replacement With Commissural Alignment: The ALIGN-ACCESS Study.

Tarantini G, Nai Fovino L, Scotti A, et al., Circ Cardiovasc Interv. 2022 Feb;15(2):e011045. doi: 10.1161/CIRCINTERVENTIONS.121.011045. Epub 2022 Feb 15. PMID: 35167332.

Coronary Cannulation After Transcatheter Aortic Valve Replacement: The RE-ACCESS Study.

Barbanti M, Costa G, Picci A, et al., JACC Cardiovasc Interv. 2020 Nov 9;13(21):2542-2555. doi: 10.1016/j.jcin.2020.07.006. Epub 2020 Oct 14. JACC Cardiovasc Interv. 2020. PMID: 33069648

Coronary Angiography After Transcatheter Aortic Valve Replacement (TAVR) to Evaluate the Risk of Coronary Access Impairment After TAVR-in-TAVR.

Nai Fovino L, Scotti A, Massussi M et al., J Am Heart Assoc. 2020 Jul 7;9(13):e016446. doi: 10.1161/JAHA.120.016446. Epub 2020 Jun 24. PMID: 32578484



Neo-Commissural Alignment Technique for Transcatheter Aortic Valve Replacement Using the ACURATE Neo Valve.

Abdel-Wahab M, Kitamura M, Fitzgerald SJ et al., Circ Cardiovasc Interv. 2022 Jul 1:101161CIRCINTERVENTIONS122011993. doi: 10.1161/CIRCINTERVENTIONS.122.011993. Epub ahead of print. PMID: 35775417.

Commissural Alignment of the ACURATE neo Valve in Transcatheter Aortic Valve Replacement.

Kitamura M, Wilde J, Gohmann R et al., JACC Cardiovasc Interv. 2021 Aug 9;14(15):1740-1742. doi: 10.1016/j.jcin.2021.05.037. PMID: 34353607.

Alignment of Transcatheter Aortic-Valve Neo-Commissures (ALIGN TAVR): Impact on Final Valve Orientation and Coronary Artery Overlap.

Tang GHL, Zaid S, Fuchs A et al., JACC Cardiovasc Interv. 2020 May 11;13(9):1030-1042. doi: 10.1016/j.jcin.2020.02.005. Epub 2020 Mar 16. PMID: 32192985

Definition of transcatheter heart Valve orientation in biCuspld aortic valve: The DA VINCI pilot study.

Tarantini G, Fabris T, Nai Fovino L et al., Front Cardiovasc Med. 2022 Dec 12;9:1056496. doi: 10.3389/fcvm.2022.1056496. PMID: 36578836; PMCID: PMC9790995.

PARAVALVULAR LEAKAGE AND AORTIC REGURGITATION

Paravalvular Aortic Regurgitation Severity Assessed by Quantitative Aortography: ACURATE neo 2 versus ACURATE neo Transcatheter Aortic Valve Implantation.

Rück A, Kim WK, Kawashima H et al., J Clin Med. 2021 Oct 9;10(20):4627. doi: 10.3390/jcm10204627. PMID: 34682750; PMCID: PMC8539505.

Haemodynamic performance and clinical outcomes of transcatheter aortic valve replacement with the self-expanding ACURATE neo2.

Scotti A, Pagnesi M, Kim WK et al., EuroIntervention. 2022 Jun 8:EIJ-D-22-00289. doi: 10.4244/EIJ-D-22-00289. Epub ahead of print. PMID: 35678222.





Transcatheter Aortic Valve Replacement with Self-Expanding ACURATE neo2: Postprocedural Hemodynamic and Short-Term Clinical Outcomes.

Buono A, Gorla R, Ielasi A et al., JACC Cardiovasc Interv. 2022 Jun 13;15(11):1101-1110. doi: 10.1016/j.jcin.2022.02.027. Epub 2022 May 17. PMID: 35595675.



Comparison of the Acurate Neo Vs Neo2 Transcatheter Heart Valves.

Kim WK, Eckel C, Renker M et al., J Invasive Cardiol. 2022 Nov;34(11):E804-E810. Epub 2022 Oct 21. PMID: 36269158.



Predictors of Prosthetic Valve Regurgitation After Transcatheter Aortic Valve Implantation with ACURATE neo in the SCOPE I Trial.

Kim WK, Walther T, Burgdorf C et al., SCOPE I Investigators. JACC Cardiovasc Imaging. 2022 Feb;15(2):367-369. doi: 10.1016/j.jcmg.2021.08.019. Epub 2021 Oct 13. PMID: 34656481.



Effectiveness and Safety of the ACURATE Neo Prosthesis in 1,000 Patients With Aortic Stenosis.

Kim WK, Möllmann H, Liebetau C, et al., Am J Cardiol. 2020 Sep 15;131:12-16. doi: 10.1016/j.amjcard.2020.05.044. Epub 2020 Jun 7. PMID: 32622500



The ACURATE neo Transcatheter Heart Valve: A Comprehensive Analysis of Predictors of Procedural Outcome.

Kim WK, Möllmann H, Liebetau C, Renker M et al., JACC Cardiovasc Interv. 2018 Sep 10;11(17):1721-1729. doi: 10.1016/j.jcin.2018.04.039. Epub 2018 May 23. PMID: 29803694



Impact of device landing zone calcification patterns on paravalvular regurgitation after transcatheter aortic valve replacement with different next-generation devices.

Mauri V, Frohn T, Deuschl F et al., Open Heart. 2020 May;7(1): e001164. doi: 10.1136/openhrt-2019-001164. PMID: 32393655



Quantitative Assessment of Acute Regurgitation Following TAVR: A Multicenter Pooled Analysis of 2,258 Valves.

Modolo R, Chang CC, Abdelghani M et al., JACC Cardiovasc Interv. 2020 Jun 8;13(11):1303-1311. doi: 10.1016/j.jcin.2020.03.002. Epub 2020 Mar 16. PMID: 32499020



PRE-DILATATION AND POST-DILATATION

Efficacy and Safety of Postdilation for a Self-Expanding Transcatheter Heart Valve.

Kim WK, Renker M, Nef H, Möllmann H et al., J Invasive Cardiol. 2022 Jun;34(6): E448-E454. Epub 2022 May 19. PMID: 35593542.

Impact of Predilatation Prior to Transcatheter Aortic Valve Implantation with the Self-Expanding ACURATE neo Device (from the Multicenter NEOPRO Registry).

Pagnesi M, Kim WK, Conradi L et al., Am J Cardiol. 2020 May 1;125(9):1369-1377. doi: 10.1016/j.amjcard.2020.02.003. Epub 2020 Feb 8. PMID: 32098656

Impact of Predilation Before Transcatheter Aortic Valve Implantation with New-Generation Devices.

Giordano A, Corcione N, Ferraro P et al., RISPEVA (Registro Italiano GISE sull'impianto di Valvola Aortica Percutanea) Study Investigators. Cardiovasc Revasc Med. 2019 Dec;20(12):1096-1099. doi: 10.1016/j.carrev.2019.01.017. Epub 2019 Jan 23. PMID: 30711476

Simplifying transfemoral ACURATE neo implantation using the TrueFlow nonocclusive balloon catheter.

Toggweiler S, Loretz L, Brinkert M et al., Catheter Cardiovasc Interv. 2020 Nov;96(6):E640-E645. doi: 10.1002/ccd.28741. Epub 2020 Jan 23. PMID: 31971346

Transfemoral aortic valve implantation using a self-expanding transcatheter heart valve without pre-dilation.

Kim WK, Liebetrau C, Renker M, et al., Int J Cardiol. 2017 Sep 15;243:156-160. doi: 10.1016/j.ijcard.2017.05.092. Epub 2017 May 27. PMID: 28579166

TAVI HEMODYNAMICS AND DURABILITY

Early Hemodynamic Outcomes in Self-Expandable Valves: Comparison of ACURATE Neo Versus ACURATE Neo2.

Miyashita H, Moriyama N, Laine M et al., Cardiovasc Revasc Med. 2022 Aug 28:S1553-8389(22)00758-8. doi: 10.1016/j.carrev.2022.08.029. Epub ahead of print. PMID: 36114124.



Mid- to long-term clinical and echocardiographic outcomes after transcatheter aortic valve replacement with a new-generation, self-expandable system.

Siqueira DA, Simonato M, Ramos AA et al., Catheter Cardiovasc Interv. 2020 May 23. doi: 10.1002/ccd.28999. Online ahead of print. PMID: 32445607.

Self-expandable transcatheter heart valves for aortic stenosis. Short-term outcome and matched hemodynamic performance.

Vera Vera S, Nombela-Franco L, Santos-Martínez S et al., Rev Esp Cardiol (Engl Ed). 2020 Nov 3:51885-5857(20)30422-9. doi: 10.1016/j.rec.2020.09.014.

Hemodynamics of self-expanding versus balloon-expandable transcatheter heart valves in relation to native aortic annulus anatomy.

Schmidt, S., Fortmeier, V., Ludwig, S et al., Clin Res Cardiol 111, 1336–1347 (2022). <https://doi.org/10.1007/s00392-022-02046-7>.

HORIZONTAL AORTA

Horizontal Aorta in Transcatheter Self-Expanding Valves: Insights from the HORSE International Multicentre Registry.

Gallo F, Gallone G, Kim WK et al., Circ Cardiovasc Interv. 2021 Sep;14(9):e010641. doi: 10.1161/CIRCINTERVENTIONS.121.010641. Epub 2021 Aug 30. PMID: 34455799.

Impact of aortic angle on transcatheter aortic valve implantation outcome with Evolut-R, Portico, and Acurate-NEO.

Gorla R, De Marco F, Garatti A et al., Catheter Cardiovasc Interv. 2020 May 13. doi: 10.1002/ccd.28957. Online ahead of print. PMID: 32400068

EX VIVO ANALYSIS

A bench study of balloon-expandable valves for the treatment of self-expanding valve failure.

Akodad M, Meier D, Sellers S et al., EuroIntervention 2023;18-online publish-ahead-of-print January 2023 DOI: 10.4244/EIJ-D-22-00769



Coronary access after ACURATE neo2 implantation for valve-in-valve TAVR: Insights from ex vivo simulations.

Khokhar AA, Giannini F, Zlahoda-Huzior A et al., Catheter Cardiovasc Interv. 2022 Sep 18. doi: 10.1002/ccd.30367. Epub ahead of print. PMID: 36116020.

Coronary access techniques following ACURATE neo2 implantation in surgical bioprosthesis.

Khokhar AA, Ponticelli F, Zlahoda-Huzior A et al., EuroIntervention. 2022 Jun 16:EIJ-D-22-00149. doi: 10.4244/EIJ-D-22-00149. Epub ahead of print. PMID: 35708304.

Long-Term Durability of the Next-Generation ACURATE neo2 Transcatheter Heart Valve: Insights From Bench Testing to 25 Years.

Sathananthan J, Landes U, Flaction L et al., JACC Cardiovasc Interv. 2021 Mar 8;14(5):586-588. doi: 10.1016/j.jcin.2020.10.030. Epub 2021 Mar 1. PMID: 33663791.

Coronary access following ACURATE neo implantation for transcatheter aortic valve-in-valve implantation: Ex vivo analysis in patient-specific anatomies.

Khokhar AA, Ponticelli F, Zlahoda-Huzior A et al., Front Cardiovasc Med. 2022 Sep 14;9:902564. doi: 10.3389/fcvm.2022.902564. PMID: 36187005; PMCID: PMC9515364.

Computed tomography analysis of coronary ostia location following valve-in-valve transcatheter aortic valve replacement with the ACURATE neo valve: Implications for coronary access.

Khokhar AA, Laricchia A, Ponticelli F et al., Catheter Cardiovasc Interv. 2021 Sep;98(3):595-604. doi: 10.1002/ccd.29503. Epub 2021 Feb 14. PMID: 33586278.

A novel 3D imaging approach to evaluate coronary access following ACURATE neo implantation to treat a degenerated surgical bioprosthesis.

Khokhar AA, Ponticelli F, Zlahoda-Huzior A et al., EuroIntervention. 2021 Aug 3:EIJ-D-21-00634. doi: 10.4244/EIJ-D-21-00634. Epub ahead of print. PMID: 34338640.

Coronary Access Following Redo TAVR: Impact of THV Design, Implant Technique, and Cell Misalignment.

David Meier, Mariama Akodad, Uri Landes, et al., JACC: Cardiovascular Interventions, 13 July 2022.



Impact of Commissural Misalignment on Hydrodynamic Function Following Valve-in-Valve Intervention With the ACURATE neo.

David Meier, Mariama Akodad, Andrew G et al., JACC: Cardiovascular Interventions, 13 July 2022.

Repeat transcatheter aortic valve implantation and implications for transcatheter heart valve performance: insights from bench testing.

Sathananthan J, Fraser R, Landes U et al., EuroIntervention. 2021 Nov 19;17(10):856-864. doi: 10.4244/EIJ-D-20-00697. PMID: 33495142.

A bench test study of bioprosthetic valve fracture performed before versus after transcatheter valve-in-valve intervention.

Sathananthan J, Fraser R, Hatoum H et al. EuroIntervention. 2020 Mar 20;15(16):1409-1416. doi: 10.4244/EIJ-D-19-00939. PMID: 31854302

Impact of implant depth on hydrodynamic function with the ACURATE neo transcatheter heart valve following valve-in-valve transcatheter aortic valve replacement in Mitroflow bioprosthetic valves: an ex vivo bench study.

Sathananthan J, Sellers SL, Fraser R et al., EuroIntervention. 2019 May 20;15(1):78-87. doi: 10.4244/EIJ-D-18-00947. PMID: 30479308

Implications of hydrodynamic testing to guide sizing of self-expanding transcatheter heart valves for valve-in-valve procedures.

Sathananthan J, Hensey M, Fraser R et al., Catheter Cardiovasc Interv. 2020 Sep 1;96(3):E332-E340. doi: 10.1002/ccd.28537. Epub 2019 Oct 24. PMID: 31647178

Valve-in-Valve Transcatheter Aortic Valve Replacement and Bioprosthetic Valve Fracture Comparing Different Transcatheter Heart Valve Designs: An Ex Viv Bench Study.

Sathananthan J, Sellers S, Barlow AM et al., JACC Cardiovasc Interv. 2019 Jan 14;12(1):65-75. doi: 10.1016/j.jcin.2018.10.043. PMID: 30621980

Neosinus and Sinus Flow After Self-Expanding and Balloon-Expandable Transcatheter Aortic Valve Replacement.

Hatoum H, Gooden SCM, Sathananthan J et al., JACC Cardiovasc Interv. 2021 Dec 27;14(24):2657-2666. doi: 10.1016/j.jcin.2021.09.013. Epub 2021 Nov 24. PMID: 34838462



Radial Force: An Underestimated Parameter in Oversizing Transcatheter Aortic Valve Replacement Prostheses: In Vitro Analysis with Five Commercialized Valves.

Egron S, Fujita B, Gullón L et al., ASAIO J. 2018 Jul/Aug;64(4):536-543. doi: 10.1097/MAT.0000000000000659. PMID: 28885378

Coronary Access and PCI After Chimney Stenting During Self-Expanding Transcatheter Aortic Valve Replacement in Surgical Bioprosthesis:

Francesco Giannini F, Khokhar A et al., JACC Cardiovasc Interv. 2022 Aug 8;15(15):1585-1589. doi: 10.1016/j.jcin.2022.05.015. Epub 2022 Jun 15. PMID: 35717397.

GENERAL TOPICS

Sizing Considerations of the ACURATE Transcatheter Heart Valve Platform: Insights From 2,865 Patients.

Charitos EI, Renker M, Eckel C, et al., JACC Cardiovasc Interv. 2022 Dec 12;15(23):2457-2459. doi: 10.1016/j.jcin.2022.08.038.

Single Access for Transfemoral Transcatheter Aortic Valve Implantation With the Acurate neo/neo 2 Self-Expanding Valve.

Bagur R, Chu MWA, Ordoñez S, et al., Can J Cardiol. 2022 Oct 31:S0828-282X(22)00992-8. doi: 10.1016/j.cjca.2022.10.026. Epub ahead of print. PMID: 36328177.

Transcatheter Aortic Valve Replacement Through a Single Femoral Access: A Multicenter Experience.

Toggweiler S, Bagur R, Agatiello CR, et al. J Invasive Cardiol. 2022 Oct;34(10):E739-E742. Epub 2022 Sep 16. PMID: 36121924.

Micro-dislodgement of a self-expanding transcatheter heart valve: Incidence, predictors, and outcomes.

Kim WK, Renker M, Dörr O et al., Int J Cardiol. 2022 Jul 1;358:77-82. doi: 10.1016/j.ijcard.2022.04.031. Epub 2022 Apr 15. PMID: 35436560.



Clinical and hemodynamic performance of the ACURATE neo bioprosthesis: Looking far and beyond randomized trials.

Franzone A, Esposito G. Catheter Cardiovasc Interv. 2021 Jan 1;97(1):175-176. doi:10.1002/ccd.29447. PMID: 33460273.

The impact of learning-curve-experience on transcatheter aortic valve replacement outcomes: Insights from the United Kingdom and Ireland all-comers second-generation ACURATE neo™ transcatheter aortic heart valve registry.

Gunarathne A, Tan J, Varcoe R et al., Catheter Cardiovasc Interv. 2021 Nov 18. doi: 10.1002/ccd.30000. Epub ahead of print. PMID: 34792842.

Procedural Outcomes of a Self-Expanding Transcatheter Heart Valve in Patients with Porcelain Aorta.

Eckel C, Blumenstein J, Grothusen C et al., J Clin Med. 2023 Jan 26;12(3):945. doi: 10.3390/jcm12030945. PMID: 36769593; PMCID: PMC9917710.

One-year clinical outcome with a novel self-expanding transcatheter heart valve.

Pellegrini C, Rheude T, Trenkwalder T et al., Catheter Cardiovasc Interv. 2019 Nov 15;94(6):783-792. doi: 10.1002/ccd.28144. Epub 2019 Feb 24. PMID: 30801906

Transfemoral implantation of the ACURATE neo prosthesis using a low-profile expandable introducer system: A multicenter registry.

Kim WK, Brinkert M, Mangner N et al. Int J Cardiol. 2019 Apr 15;281:76-81. doi: 10.1016/j.ijcard.2019.01.081. Epub 2019 Jan 25. PMID: 30711266

Accuracy of device landing zone calcium volume measurement with contrast-enhanced multidetector computed tomography.

Kim WK, Renker M, Rolf A, et al., Int J Cardiol. 2018 Jul 15;263:171-176. doi: 10.1016/j.ijcard.2018.02.042. PMID: 29754916



OFF-LABEL USE

(SAFETY AND EFFECTIVNESS HAVE NOT BEEN ESTABLISHED IN THESE PATIENT POPULATIONS)

VALVE-IN-VALVE

Comparison of two self-expanding transcatheter heart valves for degenerated surgical bioprostheses: The AVENGER multicentre registry.

Kim WK, Seiffert M, Ruck A et al., EuroIntervention. 2023 Sep. DOI: 10.4244/EIJ-D-23-00779

Valve-in-Valve Implantation Using the ACURATE Neo in Degenerated Aortic Bioprostheses: An International Multicenter Analysis.

Holzamer A, Kim WK, Rück A et al., JACC Cardiovasc Interv. 2019 Nov 25;12(22):2309-2316. doi: 10.1016/j.jcin.2019.07.042. PMID: 31753302

PURE AORTIC INSUFFICIENCY / AORTIC REGURGITATION

Transcatheter treatment of native aortic valve regurgitation: Results from an international registry using the transfemoral ACURATE neo valve.

Purita PAM, Tahoces LS, Fraccaro C, et al., Int J Cardiol Heart Vasc. 2020 Feb 12;27:100480. Doi: 10.1016/j.ijcha.2020.100480. eCollection 2020 Apr. PMID: 32083165

Transfemoral Implantation of the Acurate neo for the Treatment of Aortic Regurgitation.

Toggweiler S, Cerillo AG, Kim WK, et al., J Invasive Cardiol. 2018 Sep;30(9):329-333. Epub 2018 Jul 15. PMID: 30012890

TAVI for Pure Non-calcified Aortic Regurgitation Using a Self-Expandable Transcatheter Heart Valve.

Schneeberger Y, Seiffert M, Schaefer A et al., Front Cardiovasc Med. 2022 Jan 25;8:743579. doi: 10.3389/fcvm.2021.743579. PMID: 35146003; PMCID: PMC8823702.



BICUSPID VALVE

Transcatheter aortic valve replacement in raphe-type bicuspid valves with the ACURATE neo2 according to the LIRA method.

Beneduce A, Russo F, Ghizzoni G, et al., *AsiaIntervention*. 2022 Oct 6;8(2):145-149. doi: 10.4244/AIJ-D-22-00006. PMID: 36483274; PMCID: PMC9706757.

Outcome of transcatheter aortic valve replacement in bicuspid aortic valve stenosis with new-generation devices.

Gorla R, Casenghi M, Finotello A et al., *Interact Cardiovasc Thorac Surg*. 2021 Jan 1;32(1):20-28. doi: 10.1093/icvts/ivaa231. PMID: 33201

Transcatheter aortic valve implantation using the ACURATE neo in bicuspid and tricuspid aortic valve stenosis: a propensity-matched analysis of a European experience.

Mangieri A, Chieffo A, Kim WK et al., *EuroIntervention*. 2018 Dec 20;14(12):e1269-e1275. doi: 10.4244/EIJ-D-18-00281. PMID: 30375336

ALTERNATIVE ACCESS ROUTE

Procedural success in transaxillary transcatheter aortic valve implantation according to type of transcatheter heart valve: results from the multicenter TAXI registry.

Schaefer A, Bhadra OL, Conradi L et al., *Clin Res Cardiol*. 2023 May 4. doi: 10.1007/s00392-023-02216-1. Epub ahead of print. PMID: 37138103.

Transaxillary transcatheter ACURATE neo aortic valve implantation - The TRANSAX multicenter study.

Amat-Santos IJ, Santos-Martínez S, Conradi L, et al., *Catheter Cardiovasc Interv*. 2020 Dec 14. doi: 10.1002/ccd.29423. Online ahead of print. PMID: 33315296

Direct percutaneous transaxillary implantation of a novel self-expandable transcatheter heart valve for aortic stenosis.

Deuschl F, Schofer N, Seiffert M et al., *Catheter Cardiovasc Interv*. 2017 Dec 1;90(7):1167-1174. doi: 10.1002/ccd.26986. Epub 2017 Mar 15. PMID: 28296065





Interventional Cardiology
300 Boston Scientific Way
Marlborough, MA 01752-1234
www.bostonscientific.com

*To order product or for more information
contact customer service at 1.888.272.1001.*

© 2024 Boston Scientific Corporation
or its affiliates. All rights reserved.

CAUTION: The ACURATE Prime Aortic Valve System is an investigational device.
Not available for sale. In Europe, ACURATE neo2 Aortic Valve System is CE-marked.
In the USA, ACURATE neo2 is an investigational device and restricted under federal
law to investigational use only. Not available for sale.

SH-1376306-AD