







4G USB Cellular Adapter Setup Guide

The 4G USB cellular adapter includes two (2) items:

- USB cellular adapter
- USB extension cable

Follow the steps below to setup your 4G cellular adapter:

<p>Step 1</p>	<p>If you still have your old USB cellular adapter connected, please remove it now.</p> <p><u>Unplug</u> the power supply cord from the back of your LATITUDE™ Communicator or from the electrical outlet.</p>	
<p>Step 2</p>	<p>Connect the new 4G USB cellular adapter to the supplied USB extension cable by inserting the metal end of the cellular adapter securely into the female end of the extension cable</p> <p>NOTE: There will be a slight gap between the adapter and the cable. Always keep the two parts connected together.</p>	
<p>Step 3</p>	<p>Connect the USB extension cable into either one of the USB ports on the back of the LATITUDE Communicator.</p>	
<p>Step 4</p>	<p><u>Reconnect</u> the power supply cord to the back of your LATITUDE Communicator and/or the electrical outlet. The Communicator will take about 15 seconds to turn on. Wait until you see the green LATITUDE indicator light before proceeding to the next step.</p> <p>NOTE: The USB cellular adapter LED light will begin to flash when properly connected and may flash or steadily display various colors (example: white, red, blue, orange).</p>	

	This is normal behavior and is not related to your implanted device.	
Step 5	To test that your new 4G cellular adapter is setup correctly, press and hold the white indented status button on the back of the LATITUDE Communicator for 3 seconds to complete the setup. DO NOT press the heart button on the front of the Communicator. The wave lights on your Communicator will flash green in sequence for several minutes on the right side of the Heart Button.	
Step 6	Once the connection is complete, the waves on your Communicator will turn solid green and a blue doctor icon will appear. Your new adapter is connected and no further action is needed. After two minutes, the lights will go off except for the green light under the word LATITUDE. The USB cellular adapter's LED light will remain lit and/or flashing various colors. This is normal behavior and is not related to your implanted device.	

Questions? Call Boston Scientific RhythmCARE™ Patient Services at 1-866-484-3268; Monday-Friday 8am-5pm CST.

Download the FREE MyLATITUDE™ Patient App (Recommended)

The MyLATITUDE Patient App is an optional resource for you to learn more about your LATITUDE Communicator and implanted device.

Download the MyLATITUDE Patient App for:

- The connection status of your LATITUDE Communicator
- Notification if your Communicator is not working and how to resolve the issue
- Information about your scheduled transmissions
- The status of your implanted device battery
- Communicator setup instructions
- The above content in Spanish



Get Started:

- 1** Open the App Store within your mobile device
- 2** Search for "MyLATITUDE"
- 3** Select the "MyLATITUDE Patient App"