



**THERASPHERE™**  
Y-90 Glass Microspheres

**CONTROL FOR  
DURABLE OUTCOMES**





# THERASPHERE™

## Y-90 Glass Microspheres

Designed for control to achieve durable patient outcomes, TheraSphere provides reproducible results for you and your referrers.



**Proven, well-tolerated therapy since 1999  
with more than 100,000 patient treatments globally**

**1999**

TheraSphere receives **FDA HDE approval**

**2011**

**Radiation segmentectomy:** high-dose radiation delivered to  $\leq 2$  hepatic segments<sup>1</sup>

**2016**

**PREMIERE:** TheraSphere significantly prolongs TTP compared with cTACE<sup>2</sup>

**2020**

**LEGACY** confirms neoadjuvant or standalone treatment in HCC<sup>3</sup>

**DOSISPHERE-01** personalized dosimetry approach improves OS<sup>4</sup>

**2021**

**1st and only Y-90 therapy PMA-approved by FDA for HCC**

**BCLC Guidelines: LEGACY study led to inclusion of TARE for very early and early-stage unresectable HCC<sup>5</sup>**

**TARGET** global real-world study confirms tumor absorbed dose is critical for predictable tumor response and OS in broad population.<sup>6</sup>

**2022**

**RASER:** Prospective study assessing radiation segmentectomy using TheraSphere<sup>7</sup>

**TRACE:** TheraSphere offers superior tumor control and survival compared with DEB-TACE\*<sup>8</sup>

**2024**

Ochsner tumor dose analysis demonstrated reproducible outcomes with high dose and high RPM<sup>9</sup>



**Targeted for control.**  
Unique, optimal design for tumor control.

**Maximize  
dose to tumor**




**Minimize  
dose to  
healthy tissue**

TheraSphere is uniquely engineered to have unmatched radiation per microsphere (RPM) to maximize repetitive and cumulative radiation exposure to tumor cells.

This allows you to confidently destroy tumor cells with proven ablative dosing and drive complete pathological necrosis (CPN).<sup>3,10</sup> The flexibility of personalized dosing serves your treatment intent across various clinical scenarios.<sup>11</sup>

**Targeted delivery with unmatched RPM helps achieve CPN:<sup>10</sup>**




TheraSphere maximizes lethal hits to tumor DNA, leading to tumor cell death<sup>9</sup>

# Durability leads to quality of life.

Tumor control yields durable outcomes and reproducible results.

## Durable outcomes for your patient



### LEGACY

**88%** best overall response rate using localized mRECIST, with **76.1%** of patients exhibiting a DoR  $\geq 6$  months<sup>3</sup>


### RASER

**90%** of patients achieved a sustained complete response with a median duration of **20.9 months (635 days)**<sup>7</sup>

## Reproducible results widely published



### Proven absorbed dose-response relationship with TheraSphere



Reproducible results backed by international, multi-disciplinary Dosimetry Steering Committee recommendations.<sup>11</sup>


# Quality of life creates long-term value.

Minimal side effects in an outpatient procedure.

## Durable outcomes are cost effective



Glass Y-90 treatment showed longer time to tumor progression than cTACE<sup>2</sup>



## Fewer procedures positively impact patient lives



TheraSphere required fewer total procedures per patient on average when compared to TACE<sup>16</sup>



At 5 years, TARE is a cost-effective treatment compared to cTACE, **saving more than \$15,779 per person.**<sup>17</sup>

# Committed to your cause.

With TheraSphere resources you can:

- Reduce time between ordering and treatment with expedited delivery, arriving as early as next day
- Order personalized doses online through TheraSphere Now – a free cloud-based platform with real-time inventory visibility, delivery tracking and order management
- Get case support, dosimetry insights, product expertise and educational resources from experienced field teams
- Apply recommendations from the Dosimetry Steering Committee with Simplicit<sup>90Y</sup>™ software
- Leverage patient education materials to help navigate the treatment journey with your patients



Visit **TheraSphere.com** to find resources or request to speak with a field consultant



# To place an order for TheraSphere™ Y-90 Glass Microspheres, please contact your TheraSphere consultant at [TherasphereCustomerSupport@bsci.com](mailto:TherasphereCustomerSupport@bsci.com)

or call Customer Service (US/CANADA): 1-866-363-3330.

Visit our website at [www.therasphere.com](http://www.therasphere.com)

1. Riaz A, Gates VL, Atassi B et al. Int J Radiat Oncol Biol Phys. 2011 Jan 1;79(1):163-71.
2. Salem R, Gordon AC, Mouli S, et al. Y90 Radioembolization Significantly Prolongs Time to Progression Compared With Chemoembolization in Patients With Hepatocellular Carcinoma. Gastroenterology. 2016;151(6):1155-1163.e2. doi:10.1053/j.gastro.2016.08.029.
3. Salem R, Johnson GE, Kim E, Riaz A, Bishay V, Boucher E, Fowers K, Lewandowski R, Padia SA. Yttrium-90 Radioembolization for the Treatment of Solitary, Unresectable Hepatocellular Carcinoma: The LEGACY Study. Hepatology. 2021 Mar 19. doi:10.1002/hep.31819. Funding provided by Boston Scientific.
4. Garin et al. J Clin Oncol 38, 2020 (suppl 4; abstr 516)
5. Reig M, Forner A, Rimola J, et al. BCLC strategy for prognosis prediction and treatment recommendation: The 2022 update. J Hepatol. 2021;S0168-8278(21)02223-6. doi:10.1016/j.jhep.2021.11.018. [published online ahead of print, 2022 May 23]. Lancet Gastroenterol Hepatol. 2022;S2468-1253(22)00091-7. doi:10.1016/S2468-1253(22)00091-7. Funding provided by Boston Scientific.
6. Lam, M., Garin, E., Maccauro, M. et al. A global evaluation of advanced dosimetry in transarterial radioembolization of hepatocellular carcinoma with Yttrium-90: the TARGET study. Eur J Nucl Med Mol Imaging (2022). <https://doi.org/10.1007/s00259-022-05774-0>. Funding provided by Boston Scientific.
7. Kim E, Sher A, Abboud G, et al. Radiation segmentectomy for curative intent of unresectable very early to early stage hepatocellular carcinoma (RASER): a single-centre, single-arm study [published online ahead of print, 2022 May 23]. Lancet Gastroenterol Hepatol. 2022;S2468-1253(22)00091-7. doi:10.1016/S2468-1253(22)00091-7. Funding provided by Boston Scientific.
8. Dhondt E, Lambert B, Hermie L, et al. 90Y Radioembolization versus Drug-eluting Bead Chemoembolization for Unresectable Hepatocellular Carcinoma: Results from the TRACE Phase II Randomized Controlled Trial. Radiology. 2022;303(3):699-710. doi:10.1148/radiol.211806.
9. Sandow T, Gimenez J, Nunez K, Tramel R, Gilbert P, Oliver B, Cline M, Fowers K, Cohen A, Thevenot P. Using Voxel-based Dosimetry to evaluate sphere concentration and tumor dose in Hepatocellular Carcinoma treated with Y-90 Radiation Segmentectomy with glass microspheres. Journal of Vascular and Interventional Radiology (2024). doi: <https://doi.org/10.1016/j.jvir.2024.05.020>.
10. Gabr A, Riaz A, Johnson GE, et al. Correlation of Y90-absorbed radiation dose to pathological necrosis in hepatocellular carcinoma: confirmatory multicenter analysis in 45 explants. Eur J Nucl Med Mol Imaging. 2020;48(2):580-583. doi:10.1007/s00259-020-04976-8.
11. Salem, R., Padia, S.A., Lam, M. et al. Clinical, dosimetric, and reporting considerations for Y-90 glass microspheres in hepatocellular carcinoma: updated 2022 recommendations from an international multidisciplinary working group. Eur J Nucl Med Mol Imaging 50, 328–343 (2023). <https://doi.org/10.1007/s00259-022-05956-w>. Funding provided by Boston Scientific.
12. Vouche M, Habib A, Ward TJ, et al. Unresectable solitary hepatocellular carcinoma not amenable to radiofrequency ablation: Multicenter radiology-pathology correlation and survival of radiation segmentectomy. Hepatology. 2014;60(1):192-201. doi:10.1002/hep.27057.
13. Montazeri SA, De la Garza-Ramos C, Lewis AR, et al. Hepatocellular carcinoma radiation segmentectomy treatment intensification prior to liver transplantation increases rates of complete pathologic necrosis: an explant analysis of 75 tumors. Eur J Nucl Med Mol Imaging. 2022;49(11):3892-3897. doi:10.1007/s00259-022-05776-y.
14. De la Garza-Ramos C, Montazeri SA, Croome KP, et al. Radiation Segmentectomy for the Treatment of Solitary Hepatocellular Carcinoma: Outcomes Compared with Those of Surgical Resection. Journal of Vascular and Interventional Radiology. 2022;33(7):775-785.e2. doi:10.1016/j.jvir.2022.03.021.
15. Padia SA, Kwan SW, Roudsari B, Monsky WL, Coveler A, Harris WP. Superselective Yttrium-90 Radioembolization for Hepatocellular Carcinoma Yields High Response Rates with Minimal Toxicity. Journal of Vascular and Interventional Radiology. 2014;25(7):1067-1073. doi:10.1016/j.jvir.2014.03.030.
16. Moreno-Luna LE et al. Cardiovasc Intervent Radiol 2013;36(3):714–23.
17. Patel, M. et al (2023) Transarterial therapies in patients with hepatocellular carcinoma eligible for Transarterial embolization: a US cost-effectiveness analysis, Journal of Medical Economics, 26:1, 1061-1071, DOI: 10.1080/13696998.2023.2248840.

## TheraSphere™ Yttrium-90 Glass Microspheres

**CAUTION:** Federal (USA) law restricts this device to sale by or on order of a physician. **INDICATION FOR USE:** TheraSphere is indicated for use as selective internal radiation therapy (SIRT) for local tumor control of solitary tumors (1-8 cm in diameter), in patients with unresectable hepatocellular carcinoma (HCC), Child-Pugh Score A cirrhosis, well-compensated liver function, no macrovascular invasion, and good performance status. **CONTRAINDICATIONS:** TheraSphere is contraindicated in patients: whose Tc-99m macroaggregated albumin (MAA) hepatic arterial perfusion scintigraphy shows any deposition to the gastrointestinal tract that may not be corrected by angiographic techniques • who show shunting of blood to the lungs that could result in delivery of greater than 16.5 mCi (0.61 GBq) of Y-90 to the lungs. Radiation pneumonitis has been seen rarely in patients receiving doses to the lungs greater than 30 Gy in a single treatment. • in whom hepatic artery catheterization is contraindicated, such as patients with vascular abnormalities or bleeding diathesis • who have pulmonary insufficiency (conventionally defined by an arterial oxygen pressure (PaO<sub>2</sub>) of < 60 mmHg, or oxygen saturation (SaO<sub>2</sub>) of < 90%) or severe liver dysfunction, including hepatic encephalopathy, clinically evident ascites or treatment with diuretics for ascites • with portal vein thrombosis (PVT) Type I involvement and lack of Tc-99m MAA deposition on the PVT seen on the Tc-99m MAA imaging with >70% tumor replacement in the liver • with comorbidities or poor overall health (e.g., ECG performance status rating > 2) which may make the patient a poor candidate for locoregional radiation treatment. • who are pregnant. **WARNINGS:** The following pre-treatment, high-risk factors (disease characteristics) have been associated with serious adverse events deemed possibly related to use of the device: infiltrative tumor type • tumor nodules too numerous to count • AST or ALT > 5 times ULN • bilirubin > 2 mg/dL • tumor volume > 50% combined with albumin < 3 g/dL. Keep the TheraSphere dose vial upright and stored in its lead pot before and during patient treatment, except as required for radiation measurement. Do not open the dose vial acrylic shield prior to patient treatment. Post-treatment, waste materials require caution to prevent contamination and beta shielding due to residual glass microspheres. **PRECAUTIONS: GENERAL PRECAUTIONS:** As in any intra-arterial procedure, aseptic technique should be practiced, and care should be taken to ensure minimum patient anesthesia exposure extraneous to therapeutic objective. • Consideration of patient comorbidities should be used when determining the type and volume of fluid to infuse via catheter to avoid electrolyte imbalance, fluid shift, and hyperglycemia. • It is important to avoid any aggressive arterial procedure that may lead to arterial spasm that impairs TheraSphere distribution into the perfused liver target volume which may lead to underdosing or non-target deposition of TheraSphere. **PRECAUTION IN PATIENTS WITH IMPAIRED LIVER FUNCTION:** No efficacy or safety data from the LEGACY study are available to support the use of the device in patients with Child-Pugh score B or C cirrhosis. **PRECAUTION IN VULNERABLE PATIENTS:** No effectiveness or safety data are available to support the use of the device in children or breast-feeding women. **ENDOCRINE DISRUPTION, CARCINOGENICITY, MUTAGENICITY, TOXICITY TO REPRODUCTION:** Ideally the use of this radioactive device in women of childbearing capability should be performed during the first few (approximately 10) days following the onset of menses. **RADIATION SAFETY:** Radioactive products should be used only by healthcare professionals who are qualified by specific training in the safe use and handling of radionuclides and whose experience and training have been approved by the appropriate government agency authorized to license the use of radionuclides. • As in the use of any radioactive material, ensure minimum radiation exposure to the patient extraneous to the therapeutic objective, and to minimize radiation exposure to workers and others in contact with the patient. **RELEASE AND POST-TREATMENT PRECAUTIONS:** Post treatment patient care: use universal precautions for body fluid contact. Trace Y-90 may be detectable in blood and urine; handle with gloves and dispose as normal body fluids. The radiation field is expected to be less than 1 mrem/h (10 μSv/h) at 3 ft (1 m) from the patient's abdomen. Supplemental shielding and segregation of the patient are not required to maintain exposure to others below regulated limits. • Release instructions: The patient should follow good hygiene (e.g., proper hand washing). Caregivers, family, and others do not require restrictions on patient contact; however, they can minimize their radiation exposure by avoiding prolonged time (>12 hours per day) within 1 ft (0.3 m) of the patient's abdomen for the first week post therapy. Patients should be advised that radiation emitted from the patient may be detectable at security screening (e.g., international travel). • Special precautions post-administration: If the patient requires hospitalization, surgery, medical assessment or treatment regarding any part of their thorax or abdomen within first 2 weeks of treatment, the patient should advise the hospital and treating physician of the Y-90 TheraSphere implant. The physician should consult their radiation safety staff for handling and disposal of liver tissue. • Special liver tissue handling: Special liver tissue handling may be required for post-treatment surgery, explant, or transplant since the glass microspheres remain permanently implanted in the liver tissue. Disclosure of the treatment will be required if cremation is considered. **POTENTIAL ADVERSE EVENTS:** The use of this product leads to irradiation of both tumorous and normal liver tissue. As a result, patients with compromised liver function may be at greater risk of liver function impairment and hence could experience complications. Clinical side effects usually occur within the first 4 to 6 weeks after treatment. Based on clinical trial data, literature reviews and post market surveillance, adverse events potentially associated with treatment using Y-90 microspheres, including TheraSphere, may include the following: Allergic reaction • Altered liver function, acute or chronic • Anorexia • Anxiety • Ascites • Bile Duct injury • Bleeding/hemorrhage • Chills / rigors • Cholecystitis (inflammatory or infectious) • Colitis • Death • Dehydration • Diarrhea • Dizziness • Dyspnea • Edema (any location) • Electrolyte abnormalities • Elevated BUN/creatinine • Fall • Fatigue • Fever • Gastrointestinal bleeding / hemorrhage • Gastrointestinal ulcer or ulceration • Hepatic encephalopathy • Hepatorenal failure • Hiccups • Hypertension • Hypotension • Infection (any location) • Liver failure, acute or chronic • Lymphopenia • Malaise • Mood alteration • Muscle weakness • Nausea • Neutropenia • Pain (any location) • Pancreatitis • Platelet count abnormalities • Pleural effusion • Portal hypertension • Pre-existing chronic liver disease decompensation • Pulmonary edema • Pulmonary fibrosis • Radiation hepatitis • Radiation induced disease, acute • Radio Embolization Induced Liver Disease (REILD) • Sepsis • Supraventricular arrhythmia • Thrombosis (arterial or venous) • Tumor inflammation (including tumor edema) • Tumor-lysis syndrome • Vomiting • Weight loss. Complications related to the administration procedure itself may include: Allergic reaction: Arterial injury including vessel dissection • Aspiration pneumonia • Bruising/bleeding/hematoma at site • Constipation/abdominal distension • Fatigue • Flushing • Infection • Nausea • Nerve damage. PI-992004-AA.

TheraSphere is a registered trademark of Theragnostics Corporation used under license by Boston Scientific Medical Device Limited, a wholly owned indirect subsidiary of Boston Scientific Corporation.

**Boston Scientific**  
Advancing science for life™

**Peripheral Interventions**  
300 Boston Scientific Way  
Marlborough, MA 01752-1234  
[www.bostonscientific.com](http://www.bostonscientific.com)

© 2024 Boston Scientific Corporation or its affiliates. All rights reserved.

PI-1941301-AA