REF 6386

Magnet in shielded box
Aimant dans coffret blindé
Magnet in abgeschirmtem Gehäuse
Magnete in scatola schermata
Magneet in afgeschermdoos
Imán en caja aislante

ENGLISH

WARNING: The magnet provided with the ICM system may cause interference with devices sensitive to magnetic fields such as hearing aids, pacemakers, and other implanted devices. It can also permanently disable some magnetic strip cards. Keep the magnet at least 15 cm (6 inches) away from items sensitive to magnetic fields, including the ICM device when the magnet is not being used to initiate communication between the device and the patient or clinic app.

WARNING: The magnet is MR Unsafe and must remain outside the MRI Site Zone III (and higher) as defined by the American College of Radiology Guidance Document on MR Safe Practices¹. Under no circumstances should the magnet be brought into the MRI scanner room, the control room, or the MRI site Zone III or IV areas.

For a full list of magnet-related warnings and precautions, see the User's Manual or Patient Handbook provided with the ICM system.

Patients: Report any serious incident that involves this magnet (i.e., an event requiring medical attention) to your health care provider, to Boston Scientific, and to the relevant local regulatory authority for medical devices in your country.

Items Included in the Packaging: Magnet in shielded box, adhesive strip, product literature.

If the adhesive strip wears over time, contact Boston Scientific for a replacement.

¹Kanal E, Barkovich AJ, Bell C, et al. ACR guidance document on MR safe practices: 2013. J. Magn. Reson. Imaging 2013;37:501-530.

Intended Use

The magnet is used to initiate communication between a LUX-Dx™ insertable cardiac monitor (ICM) and the LUX-Dx Clinic Assistant or myLUX™ software applications.

Related Information

Refer to the User's Manual or Patient Handbook and app screens for instructions on how and when to use the magnet, including any warnings, instructions related to establishing communication between the apps and the ICM device, verification that the magnet is performing as intended, and circumstances in which a patient should contact their healthcare provider.

Clinical Benefits of the Device

The magnet provided with the ICM system is part of a feature that allows the patient or clinician to initiate communication between the ICM device and the patient or clinic app. The magnet as an isolated device has no clinical benefit except for the capabilities that it enables when used as part of the ICM system. For clinical benefits of the ICM system, see the instructions for use provided with the ICM device or mobile software applications.

Summary of Safety and Clinical Performance

For customers in the European Union, use the Reference Number found in the labeling to search for the device's Summary of Safety and Clinical Performance, which is available on the European database on medical devices (Eudamed) website: <https://ec.europa.eu/tools/eudamed>.

Patient Contacting Materials

100% of the coated magnet surface area is made from a plastic material called PEBAX™ Polyether block amide or PEBA.²

Care and Disposal

Use the following instructions to avoid damage to the magnet and magnet case:

- Store the magnet between -20 to 70° C (-4 to 158° F).
- Clean your magnet only with the following cleaning agents: warm tap water, Isopropyl alcohol (70% concentration), or a weak concentration of mild detergent mixed with warm water. Apply cleaning agents to a soft cloth and use gentle pressure. Do not submerge the magnet or place under running water.
- If you lose your magnet or if the plastic coating surrounding the magnet cracks, contact Boston Scientific for a replacement.

²PEBAX is a trademark of Arkema.

There are no special requirements for disposal of the magnet or its packaging.

The following are trademarks of Boston Scientific Corporation or its affiliates: LUX-Dx, myLUX.

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All rights reserved
© 2022 Boston Scientific Corporation or its affiliates