AutoCap™ RX

Integrated Biopsy Cap and
Guidewire Locking Device

REF

M00545100, M00545130

en	Instructions for Use	2
es	Instrucciones de uso	6
fr	Mode d'emploi	10
de	Gebruichsanweisung	14
it	Istruzioni per l'uso	18
nl	Instructies voor gebruik	22
pt-EU	Instruções de Utilização	27

51474204-01 Rev. A

2022-08

TABLE OF CONTENTS

REUSE WARNING	3
DEVICE DESCRIPTION	3
Contents.....	3
Operating Principle.....	3
User Information.....	3
INTENDED USE/INDICATIONS FOR USE	3
Clinical Benefit Statement.....	3
CONTRAINDICATIONS	3
WARNINGS	3
PRECAUTIONS	3
ADVERSE EVENTS	3
HOW SUPPLIED	3
Device Details.....	3
Handling and Storage.....	4
Device compatibility.....	4
PRIOR TO USE	4
OPERATIONAL INSTRUCTIONS	4
Figure 1. Device Attachment.....	4
Figure 2. Device Locking Anchor.....	4
DEVICE REMOVAL	4
Figure 3. Device Removal.....	5
Disposal.....	5
Post-Procedure.....	5
WARRANTY	5
SYMBOL DEFINITIONS	5

AutoCap™ RX

Integrated Biopsy Cap and Guidewire Locking Device

RX ONLY

Caution: Federal Law (USA) restricts this device to sale by or on the order of a physician.

REUSE WARNING

For single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

DEVICE DESCRIPTION

The AutoCap RX is an integrated biopsy cap and guidewire locking device that fits on the working channel port of endoscopes. The device is an accessory to be used with endoscopic biliary devices to facilitate device passage, maintain insufflation, and lock the guidewire(s) in place during ERCP procedures.

Contents

Each packaged pouch contains the following:

- One (1) AutoCap RX Integrated Biopsy Cap and Guidewire Locking Device

Operating Principle

The ILD will be attached to the duodenoscope at the biopsy port. Physicians will pass devices through the biopsy cap and lock the biliary guidewire on the locking arm. Devices will be exchanged (pushed in and pulled out of the cap over or adjacent to the guidewire(s) or by themselves without any guidewire) during the procedure.

User Information

The AutoCap RX Integrated Biopsy Cap and Guidewire Locking Device should only be used by or under the supervision of physicians trained in Endoscopic Retrograde Cholangiopancreatography (ERCP) or Endoscopic Sphincterotomy (ES). A thorough understanding of the technical principles, clinical applications and risks associated with ERCP/ES is necessary before using this device.

INTENDED USE/INDICATIONS FOR USE

The AutoCap RX Integrated Biopsy Cap and Guidewire Locking Device is intended to facilitate placement of the guidewire and lock it into place during ERCP procedures. The device also provides access for endoscopic device passage and exchange, helps maintain insufflation, minimizes leakage of biomaterial from the biopsy port throughout the ERCP procedure, and provides access for irrigation.

Clinical Benefit Statement

The AutoCap RX clinical benefit is to facilitate placement of the guidewire and lock it into place during ERCP procedures and allow instrument access through the endoscope and maintain insufflation for adequate distension of the GI lumen.

CONTRAINDICATIONS

None known.

WARNINGS

None known.

PRECAUTIONS

- Any use of this device, other than those indicated in these instructions, is not recommended.

ADVERSE EVENTS

None known.

HOW SUPPLIED

Device Details

Contents supplied sterile using an ethylene oxide (EO) process.

Do not use if package is opened or damaged. Do not use if labeling is incomplete or illegible.

Handling and Storage

This product has no special handling or storage requirements.

Device compatibility

The AutoCap[®] RX Integrated Biopsy Cap and Guidewire Locking Device is compatible with 0.035 in (0.90 mm) and 0.025 in (0.63 mm) guidewires.

PRIOR TO USE

The AutoCap[®] RX Integrated Biopsy Cap and Guidewire Locking Device is supplied sterile in a sealed pouch and is intended for single use only. Do not open the package until ready to use. Carefully examine the unit to verify that neither the contents nor the sterile package have been damaged in shipment. DO NOT USE if damaged. Immediately return the damaged product to Boston Scientific.

OPERATIONAL INSTRUCTIONS

The AutoCap[®] RX Integrated Biopsy Cap and Guidewire Locking Device should be attached to the endoscope prior to the ERCP procedure.

1. Push the AutoCap[®] RX Integrated Biopsy Cap and Guidewire Locking Device onto the endoscope channel inlet. See Figure 1 for orientation.

Confirm that the device is firmly attached to the accessory channel.


Figure 1. Device Attachment

2. To insert devices through the AutoCap[®] RX Integrated Biopsy Cap and Guidewire Locking Device, align the device with the scope inlet and insert slowly. If not properly done, this may cause more friction on the system and may damage the seal or the scope's working channel.

Note: If resistance is felt during device introduction apply water-soluble lubricant to top of cap.

3. Once the AutoCap[®] RX Integrated Biopsy Cap and Guidewire Locking Device is in place, you may lock up to 2 guidewires at a time by moving the guidewire behind and around the "Locking Arm" and into either of the 2 "Locking Anchors" (See Figure 2).


Figure 2. Device Locking Anchor

4. To unlock the guidewire follow the steps given in Figure 2 in reverse order.

DEVICE REMOVAL

The AutoCap[®] RX Integrated Biopsy Cap and Guidewire Locking Device is removed by squeezing the device while rotating and pulling the device off the barrel of the endoscope (See Figure 3). The AutoCap[®] RX Integrated Biopsy Cap and Guidewire Locking Device should be removed and discarded after each patient use.


Figure 3. Device Removal

Disposal

To minimize the risk of infection or microbial hazards after use, dispose device and packaging as follows:

After use, device may contain biohazardous substances. The device and packaging should be treated and disposed of as biohazardous waste or have them treated and disposed of in accordance with any applicable hospital, administrative, and/or local government regulations. Use of a biohazardous container with biological hazard symbol is recommended. Untreated biohazardous waste should not be disposed of in the municipal waste system.

Post-Procedure

Any serious incident that occurs in relation to this device should be reported to the manufacturer, and to the relevant local regulatory authority.

WARRANTY

For device warranty information, visit (bostonscientific.com/warranty).

Olympus is a trademark of Olympus Corporation.

SYMBOL DEFINITIONS

Commonly used medical device symbols that appear on the labeling are defined at bostonscientific.com/SymbolsGlossary.

Additional symbols are defined at the end of this document.

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© 2022 Boston Scientific Corporation or its affiliates. All rights reserved.
Boston Scientific (Australia) Pty Ltd
Boston Scientific Limited
Boston Scientific Corporation