# OverTube™

<b>English</b> [Endoscopic Access System] INSTRUCTIONS FOR USE	<b>page 2</b>	<b>Dansk</b> [Endoskopisk adgangssystem] BRUGSANVISNING	<b>side 23</b>	<b>Eesti</b> [Endoskoopiline juurdepääsusüsteem] KASUTUSJUHISED	<b>lk 47</b>
<b>Español</b> [Sistema de acceso endoscópico] INSTRUCCIONES DE USO	<b>página 5</b>	<b>Norsk</b> [Endoskopisk tilgangssystem] BRUKSANVISNING	<b>side 26</b>	<b>Srpski</b> [Sistem za endoskopski pristup] UPUTSTVO ZA UPOTREBU	<b>strana 50</b>
<b>Deutsch</b> [System für den endoskopischen Zugang] GEBRAUCHSANWEISUNG	<b>Seite 8</b>	<b>Suomi</b> [Endoskoopinen yhteysjärjestelmä] KÄYTTÖOHJEET	<b>sivu 29</b>	<b>Slovenčina</b> [Systém na endoskopický prístup] NÁVOD NA POUŽITIE	<b>strana 53</b>
<b>Français</b> [Système d'accès endoscopique] MODE D'EMPLOI	<b>page 11</b>	<b>Svenska</b> [Endoskopiskt åtkomstsystem] BRUKSANVISNING	<b>sida 32</b>	<b>Slovenščina</b> [Endoskopski dostopni sistem] NAVODILA ZA UPORABO	<b>stran 56</b>
<b>Italiano</b> [Sistema di accesso endoscopico] ISTRUZIONI PER L'USO	<b>pagina 14</b>	<b>Polski</b> [Endoskopowy system dostępu] INSTRUKCJA UŻYTKOWANIA	<b>strona 35</b>	<b>Română</b> [Sistem de abord endoscopic] INSTRUCȚIUNI DE UTILIZARE	<b>pagina 59</b>
<b>Português</b> [Sistema de acesso endoscópico] INSTRUÇÕES DE UTILIZAÇÃO	<b>página 17</b>	<b>Čeština</b> [Systém pro endoskopický přístup] NÁVOD K POUŽITÍ	<b>strana 38</b>	<b>Українська</b> [Система ендоскопічного доступу] ІНСТРУКЦІЯ З ВИКОРИСТАННЯ	<b>стр. 62</b>
<b>Nederlands</b> [Endoscopisch toegangssysteem] GEBRUIKSAANWIJZING	<b>pagina 20</b>	<b>Magyar</b> [Endoszkópos hozzáférési rendszer] HASZNÁLATI ÚTMUTATÓ	<b>41. oldal</b>	<b>Ελληνικά</b> [Συστήμα ενδοσκοπικής προσβασης] ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ	<b>σελίδα 65</b>
 <b>OVERTUBE</b>		<b>български</b> [Система за ендоскопски достъп] ИНСТРУКЦИИ ЗА УПОТРЕБА	<b>страница 44</b>	<b>Türkçe</b> [Endoskopik erişim sistemi] KULLANMA TALİMATI	<b>sayfa 68</b>
		<b>العربية</b> [نظام الوصول بالمنظار الداخلي] تعليمات الاستعمال		<b>الصفحة 71</b>	

**Product Code:**

Non-Sterile  
Sterile

OVT-027-160  
OVT-S27-160

**Caution:**

Please read all instructions prior to use

SINGLE PATIENT USE  
DISPOSABLE

**GRF-00081-00R07**

**Single use only. Disposable. Do not resterilize.**

**Caution: Federal law (USA) restricts this device to sale by or on the order of a physician.**

**Patent Pending.**

**DISCLAIMER OF WARRANTY AND LIMITATION OF REMEDY**

THERE IS NO EXPRESS OR IMPLIED WARRANTY, INCLUDING WITHOUT LIMITATION ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE, ON THE APOLLO ENDOSURGERY, INC. PRODUCT(S) DESCRIBED IN THIS PUBLICATION. TO THE FULLEST EXTENT PERMITTED BY APPLICABLE LAW, APOLLO ENDOSURGERY, INC. DISCLAIMS ALL LIABILITY FOR ANY INDIRECT, SPECIAL, INCIDENTAL, OR CONSEQUENTIAL DAMAGES, REGARDLESS OF WHETHER SUCH LIABILITY IS BASED ON CONTRACT, TORT, NEGLIGENCE, STRICT LIABILITY, PRODUCTS LIABILITY OR OTHERWISE. THE SOLE AND ENTIRE MAXIMUM LIABILITY OF APOLLO ENDOSURGERY, INC., FOR ANY REASON, AND BUYER'S SOLE AND EXCLUSIVE REMEDY FOR ANY CAUSE WHATSOEVER, SHALL BE LIMITED TO THE AMOUNT PAID BY THE CUSTOMER FOR THE PARTICULAR ITEMS PURCHASED. NO PERSON HAS THE AUTHORITY TO BIND APOLLO ENDOSURGERY, INC. TO ANY REPRESENTATION OR WARRANTY EXCEPT AS SPECIFICALLY SET FORTH HEREIN. DESCRIPTIONS OR SPECIFICATIONS IN APOLLO ENDOSURGERY, INC PRINTED MATTER, INCLUDING THIS PUBLICATION, ARE MEANT SOLELY TO GENERALLY DESCRIBE THE PRODUCT AT THE TIME OF MANUFACTURE AND DO NOT CONSTITUTE ANY EXPRESS WARRANTIES OR RECOMMENDATIONS FOR USE OF THE PRODUCT IN SPECIFIC CIRCUMSTANCES. APOLLO ENDOSURGERY, INC. EXPRESSLY DISCLAIMS ANY AND ALL LIABILITY, INCLUDING ALL LIABILITY FOR ANY DIRECT, INDIRECT, SPECIAL, INCIDENTAL, OR CONSEQUENTIAL DAMAGES, RESULTING FROM REUSE OF THE PRODUCT.

**Table of Symbols**

Description	Symbol	Description	Symbol	Description	Symbol
Consult Instructions for Use		Manufacturer		Use By	
Do Not Re-use		Caution: Federal law (USA) restricts this device to sale by or on the order of a physician		Lot Number	
Sterilized Using Ethylene Oxide		Caution		Do not use if package is damaged	
Reference Number		Do Not Resterilize		Non-Sterile	
Authorized Representative in the European Community					

**1. Instructions for Use**


Carefully read all instructions prior to use. Observe all contraindications and warnings in these instructions. Failure to do so may result in complications. Apollo Endosurgery recommends that the surgical team carefully read and familiarize themselves with the instructions before removing the device from its packaging.

**2. Indications for Use**

The OverTube is intended to be used in conjunction with an endoscope for foreign body removal or endoscopic procedures requiring multiple insertions of the endoscope into the lower or upper gastrointestinal tract.


**3. General System Description**

The Apollo Endosurgery OverTube™ is comprised of a polymer extrusion reinforced with a metallic coil. The Hub located at the proximal end of the device contains a Cuff Seal which can be inflated with air using a Syringe to prevent loss of pressure when insufflation is used.

## 4. Packaged System

The OverTube is packaged inside a sealed protective pouch.

## 5. Recommended Accessories

- Water soluble lubricant
- 10ml Syringe (Aspiration Source)
- 63FR bite-block or larger

## 6. Contraindications

Contraindications include those specific to use of an overtube and any endoscopic procedure, which may include, but not be limited to, the following:

### Esophageal Procedures:

- Esophageal bleeding, ulceration and/or laceration
- Laryngeal perforation and/or stenosis
- Trauma to teeth, gums and/or pharynx which could be aggravated during the procedure
- Bleeding disorders/anticoagulation (non reversible bleeding disorders with coagulopathy)
- The patient has any other condition or anatomical limitation that would contraindicate them from an upper GI endoscopic procedure

### Colonic Procedures:

- Varices
- Bleeding disorders/anticoagulation (non-reversible bleeding disorders with coagulopathy)
- The patient has any other condition or anatomical limitation that would contraindicate them from a lower GI endoscopic procedure

## 7. Warnings

- Do not use device where the integrity of the packaging has been compromised or if the device appears damaged
- Lubrication required. Generously lubricate inner and outer surfaces of the OverTube with water soluble lubricant prior to use
- During insertion or removal of the OverTube, do not advance or retract against resistance
- Do not advance, remove, or reposition the OverTube unless under endoscopic visual guidance
- Reuse or reprocessing of the OverTube could result in device malfunction or patient consequences to include:
  - Infection or the transmission of disease
  - Damage to the polymer extrusion, compromising safe access
  - Impaired sealing when using insufflation
  - Accumulation of lubricant and debris

## 8. Precautions

- The system may only be used if purchased from Apollo Endosurgery Inc. or one of its authorized agents
- Prior to use, performing a baseline EGD can lower the likelihood of patient complications with the use of the OverTube
- Endoscopic procedures and foreign body retrieval should only be performed by clinicians with adequate training and familiarity with endoscopic techniques. Use of an overtube should not be attempted unless proficiency in technique has been developed
- Once placed, the OverTube should not be advanced beyond the endoscope as tissue entrapment and mucosal injury could occur. For any repositioning required, refer to Section 13. Endoscope Removal and Reinsertion

## 9. Adverse Events

Possible complications that may result from using the OverTube include, but may not be limited to:

- Hemorrhage
- Hematoma
- Infection / Sepsis
- Allergic reaction
- Pharyngeal, esophageal or bowel perforation
- Pharyngeal, esophageal or bowel laceration

## 10. Preparations for Use

### 10.1 Endoscope Selection

- 10.1.1. The Endoscope must have an outer diameter (OD) between 10.0mm and 16.0mm.

### 10.2 Accessory Preparation


- 10.2.1. Select and prepare all required accessories (listed in Section 5. Recommended Accessories).

### 10.3 Patient Preparation


- 10.3.1. Prepare the patient for endoscopy as per institutional guidelines or physician preference.
- 10.3.2. The use of a 63FR bite block is recommended to prevent collapse of the OverTube by patient gums or dentition.

### 10.4 Device Preparation


- 10.4.1. Remove the device from its outer package. Verify the device has not been compromised by assuring package integrity has been maintained.
- 10.4.2. Inspect the device for any bends, kinks, or other signs of damage. If damaged, do not use the device.  
NOTE: Check that Inflation Line is secure and tighten if required.
- 10.4.3. Generously lubricate the entire length of the Endoscope's outer sheath with a water soluble lubricant.

- 10.4.4. Load the Hub opening of the OverTube with a generous volume of water soluble lubricant.


- 10.4.5. Insert the distal end of the Endoscope through the Hub of the OverTube and slide the OverTube fully onto the Endoscope.

- 10.4.6. Ensure that the tapered tip of the OverTube forms a gap-free interface with the outer diameter of the Endoscope, and that it slides freely over the Endoscope.

- 10.4.7. Generously lubricate the outside of the OverTube including the tapered Tip.

## 11. Insertion

- 11.1. Perform baseline EGD.
- 11.2. Gently advance the OverTube over the well lubricated Endoscope until it reaches the desired anatomical location.


**WARNING:** During insertion or removal of the OverTube, do not advance or retract against resistance.

**WARNING:** Never advance or reposition the OverTube without endoscopic visual guidance.


**NOTE:** The OverTube Hub includes a 27cm Marker intended to aid positioning.

## 12. Seal Inflation (if insufflating)


- 12.1. Fill the 10ml Syringe with air.
- 12.2. Using the Luer fittings, connect the Syringe to the Aspiration Line.


- 12.3. Ensure the Seal Valve is open.

- 12.4. Squeeze the Syringe to transfer air into the Cuff Seal until adequate insufflation is achieved.

- 12.5. With the Syringe still depressed, close the Seal Valve.

NOTE: The Seal prevents loss of insufflation pressure while allowing for endoscopic movement.


**CAUTION:** If the Seal is over-inflated, Endoscope and accessory movement can be affected. If required, reduce the Seal force by opening the Seal Valve and withdrawing a small volume of air back into the Syringe.

### 13. Endoscope Removal and Reinsertion

- 13.1. The OverTube Seal should be deflated prior to removal and reinsertion of the Endoscope.


NOTE: If using insufflation, de-activating the Seal will cause the loss of gas / pressure.

- 13.2. Open the Seal Valve.

- 13.3. Fully withdraw the Syringe plunger to remove all air from the Seal.

- 13.4. Close the Seal Valve.

- 13.5. Slowly remove the Endoscope from the patient.

- 13.6. Re-apply additional lubrication to the Endoscope sheath and Hub opening as required prior to reinsertion.

NOTE: This is the same for procedure used for specimen or foreign body retrieval.

**CAUTION:** Do not try to remove specimens or foreign bodies that are larger than the OverTube can accommodate.

### 14. OverTube Removal

- 14.1. Slowly retract the OverTube and Endoscope out of the patient.

**WARNING:** During removal of the OverTube, do not advance or retract against resistance.

### 15. Product Disposal

This product must be treated as a potential biohazard after use. Dispose of in accordance with accepted medical practice and applicable local, state and federal regulations.

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All rights reserved.