OverTube™

English [Endoscopic Access System] INSTRUCTIONS FOR USE	page 2	Svenska [Endoskopiskt åtkomstsystem] BRUKSANVISNING	sida 32	Ελληνικά [Σύστημα ενδοσκοπικής πρόσβασης] ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ	σελίδα 65
Español [Sistema de acceso endoscópico] INSTRUCCIONES DE USO	página 5	Polski [Endoskopowy system dostępu] INSTRUKCJA UŻYTKOWANIA	strona 35	Türkçe [Endoskopik erişim sistemi] KULLANMA TALİMATI	sayfa 68
Deutsch [System für den endoskopischen Zugang] GEBRAUCHSANWEISUNG	seite 8	Čeština [Systém pro endoskopický přístup] NÁVOD K POUŽITÍ	strana 38	الصفحة 71 العربية [نظام الوصول بالمنظار الداخلي] تعليمات الاستعمال	
Français [Système d'accès endoscopique] MODE D'EMPLOI	page 11	Magyar [Endoszkópos hozzáférési rendszer] HASZNÁLATI ÚTMUTATÓ	41. oldal	Bahasa Indonesia [Sistem Akses Endoskopik] PETUNJUK PENGGUNAAN	halaman 74
Italiano [Sistema di accesso endoscopico] ISTRUZIONI PER L'USO	página 14	български [Система за ендоскопски достъп] ИНСТРУКЦИИ ЗА УПОТРЕБА	страница 44	简体中文 [内窥镜接入系统] 使用说明	第 77 页
Português [Sistema de acesso endoscópico] INSTRUÇÕES DE UTILIZAÇÃO	página 17	Eesti [Endoskoopiline juurdepääsusüsteem] KASUTUSJUHISED	lk 47	日本語 [内視鏡アクセスシステム] 取扱説明書	80ページ
Nederlands [Endoscopisch toegangssysteem] GEBRUIKSAANWIJZING	pagina 20	Srpski [Sistem za endoskopski pristup] UPUTSTVO ZA UPOTREBU	strana 50	Latviski [Endoskopiskā piekļuves sistēma] LIETOŠANAS INSTRUKCIJA	83. lappuse
Dansk [Endoskopisk adgangssystem] BRUGSANVISNING	side 23	Slovenčina [Systém na endoskopický přístup] NÁVOD NA POUŽITIE	strana 53	Lietuvių k. [Endoskopinės prieigos sistema] NAUDOJIMO INSTRUKCIJOS	86 puslapis
Norsk [Endoskopisk tilgangssystem] BRUKSANVISNING	side 26	Slovenščina [Endoskopski dostopni sistem] NAVODILA ZA UPORABO	stran 56	Hrvatski [Sustav za endoskopski pristup] UPUTE ZA UPORABU	stranica 89
Suomi [Endoskooppinen yhteysjärjestelmä] KÄYTTÖOHJEET	sivu 29	Română [Sistem de abord endoscopic] INSTRUCȚIUNI DE UTILIZARE	pagina 59	Русский язык [Система эндоскопического доступа] ИНСТРУКЦИЯ ПО ПРИМЕНЕНИЮ	страница 92
		Українська [Система ендоскопічного доступу] ІНСТРУКЦІЯ З ВИКОРИСТАННЯ	стр. 62		

OVERTUBE

Non-Sterile

Product Code:

OVT-027-160

Caution:

Please read all instructions prior to use

GRF-00081-00R11

Single use only. Disposable.

Caution: Federal law (USA) restricts this device to sale by or on the order of a physician.

Patent Pending.

DISCLAIMER OF WARRANTY AND LIMITATION OF REMEDY

THERE IS NO EXPRESS OR IMPLIED WARRANTY, INCLUDING WITHOUT LIMITATION ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE, ON THE APOLLO ENDOSURGERY, INC. PRODUCT(S) DESCRIBED IN THIS PUBLICATION. TO THE FULLEST EXTENT PERMITTED BY APPLICABLE LAW, APOLLO ENDOSURGERY, INC. DISCLAIMS ALL LIABILITY FOR ANY INDIRECT, SPECIAL, INCIDENTAL, OR CONSEQUENTIAL DAMAGES, REGARDLESS OF WHETHER SUCH LIABILITY IS BASED ON CONTRACT, TORT, NEGLIGENCE, STRICT LIABILITY, PRODUCTS LIABILITY OR OTHERWISE. THE SOLE AND ENTIRE MAXIMUM LIABILITY OF APOLLO ENDOSURGERY, INC., FOR ANY REASON, AND BUYER'S SOLE AND EXCLUSIVE REMEDY FOR ANY CAUSE WHATSOEVER, SHALL BE LIMITED TO THE AMOUNT PAID BY THE CUSTOMER FOR THE PARTICULAR ITEMS PURCHASED. NO PERSON HAS THE AUTHORITY TO BIND APOLLO ENDOSURGERY, INC. TO ANY REPRESENTATION OR WARRANTY EXCEPT AS SPECIFICALLY SET FORTH HEREIN. DESCRIPTIONS OR SPECIFICATIONS IN APOLLO ENDOSURGERY, INC PRINTED MATTER, INCLUDING THIS PUBLICATION, ARE MEANT SOLELY TO GENERALLY DESCRIBE THE PRODUCT AT THE TIME OF MANUFACTURE AND DO NOT CONSTITUTE ANY EXPRESS WARRANTIES OR RECOMMENDATIONS FOR USE OF THE PRODUCT IN SPECIFIC CIRCUMSTANCES. APOLLO ENDOSURGERY, INC. EXPRESSLY DISCLAIMS ANY AND ALL LIABILITY, INCLUDING ALL LIABILITY FOR ANY DIRECT, INDIRECT, SPECIAL, INCIDENTAL, OR CONSEQUENTIAL DAMAGES, RESULTING FROM REUSE OF THE PRODUCT.

Table of Symbols

Description	Symbol	Description	Symbol	Description	Symbol
Consult Instructions for Use		Manufacturer		Use By	
Do Not Re-use		Caution: Federal law (USA) restricts this device to sale by or on the order of a physician		Lot Number	
Medical Device		Caution		Do not use if package is damaged	
Reference Number		Authorized Representative in the European Community		Non-Sterile	
Inner Diameter (mm)					

1. Instructions for Use


Carefully read all instructions prior to use. Observe all contraindications and warnings in these instructions. Failure to do so may result in complications. Apollo Endosurgery recommends that the surgical team carefully read and familiarize themselves with the instructions before removing the device from its packaging.

2. Indications for Use

The OverTube is intended to be used in conjunction with an endoscope for foreign body removal or endoscopic procedures requiring multiple insertions of the endoscope into the lower or upper gastrointestinal tract.

3. General System Description

The Apollo Endosurgery OverTube™ is comprised of a polymer extrusion reinforced with a metallic coil. The Hub located at the proximal end of the device contains a Cuff Seal which can be inflated with air using a Syringe to prevent loss of pressure when insufflation is used.


4. Clinical Benefit Statement

The clinical benefit of the Apollo Endosurgery OverTube™ is to protect the GI mucosa from trauma and limit the risk of aspiration.

5. Packaged System

The OverTube is packaged inside a sealed protective pouch.

6. Recommended Accessories

- Water soluble lubricant
- 10ml Syringe (Aspiration Source)
- 63FR bite-block or larger

7. Contraindications

Contraindications include those specific to use of an overtube and any endoscopic procedure, which may include, but not be limited to, the following:

Esophageal Procedures:

- Esophageal bleeding, ulceration and/or laceration
- Laryngeal perforation and/or stenosis
- Trauma to teeth, gums and/or pharynx which could be aggravated during the procedure
- Bleeding disorders/anticoagulation (non reversible bleeding disorders with coagulopathy)
- The patient has any other condition or anatomical limitation that would contraindicate them from an upper GI endoscopic procedure

Colonic Procedures:

- Varices
- Bleeding disorders/anticoagulation (non-reversible bleeding disorders with coagulopathy)
- The patient has any other condition or anatomical limitation that would contraindicate them from a lower GI endoscopic procedure

8. Warnings

- Do not use device where the integrity of the packaging has been compromised or if the device appears damaged
- Lubrication required. Generously lubricate inner and outer surfaces of the OverTube with water soluble lubricant prior to use
- During insertion or removal of the OverTube, do not advance or retract against resistance
- Do not advance, remove, or reposition the OverTube unless under endoscopic visual guidance
- Reuse or reprocessing of the OverTube could result in device malfunction or patient consequences to include:
 - Infection or the transmission of disease
 - Damage to the polymer extrusion, compromising safe access
 - Impaired sealing when using insufflation
 - Accumulation of lubricant and debris

9. Precautions

- The system may only be used if purchased from Apollo Endosurgery Inc. or one of its authorized agents.
- Prior to use, performing a baseline EGD can lower the likelihood of patient complications with the use of the OverTube.
- Endoscopic procedures and foreign body retrieval should only be performed by clinicians with adequate training and familiarity with endoscopic techniques. Use of an overtube should not be attempted unless proficiency in technique has been developed.
- Once placed, the OverTube should not be advanced beyond the endoscope as tissue entrapment and mucosal injury could occur. For any repositioning required, refer to Section 14. Endoscope Removal and Reinsertion.
- The OverTube is a single-use device intended to support one endoscopic procedure. The duration of the procedure can be over one hour. As the number of endoscopic passes increases, the endoscope may need to be relubricated to ensure smooth passage down the OverTube. Additionally, the cuff may need to be reinflated periodically during the procedure.

10. Adverse Events

Possible complications that may result from using the OverTube include, but may not be limited to:

- Hemorrhage
- Hematoma
- Infection / Sepsis
- Allergic reaction
- Pharyngeal, esophageal or bowel perforation
- Pharyngeal, esophageal or bowel laceration

NOTE: Any serious incident that has occurred in relation to the device should be reported to Apollo Endosurgery (see contact information at the end of this document) and any appropriate government entity.

11. Preparations for Use

11.1. Endoscope Selection

- 11.1.1 The Endoscope must have an outer diameter (OD) between 10.0mm and 16.0mm.

11.2. Accessory Preparation

- 11.2.1 Select and prepare all required accessories (listed in Section 6. Recommended Accessories).

11.3. Patient Preparation


- 11.3.1 Prepare the patient for endoscopy as per institutional guidelines or physician preference.
- 11.3.2 The use of a 63FR bite block is recommended to prevent collapse of the OverTube by patient gums or dentition.

11.4. Device Preparation


- 11.4.1 Remove the device from its outer package. Verify the device has not been compromised by assuring package integrity has been maintained.
- 11.4.2 Inspect the device for any bends, kinks, or other signs of damage. If damaged, do not use the device.
NOTE: Check that Inflation Line is secure and tighten if required.
- 11.4.3 Generously lubricate the entire length of the Endoscope's outer sheath with a water-soluble lubricant.

- 11.4.4 Load the Hub opening of the OverTube with a generous volume of water soluble lubricant.


- 11.4.5 Insert the distal end of the Endoscope through the Hub of the OverTube and slide the OverTube fully onto the Endoscope.

- 11.4.6 Ensure that the tapered tip of the OverTube forms a gap-free interface with the outer diameter of the Endoscope, and that it slides freely over the Endoscope.

- 11.4.7 Generously lubricate the outside of the OverTube including the tapered tip.

12. Insertion

- 12.1 Perform baseline EGD.
- 12.2 Gently advance the OverTube over the well lubricated Endoscope until it reaches the desired anatomical location.


WARNING: During insertion or removal of the OverTube, do not advance or retract against resistance.

WARNING: Never advance or reposition the OverTube without endoscopic visual guidance.


NOTE: The OverTube Hub includes a 27cm Marker intended to aid positioning.

13. Seal Inflation (if insufflating)


- 13.1 Fill the 10ml Syringe with air.
- 13.2 Using the Luer fittings, connect the Syringe to the Aspiration Line.


- 13.3 Ensure the Seal Valve is open.

- 13.4. Squeeze the Syringe to transfer air into the Cuff Seal until adequate insufflation is achieved.

- 13.5. With the Syringe still depressed, close the Seal Valve.


NOTE: The Seal prevents loss of insufflation pressure while allowing for endoscopic movement.


CAUTION: If the Seal is over-inflated, Endoscope and accessory movement can be affected. If required, reduce the Seal force by opening the Seal Valve and withdrawing a small volume of air back into the Syringe.

14. Endoscope Removal and Reinsertion

- 14.1. The OverTube Seal should be deflated prior to removal and reinsertion of the Endoscope.
NOTE: If using insufflation, de-activating the Seal will cause the loss of gas / pressure.
- 14.2. Open the Seal Valve.

- 14.3. Fully withdraw the Syringe plunger to remove all air from the Seal.
- 14.4. Close the Seal Valve.

- 14.5. Slowly remove the Endoscope from the patient.
- 14.6. Re-apply additional lubrication to the Endoscope sheath and Hub opening as required prior to reinsertion.
NOTE: This is the same for procedure used for specimen or foreign body retrieval.

CAUTION: Do not try to remove specimens or foreign bodies that are larger than the OverTube can accommodate.

15. OverTube Removal

- 15.1. Slowly retract the OverTube and Endoscope out of the patient.

WARNING: During removal of the OverTube, do not advance or retract against resistance.

16. Product Disposal

This product must be treated as a potential biohazard after use. Dispose of in accordance with accepted medical practice and applicable local, state and federal regulations.

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All rights reserved.