OverStitch Sx™

English **page 2**

[Endoscopic Suturing System] (ESS)

INSTRUCTIONS FOR USE

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[Endoscopic Suturing System] (ESS)

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[Endoskopisches Nahtsystem] (ESS)

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OVERSTITCH Sx™

HELIX

OverStitch Sx™ Endoscopic Suturing System
OverStitch™ 2-0 Polypropylene Suture
OverStitch™ Suture Cinch
Tissue Helix

Product Code:

ESS-G02-Sx1
PLY-G02-020-APL
CNH-G01-000
THX-165-028

Caution:

Please read all instructions prior to use

STERILE EO

**SINGLE PATIENT USE
DISPOSABLE**

GRF-00554-00R01

English

Single use only. Disposable. Do not resterilize.

Patent Pending.

DISCLAIMER OF WARRANTY AND LIMITATION OF REMEDY

Table of Symbols

| Description | Symbol | Description | Symbol | Description | Symbol |
|---------------------------------|--------|---|--------|--|--------|
| Consult Instructions for Use | | Manufacturer | | Use By | |
| Do Not Re-use | | Single Sterile Barrier System. Sterilized by Ethylene Oxide | | Lot Number | |
| Sterilized Using Ethylene Oxide | | Date of manufacture | | Do not use if package is damaged | |
| Reference Number | | Do Not Resterilize | | Authorized Representative in the European Community | |
| Medical Device | | Caution, contact accompanying documents | | Contains hazardous substance (CAS 7440-48-4, Cobalt) | |

1. Intended Use

The OverStitch Sx Endoscopic Suture System is indicated for the endoscopic placement of anchor-sutures. Anchor-sutures can be placed for defect closure (e.g. perforation, ESD/EMR and fistula/leak), stent fixation to reduce migration rates in the upper GI tract, endoscopic sleeve gastropasty, and transoral outlet reduction.

1.1 Intended Users

The OverStitch Endoscopic Suturing System is operated by the physician (e.g., doctors performing endoscopic procedures) and supported by allied health personnel (e.g. nurses, physician assistants). Apollo Endosurgery offers basic training on the use of OverStitch and supplemental training on endoscopic sleeve gastropasty and transoral outlet reduction. This training covers patient selection, potential adverse events, prophylactic techniques, how to perform the procedure, and after care for the patient. Physicians performing bariatric procedures should have this supplemental training. Contact your local Apollo Endosurgery representative to inquire about training.

1.2 Intended Patient Population

The system is designed to operate in the gastrointestinal (GI) tract. Prospective patients are from the general adult population who have a pathology in their gastrointestinal tract or have obesity, except for those from whom endoscopic procedures are contraindicated.

1.3 Clinical Benefit Statement

The OverStitch Endoscopic Suturing System is intended to endoscopically place sutures and approximate soft tissue within the GI tract. The clinical benefit can be measured by overall clinical outcomes, including, but not limited to, successful placement of sutures to close defects, to reduce the rate of stent migration and to induce weight loss through reducing the volume of the stomach or reducing a dilated gastric outlet.

1.4 Summary of Safety and Clinical Performance

The Summary of Safety and Clinical Performance document, as required by the European Medical Device Regulation, is located at <https://ec.europa.eu/tools/eudamed>.

1.5 Information Provided to the Patient

The packaging includes a Patient Implant Card and leaflet. The leaflet instructs the medical staff on how to complete the card. The card documents the date of the procedure, patient's name, physician's contact information, device tracking information, and MR Safety information. Patients should receive the completed card and leaflet after their procedure.

1.6 Contraindications

Contraindications include those specific to use of an endoscopic suturing system, and any endoscopic procedure, which may include, but not limited to, the following:

- This system is not for use where endoscopic techniques are contraindicated.
- This system is not for use with malignant tissue.

1.7 Warnings

- Do not use a device where the integrity of the sterile packaging has been compromised or if the device appears damaged.
- Only physicians possessing sufficient skill and experience in similar or the same techniques should perform endoscopic procedures.
- Contact of electrosurgical components with other components may result in injury to the patient and/or operator as well as damage to the device and/or endoscope.
- Verify compatibility of endoscopic instruments and accessories and ensure performance is not compromised.

NOTE: Refurbished scopes may no longer conform to original specifications.
- Ensure that there is sufficient space for the Needle to open.
- Ensure that the Handle Grip of the Endoscopic Suturing System is closed and locked during intubation and extubation.
- Reuse or reprocessing of the OverStitch system could result in device malfunction or patient consequences to include:
 - Infection or the transmission of disease
 - Failure of the handle mechanism causing the device to become locked on tissue that may require surgical intervention
 - Reduced retention on the Endoscope, causing the Endcap to detach during use that may require surgical intervention to retrieve
 - Reduced retention of the Anchor to the Needle Body, resulting in an inadvertent Anchor drop causing procedural delay or requiring subsequent intervention
 - Bending of the Needle Body, preventing the physician from driving the Needle correctly or performing the intended procedure
 - Failure of the Helix to extend fully, limiting the ability to acquire tissue and perform intended procedure
- If subject device is used to oversew foreign objects, such as staples, stents, clips or mesh, it is possible for the needle to become trapped in the foreign body, requiring surgical intervention.

- In situations where the operative site poses a risk of harm to adjacent anatomic structures, use of endoscopic accessories such as the OverStitch Tissue Helix is recommended to retract the tissue intended to be sutured away from these unseen structures.
- It is important to ensure the Tissue Helix is carefully deployed and correctly retracted to avoid entrapping tissue and potentially causing trauma. Avoid using excessive pressure or applying excess turns when deploying the Tissue Helix. Performing more turns than necessary to retract tissue may increase the risk of capturing and suturing an adjacent organ and the risk of the helix entrapping tissue, complicating removal of the instrument.
- For bariatric cases, Carbon Dioxide (CO₂) is required for insufflation. Room air should not be used to insufflate and could contribute to serious adverse events including pneumoperitoneum, pneumothorax, pneumomediastinum, and death.
- Avoid placing plications in the fundus. The fundus is relatively thin walled and located close to the spleen and diaphragm. Sutures placed in the fundus may increase the risks of leakage and inadvertent suturing of the adjacent organs.
- Maintain awareness of the potential to disrupt a short gastric artery along the greater curve. Post procedure pain with any hemodynamic instability should immediately raise concern for extra-gastric bleeding and/or hematoma formation. Management of this should include imaging, e.g., with CT along with serum hemoglobin measurements.
- When cinching the suture anchor to form the plications, use the minimum tension necessary to maintain the plication. Excessive tension may increase the risk of gastrointestinal bleeding or creating a leak. Excessive tension may also increase the risk of the suture-anchor breaking and compromising the gastric sleeve.
- Patients who develop significant persistent upper abdominal pain at any time after an ESG, with radiation to the back or supraclavicular area along with pleuritic symptoms or even dyspnea, may have developed a needle puncture site leak with the development of a sterile or infected fluid collection and inflammatory pleural effusion. These symptoms warrant investigation with an imaging study, e.g., CT.
- Reuse or reprocessing of the OverStitch system could result in device malfunction or patient consequences to include infection, failure of the device to function properly.

1.8 Precautions

- The System may only be used if purchased from Apollo Endosurgery, Inc. or one of its authorized agents.
- With the Endoscopic Suturing System installed, the endoscope's primary channel effectively becomes a 3.2 mm channel.

- An overtube with an internal diameter of at least 16.7 mm may be used with the system to protect the esophagus.

1.9 System Compatibility

The OverStitch Sx ESS is compatible with PLY-G02-020-APL sutures.

The system is compatible with endoscopes having an insertion tube and distal diameter between 8.8 mm and 9.8 mm, a working length up to 110 cm, and using overtubes having an inner diameter of at least 16.7 mm.

Availability of compatible devices may vary depending on geography.

1.10 Adverse Events

Possible complications that may result from using the Endoscopic Suturing System include, but may not be limited to:


- Pharyngitis / Sore throat
- Vomiting
- Nausea
- Moderate abdominal pain more than 24 hours after procedure. In some cases, abdominal pain may be severe and require medical intervention
- Constipation

- Generalized weakness after procedure
- Heartburn
- Fever
- Gastrointestinal bleeding (with or without melena or hematemesis)
- Dehydration and/or nutritional deficiency requiring hospital admission
- Perigastric fluid collection
- Leak
- Hemoperitoneum
- Hematoma
- Paresthesia
- GERD
- Peritonitis
- Pneumoperitoneum
- Pulmonary Embolism
- Perforation (gastric or esophageal)
- Pneumothorax
- Pneumomediastinum
- Gall bladder suture

- Spleen Laceration
- Deep Vein Thrombosis
- Esophageal tear
- Pleural effusion
- Persistent Vomiting
- Bowel obstruction
- Infection/sepsis
- Bloating
- Stricture
- Liver abscess
- Intra-abdominal (hollow or solid) visceral injury
- Aspiration
- Shortness of breath
- Acute inflammatory tissue reaction
- Death


NOTE: Any serious incident that has occurred in relation to the device should be reported to Apollo Endosurgery (see contact information at the end of this document) and any appropriate government entity.

2. Packaged System


Package 1

Needle Driver & Anchor Exchange


Package 2
Suture Assembly

Polypropylene (nonabsorbable) - BLUE


Package 3
Cinch


Package 4
Helix

(AVAILABLE SEPARATELY)

3. Nomenclature


Needle Driver Handle


Anchor Release Button

Helix Handle

Cinch Handle

Needle Driver


Anchor Exchange


Suture Assembly

Helix

Cinch


Endcap

Anchor Receptacle

Anchor Body


Suture


Helix Tip

Suture Loading Loop

Collar Plug

English

4. Assembly

Prepare Endoscope


WARNING: Endoscope distal end and insertion tube must be dry and free of any lubricants prior to device installation. Failure to do so may cause the device to slip off endoscope during use.

- 4.1. Connect endoscope to video processor and light source to assist device installation.
- 4.2. Ensure endoscope is straight with wheels locked prior to device installation.


Prepare Device

- 4.3. Remove Needle Driver from packaging.


CAUTION: Ensure Endcap is not dropped or otherwise damaged.

- 4.4. Close Needle Body.


- 4.5. From the Actuation Catheter side of the device, thread the endoscope between Actuation Catheter and Catheter Sheath.



Secure Handle

- 4.6. Slide Needle Driver Handle Collar firmly onto endoscope handle until the channel ports are positioned near the endoscope channel port.




- 4.7. Secure by stretching the rubber Handle Strap around the rear of endoscope and fitting end back onto attachment peg.

Position Endcap

- 4.8. Fully insert endoscope distal end through clear Endcap Strap loops until endoscope face touches the Overhang Stop feature.


CAUTION: Ensure the Catheter Sheath is not twisted along the length of the scope.

- 4.9. Use the monitor to orient the Endcap as desired.

NOTE: Final orientation on monitor may shift clockwise during assembly and can be compensated for here.


- 4.10. Confirm both ends of the Needle Body are visible on the monitor.

Secure Endcap Straps

- 4.11. With left thumb placed on the Apollo logo of Endcap clip, slowly pull Strap 1 until a watertight seal is achieved.

NOTE: The strap appears clear in colour when sufficiently tightened.


WARNING: Straps should not be stretched beyond 3 cm (1 inch) of travel as damage may occur.

WARNING: Do not attempt to twist or reposition Endcap. Doing so may damage device or endoscope.

- 4.12. Repeat for Strap 2.

- 4.13. To lock the Straps, apply further stretching of both Straps 1 and 2 to remove any slack and firmly press the Apollo logo on Clip Hinge until Hinge Step is fully recessed.


NOTE: The Hinge Step is easier to close whilst a stretch is applied to both Straps.

Fix Catheter Sheath to Endoscope

CAUTION: Ensure endoscope is flat, straight and with scope wheels locked.

- 4.14. Working from the distal end, remove protective backing from each side of the Sheath Strap and wrap Strap around endoscope, pressing back onto itself.



- 4.15. Work along the length of the Catheter Sheath and secure remaining 4 Straps.

CAUTION: The endoscope should fit within the inner surface of the Catheter Sheath along its length.

CAUTION: Do not overstretch any Sheath Strap.


WARNING: Ensure the Catheters and Sheath are flush along the length of the scope.

WARNING: Ensure all Sheath Straps are adhered to themselves and not misaligned. An exposed corner may peel if force is applied.

Remove Endcap Clip

4.16. While securing Endcap, twist Endcap Clip counter-clockwise (indicated by arrows) to remove from Endcap.


4.17. Discard Endcap Clip.

Trim Excess Endcap Straps

WARNING: Confirm the concealed ends of both Endcap Straps are fully recessed within the Endcap.


4.18. Use a surgical blade to trim the excess lengths of both Endcap Straps flush with the Endcap surface.

Install Biopsy Valves

4.19. Fit standard biopsy valves to all endoscope and device working channels.

4.20. The biopsy valve for the Anchor Exchange (outside) channel should be in the open position.


Load Anchor

4.21. Select appropriate Suture.

4.22. Remove Suture Assembly and Anchor Exchange from packaging.

4.23. Load the Anchor onto the Anchor Exchange.


4.24. Remove Suture from Suture Cartridge by holding and pulling the Suture, not the Anchor or Anchor Exchange.

NOTE: Ensure Suture is not tangled after removal from the cartridge.

4.25. Insert Anchor Exchange into the open Anchor Exchange (outside) channel of the device.

CAUTION: Do not use when valve covers are closed as Suture drag will be increased.

4.26. Using a 'pencil-grip' on the catheter for optimum control, advance Anchor Exchange until Anchor is positioned near to the distal end of the device.

CAUTION: If resistance is encountered when advancing the Anchor Exchange through the Anchor Exchange Channel, reduce the endoscope angulation until the device passes smoothly.

Lubricate Device

4.27. Fully lubricate the distal 50 cm of the installed device and endoscope, and full length of the inner walls of the overtube.

WARNING: Use of an appropriate overtube (with an inner diameter of at least 16.7 mm) is recommended for transoral procedures.


WARNING: Confirm installed device is of compatible size with overtube prior to use.

4.28. Insert endoscope into patient.

WARNING: Do not introduce the device with the Needle Body in its open position.

5. Navigating to target anatomy


5.1. Advance until target anatomy is located.


CAUTION: When intubating or extubating with the Endoscopic Suturing System, ensure that the working length of the endoscope and the external Sheath are advanced or retracted together and any slack in Actuation Catheter is removed.

CAUTION: Do not articulate endoscope without device accessories installed up to the distal end of the primary channel and secondary channel (if in use) as this may result in minor kinking of catheters.


5.2. Open Needle Body.


5.3. Advance Anchor Exchange and/or manipulate endoscope to create Suture slack.


5.4. Once sufficient slack has been created, retract Anchor Exchange into the device.

5.5. Close Needle Body.


5.6. Use a 'pencil-grip' on the **white** portion of the cable to advance Anchor Exchange until Anchor fully seats onto Needle Body.

5.7. A tactile 'click' or firm stop may be felt as the Anchor seats fully onto Needle Body.


5.8. Check the monitor image to ensure Anchor is properly installed on Needle Body.


5.9. Ensure that the proximal end of the Suture is visible beyond the Handle channel valves.

English

- 5.10. Fully depress Anchor Release Button to release Anchor.


- 5.11. With Anchor Release Button still fully depressed, slightly retract Anchor Exchange.

6. Tissue and Suture Management

- 6.1. Open Needle Body.



CAUTION: If the Needle Body does not open, ensure that the Anchor has been released from the Anchor Exchange.

WARNING: Ensure appropriate Suture slack has been created for desired Suture path and pattern. Advance Anchor Exchange and/or manipulate endoscope to create Suture slack.

- 6.2. Position tissue in appropriate location for suturing, using Apollo Helix or other 2.8 mm compatible accessory if required.

7. Use of Tissue Helix (Optional for Defect Closure and Stent Fixation but Required for Endoscopic Sleeve Gastroplasty and Transoral Outlet Reduction)


- 7.1. Advance Helix into Helix Channel of ESS in the retracted position, until distal tip is visible on the monitor.

CAUTION: If resistance is encountered when advancing the Helix through the Helix Channel of the ESS, reduce the endoscope angulation until the device passes smoothly,


and ensure that the device's secondary working channel is not obstructed.

- 7.2. Fully depress Helix Handle button to expose Helix Tip.


CAUTION: Do not depress Helix Handle button while advancing Helix through the device.

- 7.3. Acquire tissue by rotating Helix Handle clockwise until correct tissue depth is achieved.


- 7.4. Maintain gentle forward pressure during tissue acquisition.

- 7.5. Advance / retract Helix to position tissue into desired location.

8. Suturing


- 8.1. Drive Needle through tissue by closing the Needle Body.

- 8.2. Ensure that the Needle Arm does not inadvertently close on any foreign object or device.


Retrieve Anchor

- 8.3. Use a 'pencil-grip' to advance the Anchor Exchange until the Anchor is engaged and resistance is felt.

NOTE: Resistance may vary due to the position of the endoscope.

- 8.4. Retract Anchor Exchange to acquire Anchor.


- 8.5. Use a 'pencil-grip' on the white section of the cable and place the remaining fingers of the same hand on the device handle to prevent damage to Suture or tissue once the Anchor 'pops' off the Needle Body.

CAUTION: Do not press the Anchor Release Button, as this may cause inadvertent drop of the Anchor.


WARNING: If sufficient slack has not been created prior to driving the Anchor through tissue, retraction of the Anchor Exchange may be difficult and the Anchor may not release correctly from the Needle Body.

Release Tissue

- 8.6. Rotate Helix handle counter-clockwise until the device is free from tissue.


- 8.7. Return the Helix handle button.


- 8.8. Check monitor to ensure tip is fully retracted into shaft before withdrawing into Helix Channel.

- 8.9. Withdraw Helix a short distance into Helix Channel.

- 8.10. Open Needle Body

CAUTION: Do not tension the Suture with the Anchor in the Needle Body.


8.11. To **continue** placing stitches with this Anchor, repeat **sections 5-8**.

If stitching is **complete** for this Anchor, proceed to **section 9** to approximate tissue, secure and cut the Suture.


Multiple Anchors can be utilized with each Needle Driver and Anchor Exchange.

9. Securing and Cutting Suture


9.1. Ensure Anchor is in Anchor Exchange and open Needle Body.

9.2. Advance Anchor distal of the device until visible on monitor.


9.3. Fully depress Anchor Release Button and retract Anchor Exchange to release Anchor.


CAUTION: Do not release the Anchor inside the Anchor Exchange channel.

9.4. Remove Anchor Exchange from ESS.


NOTE: The Anchor Exchange can be utilized for additional sutures.

9.5. Feed proximal end of Suture into removable Suture Loading Loop.


9.6. After threading, release the proximal end of the Suture to enable loading.


9.7. Pull Suture Loading Loop parallel to the device, to pull Suture into Cinch.

9.8. Holding the proximal end of the Suture, feed the Cinch down the Anchor Exchange channel until 'Plug and Collar' can be seen on the monitor.


9.9. Pull the Suture and apply counter traction to the Cinch until the tissue is approximated and the desired Suture tension is achieved between the Anchor and the Cinch Collar.

NOTE: It is the Collar that sets the final position of the Cinch, not the Plug.


WARNING: Excessive tension may damage tissue.

9.10. While holding suture tension, remove safety spacer from Cinch handle.

CAUTION: The safety spacer must only be removed immediately prior to deploying Cinch.


9.11. Firmly squeeze Cinch Handle to deploy Cinch and cut Suture.


CAUTION: Suture tension must be maintained during cinch deployment.

NOTE: Significant force is required to pull and lock the Plug into the Collar, a 'pop' can often be felt once the Suture is cut. Do not over deploy cinch - only squeeze handle until suture is cut.

9.12. Remove the Cinch.

9.13. Close Needle Body.

10. Multiple Sutures

10.1. Anchor Exchange and Needle Driver can be utilized with multiple sutures. Removal of the endoscope after cinching is not required if further sutures are to be deployed.

WARNING: If endoscope is removed, ensure Endcap is secure to endoscope before next intubation. The Endcap cannot be repositioned or refitted to the endoscope.

10.2. To use another Suture Assembly, return to **step 4.21** for Anchor loading and follow all subsequent steps.

If suturing is complete, proceed to **section 11** for removal of the device.


11. Removal of Device

11.1. Remove any accessory devices.

11.2. Ensure the Needle Body is closed and retract endoscope from patient, making sure the external Sheath is retracted together with the endoscope.

WARNING: Ensure any Actuation Catheter slack in device is removed prior to device removal from patient.

11.3. Insert a surgical blade through the long vertical slot on side of Endcap and cut each Endcap Strap.

NOTE: The vertical slot protects the endoscope from damage during removal of Endcap.

11.4. While looking down on endoscope face, twist Endcap clockwise around endoscope to remove Endcap.

11.5. Unwrap each Sheath Strap and remove Catheter Sheath from endoscope.

11.6. Undo Handle Strap and slide Handle Bracket off endoscope.

12. Device Disposal

After use, the OverStitch instrumentation, as well as any explanted implants, may contain biohazardous substances. These components should be safely disposed as biohazardous waste, in accordance with any applicable hospital, administrative, and/or local government regulations.

English

Use of a biohazardous container with biological hazard symbol is recommended. Untreated biohazardous waste should not be disposed of in the municipal waste system.

13. MRI Information

MR Conditional

Non-clinical testing has demonstrated that the Sutures, Cinches and Anchors (collectively termed Anchoring System) deployed by the OverStitch Endoscopic Suturing System are MR Conditional.

A patient with this Anchoring System can be safely scanned immediately after placement in an MR system meeting the following conditions:

Static Magnetic Field

- Static magnetic field of 1.5T or 3.0T

- Maximum spatial field gradient of 2,000 gauss/cm (20T/m)

- Maximum MR system reported, whole body

averaged specific absorption rate (SAR) of 2 W/kg

Under the scan conditions defined above, the Anchoring System is expected to produce a maximum temperature rise of less than 2°C after 15 minutes of continuous scanning.

In non-clinical testing, the image artifact caused by the Anchoring System extends approximately 10 mm from this device when imaged with a gradient echo pulse sequence and a 3.0T MRI system.

14. Information on Materials and Substances to which Patients Can Be Exposed

Each implant construct, consisting of an Anchor-Suture and Cinch, consists of 1. a polypropylene suture (typically less than 5 cm long) with an anchor (made from 316L stainless steel (0.011g), cobalt chrome alloy (0.006 g)) and 2. a Cinch to hold it the implant in place (made from PEEK (0.020g)). A patient may receive more than one implant construct.

15. Hazardous Substance Information

Contains Hazardous Substance

CMR Statement - The stainless steel and cobalt alloy components within this device contain the following substance(s) (defined as a CMR) (carcinogenic, mutagenic or toxic to reproduction (1A/1B and/or endocrine disrupting in a concentration above 0.1% weight by weight):

Cobalt (CAS No. 7740-48-4; EC No. 231-158-0)

Current scientific evidence supports medical devices manufactured from these cobalt alloys or stainless steels containing cobalt do not cause increased risk of cancer or adverse reproductive effects.

16. Troubleshooting

16.1. Needle Body will not open:

| Cause | Resolution |
|--|--|
| 16.1.1. Needle obstructed: | <i>i. Assess the space in which you are working and maneuver the Sheath and endoscope as a system, straighten the endoscope to the non-Retroflexed position.</i> |
| 16.1.2. Suture movement restricted: | <i>i. If the Anchor is on the Needle Body, ensure the suture is not held proximally near the handle during the opening operation. ii. Transfer the Anchor to the Anchor Exchange. Open the Needle Body. Slowly retract the Anchor Exchange proximally and then advance the Needle Body distally to free the Suture.</i> |
| 16.1.3. Sheath or Actuation Catheter looped or kinked: | <i>i. Check to ensure that the Sheath and Actuation Catheter running down the outside of the endoscope are not looped or kinked. Straighten the endoscope to the non-Retroflexed position. Advance the endoscope forward and slightly pull any slack from the Sheath and Actuation Catheter proximally until minimal resistance is felt. Grasp both the Sheath and endoscope and adjust by advancing and retracting as a system.</i> |
| 16.1.4. Foreign body obstruction: | <i>i. Remove the Anchor Exchange only: a. Transfer the Anchor to the Needle Body and remove the Anchor Exchange from the device. b. Load a grasper through the primary channel and push the Needle Body open. ii. Remove the Anchor and Anchor Exchange, cutting the Suture if necessary: a. Transfer the Anchor to the Needle Body and remove the Anchor Exchange from the device. b. Through either channel use an appropriate accessory to cut the Suture. c. Use an accessory to push open the Needle Body. d. Use appropriate means to remove the cut Suture. iii. Once standard endoscopic techniques have been exhausted, utilize laparoscopic techniques to remove the device.</i> |

16.2. Needle Body will not close:

| Cause | Resolution |
|-------------------------------------|---|
| 16.2.1. General obstruction: | <i>i. Follow steps 13.1.1, 13.1.2, 13.1.3 above (Needle Body will not open). ii. Ensure the Needle Driver handle is locked closed and: a. Stretch the Actuation Catheter to change the effective length of the Needle Body drive cable. b. Remove the Anchor Exchange and use graspers (through the device primary channel) to grab the Needle Body. Cut the Suture if necessary.</i> |
| 16.2.2. Needle Driver cable broken: | <i>i. Deploy the Anchor and Cinch. Advance grasper through endoscope and secure Needle Driver. Pull Needle Driver closed while removing the device. If using an overtube, advance the overtube as far distally as possible and withdraw the endoscope and device into the overtube, using the distal tip of the overtube to close the Needle Body.</i> |

16.3. Anchor Exchange will not exchange:

| Cause | Resolution |
|---|---|
| 16.3.1. Anchor Exchange will not install Anchor onto the Needle Body: | <i>i. Ensure there is sufficient Suture slack and the Suture outside the device is not entangled. ii. Ensure Anchor Exchange is properly positioned in Endcap of Needle Driver. iii. If Anchor and Suture are through tissue, either drop Anchor and deploy Cinch per Section 9 of Instructions for Use or drop Anchor and use a suitable accessory to cut and remove the Suture. iv. If Anchor and Suture are not through tissue, close handle of Needle Driver. Remove device. Replace Anchor and/or Anchor Exchange.</i> |
| 16.3.2. Anchor Exchange will not release an Anchor: | <i>i. Ensure there is sufficient Suture slack and the Suture outside the device is not entangled. ii. Ensure the Anchor Exchange Release Button is FULLY depressed while retracting Anchor Exchange. iii. Reduce the articulation / tortuosity of the endoscope (if possible) and try to release the Anchor. iv. Use accessories compatible with device secondary channel or endoscope working channel to cut and remove the Suture. v. Replace the Anchor Exchange.</i> |
| 16.3.3 Anchor Exchange will not retrieve Anchor from Needle Body: | <i>i. Ensure there is sufficient Suture slack and Needle Driver handle is in closed position. ii. Use 'pencil-grip' to advance Anchor Exchange until the Anchor is engaged and resistance is felt. Retract Anchor Exchange to acquire Anchor. iii. If the Anchor cannot be retrieved, replace the Anchor Exchange. Alternatively, use a suitable accessory to cut and remove the Suture. iv. Replace the Anchor and resume suturing per Section 8 of Instructions for Use.</i> |

16.4. Cinch does not cut the Suture when fired:

| Cause | Resolution |
|-------------------------|---|
| 16.4.1. Suture not cut: | <i>i. Use a suitable accessory through the device secondary channel or endoscope working channel to cut the Suture and remove the Cinch. ii. Use standard endoscopic techniques to remove the cut Suture.</i> |

English

16.5. Inadvertent drop of Anchor:

| Cause | Resolution |
|---|--|
| 16.5.1. Anchor Exchange button depressed out of sequence: | <p>i. Retrieve the Anchor as a foreign body, or follow the procedure to Cinch in place. If the Anchor is dropped inside the device channel, use the Anchor Exchange or 3.2 mm compatible grasper to push the Anchor out through the device.</p> <p>Do not attempt to pull the Anchor back through the device as it may get stuck in the channel or y-junction at the handle.</p> |

16.6. Suture entanglement:

| Cause | Resolution |
|---------------------------------------|--|
| 16.6.1. Suture outside field of view: | i. Close the Needle Driver Handle and manipulate the endoscope back to release. |
| 16.6.2. Suture behind Tissue Guard: | <p>i. Slightly close the Needle Body while retracting the endoscope.</p> <p>ii. If necessary, transfer the Anchor to the Anchor Exchange.</p> <p>iii. Open the Needle Body and advance the Anchor Exchange beyond the Endcap to push the Suture free.</p> |
| 16.6.3. Suture twisted: | <p>i. If the Suture is twisted, move the endoscope and transfer the Anchor between the Needle Body and Anchor Exchange, on the opposite side of the suture thread, as needed to untwist.</p> <p>NOTE: If the twist was noted immediately after articulation of the endoscope, first try to articulate in reverse order to remove.</p> <p>ii. If the Anchor has been deployed, use the Cinch to push and guide the Suture free.</p> |

16.7. Endcap detachment from endoscope:

| Cause | Resolution |
|------------------------------|---|
| 16.7.1. Detached during use: | <p>i. Close the Needle Body, remove any slack in Actuation Catheter and slowly remove the device from the patient.</p> <p>If using an overtube, advance the overtube as far distally as possible and withdraw the endoscope and device into the overtube, maintaining tight hold of the Sheath to keep the Endcap flush with the endoscope.</p> |

16.8. Helix does not unscrew:

| Cause | Resolution |
|--------------------------------|---|
| 16.8.1. Helix stuck in tissue: | <p>i. Use a suitable accessory through the device primary channel or endoscope working channel to apply counter traction to the tissue around the Helix, and pull the Helix free.</p> <p>ii. Once endoscopic techniques have been exhausted, utilize laparoscopic techniques to remove the Helix.</p> |

16.9. Endcap Strap end is not seated in Endcap prior to cutting Endcap Strap Tails:

| Cause | Resolution |
|--|--|
| 16.9.1. Endcap Clip Hinge not fully pushed in during installation: | <p>i. Gently pull on the endoscope to fully recess non-tail end of both Endcap Straps within Endcap.</p> <p>ii. If not discarded, thread Endcap Strap tails through the corresponding channel slots 1 and 2 in Endcap Clip. Re-fit Endcap Clip with a clockwise rotation onto Endcap (opposite of section 4.16) and repeat assembly steps.</p> <p>iii. If Endcap Clip is not available, use left thumb in place of the Endcap Clip to hold the non-tail end of each Endcap Strap within Endcap and repeat Endcap Strap tightening steps.</p> |

16.10. Difficulty passing installed device through overtube:

| Cause | Resolution |
|---|---|
| 16.10.1. Insufficient lubricant applied: | i. Fully lubricate device and inside of overtube. Apply a twisting motion of the endoscope with the installed device during intubating and extubating. |
| 16.10.2. Endoscope and/or overtube sizing not compatible: | i. Exchange overtube selected with a larger sized overtube. |
| 16.10.3. Catheters are not flat against scope: | i. Maneuver the scope and device to lay flat against each other. Apply a twisting motion of the endoscope with the installed device during intubating and extubating. |
| 16.10.4. Endcap Straps not trimmed flush with Endcap: | i. Remove endoscope from overtube. Trim Endcap Straps flush with Endcap per Step 4.18. |

DISCLAIMER OF WARRANTY AND LIMITATION OF REMEDY

THERE IS NO EXPRESS OR IMPLIED WARRANTY, INCLUDING WITHOUT LIMITATION ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE, ON THE APOLLO ENDOSURGERY, INC. PRODUCT(S) DESCRIBED IN THIS PUBLICATION. TO THE FULLEST EXTENT PERMITTED BY APPLICABLE LAW, APOLLO ENDOSURGERY, INC. DISCLAIMS ALL LIABILITY FOR ANY INDIRECT, SPECIAL, INCIDENTAL, OR CONSEQUENTIAL DAMAGES, REGARDLESS OF WHETHER SUCH LIABILITY IS BASED ON CONTRACT, TORT, NEGLIGENCE, STRICT LIABILITY, PRODUCTS LIABILITY OR OTHERWISE. THE SOLE AND ENTIRE MAXIMUM LIABILITY OF APOLLO ENDOSURGERY, INC., FOR ANY REASON, AND BUYER'S SOLE AND EXCLUSIVE REMEDY FOR ANY CAUSE WHATSOEVER, SHALL BE LIMITED TO THE AMOUNT PAID BY THE CUSTOMER FOR THE PARTICULAR ITEMS PURCHASED. NO PERSON HAS THE AUTHORITY TO BIND APOLLO ENDOSURGERY, INC. TO ANY REPRESENTATION OR WARRANTY EXCEPT AS SPECIFICALLY SET FORTH HEREIN. DESCRIPTIONS OR SPECIFICATIONS IN APOLLO ENDOSURGERY, INC PRINTED MATTER, INCLUDING THIS PUBLICATION, ARE MEANT SOLELY TO GENERALLY DESCRIBE THE PRODUCT AT THE TIME OF MANUFACTURE AND DO NOT CONSTITUTE ANY EXPRESS WARRANTIES OR RECOMMENDATIONS FOR USE OF THE PRODUCT IN SPECIFIC CIRCUMSTANCES. APOLLO ENDOSURGERY, INC. EXPRESSLY DISCLAIMS ANY AND ALL LIABILITY, INCLUDING ALL LIABILITY FOR ANY DIRECT, INDIRECT, SPECIAL, INCIDENTAL, OR CONSEQUENTIAL DAMAGES, RESULTING FROM REUSE OF THE PRODUCT.

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All rights reserved.