# OverStitch 2-0 Polypropylene Suture

<b>English</b> [Endoscopic Suturing System] (ESS) INSTRUCTIONS FOR USE	<b>page 2</b>	<b>Norsk</b> [Endoskopisk sutureringssystem] (ESS) BRUKSANVISNING	<b>side 18</b>	<b>Srpski</b> [Endoskopski sistem za ušivanje] (ESS) UPUTSTVO ZA UPOTREBU	<b>strana 34</b>
<b>Español</b> [Endoscopic Suturing System] (ESS) INSTRUCCIONES PARA EL USO	<b>página 4</b>	<b>Suomi</b> [Endoskooppien ommeljärjestelmä] (ESS) KÄYTTÖOHJEET	<b> sivu 20</b>	<b>Slovenčina</b> [Endoskopijski šijalni sistem] (ESS) NAVODILA ZA UPORABO	<b>strana 36</b>
<b>Deutsch</b> [Endoskopisches Nahtsystem] (ESS) GEBRAUCHSANWEISUNG	<b>Seite 6</b>	<b>Svenska</b> [Endoskopiskt sutureringssystem] (ESS) BRUKSANVISNING	<b>sidan 22</b>	<b>Slovenščina</b> [Endoskopski sistem za šivanje] (ESS) NAVODILA ZA UPORABO	<b>stran 38</b>
<b>Français</b> [Système de suture endoscopique] (ESS) MODE D'EMPLOI	<b>page 8</b>	<b>Polski</b> [System wzziernikowego zakładania szwów] (ESS) INSTRUKCJA UŻYCIA	<b>strona 24</b>	<b>Română</b> [Sistem de suturare endoscopică] (ESS) INSTRUCȚIUNI DE UTILIZARE	<b>pagina 40</b>
<b>Italiano</b> [Sistema di sutura endoscopica] (ESS) ISTRUZIONI PER L'USO	<b>pagina 10</b>	<b>Čeština</b> [Endoskopický šicí systém] (ESS) NÁVOD K POUŽITÍ	<b>strana 26</b>	<b>Ελληνικά</b> [Σύστημα ενδοσκοπικής συρραφής] (ESS) ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ	<b>σελίδα 42</b>
<b>Português</b> [Endoscopic Suturing System] (ESS) INSTRUÇÕES PARA O USO	<b>página 12</b>	<b>Magyar</b> [Endoszkópos varrórendszer] (ESS) HASZNÁLATI UTASÍTÁS	<b>28. oldal</b>	<b>Türkçe</b> [Endoskopik Sütür Sistemi] (ESS) KULLANMA TALİMATLARI	<b>sayfa 44</b>
<b>Nederlands</b> [Endoscopisch hechtsysteem] (ESS) GEBRUIKSAANWIJZING	<b>pagina 14</b>	<b>Български</b> [Система за ендоскопско зашиване] (ESS) ИНСТРУКЦИИ ЗА УПОТРЕБА	<b>страница 30</b>	<b>Latviski</b> [Endoskopiskā šuvju sistēma] (ESS) LIETOŠANAS INSTRUKCIJA	<b>46. lappuse</b>
<b>Dansk</b> [Endoskopisk sutureringssystem] (ESS) BRUGSANVISNING	<b>side 16</b>	<b>Eesti keel</b> [Endoskoopiliste õmbluste süsteem] (ESS) KASUTUSJUHEND	<b>lehekülg 32</b>	<b>Lietuvių k.</b> [Endoskopinio siuvimo sistema] (ESS) NAUDOJIMO INSTRUKCIJA	<b>48 psl.</b>
				<b>Hrvatski</b> [Endoskopski sustav za šivanje] (ESS) UPUTE ZA UPORABU	<b>stranica 50</b>

**Product Code:**

PLY-G02-020-APL

**Caution:**

Please read all instructions prior to use

**STERILE EO**

**SINGLE PATIENT USE  
DISPOSABLE**

**GRF-00555-00R01**

# English

Single use only. Disposable. Do not resterilize.

Patent Pending. **DISCLAIMER OF WARRANTY AND LIMITATION OF REMEDY**

Table of Symbols

Description	Symbol	Description	Symbol	Description	Symbol
Consult Instructions for Use		Manufacturer		Use By	
Do Not Re-use		Single Sterile Barrier System. Sterilized by Ethylene Oxide		Lot Number	
Sterilized Using Ethylene Oxide		Date of manufacture		Do not use if package is damaged	
Reference Number		Do Not Resterilize		Authorized Representative in the European Community	
Medical Device		Caution, contact accompanying documents		Contains hazardous substance (CAS 7440-48-4, Cobalt)	

## Description

The OverStitch 2-0 Polypropylene Suture is comprised of sterile packaged non-absorbable polypropylene suture attached to an implantable anchor manufactured from cobalt chromium and stainless steel. The anchor component is intended to function with the OverStitch Endoscopic Suturing System to perform stitching operations and serves as an anchor to secure suture placement. The suture is manufactured from an isotactic crystalline stereoisomer of polypropylene, a synthetic linear polyolefin. The suture pigment uses CU-Phthalocyanine Blue Dye (below 0.5 WT %) to enhance visibility. The suture material meets the requirements established by the USP.

## Actions

The polypropylene suture material comprising the OverStitch 2-0 Polypropylene Suture elicits minimal acute inflammatory reaction in tissue, which is followed by gradual encapsulation of the suture by fibrous connective tissue. The suture material is not absorbed, nor is any significant change in strength retention known to occur in vivo.

## 1 Intended Use

The OverStitch Endoscopic Suture System is indicated for the endoscopic placement of anchor-sutures. Anchor-sutures can be placed for defect closure (e.g. perforation, ESD/EMR and fistula/leak), stent fixation to reduce migration rates in the upper GI tract, endoscopic sleeve gastropasty, and transoral outlet reduction.

### 1.1 System Operation

For detailed instructions, warnings, precautions and potential adverse events, refer to the instructions for use for the OverStitch or OverStitch Sx Endoscopic Suturing System. These are available at [www.apolloendo.com/dfus](http://www.apolloendo.com/dfus).

### 1.2 Summary of Safety and Clinical Performance

The Summary of Safety and Clinical Performance document, as required by the European Medical Device Regulation, is located at <https://ec.europa.eu/tools/eudamed>.

### 1.3 Information Provided to the Patient

The OverStitch packaging includes a Patient Implant Card and leaflet. The leaflet instructs the medical staff on how to complete the card. The card documents the date of the procedure, patient's name, physician's contact information, device tracking information, and MR Safety information. Patients should receive the completed card and leaflet after their procedure.

### 1.4 Warnings

- Do not re-sterilize. Discard unused sutures.
- Users should be familiar with surgical procedures and techniques involving nonabsorbable sutures before using the OverStitch 2-0 Polypropylene Suture for wound closure, as risk of wound dehiscence may vary with the site of application and the suture material used.
- As any foreign material in the presence of bacterial contamination may enhance bacterial infectivity, acceptable surgical practice must be followed with respect to drainage and closure of infected or contaminated wounds.
- Do not reuse. The reuse of single-use devices can cause cross-contamination and affect the device safety, performance and effectiveness, exposing patients and staff to unnecessary risk. The design and material used are not compatible with conventional cleaning and sterilization procedures.
- Do not use if package is damaged.
- Use by expiry date.
- Keep in cool dry place away from any direct heat.

### 1.5 Precautions

- In handling the suture material, care should be taken to avoid damage from handling.
- Avoid crushing or crimping damage of the anchor due to application of surgical instruments such as forceps or needle holders.
- Use of an OverStitch Suture Cinch is required to secure placement of sutures.

## 1.6 Adverse Events

Adverse effects associated with the use of this device include:

- Wound dehiscence
  - Infection
  - Minimal acute inflammatory tissue reaction, and pain
  - Edema and erythema at the wound site.
  - Broken anchors may result in extended or additional surgeries or residual foreign bodies.
  - Inadvertent anchor sticks with contaminated surgical needles may result in the transmission of bloodborne pathogens.
- NOTE: Any serious incident that has occurred in relation to the device should be reported to Apollo Endosurgery (see contact information at the end of this document) and any appropriate government entity.

### 1.7 Device Disposal

After use, the OverStitch instrumentation, as well as any explanted implants, may contain biohazardous substances. These components should be safely disposed as biohazardous waste, in accordance with any applicable hospital, administrative, and/or local government regulations.

Use of a biohazardous container with biological hazard symbol is recommended. Untreated biohazardous waste should not be disposed of in the municipal waste system.

### 1.8 MR Safety Information

MR Conditional

Non-clinical testing has demonstrated that the sutures deployed by the OverStitch Endoscopic Suturing System are MR Conditional.

A patient with this Anchoring System can be safely scanned immediately after placement in an MR system meeting the following conditions:

#### Static Magnetic Field

- Static magnetic field of 1.5T or 3.0T
  - Maximum spatial field gradient of 2,000 gauss/cm (20 T/m)
  - Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 2 W/kg

Under the scan conditions defined above, the Anchoring System is expected to produce a maximum temperature rise of less than 2° C after 15 minutes of continuous scanning.

In non-clinical testing, the image artifact caused by the Anchoring System extends approximately 10 mm from this device when imaged with a gradient echo pulse sequence and a 3.0 T MRI system.

### 1.9 Information on Materials and Substances to which Patients Can Be Exposed.

Each implant construct, consisting of an Anchor-Suture and Cinch, consists of 1. a polypropylene suture (typically less than 5cm long) with an anchor (made from 316L stainless steel (0.011g), cobalt chrome alloy (0.006 g)) and 2. a Cinch to hold it the implant in place (made from PEEK (0.020g)). A patient may receive more than one implant construct.

### 1.10 Hazardous Substance Information

Contains Hazardous Substance

CMR Statement – The stainless steel and cobalt alloy components within this device contain the following substance(s) defined as a CMR (carcinogenic, mutagenic or toxic to reproduction) 1A/1B and/or endocrine disrupting in a concentration above 0.1% weight by weight:

Cobalt (CAS No. 7740-48-4; EC No. 231-158-0)

Current scientific evidence supports medical devices manufactured from these cobalt alloys or stainless steels containing cobalt do not cause increased risk of cancer or adverse reproductive effects.

### 1.11 Sterility

The OverStitch 2-0 Polypropylene Suture is sterilized by Ethylene Oxide (EO) and is sterile unless package is opened or damaged. Do not re-sterilize or reuse

## DISCLAIMER OF WARRANTY AND LIMITATION OF REMEDY

THERE IS NO EXPRESS OR IMPLIED WARRANTY, INCLUDING WITHOUT LIMITATION ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE, ON THE APOLLO ENDOSURGERY, INC. PRODUCT(S) DESCRIBED IN THIS PUBLICATION. TO THE FULLEST EXTENT PERMITTED BY APPLICABLE LAW, APOLLO ENDOSURGERY, INC. DISCLAIMS ALL LIABILITY FOR ANY INDIRECT, SPECIAL, INCIDENTAL, OR CONSEQUENTIAL DAMAGES, REGARDLESS OF WHETHER SUCH LIABILITY IS BASED ON CONTRACT, TORT, NEGLIGENCE, STRICT LIABILITY, PRODUCTS LIABILITY OR OTHERWISE. THE SOLE AND ENTIRE MAXIMUM LIABILITY OF APOLLO ENDOSURGERY, INC., FOR ANY REASON, AND BUYER'S SOLE AND EXCLUSIVE REMEDY FOR ANY CAUSE WHATSOEVER, SHALL BE LIMITED TO THE AMOUNT PAID BY THE CUSTOMER FOR THE PARTICULAR ITEMS PURCHASED. NO PERSON HAS THE AUTHORITY TO BIND APOLLO ENDOSURGERY, INC. TO ANY REPRESENTATION OR WARRANTY EXCEPT AS SPECIFICALLY SET FORTH HEREIN. DESCRIPTIONS OR SPECIFICATIONS IN APOLLO ENDOSURGERY, INC. PRINTED MATTER, INCLUDING THIS PUBLICATION, ARE MEANT SOLELY TO GENERALLY DESCRIBE THE PRODUCT AT THE TIME OF MANUFACTURE AND DO NOT CONSTITUTE ANY EXPRESS WARRANTIES OR RECOMMENDATIONS FOR USE OF THE PRODUCT IN SPECIFIC CIRCUMSTANCES. APOLLO ENDOSURGERY, INC. EXPRESSLY DISCLAIMS ANY AND ALL LIABILITY, INCLUDING ALL LIABILITY FOR ANY DIRECT, INDIRECT, SPECIAL, INCIDENTAL, OR CONSEQUENTIAL DAMAGES, RESULTING FROM REUSE OF THE PRODUCT.

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All rights reserved.