RHYTHMIA HDx™

Mapping System

Hardware Instructions for Use	2
Instrucciones de uso del hardware	46
Mode d'emploi du matériel	87
Gebrauchsanweisung für die Hardware	129
Istruzioni per l'uso dell'hardware	171
Instructies voor gebruik voor de hardware	212
Instruções de utilização do hardware	252

CONTENTS

1. DEVICE DESCRIPTION	6
1.1 Contents	6
1.1.1 Signal station	6
1.1.2 System software	6
1.1.3 Workstation	6
1.1.4 Accessories	6
1.1.5 Intended applied parts	7
1.2 Operating Principle	7
1.2.1 Continuous mapping	7
1.2.2 Continuous catheter localization and tracking	7
1.3 User Information	7
2. INTENDED USE	8
3. INDICATIONS FOR USE	8
4. CLINICAL BENEFIT STATEMENT	8
5. CONTRAINDICATIONS	8
6. WARNINGS	8
7. PRECAUTIONS	10
7.1 General	10
7.2 Signal station	10
7.3 Intended applied parts	10
7.4 Workstation	10
7.5 Cables	11
7.6 Electrical	11
7.7 Body surface electrodes	11
7.8 Environmental	11
7.9 Magnetic localization system	11
7.10 Cleaning and disinfecting	12
7.11 Repair or replacement	13
7.12 Disposal	13
8. ADVERSE EVENTS	13
8.1 Arrhythmias	13
8.2 Misinterpretation of Data	13
8.3 Electrical Hazards	13
9. CONFORMANCE TO STANDARDS	14

10. HOW SUPPLIED	14
11. OPERATIONAL INSTRUCTIONS	15
11.1 Setup and installation	15
11.2 Major hardware components	15
11.3 Signal station	15
11.4 SiS front panel description	16
Table 1. SiS front panel items	17
Table 2. ☼ SiS LED indicator states	18
11.5 SiS front panel inputs	18
11.5.1 IntellaMap Orion family of catheters	18
11.5.2 Breakout box	19
11.5.3 ECG input cables	20
11.5.4 Ablation connection box	20
11.5.4.1 Ablation connection box for non-force sensing catheters	21
11.5.4.2 Ablation connection box for force sensing catheters	22
11.5.5 Stimulator jumper cables	23
11.5.6 Location reference back patch	23
11.6 SiS front panel outputs	24
11.6.1 Direct connect IC cables	24
11.6.2 Direct connect ECG cables	25
11.6.3 Rhythmia IC Out 72 Cable	25
11.6.4 ECG output box	26
11.7 SiS rear panel description	26
Table 3. SiS rear panel items	27
11.8 SiS rear panel connections	29
11.8.1 Localization generator input port	29
11.8.2 Workstation data connection	30
11.8.3 Foot switch	30
11.8.4 Electric power input	30
11.9 Workstation description	31
11.10 Workstation system requirements	31
11.11 Workstation setup	31
12. PREPARING FOR A MAPPING STUDY	32
12.1 Prior to patient arrival	32
12.1.1 Set up the hardware and make the cable connections	32
12.1.2 Prepare the workstation	33

12.2	Upon patient arrival	33
12.2.1	Make required cable and catheter connections.....	33
12.2.2	Additional connections when using one of the IntellaMap Orion family of catheters.....	33
12.3	Important considerations during a study.....	34
12.4	Ending the study	34
13.	INDIVIDUAL CABLES ASSOCIATED WITH SYSTEM	35
14.	DIRECT CONNECT CHANNEL MAPPING.....	37
15.	PREPARING INTELLAMAP ORION FAMILY OF CATHETERS FOR USE	38
16.	CLEANING GUIDELINES	38
16.1	Routine cleaning and disinfecting	38
16.2	Decontaminating the equipment before shipping.....	38
17.	SIS TROUBLESHOOTING.....	39
	Table 4. Status Indicator States.....	39
18.	ELECTROMAGNETIC EMISSIONS AND IMMUNITY.....	40
	Table 5. Guidelines and manufacturer's declaration – electromagnetic emissions.....	40
	Table 6. Guidelines and manufacturer's declaration – electromagnetic immunity.....	41
19.	TECHNICAL SPECIFICATIONS.....	43
20.	SERVICE INFORMATION.....	44
21.	WARRANTY.....	44
22.	CONTACTS.....	44
23.	SOFTWARE LICENSE.....	45
24.	SYMBOL DEFINITIONS.....	45

Остаряла версия. Да не се използва.
Zastaralá verze. Nepoužívat.
Forældet version. Må ikke anvendes.
Version überholt. Nicht verwenden.
Aegunud versioon. Ärge kasutage.
Παλιά έκδοση. Μην την χρησιμοποιείτε.
Outdated version. Do not use.
Version obsolete. No utilizar.
Version périmée. Ne pas utiliser.
Zastarjela verzija. Nemojte upotrebljavati.
Úreлт útгáфа. Notið ekki.
Versione obsoleta. Non utilizzare.
Pasenusi versija. Neizmantot.
Elavult verzió. Ne használja!
Dit is een verouderde versie. Niet gebruiken.
Wersja przeterminowana. Nie używać.
Versão obsoleta. Não utilize.
Versiune expirată. A nu se utiliza.
Zastarana verzija. Nepoužívat.
Vanhentunut versio. Älä käytä.
Föråldrad version. Använd ej.
Güncel olmayan sürüm. Kullanmayın.

Rx ONLY

Caution: Federal Law (USA) restricts this device to sale by or on the order of a physician.

Note: The equipment documented in this manual (signal station, accessory boxes and cables, localization generator, and workstation) is supplied non-sterile and cannot be sterilized. The equipment is intended for multi-patient reuse.

1. DEVICE DESCRIPTION

The RHYTHMIA HDx Mapping System (the system) is a 3D mapping and navigation system used in electrophysiological (EP) procedures.

1.1 Contents

The package contents differ based on geography and items ordered. The following list includes all components that comprise the full RHYTHMIA HDx Mapping System.

1.1.1 Signal station

The signal station (SiS) accepts signals from the intracardiac catheters and ECG electrodes used during EP procedures. It amplifies and digitizes these signals and transfers them to the workstation for real-time processing and display. The SiS also supports catheter localization/tracking and diagnostic stimulation.

1.1.2 System software

Rhythmia software runs on the workstation. It processes data received from the SiS and provides a user interface for system operation. It also performs the following major functions:

- ECG and intracardiac signal display
- Catheter localization and tracking
- 3D mapping and visualization
- Diagnostic stimulation routing

Refer to RHYTHMIA HDx Mapping System Software Instructions for Use for further information.

1.1.3 Workstation

The workstation comprises both computer hardware (e.g. computer, monitor, keyboard, mouse, and power cords) and system software. In addition to operating the system software, the workstation can also store, retrieve, and export study data.

1.1.4 Accessories

- SiS power supply
- ECG input and output cables, AAMI and IEC variants
- Breakout boxes and connecting cables
- Intracardiac output cables
- Ablation connection boxes
- Stimulator input jumper cables
- Patch cable
- Localization generator and connecting cable
- Equipotential cable

- Foot switch
 - Fiber optic data cable
 - Workstation computer, monitor, and power cords
 - Workstation ablation data peripherals

1.1.5 Intended applied parts

The following single-use applied parts are intended for use with the system but are not included with the system:

- EP catheters, including the IntellaMap Orion™ family of catheters
- Location reference patch kit

Note: Carefully read respective device instructions for use before using these parts in an EP study.

1.2 Operating Principle

The RHYTHMIA HDx Mapping System (the system) is a 3D mapping and navigation system used in electrophysiological (EP) procedures. The system employs two mechanisms to perform 3D mapping and navigation: (a.) continuous mapping based on patient cardiac signals acquired from intracardiac catheters and surface ECG electrodes, and (b.) continuous localization of magnetic and impedance tracked catheters. These mechanisms are described in greater detail below.

1.2.1 Continuous mapping

The continuous mapping feature employs user-determined beat acceptance criteria to determine which beats are included into the map. As the user traverses the mapping catheter around the chamber, the software will continuously add or reject beats to the map based on conformance to the beat acceptance criteria. The maps are displayed as color coded maps.

1.2.2 Continuous catheter localization and tracking

The system enables the tracking of catheters through the use of magnetic- and impedance-based localization technologies.

Magnetic-based localization uses a magnetic location sensor embedded in magnetically tracked catheters to measure magnetic fields generated by a localization generator located beneath the patient table. These magnetic field sensor readings are passed to the system software which in turn uses the data to determine and display the position of the catheter.

Impedance-based localization operates by driving small currents between multiple surface electrodes and measuring the voltage on each electrode of a tracked impedance catheter. These voltage values are passed to the system software which in turn uses the data to determine and display the position of the catheter.

1.3 User Information

Operation of the system is restricted to licensed medical practitioners and Boston Scientific mapping specialists. Only users who meet these criteria should use the system.

2. INTENDED USE

The RHYTHMIA HDx Mapping System (the system) is a 3D mapping and navigation system used in EP procedures. The SiS and related accessories provide data connection pathways for external input/output devices (e.g. catheters and recording systems) and serve as the data conduit to the system workstation and software.

3. INDICATIONS FOR USE

The RHYTHMIA HDx Mapping System and accessories are indicated for catheter-based atrial and ventricular mapping. The mapping system allows real-time visualization of intracardiac catheters as well as display of cardiac maps in a number of different formats. The acquired patient signals, including body surface ECG and intracardiac electrograms, may also be recorded and displayed on the system's display screen.

4. CLINICAL BENEFIT STATEMENT

The RHYTHMIA HDx Mapping System is an effective diagnostic tool for cardiac electrophysiology (EP) procedures. When used with the IntellaMap Orion mapping catheters or the IntellaNav™ ablation catheters, the RHYTHMIA HDx Mapping System provides real-time visualization of intracardiac catheters and a display of cardiac maps in selected formats in minimally invasive procedures to assist the physician in identifying the origin of the arrhythmia within the heart chamber. Acquisition of the 3-D electroanatomical maps and other patient information such as surface ECG and intracardiac electrograms on a display screen provides diagnostic information to the physician for an overall clinical benefit to identify and treat cardiac arrhythmias. Conversely, failure to treat a cardiac arrhythmia could result in symptoms such as shortness of breath, palpitations, dizziness, syncope, chest pain, stroke or sudden cardiac death.

5. CONTRAINDICATIONS

None known.

6. WARNINGS

The RHYTHMIA HDx Mapping System (the system) is intended for use with other medical devices in an EP laboratory. Carefully read the instructions for use (IFU) documents for every medical device that will be used during a study, prior to each study. Observe all contraindications, warnings, and cautions. Failure to do so may result in user harm, patient illness, injury, or death.

- Carefully read this entire document and all other product IFUs before beginning the mapping study. Be sure to fully understand and consistently follow all warnings, cautions, and instructions. Failure to properly follow the instructions may cause equipment damage, system malfunction, or harm to the patient or user.
- Diagnosis and treatment of cardiac arrhythmias using the system in conjunction with radio frequency (RF) ablation and other medical devices may pose a risk of adverse events. Adverse events (e.g. cardiac perforation, new arrhythmias, exacerbation of existing arrhythmias) may require additional intervention.
- Do not use the system to route life-sustaining pacing signals. Only diagnostic stimulation signals (e.g. induction) may be routed through the system.
- Always use the controls on the external stimulator to start or stop stimulation. The system only routes externally generated-and-controlled stimulation signals to the selected electrode and channel.
- In the event that routing of the stimulation signal through the mapping system software fails, direct stimulation may be required. Connect the stimulator jumper cables to any single pair of the direct stimulation ports located above input ports **M**, **A**, **B** or **ABL**. The direct stimulation ports above input ports **M**, **A**, and **B** connect an external stimulator to channels 61 and 62 of the connected breakout box. The direct stimulation ports above the **ABL** port connect to channels 1 and 2 of the ablation catheter.
- Do not condition an IntellaMap Orion mapping catheter when it is in contact with the patient, either external or indwelling.

- Do not operate the system near flammable anesthetics.
- All devices that are connected to system hardware must independently meet IEC 60601-1 requirements as well as any other relevant safety standards. The combined hardware configuration must also meet IEC 60601-1 safety standards. The use of system hardware with accessories and devices that do not comply with relevant standards may reduce the safety of the system, cause equipment damage or system malfunction, or harm to the patient or user.
- System hardware must be connected solely to a functional, properly-tested supply main with protective ground (earth). Do not use extension cords or adapters for ungrounded outlets. The use of a faulty or ungrounded supply main increases the risk of electrical shock and system malfunction.
- Use only Maestro™, IBI™, Ampere™, SMARTABLATE™, EP-Shuttle™, or INTELLAGEN™ RF ablation generators with the system. Do not use the system with other RF ablation generators. Compatibility with other RF ablation generators has not been demonstrated.
- Carefully read respective RF ablation generator IFUs before beginning the mapping study. Do not exceed power limits set by the manufacturer.
- Do not simultaneously connect more than one ablation catheter to the ablation connection boxes. Doing so may lead to harm to the patient.
- To reduce the risk of electric shock or equipment damage, do not clean system hardware when it is plugged in, turned on, or connected to a patient. Cleaning the system while it is in use and connected to a power source may cause an electrical shock that could cause injury or death to the patient or user.
- To reduce the risk of electric shock, assure that any ECG cables and electrodes are not in contact with any other conductive parts, including ground.
- To reduce the risk of electric shock during defibrillation, assure that the exposed connector tips on the ECG output box are enclosed at all times by the non-conductive protective cover that is built onto the ECG output boxes. Do not use the ECG output box if the protective cover is damaged.
- The system generates electrical impedance fields as part of its normal operation. Do not use other systems that also generate electrical impedance fields in the same procedure, as this may interfere with the system's normal operation and reduce the quality of catheter localization, and signals.
- Do not operate the localization generator within 200 mm of an implanted CIED (cardiac implantable electronic device.) Doing so may affect CIED pacing, temporarily suspend tachycardia therapy delivery, or lead to patient discomfort.

7. PRECAUTIONS

7.1 General

- Carefully inspect all system components prior to using system hardware. Do not use any component that shows evidence of being damaged or defective.
- Do not drop system components or subject them to extreme shock. Dropping components or forcefully hitting them against hard objects may damage components and cause system malfunction. Contact Boston Scientific Support for device repair or replacement.
- To minimize the risk of dropping and damaging equipment, use care when inserting/ removing components. If needed, use two people to attach or remove devices such as the localization generator.
- Refer to Table 5 and Table 6 for EMC interference information.

7.2 Signal station

- Do not position the signal station (SiS) in a manner that would make it difficult to disconnect the power supply from the power mains. Disconnect the power cord if the SiS needs to be isolated from the power mains.
- Connect the equipotential port on the rear panel of the SiS to an equipotential receipt box prior to using system hardware to minimize the risk of electric shock. This connection should remain connected at all times.
- Use only the SiS power supply and power supply cable provided by Boston Scientific with the system. Using another power supply or power supply cable may damage the SiS.
- Do not connect or disconnect the SiS power supply while the unit is turned on. This will minimize the risk of equipment damage.
- Disconnect SiS from its power supply prior to cleaning in order to minimize the risk of electric shock.
- Do not place containers of water or other liquids directly on or near the SiS or other system components. This reduces the risk of electric shock and/or damage to the SiS.
- Do not block the air vent on the SiS during use. Blocking the air vent can cause the SiS to overheat, which may affect system operation.
- Do not stack other devices or equipment on top of the SiS.
- Always use two hands to carry the SiS.
- Only place the SiS and SiS-related accessories on flat, stable surfaces. This will minimize the risk of dropped or toppled equipment.

7.3 Intended applied parts

- Do not use the IntellaNav XP or IntellaNav MiFi XP ablation catheters with the EP-Shuttle ablation connection box.

7.4 Workstation

- Do not place containers of water or other liquids directly on or near the workstation or other system components. This reduces the risk of electric shock and/or damage to the workstation.
- Use only a flat, stable surface to hold or transport the workstation and workstation-related accessories. This will minimize the risk of dropped or toppled equipment.

7.5 Cables

- Use only the ECG cables supplied by Boston Scientific for use with the system. ECG cables provided by Boston Scientific are designed and tested to protect the SiS from defibrillation energy. Using other ECG cables may cause damage to the system hardware.
- Inspect all external connections and cables before using the SiS and secure any connections that are loose. Loose connections may impact the accuracy of mapping results.
- Do not use excessive force when connecting or disconnecting cable connectors. Excessive force can damage the connectors, which may cause system malfunction.
- Do not kink or bend cables. Kinks and sharp bends can damage the cables, which may cause system malfunction.
- Store unused system cables in a clean, dry, and secure location consistent with storage guidelines to minimize the risk of damage. For specific storage guidelines, see section 19.

7.6 Electrical

- Do not use ungrounded electrical outlets to power any system components. Do not use extension cords or adapters for ungrounded outlets. Using ungrounded outlets, extension cords, or adapters may cause equipment damage, system failure or malfunction.

7.7 Body surface electrodes

- Use care when attaching the body surface electrodes to lead connectors. To minimize the risk of electric shock, make sure that electrodes and lead connectors do not contact one another or contact ground.
- Properly prepare the skin prior to attaching the electrodes to prevent receiving low quality signals from body surface electrodes. Do not use excessive gel as this may lead to signal crossover between electrodes.
- To minimize signal interference, route the surface ECG cables across the torso instead of alongside it.

7.8 Environmental

- Avoid exposing system hardware to excessive moisture, heat or cold. Using the system in environmental conditions that exceed recommendations may affect system operation.
- When connecting or disconnecting system cables, protect the cable connectors from water or moisture. Wet connectors may affect system operation.
- Do not immerse any cable connectors in water or liquid. Immersion in water or liquid may damage connectors, which may cause system malfunction.
- Always follow guidelines for equipment storage and transport. Storage or transport in extreme environmental conditions can damage system components. For specific storage and handling guidelines, see "Technical Specifications"

7.9 Magnetic localization system

- Manually disabling the localization generator disables all catheter visualization and localization capabilities, including impedance tracking.
- Do not place cables used with the system hardware within 30 mm of the localization generator cable. Inaccurate tracking or "noisy" signals may occur if these cables are within 30 mm distance between each other, particularly if they are parallel.
- Do not coil the localization generator cable. Doing so can disturb the magnetic field of the localization generator, which may lead to inaccurate tracking.

- Do not use the magnetic localization system in the presence of other magnetic fields or large ferrous metal objects. Doing so may lead to inaccurate tracking.
- Do not connect or disconnect the localization generator from the SiS while the unit is turned on. This will minimize the risk of equipment damage.
- The localization generator may interfere with other systems that use magnetic-field-based technology. Consult the supplier of these systems before using them in the presence of the magnetic localization system.
- The localization generator may interfere with fluoroscopy or other imaging modalities. Consult the supplier of these systems before using them in the presence of the magnetic localization system.
- The localization generator may interfere with implanted cardiac implantable electronic devices (CIEDs). When mapping a patient with such a device, consider interrogating the device pre- and post-procedure. This will identify any changes in programmed parameters which could then be corrected before transferring the patient from the procedure room. Consult the CIED manufacturer instructions for additional information.
- If it becomes necessary to interrogate or program an implanted CIED while using the system, temporarily turn off the localization generator by using the on-screen button located on the annotating and editing maps toolbar.

7.10 Cleaning and disinfecting

- Do not immerse system components in water, cleaning solutions, or liquid. Prevent fluids from entering air vents. Ensure connectors stay dry. Failure to follow cleaning guidelines may cause equipment damage or system malfunction, and may also void any warranties or service agreements.
- To avoid equipment damage and malfunction, do not attempt to sterilize equipment that is provided non-sterile.
- To avoid equipment damage and malfunction, do not insert anything (e.g., cotton swabs or pins) into cable connectors or equipment ports or openings.
- To minimize the risk of electrical shock, prior to cleaning any electrically powered equipment, be sure to first turn off the device and then unplug it from power.
- Do not attempt to clean system components during system operation. Cleaning equipment during use increases the risk of electrical shock, system malfunction, and device drop.
- Do not use the ECG cables or any other cables or system components if they are soiled or contaminated with infectious, or potentially infectious, materials. Using soiled or contaminated items increases the risk of patients acquiring serious infections or contaminating other patients or users. Soiled or contaminated cables and equipment must be removed from use and either cleaned according to established facility protocol procedures or replaced.
- Always clean multiple-use equipment according to established facility protocol prior to each reuse.
- Do not use disinfectants such as glutaraldehyde or hydrogen peroxide to clean system components.
- Do not use solvents such as acetone to clean system components.

7.11 Repair or replacement

- Use only the equipment, supplies, and accessories supplied or recommended by Boston Scientific. The use of the other equipment, supplies, and accessories may cause equipment damage or system malfunction.
- Do not attempt to repair, modify, or open any part of system hardware. Repair attempts by untrained, unauthorized individuals may cause user injury, equipment damage, or system malfunction. Contact Boston Scientific Support for device repair or replacement.

7.12 Disposal

All external and accessible surfaces of this system should be cleaned and disinfected per the disinfection instructions in 16.2. Include any common detachable cables (power cord, video cables, accessories, etc.). Do not dispose of by incineration, burial or placement into common waste stream. System should be safely disposed of in accordance with hospital, administrative, and/or local government policy.

8. ADVERSE EVENTS

Any potential clinical complications are in large part expected to be related to the accessory diagnostic or ablation catheters that are used with the system, rather than the system itself. In order to identify potential adverse events, the user is instructed to read pertinent instructions for use documents associated with the catheters and ablation generators that will be employed during a mapping session.

As with other mapping systems, the RHYTHMIA HDx Mapping System can be incidentally associated with minor or major clinical complications intrinsic to intracardiac procedures. Potential adverse events associated with the use of the system include, but are not limited to, the following:

8.1 Arrhythmias

Due to the programmed electrical stimulation performed during EP diagnostic procedures and catheter manipulations, patients undergoing EP procedures are at potential risk of arrhythmias. The patient may experience discomfort from the rapid pacing and/or the initiation of an arrhythmia. While the system has no active role in RF ablation, a risk does exist that the effectiveness of an RF ablation procedure could be suboptimal and cause the targeted arrhythmia to reoccur.

8.2 Misinterpretation of Data

Localization

Poor catheter localization may lead to clinical data misinterpretation and the potential of resultant patient injury. To ensure correct clinical decisions, the physician must use fluoroscopy, ultrasound, pace mapping or other visualization techniques to verify 3-D mapping results and catheter position.

Incorrect Force Measurements

Incorrect force measurement displayed or misinterpretation of the force displayed may lead to the user to apply more force than desired during mapping or ablating. User must observe any system messages displayed. If User Operator applies more force than desired during mapping or ablating, myocardial perforation, myocardial contusion, or myocardial injury could result.

8.3 Electrical Hazards

With any electrical system there is a potential risk of electrical shock to the user, patient, and service representative.

Note: In the event that a serious incident occurred in relation to the device, including all patient deaths for procedures where the BSC product was used, the event should be reported to BSC and the competent authority of the Member State in which the user and/or patient is established.

9. CONFORMANCE TO STANDARDS

The system hardware components meet the following standards

- IEC 60601-1:2005+A1:2012
- IEC 60601-1-6:2010+A1:2013
- IEC 62366:2007+A1:2014
- CISPR 11:2009+A1:2010, Group 1, Class A
- IEC 60601-1-2:2014
- IEC 60601-2-27:2011 CORR 1:2012
- EN 55011:2009+A1:2010, Group 1, Class A

10. HOW SUPPLIED

The system is packaged and supplied in a series of shipping cartons handled by common carrier.

Do not use if any of the cartons have been damaged or unintentionally opened before use.

Do not use if labeling is incomplete or illegible.

Остаряла версия. Да не се използва.
Zastaralá verze. Nepoužívat.
Forældet version. Må ikke anvendes.
Version überholt. Nicht verwenden.
Aegunud versioon. Ärge kasutage.
Παλιά έκδοση. Μην την χρησιμοποιείτε.
Outdated version. Do not use.
Version obsolete. No utilizar.
Zastarjela verzija. Nemojte upotrebljavati.
Úreлт útgáfa. Notið ekki.
Versione obsoleta. Ne pas utiliser.
Novecojsi versija. Non utilizzate.
Zastarjela verzija. Neizmantot.
Úreлт útgáfa. Notið ekki.
Pasenusi versija. Nenaudokite.
Elavult versio. Ne használja!
Dit is een verouderde versie. Niet gebruiken.
Utdatert versjon. Skal ikke brukes.
Wersja przeterminowana. Nie używać.
Versão obsoleta. Não utilize.
Zastarana verzija. Nepoužívat.
Zastarela različica. A nu se utiliza.
Vanhentunut versio. Älä käytä.
Föråldrad version. Använd ej.
Güncel olmayan sürüm. Kullanmayın.

11. OPERATIONAL INSTRUCTIONS

11.1 Setup and installation

Setup and installation of the RHYTHMIA HDx Mapping System (the system) hardware must be performed by an authorized Boston Scientific representative. Only such authorized personnel should open the sealed product packaging.

11.2 Major hardware components

The system is comprised of two major hardware groups:

- The SiS and related accessories, typically placed near the examination table in the EP lab, and
- The workstation and related accessories, typically used in the EP lab control area.

11.3 Signal station

Function

The signal station (SiS) (Figure 1) acquires signals from intracardiac catheters and ECG electrodes used during EP procedures. Intracardiac signals are collected from Boston Scientific magnetically tracked catheters and from other catheters that are connected to the system through breakout boxes (Figure 1). The SiS amplifies and digitizes catheter and ECG signals and provides them to the workstation for real-time processing and display. The SiS can simultaneously generate and collect signals that are used for catheter localization and tracking (Figure 1).

Figure 1. SiS overview

11.4 SiS front panel description

Specific input and output cables connect to dedicated ports on the front panel of the SiS (Figure 2 and Table 1). Ports 1 through 5 in Figure 2 are inputs. Ports 6, 7, 8, and 9 are outputs. The three ports labeled **M**, **A**, and **B** in the upper row are intracardiac signal inputs. The **ECG** input port lies immediately to the right of port B. Each input port has a corresponding output port immediately below it. The eight input and output ports are color-keyed to identifier rings on the cable connectors. The light-blue-ringed umbilical cable used with an IntellaMap Orion mapping catheter can only be connected to input port **M**. The dark-blue-ringed breakout box cable (Figure 3) may be connected to input ports **M**, **A**, and **B**. The gray-ringed ECG input cable (Figure 4) is connected to the gray-ringed **ECG** input port.


Figure 2. SiS front panel

Direct connect IC cables (Figure 10) and the Rhythmia IC Out 72 cable (Figure 12) may be connected to the lower-row output ports **M**, **A**, and **B** in order to transfer the intracardiac input signals from the SiS to a recording system. Either an ECG direct connect cable (Figure 11) or an ECG output box (Figure 13) may be connected to the lower **ECG** port in order to transfer surface ECG data to a recording system.

Other SiS front panel items are described in Table 1 and Table 2.

Table 1. SiS front panel items


	External stimulator direct input	External stimulator signal direct inputs connected to breakout box ports 61 and 62
	Status LED	Indicates SiS readiness See (Table 2) for a description of status LED function. Do not use the SiS if the light does not perform as described—contact Boston Scientific Support.
	External stimulator direct input	External stimulator signal direct input to ablation catheter electrodes 1 and 2.
1	Input port M	<ul style="list-style-type: none"> Collects mapping signals directly from an IntellaMap Orion mapping catheter via the umbilical cable Collects signals from catheters connected to breakout box
2	Input port A	Collects signals from catheters connected to breakout box
3	Input port B	Collects signals from catheters connected to breakout box
4	ECG input	Collects signals from body-surface electrodes connected to ECG cable CAUTION: Use only ECG cables provided with the SiS by Boston Scientific. The ECG trunk is part of the SiS's defibrillation protection.
5	Ablation connection box input	Connection point for ablation connection box attached cable
6	Output port M	Signal output to recording system from port M inputs: IntellaMap Orion mapping catheter or breakout-box-connected catheters
7	Output port A + ABL	Signal output to recording system from catheters connected to port A breakout box as well as the ablation catheter
8	Output port B	Signal output to recording system from catheters connected to port B breakout box
9	ECG output	Provides ECG signals to a recording system through either a direct connect cable or the ECG output box (Figure 13)
10	Patch input	Connection point for patch cable
11	Condition Button	Starts IntellaMap Orion mapping catheter conditioning process
12	Reserved for future use	Not currently used
13	Stimulator input	Connects to external stimulator signal output WARNING: Only IEC 60601 compliance-certified stimulators should be used with the RHYTHMIA HDx Mapping System.
14	Reserved for future use	Not currently used

Table 2. ☀ SiS LED indicator states

Indicator state	Meaning	Action required
 DARK	Power off	None
 FLASHING	SiS hardware initialization in progress	Contact Boston Scientific Support if SiS initialization period takes more than two minutes.
 FLASHING	Hardware initialization complete, SiS ready to communicate with workstation	<ol style="list-style-type: none"> 1. Connect the fiber optic data cable to the SiS and workstation 2. Start a LIVE mapping session on the workstation
 SOLID	RHYTHMIA HDx Mapping System fully operational	None
 SOLID	Error identified	Do not use the SiS. Contact Boston Scientific Support.

Note: Connecting the fiber optic data cable and starting the LIVE mapping session may be performed in any order. The fiber optic data connection may also be made before any of the hardware is powered on.

11.5 SiS front panel inputs

11.5.1 IntellaMap Orion family of catheters

An IntellaMap Orion mapping catheter may be used with the system by connecting its umbilical cable to **IN-M**. This catheter must be conditioned before use, see section 15.

11.5.2 Breakout box

The breakout box and its connecting cable (Figure 3) provide the physical interface between the SiS and diagnostic catheters. One of the connecting cable end fittings is color-keyed with a dark-blue ring; the other is without color key. The color-keyed end connects to the SiS, and the other end connects to the breakout box. The breakout box is comprised of several columns of standard 2 mm touch-proof pin ports as inputs (64 in total) and a single multi-pin output port. Pin 1 (green) on the breakout box connected to port A becomes the input for the system reference electrode.

Multiple catheters or devices can be connected simultaneously through a single breakout box, depending on the procedure and user needs. Each breakout box indicates its connected SiS input port by illuminating the appropriate **M**, **A**, or **B** indicator on the breakout box.

Two of the 64 corresponding connectors may also be used for external stimulator use. Connectors 61 (red) and 62 (black) are hard-wired to the SiS direct stimulator connectors located above (and associated with) input ports **M**, **A**, and **B**. Connectors 61 and 62 are used for stimulation when the breakout box is connected to the system and an external stimulator is connected to the external stimulator direct input ports associated with the particular **M**, **A**, or **B** port in use. These connections can be used to manually connect an intracardiac catheter to the stimulation signal, if desired.

Figure 3. Breakout box and cable

11.5.3 ECG input cables

The ECG cable (Figure 4) is comprised of two cable subsets (the limb and the chest cables) that connect to a third cable (the trunk cable). The complete ECG cable set gathers body surface electrode signals for the SiS. Signals are communicated to the SiS via the **IN-ECG** port on the SiS front panel. ECG cables are supplied in either AAMI or IEC versions.


Figure 4. ECG input cable (AAMI)

11.5.4 Ablation connection box

The system itself does not provide RF ablation energy; external RF generators are used for this purpose.

WARNING: Use only Maestro, IBI, Ampere, SMARTABLATE, EP-Shuttle, or INTELLAGEN RF ablation generators with the system. Do not use the system with other RF ablation generators. Compatibility with other RF ablation generators has not been demonstrated.

An ablation connection box is required when connecting an ablation catheter to the system. The ablation connection box routes intracardiac signals and location information sensed by the ablation catheter to the mapping system and prevents RF energy from affecting catheter localization and other mapping system features. The ablation connection box also passes catheter tip temperature and catheter tip impedance information, as well as RF energy between the RF generator and ablation catheter.

The Ablation connection box connects to the

- RF generator
- Ablation catheter extension cable
- Ablation catheter input port on SiS (item 5 in both Figure 2 and Table 1).

11.5.4.1 Ablation connection box for non-force sensing catheters

There are five connection box variants to connect to different-brand RF generators. See Figure 5 for the generator brand that is compatible with each connection box. The catheter connector ports are detailed in Figure 6. IntellaNav catheters connect to the connector labeled **IntellaNav**. Third party ablation catheters connect to the connector labeled **STANDARD CATHETER**. When using an IntellaTip MiFi™ catheter, the bifurcated catheter extension cable is connected to the connectors labeled **STANDARD CATHETER** and **IntellaTip MiFi** on the ablation connection box. The ends of the bifurcated cable are colored-coded to match the connectors on the connection box.

CAUTION: Do not use the IntellaNav XP, IntellaNav MiFi XP or IntellaNav ST ablation catheters with the EP Shuttle ablation connection box.

WARNING: Do not simultaneously connect more than one ablation catheter to the ablation connection boxes. Doing so may lead to harm to the patient.

Note: The SMARTABLATE Ablation Connection Box utilizes the *SMARTABLATE Ablation Cable* from *SMARTABLATE System* to *CELSIUS Catheter (TC/THR)* to connect to the SMARTABLATE RF Generator.

Figure 5. Ablation connection box variant for non-force sensing catheters


Figure 6. Ablation connection box connectors

11.5.4.2 Ablation connection box for force sensing catheters

The RHYTHMIA HDx Ablation Connection Box - IntellaNav StablePoint and Ablation Connection Box - IntellaNav StablePoint (for INTELLAGEN) enable use of the RF generator and IntellaNav StablePoint ablation catheter with the mapping system. This ablation connection box routes force information sensed by the ablation catheter to the mapping system in addition to the other functions of an ablation connection box.

The ablation connection box (Figure 7) connects to the

- RF generator through the tethered cabled labeled "RF GENERATOR"
- Ablation input port on the mapping system through the tethered cable labeled "SIGNAL STATION"
- IntellaNav StablePoint ablation catheter through the IntellaNav StablePoint catheter cable

Figure 7. Ablation connection box for force sensing catheters

11.5.5 Stimulator jumper cables

The stimulator jumper cables (Figure 8) are used to connect an external stimulator to either the standard stimulator input ports (item 13 in both Figure 2 and Table 1) or to the direct stimulator input ports (indicated by the \square symbol in both Figure 2 and Table 1). The difference between these inputs is as follows:

- The standard stimulator input ports are located at the lower right corner of the SiS. The Rhythmia software provides the means for routing input stimulation signals to specific intracardiac electrodes.
- The direct input ports are located above input ports **M**, **A**, **B** and **ABL**. Signals received at any of the **M A B** direct input ports are routed to ports 61 and 62 on associated breakout boxes. Signals received at the **ABL** direct input port are routed to electrodes 1 and 2 on the ablation catheter via the ablation connection box.


Figure 8. Stimulator input jumper cables

11.5.6 Location reference back patch

The location reference back patch (back patch) is applied in a medial, stable position on the patient's back and serves as a location reference for the tracking system. The back patch has a cable and connector that connect to the patch cable (Figure 9). The patch cable connects the **PATCH** port on the SiS front panel (item 10 in both Figure 2 and Table 1) to the back patch using a dedicated connector.

Figure 9. Patch cable and connection

11.6 SiS front panel outputs

Green color-coding identifies three of the lower-bank connectors as IC output ports for external recording systems. Moving from left-to-right, the three ports are labeled **OUT-M**, **OUT-A+ABL**, and **OUT-B**. Data signals are sent to recording systems using either direct connect cables or the Rhythmia IC Out 72 cable.

11.6.1 Direct connect IC cables

Direct connect IC output cables (Figure 10) are identified by green-colored rings on the connectors that attach to the SiS. These cables are used with the Clearsign™ amplifier and the CardioLab™ external recording system. The two systems employ different cables and connection modalities see section 14.

- **CLEARSIGN amplifier**—The gray Direct Connect Clearsign IC Orion cable is used at the **OUT-M** port when an IntellaMap Orion mapping catheter is used. The black Direct Connect Clearsign IC A/B cable is connected to any of the three output ports when breakout box-connected catheters are used. When the black direct connect cable is connected to the **OUT-A+ABL** port, ablation data is carried on output ports 65-72.
- **CardioLab recording system**—The Direct Connect CardioLab IC cable may be used at any of the three output ports. Note that whenever this cable is used at port **OUT-A+ABL**, any data input on breakout box ports 33-40 is supplanted by ablation data on output ports 33-40.

Figure 10. Direct connect IC output cables

11.6.2 Direct connect ECG cables

Direct connect ECG output cables (Figure 11) are identified by gray-colored rings on the connectors that attach to the SiS and are provided for both the Clearsign amplifier and the CardioLab recording system. Both cables connect to the **OUT-ECG** port (item 9 in both Figure 2 and Table 1).

Figure 11. Direct connect ECG output cables

11.6.3 Rhythmia IC Out 72 Cable

The Rhythmia IC Out 72 cable (Figure 12) is used to provide output signals from any of the three IC output ports to an external recording system that employs pin-block input.


Figure 12. Rhythmia IC Out 72 cable

11.6.4 ECG output box

The ECG output box (Figure 13) routes ECG signals from the SiS to a recording system. The ECG output box cable connects to the **OUT-ECG** port on the SiS front panel (item 9 in both Figure 2 and Table 1). The ECG output box includes a protective, non-conductive cover that protects the user from high voltage during defibrillation. Two versions of the box are available-AAMI and IEC-the difference between them being connector labels and colors.

WARNING: To reduce the risk of electric shock during defibrillation, ensure that the exposed connector tips on the ECG output box are covered at all times with the protective, non-conductive material provided with the ECG output boxes. Do not use the ECG output box if the protective cover is damaged.

Figure 13. ECG output box variants

11.7 SiS rear panel description

SiS rear panel connectors (Figure 14) attach to the SiS's dedicated power supply, the workstation, and several SiS accessories. Rear panel items are described in Table 3.

Rear panel connections are diagrammed in Figure 15.


Figure 14. Rear panel ports

Table 3. SiS rear panel items

1	Localization generator input	Connection point for localization generator cable
2	Workstation data connection	Fiber optic data link to workstation
3	Equipotential terminal	Connection point for equipotential cable used to reduce risk of electric shock
4	Foot switch	Connection point for foot switch used to start/stop mapping
5	Power switch	Turn the SiS on or off; illuminated blue ring indicates power on
6	Electric power input	Connection point for the external power supply. CAUTION: Use only the SiS power supply and power supply cable provided by Boston Scientific with the system. Using another power supply or power supply cable may damage the SiS.
7	Reserved for future use	Not currently used


Figure 15. Rear panel connections

11.8 SiS rear panel connections

11.8.1 Localization generator input port

The localization generator produces the magnetic field used for the magnetic localization technology. The connection cable (Figure 16) connects the generator to the input port (item 1 in both Figure 14 and Table 3).

Figure 16. Localization generator cable

Note: The potential exists for localization generator interference with other medical equipment. See the following cautions for important information on mitigating potential localization generator interference.

CAUTION: The localization generator may interfere with other systems that use magnetic-field-based technology. Consult the supplier of these systems before using them in the presence of the magnetic localization system.

CAUTION: The localization generator may interfere with implanted CIEDs (cardiac implantable electronic device). When mapping a patient with such a device, consider interrogating the device pre – and post-procedure. This will identify any changes in programmed parameters which could then be corrected before transferring the patient from the procedure room. Consult the CIED manufacturer instructions for additional information.

CAUTION: If it becomes necessary to interrogate or program a CIED while using the system, turn off the localization generator using the on-screen button located on the annotating and editing maps toolbar.

11.8.2 Workstation data connection

A fiber optic data cable (Figure 17) links the workstation to the SiS and carries the signals and data required for localization and mapping.

Figure 17. Fiber optic data cable

11.8.3 Foot switch

The foot switch (Figure 18) enables the EP to directly and immediately control the mapping process from within the procedure room.

Figure 18. Foot switch

11.8.4 Electric power input

The dedicated power supply (Figure 19) provides 24 V DC power for the SiS. Maximum power consumption for the power supply is 250 W.

The male end of the power supply's attached cable connects to the power supply port on the rear panel of the SiS. The detachable power cord plugs into a grounded AC mains electrical outlet or, preferably, an Uninterrupted Power Supply (UPS).

Figure 19. SiS external power supply

11.9 Workstation description

Function

The workstation (Figure 20) receives signals from the SiS. It processes the signals and generates data for real-time display on the workstation monitor.


Figure 20. Workstation

The workstation includes the computer hardware, accessories and software that receives, interprets, and displays data received from the SiS. It provides options that permit study archiving and retrieval. The workstation can be stored on a cart for mobility purposes.

11.10 Workstation system requirements

Power input: 100-240 V, 50/60Hz, 8A

11.11 Workstation setup

Initial workstation setup and cable connections are performed by an authorized Boston Scientific representative.

CAUTION: Use only the equipment, supplies, and accessories supplied or recommended by Boston Scientific for use with the RHYTHMIA HDx Mapping System. The use of other equipment, supplies, and accessories may cause equipment damage or system malfunction.

12. PREPARING FOR A MAPPING STUDY

12.1 Prior to patient arrival

12.1.1 Set up the hardware and make the cable connections

- Place the SiS on a dedicated table or cart that is alongside and slightly higher than the patient table. The SiS front panel should face the table and be positioned so that the cables can be attached without pulling, bending, kinking or otherwise straining them.
- Attach the localization generator to the patient table using the localization generator mount or the supplied straps. Make certain that the localization generator is located immediately below the patient's torso.
- Verify that the localization generator cable is properly connected to both the localization generator and the SiS rear panel.
- Verify that the power supply connector and fiber optic connector on the rear panel of the SiS are properly connected and secure.
- Power on the SiS and workstation.
- Observe the SiS status LED to confirm that the SiS is working properly.
- Connect breakout box(es) to ports **IN-M**, **IN-A**, or **IN-B** as needed.
- Verify that the illuminated port indicator (**M**, **A**, or **B**) on the breakout box correctly identifies the connected port.
- Connect the ECG input cable to the SiS.
- Connect the ablation connection box (specific to the RF generator being used) to the SiS and RF generator.
- Connect the external recording system to an output port (**OUT-M**, **OUT-A+ABL**, or **OUT-B**):
 - **Direct connect to Boston Scientific Clearsign amplifier:** The gray Direct Connect Clearsign IC Orion cable is used at the OUT-M port when an IntellaMap Orion mapping catheter is used. The black Direct Connect Clearsign IC A/B cable is connected to any of the three output ports when breakout box-connected catheters are used.
 - **Direct connect to GE CardioLab recording system:** The Direct Connect CardioLab IC cable may be used at any of the three output ports. Note that whenever this cable is used at port OUT-A+ABL, any data input on breakout box ports 33-40 is supplanted by ablation data on output ports 33-40.
 - **The Rhythmia IC Out 72 cable is universal** and may be used with any combination of catheters, ports, and recording devices.
- Connect the proper ECG direct connect cables or the ECG output box to the ECG output connector.
- Verify the patch cable is connected to the SiS **PATCH** input connector.
- If stimulator signals are to be routed through the SiS, connect the stimulator input jumper cable to the **STIM IN** ports and to the external stimulator.
- Verify that additional rear panel connectors (e.g. equipotential and, if applicable, foot switch) are properly connected.

12.1.2 Prepare the workstation

- Verify that the workstation fiber optic cable is properly connected and secure.
- If display of ablation information on the Rhythmia software is desired, verify cable connectivity between the workstation and the data connection port on the ablation generator.

12.2 Upon patient arrival

12.2.1 Make required cable and catheter connections

- Apply the back patch to the patient's back according to its instructions for use.
- Connect the back patch cable to the patch cable.
- Position the ECG electrodes using established facility protocol and apply according to the instructions for use provided with the Location Reference Patch Kit.
- Connect the limb and chest leads from the ECG input cable to the ECG electrodes.
- Verify the ECG electrode and limb and chest leads connectors are securely connected to reduce noise.
- Verify the limb and chest leads are properly seated in the trunk cable.
- Connect the diagnostic catheters to breakout boxes.
- Verify that the system reference electrode is connected to channel 1 of the breakout box that is connected to the **IN-A** port.
- Verify the electrode pair that may be used for stimulation are in channels 61 and 62 of a connected breakout box.
- Connect the ablation catheter to the ablation connection box via the ablation catheter extension cable.

12.2.2 Additional connections when using one of the IntellaMap Orion family of catheters

- Connect the umbilical cable to the **IN-M** port connection.
- Connect the catheter to the umbilical cable.
- Prepare the catheter according to its instructions for use.
- Condition the catheter, see section 15.

12.3 Important considerations during a study

During a study, all system operation is controlled by the workstation software and catheter manipulation. In case of electrical power interruption or system lock-up, restart the SiS and workstation and use the software 'resume' option to resume the same procedure. Verify that all catheters, signals and references are configured appropriately before continuing the procedure.

There is typically no need for user interaction with the SiS during an active study, with the exception of the following:

- Connecting and disconnecting catheters from the ablation connection box and breakout boxes according to the procedure needs.
- Monitoring the SiS Status LED – The LED should remain solid green during normal operation. Contact Boston Scientific support if this is not the case.

CAUTION: Do not connect or disconnect the SiS power supply while the unit is turned on. This will minimize the risk of equipment damage.

CAUTION: Do not connect or disconnect the localization generator from the SiS while the unit is turned on. This will minimize the risk of equipment damage.

CAUTION: Manually disabling the localization generator disables all catheter visualization and localization capabilities, including impedance tracking.

12.4 Ending the study

Disconnect all catheters and patient sensors from the SiS and discard all single use items in accordance with local procedures and regulations.

1. Power off the SiS: Push the on/off button such that the blue light is no longer illuminated.
2. Close the study and power off the workstation using the power off button on the main screen.
3. Clean the SiS and accessories, see section 16.

13. INDIVIDUAL CABLES ASSOCIATED WITH SYSTEM

Cable Type/ Label Name	Connection Scheme	Device Description
INTRACARDIAC INPUTS		
Umbilical cable	IntellaMap Orion mapping catheter to SiS front panel port IN-M	Umbilical cable is used to connect an IntellaMap Orion mapping catheter to the SiS and carry signals for 64 intracardiac electrodes, magnetic tracking sensors and memory chip.
Breakout box cable	Breakout box to SiS front panel IN-M, IN-A, and IN-B ports	Breakout box cable is an accessory cable used to connect the intracardiac inputs on the SiS (IN-M, IN-A, and IN-B) to the breakout box.
SURFACE ECG INPUTS ("AAMI" OR "IEC" APPENDED TO LABEL NAMES)		
ECG IN Chest	ECG patches to SiS front panel IN-ECG port	ECG chest and ECG limb cables are accessory cables used to connect ECG electrodes to the ECG trunk cable
ECG IN Limb		
ECG IN Trunk		ECG trunk cable connects ECG chest and ECG limb cables to the SiS and is the primary component of SiS defibrillation protection.
STIMULATOR SIGNAL INPUT		
Male to Male Tip-Jacks 2-2 Pin (Stimulator Jumper Cable)	Stimulator to front panel of SiS	Cable is used to connect the Stimulator to the SiS
RECORDING SYSTEM OUTPUTS—Direct Connect Clearsign		
Direct Connect Clearsign ECG (3m or 6m)	SiS front panel OUT-ECG port to recording system input module	Clearsign direct connect cables are accessory cables used to connect the SiS to the Boston Scientific Clearsign amplifier. Direct connect cable assemblies are fitted with a keyed connector with colored features that match the keyed receptacle and color on the SiS.
Direct Connect Clearsign IC Orion (3m or 6m)	SiS front panel OUT-M port to recording system input module—gray cable	
Direct Connect Clearsign IC A/B (3m or 6m)	SiS front panel OUT-A+ABL or OUT-B port to recording system input module—black cable	
RECORDING SYSTEM OUTPUTS—Direct Connect CardioLab		
Direct Connect CardioLab ECG (3m or 6m)	SiS front panel OUT-ECG port to recording system input module	CardioLab direct connect cables are accessory cables used to connect the SiS to the GE CardioLab recording system. Direct connect cable assemblies are fitted with a keyed connector with colored features that match the keyed receptacle and color on the SiS.
Direct Connect CardioLab IC (3m or 6m)	SiS front panel OUT-M, OUT-A+ABL, OUT-B ports to recording system input module	
RECORDING SYSTEM OUTPUTS—Other		
Rhythmia IC Out 72 Cable	SiS front panel OUT-M, OUT-A+ABL, OUT-B ports to recording system input module	Accessory cable used to connect the intracardiac outputs on the SiS (OUT-M, OUT-A+ABL, OUT-B) to the recording system pin blocks.

Cable Type/ Label Name	Connection Scheme	Device Description
WORKSTATION ABLATION DATA CABLES		
USB to RS232 serial adapter cable	RF generator to workstation	Serial adapter is used to connect the RF generators to the workstation
DB9 RS232 serial Cable Strt Thr M/F 10ft	RF generator to USB to RS232 Serial adapter to workstation	Serial cable is used to connect the Maestro, IBI, Ampere RF generators to the work station.
DB9 RS232 serial cable null F/F 10ft	RF generator to USB to RS232 serial adapter to workstation	Serial cable is used to connect the EP-Shuttle, SMARTABLATE, and INTELLAGEN RF generators to the workstation.
Serial splitter	RF generator to the workstation and recording system	The serial splitter is used to share ablation information from the RF generator to the workstation and the recording system.
OTHER FRONT PANEL CABLES		
Rhythmia patch cable	Back patch to front panel of SiS	The patch cable is an accessory cable used to connect the back patch to the SiS PATCH input.
REAR PANEL CABLES		
RHYTHMIA HDx power supply	Input: AC mains Output: SiS rear panel	Power supply is used to convert mains AC voltage to low level DC voltage for the SiS. It includes the power converter itself as well as a connection cable to the SiS.
Equipotential cable	Connects to common facility equipotential terminal	Equipotential cable prevents touch voltages by bringing the SiS to the same potential as other conductive surfaces in the room.
Localization generator cable	Localization generator to the rear panel of SiS	Localization generator generates the magnetic field used for the magnetic localization technology. The localization generator cable connects the localization generator to the SiS.
Fiber Ethernet cable	Rear panel of SiS to workstation	Fiber optic data link to the workstation from the SiS
OTHER CABLES		
Power cord	AC mains to workstation	Power cord used to provide power to the workstation
	AC mains to workstation monitor	Power cord used to provide power to the workstation monitor
	AC mains to SiS power supply	Power cord used to provide power to the SiS power supply. Length: 1.8 m to 3.05 m

14. DIRECT CONNECT CHANNEL MAPPING

Clearsign IC Orion cable (Gray)		Clearsign IC A/B cable (Black)				CardioLab cable (3 end connectors)					
OUT-M		OUT-A+ABL		OUT-M and OUT-B		OUT-M		OUT-A+ABL		OUT-B	
Orion spline	Output channels	Input channels	Output channels	Input channels	Output channels	Input channels	Output channels	Input channels	Output channels	Input channels	Output channels
A1-A8	1-8	BB 1-64	1-64	BB 1-64	1-64	Orion A1-H8	1-64	BB 1-32	1-32	BB 1-64	1-64
B1-B8	11-18	ABL 1-8	65-72			BB 1-64		ABL 1-8	33-40		
C1-C8	21-28							BB 33-40	N/A		
D1-D8	31-38							BB 41-64	41-64		
E1-E8	41-48										
F1-F8	51-58										
G1-G8	61-68										
H1-H8	71-78										

Input channels: BB = Breakout Box ABL = Ablation Connection Box

15. PREPARING INTELLAMAP ORION FAMILY OF CATHETERS FOR USE

When used with the RHYTHMIA HDx Mapping System, an IntellaMap Orion mapping catheter must be conditioned before use.

1. Connect an IntellaMap Orion mapping catheter to either end of the umbilical cable.
2. Connect the other end of the umbilical cable to the **M** input connection on the front panel of the SiS.
3. Verify that the SiS is turned on.
4. Confirm that the localization generator and back patch are connected to the SiS.
5. Confirm that the back patch is properly positioned on the patient's back and is within the localization generator magnetic field.
6. Make certain that the catheter is *outside* the patient's body and the localization generator's magnetic field.
7. Completely submerge an IntellaMap Orion electrode array in sterile saline. It is essential that all electrodes remain immersed in the saline solution during the entire conditioning process.
8. Wait for the blue circle on the SiS front panel **Condition** button to begin flashing and then press the button to begin the conditioning process. The circle will remain illuminated for approximately 10 seconds and a status message will be shown on the workstation screen.
9. Conditioning process completion is indicated when the blue circle is no longer illuminated and a confirmation message is shown on the workstation screen.

Note: Consult RHYTHMIA HDx Mapping System Software Instructions for Use for guidelines on overriding pre-conditions.

16. CLEANING GUIDELINES

16.1 Routine cleaning and disinfecting

Observe facility cleaning protocol while following the instructions below:

- Turn off and unplug all system components that are connected to a power supply.
- Wipe the external system components with a cloth slightly dampened with a mild cleaning solution. Do not allow any amount of cleaning solution or moisture to come in contact with the cable connection ports and do not spray water or other liquids onto system components.
- The system components do not require sterilization or disinfection.
- Allow the cleaned surfaces to air dry prior to reconnecting or using the system components.

16.2 Decontaminating the equipment before shipping

To comply with United States postal and transportation laws, equipment returned to Boston Scientific must be properly decontaminated with chemical germicide that is approved for use as a hospital disinfectant. If equipment is received without having been decontaminated, Boston Scientific will charge a cleaning surcharge. Any equipment return to Boston Scientific without having been properly decontaminated must have an appropriate biohazard warning.

CAUTION: Do not immerse system components in water, cleaning solutions, or liquid. Prevent fluids from entering air vents. Ensure connectors stay dry. Failure to follow cleaning guidelines may cause equipment damage or system malfunction, and may also void any warranties or service agreements.

CAUTION: To avoid equipment damage and malfunction, do not attempt to sterilize equipment that is provided non-sterile.

CAUTION: To avoid equipment damage and malfunction, do not insert anything (e.g., cotton swabs or pins) into cable connectors or equipment ports or openings.

CAUTION: To minimize the risk of electrical shock, prior to cleaning any electrically powered equipment, be sure to first turn off the device and then unplug it from power.

CAUTION: Do not attempt to clean system components during system operation. Cleaning equipment during use increases the risk of electrical shock, system malfunction, and device drop.

17. SiS TROUBLESHOOTING

The front panel status indicator LED displays system hardware and communication status. It should be monitored during both SiS initial start-up and live mapping sessions to assure correct system operation.

Table 4. Status Indicator States

Indicator state	Meaning	Action required
 DARK	Power off	None
 FLASHING	SiS hardware initialization in progress	Contact Boston Scientific Support if SiS initialization period takes more than two minutes.
 FLASHING	Hardware initialization complete, SiS ready to communicate with workstation	<ol style="list-style-type: none"> 1. Connect the fiber optic data cable to the SiS and workstation 2. Start a LIVE mapping session on the workstation
 SOLID	RHYTHMIA HDx Mapping System fully operational	None
 SOLID	Error identified	<p>Do not use the SiS.</p> <p>Contact Boston Scientific Support.</p>

Note: Connecting the fiber optic data cable and starting the LIVE mapping session may be performed in any order. The fiber optic data connection may also be made before any of the hardware is powered on.

18. ELECTROMAGNETIC EMISSIONS AND IMMUNITY

Table 5. Guidelines and manufacturer's declaration – electromagnetic emissions

The RHYTHMIA HDx Mapping System is intended for use in the electromagnetic environment specified below. The customer or the user of the RHYTHMIA HDx Mapping System should assure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment – guidance
RF Emissions EN 55011 CISPR 11	Group 1	The RHYTHMIA HDx Mapping System uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment. The RHYTHMIA HDx Mapping System is suitable for use in all establishments other than domestic, and may be used connected to the public low voltage power supply network that supplies buildings used for domestic purposes, provided the following warning is heeded: WARNING: The RHYTHMIA HDx Mapping System is intended for use by healthcare professionals only. This system may cause radio interference or may disrupt the operation of nearby equipment. It may be necessary to take mitigation measure such as re-orienting or relocating the RHYTHMIA HDx Mapping System or shielding the location.
RF Emissions EN 55011 CISPR 11	Class A	
Harmonic Emission EN 61000-3-2	Class A	
Voltage Fluctuations/ flicker Emission EN 61000-3-3	Class A	

Note: The RHYTHMIA HDx Mapping System is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the RHYTHMIA HDx Mapping System can help prevent electromagnetic interference by maintaining a minimum distance of 30 cm (12 inches) between portable and mobile RF communications equipment (transmitters) and the RHYTHMIA HDx Mapping System.

Table 6. Guidelines and manufacturer's declaration – electromagnetic immunity

The RHYTHMIA HDx Mapping System is intended for use in the electromagnetic environment specified below. The customer or the user of the RHYTHMIA HDx Mapping System should assure that it is used in such an environment.

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment
Electrostatic Discharge (ESD) IEC 61000-4-2	±8 kV contact discharge ±15 kV air discharge	Pass	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical Fast Transient/Burst IEC 61000-4-4	±2 kV AC Mains ±1 kV I/O Lines	Pass	Mains power quality should be that of a typical commercial or hospital environment.
Surge Line to Line (AC Power) IEC 61000-4-5	±1 kV Line to Line ±2 kV Line to Ground	Pass	Mains and power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	100% dip in 1 cycle 30% dip in 25/30 cycles 100% dip in 250/300 cycles	Pass	Mains power quality should be that of a typical commercial or hospital environment. If the user of the RHYTHMIA HDx Mapping System requires continued operation during power mains interruptions, it is recommended that the RHYTHMIA HDx Mapping System be powered by an uninterruptable power supply (UPS).
Power frequency (50/60Hz) magnetic field IEC 61000-4-8	30 A/m	Pass	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
Conducted RF IEC 61000-4-6	6 V/m in ISM bands from 0.15 MHz to 80 MHz	Pass	Portable and mobile RF communications equipment should be used no closer to any part of the RHYTHMIA HDx Mapping System, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.

Table 6. Guidelines and manufacturer's declaration – electromagnetic immunity (continued)

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment	
Radiated RF IEC 61000-4-3	3 V/m	Pass	Recommended separation distance where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m)	
	80 MHz to 2.7 GHz: 27 V/m			
	380 - 390 MHz: 28 V/m			
	430 - 470 MHz: 9 V/m			
	704 - 787 MHz: 28 V/m			
	800 - 960 MHz: 28 V/m		$d = \left[\frac{1.17}{V_i} \right] \sqrt{P}$	150 KHz to 80 MHz
	1700 - 1900 MHz: 28 V/m		$d = \left[\frac{1.17}{E_i} \right] \sqrt{P}$	80 MHz to 800 MHz
	2400 - 2570 MHz: 9 V/m		$d = \left[\frac{2.33}{E_i} \right] \sqrt{P}$	800 MHz to 2.5 GHz
	5100 - 5800 MHz:			
<p>Note 1: UT is the AC mains voltage prior to application of the test level.</p> <p>Note 2: At 80 MHz and 800 MHz, the higher frequency range applies.</p> <p>Note 3: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.</p> <p>Note 4: The ISM bands between 0.15 MHz and 80 MHz are 6.765 MHz to 6.795 MHz, 13.553 MHz to 13.567 MHz, 26.957 MHz to 27.283 MHz, and 40.66 MHz to 40.70 MHz.</p> <p>The amateur radio bands between 0.15 MHz and 80 MHz are: 1.8 MHz to 2.0 MHz, 3.5 MHz to 4.0 MHz, 5.3 MHz to 5.4 MHz, 7.0 MHz to 7.3 MHz, 10.1 MHz to 10.15 MHz, 14.0 MHz to 14.2 MHz, 18.07 MHz to 18.17 MHz, 21.0 MHz to 21.4 MHz, 24.89 MHz to 24.99 MHz, 28.0 MHz to 29.7 MHz, and 50.0 MHz to 54.0 MHz.</p>				

19. TECHNICAL SPECIFICATIONS

Device	Category	
SiS	Position	The SiS should be placed on a dedicated table or cart that is adjacent to the foot of the patient table or on the foot of the patient table.
	Patient Inputs	12 lead ECG
		64 intracardiac channels of an IntellaMap Orion mapping catheter
		128 additional diagnostic channels plus eight dedicated ablation signal channels
	Connectivity with other lab equipment	Stimulator: up to two stimulator channels
		Recording system: Output all signals to a third-party recording system through direct connect cables or Rhythmia IC Out Cable 72
		RF generator: dedicated connection to Maestro, IBI, Ampere, SMARTABLATE, EP-Shuttle, and INTELLAGEN RF generators
	Power Inputs	Power supply manufacturer: XP Power
		Model: AHM250PS24-XD0653
		Input voltage: 100-240 VAC/50 Hz; 100-220 VAC/60 Hz
		Input current: 2.3 A at 115 VAC; 1.2 A at 230 VAC
		Output voltage: 24 VDC
		Output power: 250 W; Output current: 10.41 A
	Safety	IEC 60601-1:2005+A1:2012
		IEC 60601-1-2:2014
		IEC 60601-1-6:2010+A1:2013
		IEC 60601-2-27:2011 CORR 1:2012
IEC 62366:2007+A1:2014		
Protection against electric shock: Class I, Type CF	EN 55011:2009+A1:2010, Group 1, Class A;	
	CISPR 11:2009+A1:2010, Group 1, Class A	
	Protection against fluid ingress: IPX1 (protection against vertically falling water drops)	
Defibrillation recovery time: 5 seconds		
Transmitter	Frequency: 800 Hz	
	Output Power: 72 dB μ A/m @ 10 m	
Input Accuracy	Intracardiac: 5%, 1-300Hz	
	ECG: 5%, 0.67-150Hz	

Device	Category	
Workstation	Position	The workstation is typically positioned in the control room.
	Operating System	Linux
	Connection to SiS	Fiber optic cable
	Power Inputs	Input voltage: 100-240 VAC
		Input current: 8.0 A
		Input frequency: 50-60Hz
Accessories	Monitor, keyboard, and mouse	
Environmental Operating Conditions		Equipment Storage and Transport
<ul style="list-style-type: none"> • Temperature: 10°C (50°F) – 30°C (86°F) • Humidity: 30% – 75% non-condensing • Atmospheric pressure: 580 –760 mmHg 		<ul style="list-style-type: none"> • Temperature -29°C to +60°C • Humidity 10% -90% non-condensing

20. SERVICE INFORMATION

Equipment Service and Maintenance

The RHYTHMIA HDx Mapping System does not have any parts that require periodic service from users. In case of malfunction, contact Boston Scientific Support.

21. WARRANTY

For device warranty information, visit (www.bostonscientific.com/warranty).

EU Importer: Boston Scientific International B.V., Vestastraat 6, 6468 EX Kerkrade, The Netherlands
RHYTHMIA HDx, IntellaNav, IntellaMap Orion, IntellaTip MiFi, Maestro, Clearsign, EP-Shuttle, and INTELLAGEN are trademarks of Boston Scientific Corporation or its affiliates.

CardioLab is a trademark of GE Medical Systems.

IBI and Ampere are trademarks of St. Jude Medical, Atrial Fibrillation Division, Inc.

SMARTABLATE is a trademark of Biosense Webster, Inc. Biosense Webster has no association with the RHYTHMIA HDx Ablation Connection Box by Boston Scientific.

While testing was done by Boston Scientific, Biosense Webster has not tested or validated the use of this device with its SMARTABLATE System.

All other trademarks are the property of their respective owners.

22. CONTACTS

For service and support in using this system please contact Boston Scientific Support using the resources given below. Do not send any parts or equipment for service to Boston Scientific without prior authorization.

Technical Support (North America)

Tel +1 800 949 6708

Fax +1 510 624 2493

CETechSupportUSA@bsci.com

Technical Support (Europe, Middle East, Africa)

Tel 0031 (0)45 5467707

Fax 0031 (0)45 5467805

CETechSupportEMEA@bsci.com

Technical Support (Japan)

Tel +81 03 6853 1000

Fax +81 45 444 2799

japantsc@bsci.com

23. SOFTWARE LICENSE

You have acquired the RHYTHMIA HDx Mapping System which includes software developed and owned by Boston Scientific Corporation as well as software licensed by Boston Scientific from various software licensors. Refer to RHYTHMIA HDx Mapping System Software Instructions for Use for further information.

24. SYMBOL DEFINITIONS

Commonly used medical device symbols that appear on the labeling are defined at www.bostonscientific.com/SymbolsGlossary.

Additional symbols are defined at the end of this document.

Остаряла версия. Да не се използва.
Zastaralá verze. Nepoužívat.
Forældet version. Må ikke anvendes.
Version überholt. Ärge kasutage.
Aegunud versioon. Μην την χρησιμοποιείτε.
Παλιά έκδοση. Μην την χρησιμοποιείτε.
Outdated version. Do not use.
Version obsolete. Ne pas utiliser.
Zastarjela verzija. Nemojte upotrebljavati.
Úrejt útgáfa. Ne pas utilizar.
Versione obsolete. Ne pas utiliser.
Zastarjela verzija. Neizmantot.
Úrejt útgáfa. Notið ekki.
Novecojsi versija. Non utilizzate.
Elavult verzió. Ne használjate.
Dit is een verouderde versie. Niet gebruiken.
Wersja przeterminowana. Nie używać.
Versão obsoleta. Não utilize.
Utdatert versjon. Skal ikke brukes.
Wersja przeterminowana. Nie używać.
Zastaraná verzia. Nepoužívať.
Zastarela različica. Ne uporabite.
Vanhentunut versio. Älä käytä.
Föråldrad version. Använd ej.
Güncel olmayan sürüm. Kullanmayın.

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All rights reserved.
Boston Scientific AG
Boston Scientific (Australia) Pty Ltd
© 2022 Boston Scientific Corporation or its affiliates.