# ACURATE neo2™

## Transfemoral Delivery System Loading Kit

**REF**

SYM-DS-010

SYM-AC-010

SYM-SV23-004

SYM-SV25-004

SYM-SV27-004

---

**en Instructions for Use 3**

---

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---

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---

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---

51081295-01

2022-10

## TABLE OF CONTENTS

<b>1. REUSE WARNING</b> .....	<b>3</b>
<b>2. DEVICE DESCRIPTION</b> .....	<b>3</b>
2.1 ACURATE neo2 Transfemoral Delivery System.....	3
2.1.1 Delivery System Contents .....	3
2.1.2 Delivery System Operating Principle.....	3
2.1.3 Delivery System Materials .....	5
2.1.4 Delivery System Non-Pyrogenic.....	5
2.2 ACURATE neo2 Loading Kit.....	5
2.2.1 Loading Kit Contents.....	5
2.2.2 Loading Kit Operating Principle.....	5
Table 1: Loading Kit Components.....	6
Table 2: Device Catalog Numbers.....	7
2.3 User Information.....	7
<b>3. INTENDED USE/INDICATIONS FOR USE</b> .....	<b>7</b>
<b>4. CLINICAL BENEFIT STATEMENT</b> .....	<b>7</b>
<b>5. SUMMARY OF SAFETY AND CLINICAL PERFORMANCE</b> .....	<b>8</b>
<b>6. CONTRAINDICATIONS</b> .....	<b>8</b>
<b>7. WARNINGS</b> .....	<b>8</b>
<b>8. PRECAUTIONS</b> .....	<b>9</b>
8.1 Precautions Prior to Use.....	9
8.2 Precautions During Use.....	9
<b>9. ADVERSE EVENTS</b> .....	<b>9</b>
<b>10. HOW SUPPLIED</b> .....	<b>10</b>
10.1 Device Details.....	10
10.2 Handling and Storage.....	11
<b>11. OPERATIONAL INSTRUCTIONS</b> .....	<b>11</b>
11.1 Additional Items for Safe Use.....	11
11.1.1 Sterile Surgical Tools.....	11
11.1.2 Sterile Rinsing Material.....	11
11.2 Device Preparation.....	11
11.2.1 Device Selection and Pre-Implantation Requirements.....	11
11.2.2 Valve Preparation.....	12
11.2.3 Inspection of Delivery System and Loading Kit.....	12
11.2.4 Valve Rinsing.....	12
11.3 Loading Procedure.....	13
11.3.1 Loading of the Valve onto the Delivery System.....	13
11.3.2 Flushing the Delivery System.....	13
11.4 Implantation Procedure.....	13
11.4.1 Coronary Occlusion Risk Mitigation.....	13
Table 3: Occlusion Risk Mitigation.....	14
11.4.2 Transfemoral Access.....	14

11.4.3 Pre-Dilation of Native Valve .....	14
11.4.4 Valve Delivery .....	14
11.4.5 Verification of the Valve Position and Post-Implant Monitoring .....	16
Table 4: Maximum Inflated Balloon Diameter .....	16
11.5 Disposal .....	17
11.6 Post-Procedure .....	17
<b>12. WARRANTY .....</b>	<b>17</b>
<b>13. SYMBOL DEFINITIONS .....</b>	<b>17</b>

# ACURATE neo2™

## Transfemoral Delivery System

### Loading Kit

#### Rx ONLY

**Caution:** Federal Law (USA) restricts this device to sale by or on the order of a physician.

#### 1. REUSE WARNING

Contents supplied STERILE using a radiation process. Do not use if sterile barrier is damaged. If damage is found, call your Boston Scientific representative.

For single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

After use, dispose of product and packaging in accordance with hospital, administrative and/or local government policy.

#### 2. DEVICE DESCRIPTION

##### 2.1 ACURATE neo2 Transfemoral Delivery System


###### 2.1.1 Delivery System Contents

One (1) ACURATE neo2 Transfemoral Delivery System.

**Table 2** indicates the available catalog number.

###### 2.1.2 Delivery System Operating Principle

The ACURATE neo2 Transfemoral Delivery System (or Delivery System) is used to position and deploy the ACURATE neo2 Aortic Valve (or Valve) at the intended position over the patient's native, calcified aortic valve via transfemoral access. The Delivery System is comprised of a catheter, a release handle, and an insertion aid. Key components of the Delivery System are shown in **Figure 1: Transfemoral Delivery System**.

**Figure 1: Transfemoral Delivery System**

1. The catheter is comprised of three components with the following features:
  - A flexible inner member which contains a guidewire lumen bonded proximally onto a female luer lock and distally to a radio-opaque atraumatic tip and capsule which contains the distal section of the Valve when loaded onto the Delivery System. The inner member is fixed proximally to the release handle.
  - A flexible middle member with a stent holder to avoid premature delivery of the Valve during the release procedure and a sliding shuttle with an atraumatic soft tip to facilitate the Delivery System retrieval through the introducer sheath and full closure of the Delivery System if required. A self-expanding cage positioned over the stent holder guarantees a smooth transition to the capsule of the inner member. In particular, the middle member of the Delivery System presents an additional marker band for easier Valve positioning within the native aortic annulus. The middle member encompasses the inner member and is fixed proximally to the release handle.
  - A flexible outer member which contains distally the proximal section of the Valve when loaded onto the Delivery System. The outer member encompasses the middle member and is fixed proximally to the release handle.
2. The release handle provides an ergonomic fit to the physician's hand to facilitate the precise two step deployment of the Valve. The release handle comprises the following features:
  - A first rotating knob acting on the outer member for the controlled release of the proximal section of the Valve from the Delivery System
  - A second rotating knob acting on the inner member & capsule for the controlled release of the distal section of the Valve from the Delivery System
  - A safety button to avoid premature release of the distal section of the Valve
  - A female luer lock / flushing line for flushing the annular space between the inner, middle and outer member
  - A positioning sheath proximally connected to the release handle and covering the flexible outer member
3. The insertion aid slides over the outer member and facilitates the insertion of the Delivery System into the compatible introducer sheath.

All three flexible members of the catheter are coaxially arranged. The inner and outer members slide longitudinally relative to the middle member during the two-step release procedure of the Valve.

The ACURATE neo2 Transfemoral Delivery System has a nominal outer diameter of 18F (6.0 mm) and a usable length of 113 cm and is compatible with a 14F iSLEEVE™ Expandable Introducer Set and with 0.035 in (0.89 mm) Safari<sup>2</sup>™ Pre-Shaped TAVI Guidewires or 0.035 in (0.89 mm) guidewires.

Refer to the ACURATE neo2 Aortic Valve Instructions for Use for a detailed description of the valves.

---

**NOTE:** Refer to Introducer Sheath Instructions for Use for femoral artery diameter compatibility.

---

### 2.1.3 Delivery System Materials

The Delivery System is comprised of a catheter<sup>1</sup>, a release handle, and an insertion aid consisting of multiple polymers and metal alloys.

<sup>1</sup>Contains cobalt: CAS No. 7440-48-4; EC No. 231-158-0. Defined as a 1B carcinogen and reproductive toxicant according to the European Commission in a concentration above 0.1% weight by weight.

---

**Note:** This device is made with a metal alloy which contains cobalt. Current scientific evidence supports that metal alloys containing cobalt used in medical devices, do not cause an increased risk of cancer or adverse reproductive effects.

---

### 2.1.4 Delivery System Non-Pyrogenic

This device meets pyrogen limit specifications.

## 2.2 ACURATE neo2 Loading Kit

### 2.2.1 Loading Kit Contents

One (1) Crimping Spatula

One (1) Flushing Line

One (1) Loading Tray and Lid

One (1) Loading Tube

One (1) Safety Button

One (1) Shuttle Pusher

One (1) Stylet

Three (3) Rinsing Bowls

One (1) Distal Loader

One (1) Compressor Ring

One (1) Loading Funnel

One (1) Loading Nut

One (1) Silicon Ring

One (1) Split Ring



One (1) Star Pusher

One (1) Proximal Loader


### 2.2.2 Loading Kit Operating Principle

The Loading Kit is used to load the ACURATE neo2 Aortic Valve onto the ACURATE neo2 Transfemoral Delivery System. See **Table 1** for visual representation of Loading Kit components and **Table 2** indicates the available catalog number.

**Table 1: Loading Kit Components**

Component Name	Visual Representation
Crimping Spatula	
Flushing Line	
Loading Tube	
Loading Tray and Lid	
Safety Button	
Shuttle Pusher	
Stylet	
Three (3) Rinsing Bowls	
<b>Distal Loader Components</b>	
 Compressor Ring	 Loading Funnel
 Loading Nut	 Silicon Ring

**Table 2: Device Catalog Numbers**

Catalog Numbers	Product Name
SYM-DS-010	ACURATE neo2 Transfemoral Delivery System
SYM-AC-010	ACURATE neo2 Loading Kit

### 2.3 User Information

Intended Users of the ACURATE neo2 Transfemoral Delivery System are physicians who have completed the Boston Scientific Corporation (BSC) Physician Training Plan. Use of the ACURATE neo2 Transfemoral Delivery System should only be performed by these Intended Users.

The ACURATE neo2 Aortic Valve loading procedure must only be performed by personnel trained by BSC. Other personnel are not authorized to execute the procedure.

To obtain optimal clinical results, proper patient selection, careful sizing and attention to procedural execution is warranted. A more detailed guidance can be found in the following publications: Kim W-K et al. Am J Cardiol 2020;131:12-16 and Okuno T et al. Current Cardiology Reports (2020) 22: 107. ACURATE neo2 is subject to continuous clinical monitoring.

### 3. INTENDED USE/INDICATIONS FOR USE

The ACURATE neo2 Aortic Valve, in combination with the ACURATE neo2 Transfemoral Delivery System and the ACURATE neo2 Loading Kit, is indicated for relief of aortic stenosis in patients with symptomatic heart disease due to severe native calcific aortic stenosis who are judged by a Heart Team, including a cardiac surgeon, to be appropriate for the transcatheter heart valve replacement therapy.

**FOR CANADA ONLY:** The Indications for Use in this document are not valid in Canada. See Addendum attached to product for Indications for Use

### 4. CLINICAL BENEFIT STATEMENT

The clinical benefit of ACURATE neo2 Aortic Valve System is to improve aortic valve function with an intent to relieve symptoms and improve survival in patients with symptomatic heart disease due to

severe native calcific aortic stenosis who are judged by a heart team to be appropriate for transcatheter heart valve replacement therapy.

## 5. SUMMARY OF SAFETY AND CLINICAL PERFORMANCE

For customers in the European Union, use the device name found in the labeling to search for the device's Summary of Safety and Clinical Performance, which is available on the European database on medical devices (EUDAMED) website: (<https://ec.europa.eu/tools/eudamed>).

## 6. CONTRAINDICATIONS

The ACURATE neo2 Aortic Valve, in combination with the ACURATE neo2 Transfemoral Delivery System and the ACURATE neo2 Loading Kit is contraindicated in patients with:

- Non-valvular aortic stenosis
- Congenital aortic stenosis or unicuspid or bicuspid aortic valve
- Presence of a prosthetic mitral valve
- Presence of previously implanted aortic bioprosthesis
- Non-calcific acquired aortic stenosis
- Non-stenotic aortic insufficiency
- Severe eccentricity of calcification
- Evidence of intracardiac mass, thrombus or vegetation
- Severe coagulation problems
- Active bacterial endocarditis or other active infections
- Severe ventricular dysfunction with ejection fraction < 20%
- Inability to tolerate anticoagulation therapy
- Hypertrophic cardiomyopathy with or without obstruction (HOCM)
- Known allergy to nickel, aspirin or contrast media
- Distance between the base of the coronary cusp and the respective coronary ostium less than 8 mm
- Leaflet geometry with respect to coronary ostium location presenting a risk of overlapping
- Anatomy NOT appropriate for transfemoral implant due to the size, disease and degree of calcification or tortuosity of the aorta or ilio-femoral arteries

The ACURATE neo2 Aortic Valve, in combination with the ACURATE neo2 Transfemoral Delivery System and the ACURATE neo2 Loading Kit should not be used if the implanting physician believes its implantation would be contrary to the best interest of the patient.

The ACURATE neo2 Aortic Valve shall not be placed in positions other than in the native aortic valve.

---

## 7. WARNINGS

- Implantation of the ACURATE neo2 Aortic Valve with the ACURATE neo2 Transfemoral Delivery System and Loading Kit should only be performed by physicians trained according to the BSC Physician Training Plan. Other physicians are not authorized to perform valve implantation.
- The Valve loading procedure must only be performed by personnel trained by BSC. Other personnel are not authorized to execute the procedure.

- The outer surface of the peelable pouches of the Delivery System and Loading Kit are not sterile and must not be placed in the sterile field.
- If the Delivery System or Loading Kit remains unused after being removed from its primary packaging it shall be considered non-sterile and must not be used.
- Handle the Delivery System with care. Catheter damage or malfunction may result from kinking, stretching or torquing/rotating the Delivery System.
- Ensure full detachment of the Valve from the stent holder prior to retraction of the Delivery System.
- Proceed carefully if resistance is felt during retraction of the Delivery System.

Please refer to the Instructions for Use of the ACURATE neo2 Aortic Valve for further warnings related to the Valve.

## 8. PRECAUTIONS

### 8.1 Precautions Prior to Use

- The Delivery System contains cobalt. Current scientific evidence supports that metal alloys containing cobalt used in medical devices do not cause an increased risk of cancer or adverse reproductive effects. However, the ACURATE neo2 Aortic Valve, in combination with the ACURATE neo2 Transfemoral Delivery System and the ACURATE neo2 Loading Kit should not be used if the implanting physician believes its implantation would be contrary to the best interest of the patient, including treatment of children or treatment of pregnant or breast-feeding women, or other patient groups considered particularly vulnerable to cobalt.
- Do not use the Delivery System if it is observed to be damaged or malfunctioning.
- Do not use the Delivery System if unable to flush.

### 8.2 Precautions During Use

- The implant procedure should be conducted under fluoroscopic guidance.

Please refer to the Instructions for Use of the ACURATE neo2 Aortic Valve for further PRECAUTIONS related to the Valve.

## 9. ADVERSE EVENTS

Risks, complications and adverse events that may be associated with the use of the ACURATE neo2 Aortic Valve, in combination with the ACURATE neo2 Transfemoral Delivery System and the ACURATE neo2 Loading Kit include those risks related to conventional surgical aortic valve replacement (SAVR), as well as those related to transcatheter aortic valve implantation (TAVI).

Known or anticipated risks are listed below in alphabetical order:

- Allergic or adverse reaction (including medications, anesthesia, contrast, or device materials)
- Angina
- Arrhythmia or new conduction system injury (including need for pacemaker insertion)
- Bleeding including hemorrhage or hematoma (possibly requiring transfusion or additional intervention)

- Cardiac failure leading to low cardiac output (cardiogenic shock) or pulmonary edema
- Cerebrovascular accident, stroke, transient ischemic attack or cerebral infarction including asymptomatic neuroimaging findings
- Coronary obstruction
- Death
- Device misplacement, migration or embolization
- Emboli (including air, tissue, thrombus or device materials)
- Fever and pyrogen reaction
- Heart failure
- Hemolysis and/or hemolytic anemia
- Hypertension/hypotension
- Infection (local or systemic including endocarditis)
- Mitral valve insufficiency or injury
- Myocardial infarction
- Myocardial or valvular injury (including perforation or rupture)
- Nerve injury or neurologic deficits (including encephalopathy)
- Pain or inflammation
- Pericarditis, pericardial effusion or tamponade
- Peripheral ischemia or infarction
- Pleural effusion
- Radiation injury
- Renal insufficiency or failure
- Respiratory insufficiency or failure
- Valve dysfunction, deterioration or failure
- Valve or device thrombosis
- Valvular stenosis or regurgitation (central or paravalvular)
- Vessel injury (including access site) such as spasm, lymphatic problems, pseudoaneurysm, arteriovenous fistula, trauma, dissection, perforation, and rupture

As a result of these adverse events, the subject may require medical, percutaneous or surgical intervention, including implantation of a second valve, re-operation and replacement of the valve. These events may lead to fatal outcomes.

## 10. HOW SUPPLIED

### 10.1 Device Details

The ACURATE neo2 Transfemoral Delivery System and the ACURATE neo2 Loading Kit are supplied sterile and non-pyrogenic.

The primary packaging of the Delivery System consists of a tray containing the device which is placed in a sealed peelable pouch. The primary packaging of the Loading Kit consists of a tray containing the loading components which are placed in a sealed peelable pouch.

The secondary packaging for the Delivery System and the Loading Kit consists of a carton.

Do not use if package is damaged or unintentionally opened before use.

Do not use if labeling is incomplete or illegible.

## 10.2 Handling and Storage

The ACURATE neo2 Transfemoral Delivery System and ACURATE neo2 Loading Kit shall be stored in a dry area at room temperature thereby avoiding direct contact with sunlight. Storage at elevated temperatures may potentially damage the polymeric components and adhesives thus compromising the product performances.

## 11. OPERATIONAL INSTRUCTIONS

Physicians using the ACURATE neo2 Transfemoral Delivery System and Loading Kit with the ACURATE neo2 Aortic Valve must be experienced in:

- Balloon aortic valvuloplasty
- Transfemoral access and catheterization
- Aortic bioprosthesis implantation using transcatheter / transfemoral procedures

### 11.1 Additional Items for Safe Use

Specific equipment required for the preparation of the ACURATE neo2 Transfemoral Delivery System and the ACURATE neo2 Aortic Valve are as follows:

#### 11.1.1 Sterile Surgical Tools

- Scalpel
- Large Clamp or Tweezers

#### 11.1.2 Sterile Rinsing Material

- 3 bowls (included)
- Loading Tray (included)
- 3.0 L refrigerated saline  $\leq 5^{\circ}\text{C}$
- 1.5 L saline, room temperature
- 5,000 units Heparin
- 2 x 20 ml luer syringe (with luer lock)

General equipment required for the use of the ACURATE neo2 Transfemoral Delivery System and the Loading Kit with the ACURATE neo2 Aortic Valve are as follows:

- Standard catheter laboratory equipment such as fluoroscopy, echocardiography, pressure measuring equipment.
- Standard catheterization devices such as guidewires, introducer sheaths, balloon valvuloplasty catheters, pacemakers.

### 11.2 Device Preparation

#### 11.2.1 Device Selection and Pre-Implantation Requirements

Native aortic valve pre-dilation and transfemoral aortic valve implantation shall be performed in a catheter laboratory or hybrid room with hemodynamic monitoring, high-quality fluoroscopic and echocardiographic imaging capabilities.

The procedure for the transfemoral Valve implantation consists of a standard catheterization of the right or left femoral artery followed by the introduction of the Delivery System in the vasculature for retrograde access to the aortic root and positioning and deploying the Valve over the stenotic, native aortic valve.

Different techniques can be used for the measurement of the aortic annulus diameter and consequently for the selection of the size of the Valve to be implanted. Among these are transesophageal echocardiography (TEE), transthoracic echocardiography

(TTE), computerized tomography (CT) scan, and angiography (ventriculography).

3D CT-scan reconstruction shall be considered as the best choice to measure the annulus diameter.

The technique for the dimensional assessment of the aortic valve and the supra-annular apparatus (e.g. coronary ostia positioning) is left at the physician's discretion.

Finally, the implantation shall be performed using general or regional anesthesia and without cardiopulmonary bypass.

#### **11.2.2 Valve Preparation**

The Valve and its preserving solution are sterile and do not require further treatment. The Valve is packaged sterile in a bottle with a sealed screw cap and is identified with a serial number tag attached to a valve holder.

Before opening the primary packaging, the following steps shall be performed:

- Confirm no signs of broken seal and/or signs of moisture
- Confirm no solution residue
- Confirm no leakage
- Ensure level of preserving solution totally covers the Valve
- Confirm integrity of the sealed screw cap to ensure the bottle has not been opened previously
- Confirm Use By date has not passed
- Ensure temperature indicator displays "OK" on screen

#### **11.2.3 Inspection of Delivery System and Loading Kit**

Before opening the primary packaging of the Delivery System and Loading Kit, the following controls shall be performed:

- Confirm integrity of the peelable pouch with no pin holes or other damage
- Confirm the Use By date has not passed

Then:

- Peel open the pouch of the Delivery System and remove the tray containing the device
- Peel open the pouch of the Loading Kit and remove the tray containing the components
- Unpack the Delivery System and Loading Kit from their trays
- Perform a visual inspection of all components

#### **11.2.4 Valve Rinsing**

The Valve rinsing shall be performed as follows:

- Set up three (3) sterile bowls with 500 ml of fresh sterile saline solution.
- Open the bottle by unscrewing the cap and remove the ACURATE neo2 Aortic Valve from the bottle via the valve holder using tweezers and transfer it to the first rinsing bowl.
- Ensure that the saline solution completely covers the ACURATE neo2 Aortic Valve and the valve holder. Gently stir the Valve for a minimum of 2 minutes, holding it by the valve holder.

- Repeat the rinsing procedure two (2) additional times, always in 500 ml of fresh saline solution and stirring for at least 2 minutes. The last rinse should include 5,000 units of heparin added to the saline solution.
- Carefully remove the valve holder / serial number tag from the Valve by cutting the anchoring thread with a scalpel. Check that no free thread loops remain attached to the Valve.
- Inspect the Valve for any signs of damage to the metallic stent or to the biological component.
- Leave the Valve in the final rinse solution until the time of loading onto the Delivery System in order to prevent the biological tissue from drying.

---

**NOTE:** The recommended rinsing is to remove residual glutaraldehyde-based preserving solution.

---

**CAUTION:** Do not place other objects in the rinse bowls.

---

### 11.3 Loading Procedure

#### 11.3.1 Loading of the Valve onto the Delivery System

---

**WARNING:** The ACURATE neo2 Aortic Valve loading procedure must only be performed by personnel trained by BSC. Other personnel are not authorized to execute the procedure.

---

- Flush the Delivery System
- Place the Valve into the loader within the loading tray filled with refrigerated saline solution
- Align the Delivery System with the Valve / loader
- Insert the tip into the Valve until the stent-holder pins are aligned with the stent fixation loops/hooks
- Load the Valve in refrigerated saline solution into the Delivery System using the rotating handle
- Verify if the pins are engaged

Refer to the training material and Loading Booklet for detailed instructions on how to load the device.

#### 11.3.2 Flushing the Delivery System

- Flush the guidewire lumen with heparinized saline solution
- Flush the Delivery System two to three additional times with heparinized saline solution, tapping gently its shaft to de-bubble
- Verify that the fixation loops/hooks are still engaged
- Inspect for any bubbles

### 11.4 Implantation Procedure

#### 11.4.1 Coronary Occlusion Risk Mitigation

The distance between the base of the coronary cusp and the respective coronary ostium shall be measured pre-operatively by TEE and if available by Angio CT to evaluate the risk of coronary ostia occlusion due to overlapping of a leaflet. **Table 3** indicates the occlusion risk mitigation as a function of that distance:

**Table 3: Occlusion Risk Mitigation**

< 8 mm	8 mm-14 mm	14 mm-23 mm	> 23 mm
Very High	High	Moderate	Low

The presence of large and/or inhomogeneous blocks of calcium is a factor potentially increasing the risk level.

**CAUTION:** Performing a percutaneous coronary intervention (PCI) following Valve implantation may be difficult as a portion of the Valve, if positioned in the vicinity of the coronary ostia, may render access challenging for PCI devices.

#### 11.4.2 Transfemoral Access

- Access and puncture of the femoral arteries shall be performed in a routine fashion.
- Insert and place a pacing wire according to standard practice and test it for pacemaker capture.
- Advance a compatible introducer sheath under fluoroscopic control in a routine fashion. Secure the final position of the introducer to avoid any displacement.

#### 11.4.3 Pre-Dilation of Native Valve

All steps described here below shall be performed under fluoroscopic control and with optional TEE guidance:

- Advance a pigtail catheter into the ascending aorta to perform a supra-aortic angiogram with the projection of the native aortic valve in a perpendicular angle to the screen. The fluoroscopy projection is correctly positioned if all the three native aortic valve cusps are in the same plane.
- Prepare the appropriate Balloon Valvuloplasty Catheter (BVC) according to its Instructions for Use.
- Advance a 275 cm x 0.035 inch (0.89 mm) guidewire retrograde across the native stenotic valve into the left ventricle.
- Perform balloon valvuloplasty during rapid ventricular pacing:
  - Begin rapid pacing between 170 bpm to 220 bpm. Balloon inflation can start once there is no significant pulse pressure (temporary cessation of Left Ventricle [LV] ejection).
  - Verify the desired axial position and rapidly inflate the BVC with the entire content of the inflation syringe, aiming at homogenous inflation of the BVC without dent.
  - Rapidly deflate the BVC until it is completely deflated.
  - After confirming the BVC is completely deflated, turn off the pacemaker and remove the BVC, leaving the 0.035 inch (0.89 mm) guidewire in position.

#### 11.4.4 Valve Delivery

All steps described here below shall be performed under fluoroscopic control and with optional TEE guidance:

- Place the pigtail catheter in at the deepest point in the non-coronary cusp (NCC) as reference for positioning.
- Advance the Valve loaded onto the Delivery System over the 0.035 inch (0.89 mm) guidewire while holding the guidewire. Continue to hold the guidewire throughout the procedure.

- To assist commissural alignment, the orientation may be adjusted by rotating the Delivery System handle. Where torque is built up and/or held in the Delivery System handle, it should be released in advance of proximal deployment.
- Cross the native aortic valve with the Delivery System, until the radio-opaque marker band is positioned at the level of the virtual basal ring of the native stenotic aortic valve.

---

**NOTE:** If the valve cannot be placed in the native aortic valve, perform assessment under fluoroscopy before proceeding. Possible options to consider include, but are not limited to, removal of the valve and introducer sheath as one unit, ectopic placement of the valve, or surgical intervention.

---

**WARNING:** If retraction of the Delivery System with a partially deployed valve is required use caution when attempting to retract the Delivery System with the partially deployed valve around the aortic arch. If resistance is felt, assess the structural integrity of the valve, delivery system, introducer sheath or guidewire before proceeding. Do not attempt to retract the undeployed valve into the introducer sheath.

---

- Under fluoroscopic control, turn the first rotating knob of the release handle counter-clockwise, until full stop. Holding the Delivery System at the level of the positioning sheath will allow a more controlled placement of the Valve. At this stage the stabilization arches are fully deployed. The upper crown is partly deployed as well. The Valve is still anchored to the Delivery System through the fixation loops/hooks of the stent.
- Verify the axial positioning of the Valve with respect to the native valve. The radio-opaque marker band (embedded into the middle member of the Delivery System and located in correspondence of the lower crown) shall be positioned at the level of the virtual basal ring of the native stenotic aortic valve. An extra episode of rapid ventricular pacing may facilitate the final confirmation of the optimal positioning.
- Remove the safety button.
- Performing final Valve deployment under rapid ventricular pacing is optional. In case rapid ventricular pacing is performed, follow the guideline here below:
  - Begin rapid pacing between 170 bpm to 220 bpm. Once there is no significant pulse pressure (temporary cessation of LV ejection) final Valve deployment can start.
- Retrieve the pigtail catheter into the ascending aorta.
- Under fluoroscopic guidance, turn the second rotating knob of the release handle counter-clockwise until full stop, while maintaining forward pressure on the Delivery System.
- Once the Valve is totally deployed, the pacemaker shall be turned off.
- Ensure that the Valve is detached from the Delivery System.

---

**CAUTION:** If not fully detached, ensure that the Delivery System is fully opened by rotating the second knob counter-clockwise until full stop.

---

**CAUTION:** In case the Delivery System comes in contact with the Valve during retraction, partial or full closure of the Delivery System shall be performed.

- Carefully retract the tip and capsule of the Delivery System through the functional Valve, leaving the guidewire in position across the Valve.
- Once crossing the functional Valve, carefully retract the Delivery System into the descending aorta. Under fluoroscopic control, turn the second rotating knob clockwise until it comes to a hard stop and turn the first rotating knob clockwise until the shuttle makes visual contact with the stent holder.
- Carefully retract the Delivery System into the introducer sheath under fluoroscopy, leaving the guidewire in position across the Valve.

#### 11.4.5 Verification of the Valve Position and Post-Implant Monitoring

Leaving the guidewire in position across the Valve, measure both invasive and non-invasive hemodynamic parameters to check positioning and function of the Valve. Perform an angiogram to evaluate device performance and coronary patency after Valve deployment. The use of echocardiographic imaging supports the assessment of the position of the Valve and an evaluation of the para-valvular and intravalvular leakage.

**NOTE:** The effectiveness of the new porcine pericardium sealing skirt in reducing the frequency and severity of paravalvular leakage has been evaluated in a limited clinical data set.

Post-dilation of ACURATE neo2 Aortic Valve is recommended in the presence of significant paravalvular leak.

**CAUTION:** Proceed carefully if a balloon post-dilation is deemed necessary after Valve deployment. The use of a new balloon catheter is recommended. The use of a balloon catheter with a straight shape is required. If a guidewire needs to re-cross the implanted Valve, confirm it has passed between the stabilization arches (not through a stabilization arch) and through the center of the implanted Valve. Failure to confirm correct guidewire crossing can result in migration or embolization of the Valve. Do not use balloon sizes with a maximum diameter (including the tolerance) larger than those specified in **Table 4**.

**Table 4: Maximum Inflated Balloon Diameter**

Valve Catalog Numbers	Maximum Inflated Balloon Diameter Including Tolerance
SYM-SV23-004	22 mm
SYM-SV25-004	24 mm
SYM-SV27-004	26 mm

**NOTE:** Refer to Manufacturer's Instructions for Use for proper sizing and pressure of the balloon.

**NOTE:** The presence of the guidewire crossing the Valve may create transient central leak.

---

- Upon satisfactory results remove the guidewire.
- Measure the transvalvular pressure gradients.
- Remove all catheters and sheaths when the ACT level is appropriate.
- Apply local hemostatic compression on the catheterization puncture sites, or close surgically if clinically indicated.

Administer periprocedural antiplatelet and/or anticoagulant therapy at the discretion of the physician consistent with the local standard-of-care.

### 11.5 Disposal

After use, device and primary packaging may contain biohazardous substances. The device and primary packaging should be treated and disposed of as biohazardous waste or have them treated and disposed of in accordance with any applicable hospital, administrative, and/or local government regulations. Use of a biohazardous container with biological hazard symbol is recommended. Untreated biohazardous waste should not be disposed of in the municipal waste system.

All other packaging should be safely disposed of in accordance with hospital, administrative, and/or local government policy.

### 11.6 Post-Procedure

See ACURATE neo2 Aortic Valve instructions for use for post procedure follow-up recommendations.

Any serious incident that occurs in relation to this device should be reported to the manufacturer and relevant local regulatory authority.

## 12. WARRANTY

For device warranty information, visit ([www.bostonscientific.com/warranty](http://www.bostonscientific.com/warranty)).

The following are trademarks of Boston Scientific Corporation or its affiliates: ACURATE neo2, iSLEEVE and Safari.

All other trademarks are the property of their respective owners.

## 13. SYMBOL DEFINITIONS

Commonly used medical device symbols that appear on the labeling are defined at [www.bostonscientific.com/SymbolsGlossary](http://www.bostonscientific.com/SymbolsGlossary). Additional symbols are defined at the end of this document.

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## INSTRUCTIONS FOR USE

en

The Instructions for Use (IFU) for this product are supplied in electronic form over the Internet.

Visit [www.IFU-BSCI.com](http://www.IFU-BSCI.com) to access the IFU using your internet browser. Information from the product labeling such as Product Code, UDI-DI or Product Name may be used to search for the specific IFU.

If you have difficulty accessing the IFU online, or would prefer to receive a paper copy, please contact Boston Scientific Customer Service or your local country contact. A copy will be sent to you at no charge and should arrive within seven days.

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All rights reserved.
Boston Scientific Limited
Boston Scientific (Australia) Pty Ltd
© 2022 Boston Scientific Corporation or