# ACURATE neo2™

## Aortic Valve

**REF**

SYM-AC-010

SYM-DS-010

SYM-SV23-004

SYM-SV25-004

SYM-SV27-004

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2022-11

## TABLE OF CONTENTS

<b>1. REUSE WARNING</b>	<b>3</b>
<b>2. DEVICE DESCRIPTION</b>	<b>3</b>
2.1 Contents	3
2.2 Operating Principle	3
Figure 1: ACURATE neo2 Aortic Valve	4
2.3 Materials	5
Table 1: Patient Contacting Materials	5
2.4 Non-Pyrogenic	5
Table 2: Device Catalog Numbers	5
Table 3: Compatible Delivery System and Loading Kit	6
2.5 User Information	6
<b>3. INTENDED USE/INDICATIONS FOR USE</b>	<b>7</b>
<b>4. CLINICAL BENEFIT STATEMENT</b>	<b>7</b>
<b>5. SUMMARY OF SAFETY AND CLINICAL PERFORMANCE</b>	<b>7</b>
<b>6. CONTRAINDICATIONS</b>	<b>8</b>
<b>7. WARNINGS</b>	<b>9</b>
<b>8. PRECAUTIONS</b>	<b>10</b>
8.1 Precautions Prior to Use	10
8.2 Precautions During Use	11
<b>9. ADVERSE EVENTS</b>	<b>11</b>
<b>10. MRI SAFETY INFORMATION</b>	<b>13</b>
<b>11. HOW SUPPLIED</b>	<b>14</b>
11.1 Device Details	14
11.2 Handling and Storage	14
Figure 2: Interpreting the Temperature Indicator Display	15

<b>12. OPERATIONAL INSTRUCTIONS</b> .....	<b>15</b>
12.1 Disposal.....	15
12.2 Post-Procedure.....	16
12.3 Implantable Device Patient Information.....	16
12.4 Implant Card Instructions.....	16
Table 4: Understanding symbols on the patient implant card.....	17
<b>13. INFORMATION TO BRIEF THE PATIENT</b> .....	<b>17</b>
<b>14. WARRANTY</b> .....	<b>18</b>
<b>15. SYMBOL DEFINITIONS</b> .....	<b>18</b>

# ACURATE neo2™

## Aortic Valve

### Rx ONLY

**Caution:** Federal Law (USA) restricts this device to sale by or on the order of a physician.

---

#### 1. REUSE WARNING

Contents supplied STERILE using a chemical solution. Do not use if sterile barrier is damaged. If damage is found, call your Boston Scientific representative.

For single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

After use, dispose of packaging in accordance with hospital, administrative and/or local government policy.

---

#### 2. DEVICE DESCRIPTION

##### 2.1 Contents


One (1) ACURATE neo2 Aortic Valve.

##### 2.2 Operating Principle

The ACURATE neo2 Aortic Valve (or Valve) is a bioprosthesis for minimally invasive, transcatheter aortic valve replacement in patients presenting with severe aortic valve stenosis. The Valve is comprised of a radio-opaque, Nitinol self-expanding stent, with

integrated porcine pericardium leaflets and an inner and outer skirt. It is available in three different sizes: S (Small), M (Medium) and L (Large). Key components of the Valve are shown in **Figure 1: ACURATE neo2 Aortic Valve**. **Table 2** indicates the available catalog numbers and the size of the Valve to use in regard to the native annulus diameter of the patient.

Prior to implantation, the Valve is carefully loaded onto the ACURATE neo2 Transfemoral Delivery System using the ACURATE neo2 Loading Kit.

**Figure 1: ACURATE neo2 Aortic Valve**

The biological valve is manufactured from tissue of animal origin (porcine pericardium tissue) that has been preserved in low concentration solutions of glutaraldehyde, preserving its flexibility and strength. The biological tissue for the ACURATE neo2 Aortic Valve has been treated with the BioFix™ calcification-mitigation process, which has demonstrated significant reduction in calcification in small animal models. No clinical data is available showing effectiveness in humans. The Valve is chemically sterilized according to a validated process, involving heat treatment of the tissue in sterilant solution containing glutaraldehyde. The Valve is packaged and stored in low concentration solution of glutaraldehyde. Glutaraldehyde is shown to reduce the

antigenicity of tissue valves and increase tissue stability. However, glutaraldehyde alone has not been shown to affect or reduce the calcification rate of the Valve.

### 2.3 Materials

The device is composed of a self-expanding nitinol stent and a porcine pericardium valve.

Patient contacting materials present in this device are listed in **Table 1**.

**Table 1: Patient Contacting Materials**

Material Name	Weight per component (g)
Nitinol	1.080
Porcine Pericardium	1.826
Polyethylene Terephthalate (PET)	0.039

**WARNING:** Patient contacting materials may cause allergic reaction in some patients.

### 2.4 Non-Pyrogenic

This device meets pyrogen limit specifications.

**Table 2: Device Catalog Numbers**

Catalog Numbers	Product Name	Native Aortic Annulus Diameter
SYM-SV23-004	ACURATE neo2 Aortic Valve S	20.5 mm $\leq$ annulus diameter $\leq$ 23 mm
SYM-SV25-004	ACURATE neo2 Aortic Valve M	22.5 mm $<$ annulus diameter $\leq$ 25 mm
SYM-SV27-004	ACURATE neo2 Aortic Valve L	24.5 mm $<$ annulus diameter $\leq$ 27 mm

**Table 3** indicates the current compatible Delivery System and Loading Kit catalog numbers.

**Table 3: Compatible Delivery System and Loading Kit**

Catalog Numbers	Product Name
SYM-DS-010	ACURATE neo2 Transfemoral Delivery System
SYM-AC-010	ACURATE neo2 Loading Kit

## 2.5 User Information

Implantation of the ACURATE neo2 Aortic Valve should only be performed by physicians trained according to the Boston Scientific Corporation (BSC) Physician Training Plan. Other physicians are not authorized to perform valve implantation.

The ACURATE neo2 Aortic Valve loading procedure must only be performed by personnel trained by BSC. Other personnel are not authorized to execute the procedure.

ACURATE neo2 Aortic Valve System with the Active PVseal™ was designed to minimize paravalvular leakage (PVL), provide excellent hemodynamics, minimize the risk of conduction disturbances and preserve coronary access. The safety and performance of BSC ACURATE neo2 Aortic Valve System and its predecessor ACURATE neo™ Aortic Valve System has been established in multiple BSC-sponsored pre-market and post-market studies in over 1,800 patients. Procedural success was achieved in 97.5%-98.7% of the patients, demonstrating a low rate of procedural complications. Key safety data includes rates of death (30-day: 1.4% - 3.3%; 1-year: 7.96% - 11.9%), stroke (30-day: 1.9% - 2.6% and 1-year: 2.5% - 3.6%), coronary obstruction (30-day: 0.0% - 0.8%) and permanent pacemaker implantation (30-day: 8.3%-15.0%; 1-year: 9.94% -17.8%). Key performance data includes mean aortic gradient (30-day: 6.7 mm Hg - 7.9 mm Hg; 1-year: 6.9 mm Hg - 7.6 mm Hg), paravalvular leak (Moderate/Severe: 30-day: 2.9%; 1-year: 2.3%) and New York Heart Association functional Class (I/II: 30-day:

70.1% - 90.0%; 1-year: 71.6% - 86.9%). This data supports the safety and performance of BSC ACURATE neo2 Aortic Valve System.

To obtain optimal clinical results, proper patient selection, careful sizing and attention to procedural execution is warranted. A more detailed guidance can be found in the following publications: Kim W-K et al. Am J Cardiol 2020;131:12-16 and Okuno T et al. Current Cardiology Reports (2020) 22: 107. ACURATE neo2 is subject to continuous clinical monitoring.

### **3. INTENDED USE/INDICATIONS FOR USE**

The ACURATE neo2 Aortic Valve, in combination with the ACURATE neo2 Transfemoral Delivery System and the ACURATE neo2 Loading Kit, is indicated for relief of aortic stenosis in patients with symptomatic heart disease due to severe native calcific aortic stenosis who are judged by a Heart Team, including a cardiac surgeon, to be appropriate for the transcatheter heart valve replacement therapy.

**FOR CANADA ONLY:** The Indications for Use in this document are not valid in Canada. See Addendum attached to product for Indications for Use

### **4. CLINICAL BENEFIT STATEMENT**

The clinical benefit of ACURATE neo2 Aortic Valve System is to improve aortic valve function with an intent to relieve symptoms and improve survival in patients with symptomatic heart disease due to severe native calcific aortic stenosis who are judged by a heart team to be appropriate for transcatheter heart valve replacement therapy.

### **5. SUMMARY OF SAFETY AND CLINICAL PERFORMANCE**

For customers in the European Union, use the device name found in the labeling to search for the device's Summary of Safety and Clinical Performance, which is available on the European database on medical devices (EUDAMED) website: (<https://ec.europa.eu/tools/eudamed>).

## 6. CONTRAINDICATIONS

The ACURATE neo2 Aortic Valve, in combination with the ACURATE neo2 Transfemoral Delivery System and the ACURATE neo2 Loading Kit, is contraindicated in patients with:

- Non-valvular aortic stenosis
- Congenital aortic stenosis or unicuspid or bicuspid aortic valve
- Presence of a prosthetic mitral valve
- Presence of previously implanted aortic bioprosthesis
- Non-calcific acquired aortic stenosis
- Non-stenotic aortic insufficiency
- Severe eccentricity of calcification
- Evidence of intracardiac mass, thrombus or vegetation
- Severe coagulation problems
- Active bacterial endocarditis or other active infections
- Severe ventricular dysfunction with ejection fraction < 20%
- Inability to tolerate anticoagulation therapy
- Hypertrophic cardiomyopathy with or without obstruction (HOCM)
- Known allergy to patient contacting materials (Nitinol, porcine pericardium, and Polyethylene Terephthalate), aspirin or contrast media
- Distance between the base of the coronary cusp and the respective coronary ostium less than 8 mm
- Leaflet geometry with respect to coronary ostium location presenting a risk of overlapping
- Anatomy NOT appropriate for transfemoral implant due to the size, disease and degree of calcification or tortuosity of the aorta or ilio-femoral arteries

The ACURATE neo2 Aortic Valve, in combination with the ACURATE neo2 Transfemoral Delivery System and the ACURATE neo2 Loading

Kit should not be used if the implanting physician believes its implantation would be contrary to the best interest of the patient.

The ACURATE neo2 Aortic Valve shall not be placed in positions other than in the native aortic valve.

---

## 7. WARNINGS

- Implantation of the Valve should only be performed by physicians trained according to the BSC Physician Training Plan. Other physicians are not authorized to perform valve implantation.
- The Valve loading procedure must only be performed by personnel trained by BSC. Other personnel are not authorized to execute the procedure.
- Correct sizing of the native aortic annulus before use of the Valve is essential to prevent para-valvular leak or migration. The Valve is intended for use in patients with a native aortic annulus size ranging from 20.5 mm to 27 mm.
- Do not use the Valve if the expiry date has passed or the temperature indicator has tripped. (See **Handling and Storage**)
- Do not use the Valve if the glutaraldehyde-based preserving solution does not completely cover the Valve.
- The outer surface of the Valve bottle is not sterile and must not be placed in the sterile field.
- Do not add or apply antibiotics, drugs or chemical substances to the preserving solution, to the rinse solution (other than heparin in the final rinse) or to the Valve as these could modify the biological tissue status.
- If the Valve remains unused after being removed from its primary packaging it shall be considered non-sterile and must not be used.

- Handle the Valve with maximum care and in accordance with the ACURATE neo2 Transfemoral Delivery System/Loading Kit Instructions for Use, to avoid damaging the stent and/or biological tissue leaflets.
- Do not load the Valve on the Delivery System more than twice as this may cause damage to the stent and/or the biological tissue.
- The Valve shall be kept moist through all the preparation procedure to prevent the biological tissue from drying. If dehydrated, the biological tissue suffers irreversible damage and the valve must not be used. If necessary, keep the Valve moist with a sterile saline solution during the procedure.
- Post-dilation of the Valve could damage the device integrity or cause migration of the Valve. Proceed with caution if it is necessary to post-dilate the Valve. Ensure post-dilation balloon shape, dimensions and tolerances are suitable for the Valve.
- Valve recipients should be maintained on antiplatelet therapy except where contraindicated, as determined by their physician.
- Accelerated deterioration due to calcific degeneration of the Valve (as with any glutaraldehyde crosslinked bioprosthesis) may occur in patients with altered calcium metabolism.

Please refer to the Instructions for Use of the ACURATE neo2 Transfemoral Delivery System/Loading Kit for further warnings related to the Delivery System.

---

## 8. PRECAUTIONS

### 8.1 Precautions Prior to Use

- It is recommended that all ACURATE neo2 Aortic Valve recipients are prophylactically treated for endocarditis to minimize the possibility of prosthetic valve infection.
- Do not use the ACURATE neo2 Aortic Valve if it is observed to be damaged or dehydrated.

## 8.2 Precautions During Use

- The Valve shall be handled in such a way as to avoid contact with fibers or foreign bodies in which adherence to the Valve could induce emboli or undesired reactions with the blood.
- The Valve shall be handled using aseptic technique. The preserving solution is not suitable for Valve re-sterilization.
- Glutaraldehyde may cause irritation of the skin, eyes, nose and throat. Avoid prolonged or repeated exposure or breathing of the solution. Use only with adequate ventilation. If skin contact occurs, immediately flush the affected area with water. In the event of contact with eyes, seek immediate medical attention.
- Implantation of the Valve shall be preceded by dilation of the stenotic native aortic valve by means of balloon aortic valvuloplasty.
- The long-term durability of transcatheter heart valves is not known at this time. Regular medical follow-up is recommended to evaluate how well a patient's heart valve is performing.

Please refer to the Instructions for Use of the ACURATE neo2 Transfemoral Delivery System/Loading Kit for further precautions related to the Delivery System.

## 9. ADVERSE EVENTS

Risks, complications and adverse events that may be associated with the use of the ACURATE neo2 Aortic Valve, in combination with the ACURATE neo2 Transfemoral Delivery System and the ACURATE neo2 Loading Kit include those risks related to conventional surgical aortic valve replacement (SAVR), as well as those related to transcatheter aortic valve implantation (TAVI).

Known or anticipated risks are listed below in alphabetical order:

- Allergic or adverse reaction (including medications, anesthesia, contrast, or device materials)
- Angina
- Arrhythmia or new conduction system injury (including need for pacemaker insertion)

- Bleeding including hemorrhage or hematoma (possibly requiring transfusion or additional intervention)
- Cardiac failure leading to low cardiac output (cardiogenic shock) or pulmonary edema
- Cerebrovascular accident, stroke, transient ischemic attack or cerebral infarction including asymptomatic neuroimaging findings
- Coronary obstruction
- Death
- Device misplacement, migration or embolization
- Emboli (including air, tissue, thrombus or device materials)
- Fever and pyrogen reaction
- Heart failure
- Hemolysis and/or hemolytic anemia
- Hypertension/hypotension
- Infection (local or systemic including endocarditis)
- Mitral valve insufficiency or injury
- Myocardial infarction
- Myocardial or valvular injury (including perforation or rupture)
- Nerve injury or neurologic deficits (including encephalopathy)
- Pain or inflammation
- Pericarditis, pericardial effusion or tamponade
- Peripheral ischemia or infarction
- Pleural effusion
- Radiation injury
- Renal insufficiency or failure
- Respiratory insufficiency or failure
- Valve dysfunction, deterioration or failure
- Valve or device thrombosis
- Valvular stenosis or regurgitation (central or paravalvular)

- Vessel injury (including access site) such as spasm, lymphatic problems, pseudoaneurysm, arteriovenous fistula, trauma, dissection, perforation, and rupture

As a result of these adverse events, the subject may require medical, percutaneous or surgical intervention, including implantation of a second valve, re-operation and replacement of the valve. These events may lead to fatal outcomes.

## 10. MRI SAFETY INFORMATION

Non-clinical testing has demonstrated that the ACURATE Aortic Valve is Magnetic Resonance (MR) Conditional. A patient with this device can be safely scanned in an MR system meeting the following conditions:

- Static magnetic field of 1.5 Tesla and 3 Tesla, **with**
- Maximum spatial field gradient of 13,200 G/cm (132 T/m)
- **Maximum force product of 238,000,000 G<sup>2</sup>/cm (238 T<sup>2</sup>/m)**
- **Theoretically estimated** maximum whole body averaged (WBA) specific absorption rate (SAR) of 2 W/kg (Normal Operating Mode)

If information about a specific parameter is not included, there are no conditions associated with that parameter.

Under the scan conditions defined above, the ACURATE Aortic Valve is expected to produce a maximum temperature rise of less than 1.8°C (**2 W/kg, 1.5 Tesla**) RF-related temperature increase and 1.5°C (**2 W/kg, 3 Tesla**) RF-related temperature increase after 15 minutes of continuous scanning.

In non-clinical testing, the image artifact caused by the device extends approximately 7.1 mm from the ACURATE Aortic Valve when imaged with a gradient echo pulse sequence and a 3 Tesla MR system.

MR image quality may be compromised if the area of interest is in the exact same area or relatively close to the position of the Valve.

## 11. HOW SUPPLIED

### 11.1 Device Details

The ACURATE neo2 Aortic Valve is delivered sterile and non-pyrogenic with a serial number tag and valve holder attached to the stabilization arches.

The primary packaging consists of a bottle closed with a sealed screw cap and contains the Valve immersed in a glutaraldehyde based preserving solution with a temperature indicator visible on the top of the carton. The temperature indicator serves to detect out-of-range temperature conditions during transportation.

The secondary packaging consists of two foam protective shells, to preserve the primary packaging from mechanical shocks, and a carton.

Do not use if package is damaged or unintentionally opened before use.

Do not use if labeling is incomplete or illegible.

### 11.2 Handling and Storage

The ACURATE neo2 Aortic Valve shall be stored at a temperature in the range 5 °C (41 °F) to 25 °C (77 °F). The lower temperature limit shall be observed carefully since at temperatures close to 0 °C (32 °F) the preserving solution begins to freeze, and this could cause irreversible damage to the biological tissue.

Do not leave the Valve in areas exposed to direct sunlight or close to heating or air conditioning equipment. Exposure to sunlight encourages changes in the preserving solution properties.

An activated temperature indicator shows that the Valve has been exposed to an out-of-range temperature during transportation. Only use the Valve if the indicator clearly displays OK. See **Figure 2: Interpreting the Temperature Indicator Display.**

Out-of-Range Temperature: Do not use product	OK to use product
 <p>DO NOT use product</p>	 <p>OK to use product</p>

**Figure 2: Interpreting the Temperature Indicator Display**

**NOTE:** The purpose of the temperature indicator is to monitor out-of-range temperatures during transportation. It is not intended for monitoring temperature conditions during the shelf life of the ACURATE neo2 Aortic Valve.

## 12. OPERATIONAL INSTRUCTIONS

Please refer to the Instructions for Use for the ACURATE neo2 Transfemoral Delivery System/Loading Kit for details on additional items for safe use, preparation, loading and implantation of the device.

### 12.1 Disposal

After use, device and primary packaging may contain biohazardous substances. The device and primary packaging should be treated and disposed of as biohazardous waste or have them treated and disposed of in accordance with any applicable hospital, administrative, and/or local government regulations. Use of a biohazardous container with biological hazard symbol is recommended. Untreated biohazardous waste should not be disposed of in the municipal waste system.

Do not incinerate any component containing a battery and/or electronics. Improper disposal could result in an explosion.

All other packaging should be safely disposed of in accordance with hospital, administrative, and/or local government policy.

### **12.2 Post-Procedure**

Prescribe antiplatelet therapy consistent with local standard-of-care. Post-procedure follow-up is recommended at 30 days, 1 year, and annually thereafter, or per the local standard-of-care.

Please keep a patient record of the product number (UPN) and serial number of the device implanted. Complete the Implant Card as per the instructions in Section 12.4 and provide the completed implant card to the patient.

Any serious incident that occurs in relation to this device should be reported to the manufacturer and relevant local regulatory authority.

### **12.3 Implantable Device Patient Information**

Advise the patient that additional information may be available to them on the Boston Scientific website ([www.bostonscientific.com/patientlabeling](http://www.bostonscientific.com/patientlabeling)).

### **12.4 Implant Card Instructions**

Remove the perforated hang tab from the implant card before giving to the patient. Apply a peel-off label from the product onto the supplied patient implant card and complete the required content on the card. See **Table 4**.

**Table 4: Understanding symbols on the patient implant card.**

Symbol	Action
	Provide date of device implantation
	Provide healthcare facility and/or physician information.
	Provide patient name

Give the completed implant card to the patient.

### **13. INFORMATION TO BRIEF THE PATIENT**

Ensure the patient is informed of the following important points:

- Contraindications, warnings, precautions and potential adverse events that are directly relevant to the patient.
- Post-procedure medication regimen and necessary follow-up visits as per the local standard of care.
- Complying with the follow-up visits and notifying their health care provider of any possible symptoms and signs that could potentially indicate aortic valve deterioration or malfunction (e.g. heart failure, myocardial infarction, stroke).
- The increased risk of infection with the implanted aortic valve and the need to keep the patient implant card with them at all times to inform their health care provider prior to receiving any invasive procedure or MRI examination.
- The device is a permanent implant and has been tested for structural performance (device integrity and component fatigue) for up to 5 years. The long-term durability of transcatheter heart valves is not known at this time. Regular

medical follow-up is recommended to evaluate how well a patient's heart valve is performing.

#### **14. WARRANTY**

For device warranty information, visit ([www.bostonscientific.com/warranty](http://www.bostonscientific.com/warranty)).

ACURATE neo2, ACURATE neo, Active PVseal and BioFix are trademarks of Boston Scientific Corporation or its affiliates.

All other trademarks are the property of their respective owners.

#### **15. SYMBOL DEFINITIONS**

Commonly used medical device symbols that appear on the labeling are defined at [www.bostonscientific.com/SymbolsGlossary](http://www.bostonscientific.com/SymbolsGlossary). Additional symbols are defined at the end of this document.

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## INSTRUCTIONS FOR USE

en

The Instructions for Use (IFU) for this product are supplied in electronic form over the Internet.

Visit [www.IFU-BSCI.com](http://www.IFU-BSCI.com) to access the IFU using your internet browser. Information from the product labeling such as Product Code, UDI-DI or Product Name may be used to search for the specific IFU.

If you have difficulty accessing the IFU online, or would prefer to receive a paper copy, please contact Boston Scientific Customer Service or your local country contact. A copy will be sent to you at no charge and should arrive within seven days.

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All rights reserved.
Boston Scientific Limited
Boston Scientific (Australia) Pty Ltd
© 2022 Boston Scientific Corporation or its affiliates.