# ACURATE Prime™

## Aortic Valve

**REF**

H74939690230

H749396822325

H74939690250

H749396822729

H74939690270

H749396942325

H74939690290

H749396942729

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2024-05

## TABLE OF CONTENTS

<b>1. REUSE WARNING</b> .....	<b>3</b>
<b>2. DEVICE DESCRIPTION</b> .....	<b>3</b>
2.1 Contents.....	3
2.2 Operating Principle .....	3
Figure 1: ACURATE Prime Aortic Valve .....	4
Table 1. Valve Sizing Chart with Compatible Delivery System and Loading Kit.....	5
Table 2. 23 and 25 mm Valve Sizes and Annulus Compatibility Information .....	5
Table 3. 27 and 29 mm Valve Sizes and Annulus Compatibility Information .....	6
2.3 Materials.....	6
Table 4: Patient Contacting Materials.....	6
2.4 Non-Pyrogenic .....	6
2.5 User Information .....	6
<b>3. INTENDED USE/INDICATIONS FOR USE</b> .....	<b>7</b>
<b>4. CLINICAL BENEFIT STATEMENT</b> .....	<b>7</b>
<b>5. SUMMARY OF SAFETY AND CLINICAL PERFORMANCE</b> .....	<b>7</b>
<b>6. CONTRAINDICATIONS</b> .....	<b>8</b>
<b>7. WARNINGS</b> .....	<b>8</b>
<b>8. PRECAUTIONS</b> .....	<b>10</b>
8.1 General Precautions.....	10
8.2 Precautions Prior to Use.....	11
8.3 Precautions During Use.....	11
<b>9. ADVERSE EVENTS</b> .....	<b>12</b>
<b>10. MRI SAFETY INFORMATION</b> .....	<b>14</b>
<b>11. HOW SUPPLIED</b> .....	<b>15</b>
11.1 Interpreting the Activated Temperature Indicator .....	15

Figure 2: Interpreting the Temperature Indicator Display.....	16
11.2 Handling and Storage.....	16
<b>12. OPERATIONAL INSTRUCTIONS.....</b>	<b>16</b>
12.1 Device Disposal.....	16
12.2 Post-Procedure.....	17
12.3 Implantable Device Patient Information.....	17
12.4 Implant Card Instructions.....	17
Table 5: Understanding Symbols on the Patient Implant Card.....	18
<b>13. INFORMATION TO BRIEF THE PATIENT.....</b>	<b>18</b>
<b>14. WARRANTY.....</b>	<b>19</b>
<b>15. SYMBOL DEFINITIONS.....</b>	<b>19</b>

# ACURATE Prime™

## Aortic Valve

### Rx ONLY

**CAUTION:** Federal Law (USA) restricts this device to sale by or on the order of a physician.

---

#### 1. REUSE WARNING

Contents supplied STERILE using a chemical solution. Do not use if sterile barrier is damaged. If damage is found, call your Boston Scientific representative.

For single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

After use, dispose of packaging in accordance with hospital, administrative and/or local government policy.

---

#### 2. DEVICE DESCRIPTION


##### 2.1 Contents

One (1) ACURATE Prime Aortic Valve.

##### 2.2 Operating Principle

The ACURATE Prime Aortic Valve (or Valve) is a bioprosthesis for minimally invasive, transcatheter aortic valve replacement in adult patients presenting with severe aortic valve stenosis. The Valve is comprised of a radio-opaque, self-expanding Nitinol stent, with

integrated porcine pericardium leaflets and an inner and outer skirt. Key components of the Valve are shown in **Figure 1**.

- [A] Lower Crown [B] Upper Crown [C] Commissure Post
- [D] Stabilization Arches [E] Self-Expanding Nitinol Stent
- [F] Supra-Annular Porcine Pericardial Leaflets
- [G] Porcine Pericardial Sealing Skirt [H] Fixation Loops/Hooks

### Figure 1: ACURATE Prime Aortic Valve

The biological Valve is manufactured from tissue of animal origin (porcine pericardium) that has been preserved in low concentration solutions of glutaraldehyde, preserving its flexibility and strength. The biological tissue for the ACURATE Prime Aortic Valve has been treated with the BioFix™ calcification-mitigation process, which has demonstrated significant reduction in calcification in small animal models. No clinical data is available showing effectiveness in humans. The Valve is chemically sterilized according to a validated process, involving heat treatment of the tissue in sterilant solution containing glutaraldehyde. The Valve is packaged and stored in a low concentration solution of glutaraldehyde. Glutaraldehyde is shown to reduce the antigenicity of tissue valves and increase tissue stability. However, glutaraldehyde alone has not been shown to affect or reduce the calcification rate of the Valve.

The ACURATE Prime Aortic Valve is available in four sizes: 23 mm, 25 mm, 27 mm and 29 mm. **Table 1** details Valve compatibility with the Delivery System and Loading Kit.

**Table 1. Valve Sizing Chart with Compatible Delivery System and Loading Kit**

	ACURATE Prime Aortic Valve Size			
	23 mm	25 mm	27 mm	29 mm
<b>Delivery System and Loading Kit Carton Dimensions and Flag Color*</b>	23 mm   25 mm Green		27 mm   29 mm Yellow	

\*See product carton for dimensions and flag color.

**Table 2** and **Table 3** indicate available catalog numbers, product name, and the corresponding valve size in CT-based measurements for both native aortic annulus perimeter and perimeter-derived diameter. Prior to implantation, the Valve is carefully loaded onto ACURATE Prime Delivery System using ACURATE Prime Loading Kit.

**Table 2. 23 and 25 mm Valve Sizes and Annulus Compatibility Information**

	ACURATE Prime Aortic Valve Size	
	23 mm	25 mm
<b>Catalog Number</b>	H74939690230	H74939690250
<b>Native Aortic Valve Perimeter</b>	64 mm - 72 mm	71 mm - 79 mm
<b>Perimeter-Derived Diameter</b>	20.5 mm - 23 mm	22.5 mm - 25 mm

**Table 3. 27 and 29 mm Valve Sizes and Annulus Compatibility Information**

	ACURATE Prime Aortic Valve Size	
	27 mm	29 mm
Catalog Number	H74939690270	H74939690290
Native Aortic Valve Perimeter	77 mm - 85 mm	83 mm - 91 mm
Perimeter-Derived Diameter	24.5 mm - 27 mm	26.5 mm - 29 mm

### 2.3 Materials

The device is composed of a self-expanding nitinol stent and a porcine pericardium valve comprised of three (3) porcine pericardium leaflets, and a double porcine pericardium skirt. Patient contacting materials present in this device are listed in **Table 4**.

**Table 4: Patient Contacting Materials**

Material Name	Weight per component (g)
Nitinol	≤ 1,204
Porcine Pericardium	≤ 1.413
Polyethylene Terephthalate (PET)	≤ 0.041

**WARNING:** Patient contacting materials may cause allergic reaction in some patients.

### 2.4 Non-Pyrogenic

This device meets pyrogen limit specifications.

### 2.5 User Information

Implantation of the ACURATE Prime Aortic Valve with the ACURATE Prime Delivery System shall only be performed by physicians trained according to the Boston Scientific Corporation (BSC)

Physician Training Plan. Other physicians are not authorized to perform Valve implantation.

The ACURATE Prime Aortic Valve loading procedure shall only be performed by personnel trained and certified according to the BSC Training Plan. Other personnel are not authorized to execute the Valve loading procedure.

### **3. INTENDED USE/INDICATIONS FOR USE**

The ACURATE Prime Aortic Valve, in combination with the ACURATE Prime Delivery System and the ACURATE Prime Loading Kit, is indicated for relief of aortic stenosis in patients with symptomatic heart disease due to severe native calcific aortic stenosis who are judged by a Heart Team, including a cardiac surgeon, to be appropriate for the transcatheter heart valve replacement therapy.

**FOR CANADA ONLY:** The Indications for Use in this document are not valid in Canada. See Addendum attached to product for Indications for Use

### **4. CLINICAL BENEFIT STATEMENT**

The clinical benefit of ACURATE Prime Aortic Valve System is to improve aortic valve function with an intent to relieve symptoms and improve survival in patients with symptomatic heart disease due to severe native calcific aortic stenosis who are judged by a heart team to be appropriate for transcatheter heart valve replacement therapy.

### **5. SUMMARY OF SAFETY AND CLINICAL PERFORMANCE**

For customers in the European Union, use the device name found in the labeling to search for the device's Summary of Safety and Clinical Performance, which is available on the European database on medical devices (EUDAMED) website: (<https://ec.europa.eu/tools/eudamed>).

## 6. CONTRAINDICATIONS

The ACURATE Prime Aortic Valve, in combination with the ACURATE Prime Delivery System and the ACURATE Prime Loading Kit, is contraindicated in patients with:

- Anatomy that would preclude safe transfemoral access of the bioprosthesis due to the size, disease and degree of calcification or tortuosity of the aorta or ilio-femoral arteries.
- Evidence of intracardiac mass, thrombus, vegetation, active systemic infection, or endocarditis.
- Intolerance to anticoagulation/antiplatelet regimen; known hypersensitivity or contraindication to patient contacting materials (See **Table 4** Patient Contacting Materials), or sensitivity to contrast media which cannot be adequately premedicated.

The ACURATE Prime Aortic Valve, in combination with the ACURATE Prime Delivery System and the ACURATE Prime Loading Kit should not be used if the implanting physician believes its implantation would be contrary to the best interest of the patient.

---

## 7. WARNINGS

- Implantation of the ACURATE Prime Aortic Valve shall only be performed by physicians trained according to the BSC Training Plan. Other physicians are not authorized to perform valve implantation.
- The Valve loading procedure shall only be performed by personnel trained according to the BSC Training plan. Other personnel are not authorized to perform valve loading.
- Do not attempt to load an ACURATE neo2 Aortic Valve onto an ACURATE Prime Delivery System.
- Correct sizing of the native aortic annulus before use of the Valve is essential to prevent paravalvular regurgitation or device misplacement (migration/embolization). The Valve is

intended for use in adult patients with a native aortic annulus diameter size ranging from 20.5 mm to 29 mm.

- Do not use the Valve if the expiry date is passed or the temperature indicator is activated (See Handling and Storage).
- Do not use the Valve if the glutaraldehyde-based preserving solution does not completely cover the Valve.
- The outer surface of the Valve bottle is not sterile and must not be placed in the sterile field.
- Do not add or apply antibiotics, drugs or chemical substances to the preserving solution, to the rinse solution (other than heparin in the final rinse) or to the Valve as these could modify the biological tissue status.
- If Valve remains unused after being removed from primary packaging it shall be considered non-sterile and must be disposed of (See Device Disposal).
- Handle the Valve with care and in accordance with the ACURATE Instructions for Use to avoid damaging the stent and/or biological tissue.
- Do not load the Valve on the Delivery System more than twice as this may cause damage to the stent and/or the biological tissue. Replace Valve, Delivery System and Loading Kit after two (2) loading attempts.
- The Valve shall be kept cool and moist through all the preparation and loading procedure. If dehydrated, the biological tissue suffers irreversible damage, and the valve must not be used.
- Incorrect post-dilatation technique of the ACURATE Prime Aortic Valve could damage the device integrity or cause migration of the Valve. Proceed with caution if it is necessary to post-dilate the ACURATE Prime Aortic Valve. The use of a new balloon catheter is recommended. The use of the balloon catheter with a straight shape is required.

- If performing post-dilation of the Valve, ensure post-dilation balloon shape, dimensions and tolerances are suitable for the Valve and patient anatomy. Post-dilation of the Valve can damage the device integrity or cause migration of the Valve.
- Ensure full detachment of the Valve from the stent holder prior to retraction of the ACURATE Prime Delivery System to prevent device misplacement (migration/embolization).
- Anti-platelet therapy is recommended after TAVI to decrease the risk of thrombotic or thromboembolic complications if there are no contraindications to these medications, per local standard of care.
- Accelerated deterioration due to calcific degeneration of the Valve (as with any glutaraldehyde crosslinked bioprosthesis) may occur in patients with altered calcium metabolism.

---

Please refer to the Instructions for Use of the ACURATE Prime Delivery System/Loading Kit for further warnings related to the Delivery System and Loading Kit.

---

## 8. PRECAUTIONS

### 8.1 General Precautions

The safety and effectiveness of the bioprosthesis for aortic valve replacement have not been evaluated in patient populations presenting with the following conditions:

- Unicuspid or bicuspid aortic valve.
- Pre-existing prosthetic aortic or mitral valve.
- Pure aortic regurgitation.
- Severe ventricular dysfunction with ejection fraction < 20%.
- Mixed native aortic valve disease (aortic stenosis and severe aortic regurgitation).
- Moderate or severe mitral stenosis.

- Severe mitral or tricuspid regurgitation.
- Hypertrophic cardiomyopathy.

### 8.2 Precautions Prior to Use

- It is recommended that all ACURATE Prime Aortic Valve recipients are prophylactically treated for endocarditis to minimize the possibility of prosthetic valve infection.
- Do not use the ACURATE Prime Aortic Valve if it is observed to be damaged or dehydrated.
- Device materials may cause sensitivity in some patients.

### 8.3 Precautions During Use

- The implant procedure shall be conducted under fluoroscopic guidance. The correct fluoroscopy projection for implantation of the valve is when all three (3) native aortic valve cusps are in the same plane.
- The Valve shall be handled in such a way as to avoid contact with fibers or foreign bodies in which adherence to the Valve could induce emboli or undesired reactions with the blood.
- The Valve shall be handled using aseptic technique. The preserving solution is not suitable for Valve re-sterilization.
- Glutaraldehyde may cause irritation of the skin, eyes, nose, and throat. Avoid prolonged or repeated exposure or breathing of the solution. Use only with adequate ventilation. If skin contact occurs, immediately flush the affected area with water. In the event of contact with eyes, seek immediate medical attention.
- Predilatation of the stenotic native aortic valve by means of balloon aortic valvuloplasty is recommended prior to implantation of the prosthesis.
- ACURATE Prime Aortic Valve recipients shall receive during the implant procedure anticoagulant therapy (e.g., unfractionated heparin) per local standard of care, with a recommended target activated clotting time of  $\geq 250$  seconds during the index procedure as determined by their physician.

- The long-term durability of transcatheter heart valves is not known at this time. Regular medical follow-up is recommended to evaluate how well a patient's heart valve is performing.

Please refer to the Instructions for Use of the ACURATE Prime Delivery System/Loading Kit for further precautions related to the Delivery System.

## 9. ADVERSE EVENTS

Risks, complications and adverse events that may be associated with the use of the ACURATE Prime Aortic Valve, in combination with the ACURATE Prime Delivery System and the ACURATE Prime Loading Kit may include but are not limited to the following:

- Allergic or adverse reaction (including medications, anesthesia, contrast, or device materials).
- Angina.
- Arrhythmia or new conduction system injury (including need for pacemaker insertion).
- Bleeding including hemorrhage or hematoma (possibly requiring transfusion or additional intervention).
- Cardiac failure leading to low cardiac output (cardiogenic shock) or pulmonary edema.
- Cerebrovascular accident, stroke, transient ischemic attack or cerebral infarction including asymptomatic neuroimaging findings.
- Coronary obstruction.
- Death.
- Device misplacement, migration or embolization.
- Emboli (including air, tissue, thrombus or device materials).
- Fever and pyrogen reaction.
- Heart failure.
- Hemolysis and/or hemolytic anemia.

- Hypertension/hypotension.
- Infection (local or systemic including endocarditis).
- Mitral valve insufficiency or injury.
- Myocardial infarction.
- Myocardial or valvular injury (including perforation or rupture).
- Nerve injury or neurologic deficits (including encephalopathy).
- Pain or inflammation.
- Pericarditis, pericardial effusion or tamponade.
- Peripheral ischemia or infarction.
- Pleural effusion.
- Radiation injury.
- Renal insufficiency or failure.
- Respiratory insufficiency or failure.
- Valve dysfunction, deterioration or failure.
- Valve or device thrombosis.
- Valvular stenosis or regurgitation (central or paravalvular).
- Vessel injury (including access site) such as spasm, lymphatic problems, pseudoaneurysm, arteriovenous fistula, trauma, dissection, perforation, and rupture.

As a result of these adverse events, the subject may require medical, percutaneous or surgical intervention, including implantation of a second Valve, re-operation and replacement of the Valve. These events may lead to fatal outcomes.

## 10. MRI SAFETY INFORMATION

Non-clinical testing has demonstrated that the ACURATE Prime Aortic Valve is Magnetic Resonance (MR) Conditional. A patient with this device can be safely scanned in an MR system meeting the following conditions:

- Static magnetic field of 1.5 Tesla and 3.0 Tesla, **with**
- Maximum spatial field gradient of 4,000 G/cm (40 T/m)
- **Maximum force product of 238,000,000 G<sup>2</sup>/cm (238 T<sup>2</sup>/m)**
- **Theoretically estimated** maximum whole body averaged (WBA) specific absorption rate (SAR) of 2 W/kg (Normal Operating Mode)

If information about a specific parameter is not included, there are no conditions associated with that parameter.

Under the scan conditions defined above, the ACURATE Prime Aortic Valve is expected to produce a maximum temperature rise of less than 1.8°C (**2 W/kg, 1.5 Tesla**) **RF-related temperature increase** and 1.5°C (**2 W/kg, 3.0 Tesla**) **RF-related temperature increase** after 15 minutes of continuous scanning.

In non-clinical testing, the image artifact caused by the device extends approximately 7.1 mm from the ACURATE Prime Aortic Valve when imaged with a gradient echo pulse sequence and a 3.0 Tesla MR system.

MR image quality may be compromised if the area of interest is in the exact same area or relatively close to the position of the Valve.

## 11. HOW SUPPLIED

The ACURATE Prime Aortic Valve is delivered sterile and non-pyrogenic with a serial number tag and valve holder attached to the stabilization arches.

The primary packaging consists of a labelled bottle closed with a sealed screw cap and contains the Valve immersed in a glutaraldehyde-based preserving solution.

The secondary packaging consists of two (2) foam protective shells, to preserve the primary packaging from mechanical shocks, and a labelled carton with a color-coded flag that includes valve diameter size information to facilitate device selection and has a temperature indicator visible on the top of the carton. The temperature indicator serves to detect out-of-range temperature conditions during transportation.

Do not use if package is damaged or unintentionally opened before use.

Do not use if labeling is incomplete or illegible.

### 11.1 Interpreting the Activated Temperature Indicator

Use the Valve if the temperature indicator has not been activated. Check display for "OK" in black lettering on a plain grey background. See **Figure 2 "OK to use product"** for a temperature indicator that has not been activated. This is the only acceptable display for the temperature indicator. All other displays are not acceptable, and the valve must not be used.

---

**WARNING:** Do not use the Valve if the temperature indicator has been activated. Activated temperature indicators can have many different displays. Check display for "OK" in black and grey lettering against a black background, with or without a combination of numbers 1 through 4. Also check for a blank display. See **Figure 2 "Out-of-Range Temperature: Do not use product"** for examples of temperature indicators that have been activated and are not acceptable to use.

---

**NOTE:** Additional activated temperature indicator displays can be observed that are not shown in **Figure 2**, i.e., where the numbers displayed differ from what is shown in "Out-of-Range Temperature: Do not use product".

Out-of-Range Temperature: Do not use product	OK to use product
 <p>DO NOT use product</p>	 <p>OK to use product</p>

**Figure 2: Interpreting the Temperature Indicator Display**

## 11.2 Handling and Storage

The ACURATE Prime Aortic Valve shall be stored in its product carton at a temperature in the range 5 °C (41 °F) to 25 °C (77 °F). Do not leave the Valve in areas exposed to direct sunlight or close to heating or air conditioning equipment. Exposure to sunlight encourages changes in the preserving solution properties.

## 12. OPERATIONAL INSTRUCTIONS

Please refer to the Instructions for Use for the ACURATE Prime Delivery System/Loading Kit for details on additional items for safe use, preparation, loading and implantation of the device.

### 12.1 Device Disposal

After use, device and primary packaging may contain biohazardous substances. The device and primary packaging should be treated and disposed of as biohazardous waste or have them treated

and disposed of in accordance with any applicable hospital, administrative, and/or local government regulations. Use of a biohazardous container with biological hazard symbol is recommended. Untreated biohazardous waste should not be disposed of in the municipal waste system.

Do not incinerate any component containing a battery and/or electronics. Improper disposal could result in an explosion.

All other packaging should be safely disposed of in accordance with hospital, administrative, and/or local government policy.

### 12.2 Post-Procedure

Prescribe antiplatelet therapy consistent with local standard-of-care. Post-procedure follow-up is recommended at 30 days, 1 year, and annually thereafter, or per the local standard-of-care.

Please keep a patient record of the product number (UPN) and serial number of the device implanted. Complete the Implant Card as per the instructions in **Section 12.4** and provide the completed implant card to the patient.

Any serious incident that occurs in relation to this device should be reported to the manufacturer and relevant local regulatory authority.

### 12.3 Implantable Device Patient Information

Advise the patient that additional information may be available to them on the Boston Scientific website ([www.bostonscientific.com/patientlabeling](http://www.bostonscientific.com/patientlabeling)).

### 12.4 Implant Card Instructions

Remove perforated hang tab from implant card before giving to patient. Apply a peel-off label from the product packaging onto the supplied patient implant card and complete the required content on the card. See **Table 5**.

**Table 5: Understanding Symbols on the Patient Implant Card**

Symbol	Action
	Provide date of device implantation.
	Provide healthcare facility and/or physician information.
	Provide patient name.

Give completed implant card to the patient.

### 13. INFORMATION TO BRIEF THE PATIENT

Ensure patient is informed of the following important points:

- Contraindications, warnings, precautions, and potential adverse events that are directly relevant to the patient.
- Post-procedure medication regimen and necessary follow-up visits as per the local standard of care.
- Complying with the follow-up visits and notifying their health care provider of any possible symptoms and signs that could potentially indicate aortic valve deterioration or malfunction (e.g. heart failure, myocardial infarction, stroke).
- The increased risk of infection with the implanted aortic valve and the need to always keep the patient implant card with them to inform their health care provider prior to receiving any invasive procedure or MRI examination.
- The device is a permanent implant and has been tested for structural performance (device durability and component fatigue) for up to five (5) years. The long-term durability of transcatheter heart valves beyond five (5) years is not known.

at this time. Regular medical follow-up is recommended to evaluate how well a patient's heart valve is performing.

#### **14. WARRANTY**

For device warranty information, visit ([www.bostonscientific.com/warranty](http://www.bostonscientific.com/warranty)).

ACURATE neo2, ACURATE Prime and BioFix are trademarks of Boston Scientific Corporation or its affiliates.

All other trademarks are the property of their respective owners.

#### **15. SYMBOL DEFINITIONS**

Commonly used medical device symbols that appear on the labeling are defined at [www.bostonscientific.com/SymbolsGlossary](http://www.bostonscientific.com/SymbolsGlossary). Additional symbols are defined at the end of this document.

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## INSTRUCTIONS FOR USE

en

The Instructions for Use (IFU) for this product are supplied in electronic form over the Internet.

Visit [www.IFU-BSCI.com](http://www.IFU-BSCI.com) to access the IFU using your internet browser. Information from the product labeling such as Product Code, UDI-DI or Product Name may be used to search for the specific IFU.

If you have difficulty accessing the IFU online, or would prefer to receive a paper copy, please contact Boston Scientific Customer Service or your local country contact. A copy will be sent to you at no charge and should arrive within seven days.

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All rights reserved.
Boston Scientific Limited
Boston Scientific (Australia) Pty Ltd
© 2024 Boston Scientific Corporation or its affiliates.