# IceRod™

1.5 CX 90°

# IceSphere™

1.5 CX 90°

# IcePearl™

2.1 CX STRAIGHT  
2.1 CX 90°  
2.1 CX L 90°

# IceForce™

2.1 CX STRAIGHT  
2.1 CX 90°  
2.1 CX L 90°

## Cryoablation Needle

|       |                          |    |
|-------|--------------------------|----|
| en    | Instructions for Use     | 2  |
| es    | Instrucciones de uso     | 12 |
| fr    | Mode d'emploi            | 21 |
| de    | Gebrauchsanweisung       | 31 |
| it    | Istruzioni per l'uso     | 41 |
| nl    | Instructies voor gebruik | 52 |
| pt-EU | Instruções de Utilização | 62 |

## TABLE OF CONTENTS

|  |   |
|--|---|
| <b>REUSE WARNING</b> .....   | 3 |
| <b>DEVICE DESCRIPTION</b> .....  | 3 |
| Contents .....   | 3 |
| Figure 1. 1.5 CX Cryoablation Needle – Components .....                    | 3 |
| Figure 2. 1.5 CX Cryoablation Needle – Needle Features / Shaft Marks ..... | 4 |
| Figure 3. 2.1 CX Cryoablation Needle – Components .....                    | 4 |
| Figure 4. 2.1 CX Cryoablation Needle – Needle Features / Shaft Marks ..... | 5 |
| Table 1. CX Cryoablation Needle - Handle Colors .....                      | 5 |
| Materials .....  | 5 |
| Operating Principle .....  | 5 |
| Needle Performance - Cryoablation .....                                    | 5 |
| 2.1 CX L 90° Porcine Testing/Lung .....                                    | 6 |
| Figure 5. 1.5 CX Cryoablation Needle Isotherm Data .....                   | 6 |
| Figure 6. IcePearl Cryoablation Needle 2.1 CX Isotherm Data .....          | 6 |
| Figure 7. IceForce Cryoablation Needle 2.1 CX Isotherm Data .....          | 6 |
| Needle Performance - Active Thaw .....                                     | 6 |
| Needle Performance - Track Ablation (Cautery function) .....               | 6 |
| Table 2. 1.5 CX Cryoablation Needle – Track Ablation .....                 | 6 |
| Table 3. 2.1 CX Cryoablation Needle – Track Ablation .....                 | 7 |
| CX Needles Porcine Testing/Liver .....                                     | 7 |
| User Information .....   | 7 |
| <b>INTENDED USE</b> .....  | 7 |
| Patient Groups .....   | 7 |
| <b>INDICATIONS FOR USE</b> .....   | 7 |
| Clinical Benefit Statement .....   | 7 |
| <b>CONTRAINDICATIONS</b> .....   | 7 |
| <b>WARNINGS</b> .....  | 7 |
| <b>PRECAUTIONS</b> .....   | 8 |
| General .....  | 8 |
| Handling .....   | 8 |
| Procedural .....   | 8 |

|   |    |
|---|----|
| <b>ADVERSE EVENTS</b> .....                                     | 8  |
| <b>CONFORMANCE TO STANDARDS</b> .....                           | 9  |
| <b>HOW SUPPLIED</b> .....                                       | 9  |
| Device Details .....  | 9  |
| Handling and Storage .....                                      | 9  |
| <b>OPERATIONAL INSTRUCTIONS</b> .....                           | 9  |
| Needle Preparation .....  | 9  |
| Needle Use .....  | 9  |
| Needle Handling and Insertion .....                             | 9  |
| Table 4. 1.5 CX Cryoablation Needle – Withdrawal Distance ..... | 10 |
| Table 5. 2.1 CX Cryoablation Needle – Withdrawal Distance ..... | 10 |
| Needle Removal .....  | 10 |
| Disposal .....  | 11 |
| Post-Procedure .....  | 11 |
| <b>INFORMATION TO BRIEF THE PATIENT</b> .....                   | 11 |
| <b>WARRANTY</b> .....   | 11 |
| <b>SYMBOL DEFINITIONS</b> .....                                 | 11 |

# IceRod™

1.5 CX 90°

# IceSphere™

1.5 CX 90°

# IcePearl™

2.1 CX STRAIGHT  
2.1 CX 90°  
2.1 CX L 90°

# IceForce™

2.1 CX STRAIGHT  
2.1 CX 90°  
2.1 CX L 90°

## Cryoablation Needle

### Rx ONLY

**Caution:** Federal Law (USA) restricts this device to sale by or on the order of a physician.

### REUSE WARNING

The Boston Scientific IceRod Cryoablation Needle, IceSphere Cryoablation Needle, IcePearl Cryoablation Needle, and IceForce Cryoablation Needle (CX Needles) are disposable, single-use products. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness, or death of the patient.

### DEVICE DESCRIPTION


#### Contents

One (1) of the following CX Needles:

- IceRod Cryoablation Needle, 1.5 CX 90°
- IceSphere Cryoablation Needle, 1.5 CX 90°
- IcePearl Cryoablation Needle, 2.1 CX Straight
- IcePearl Cryoablation Needle, 2.1 CX 90°
- IcePearl Cryoablation Needle, 2.1 CX L 90°
- IceForce Cryoablation Needle, 2.1 CX Straight

- IceForce Cryoablation Needle, 2.1 CX 90°
- IceForce Cryoablation Needle, 2.1 CX L 90°


The Boston Scientific CX Needles (straight or angled 90°) have a sharp cutting tip, a shaft with distal coating, a color-coded handle, gas tubing, and a connector. All components are illustrated in Figure 1 and Figure 3. Needle features, shaft marks, and coating are shown in Figure 2 and Figure 4. Needle handle colors are listed in Table 1.

| # | Component                        | # | Component                   |
|---|----------------------------------|---|-----------------------------|
| 1 | Cutting Tip                      | 5 | Needle Handle (color coded) |
| 2 | Low-friction (non-stick) Coating | 6 | Gas Tubing                  |
| 3 | Active Zone Indicator            | 7 | Electrical Connector        |
| 4 | Needle Shaft (diameter 1.5 mm)   | 8 | Gas Connector               |


Figure 1. 1.5 CX Cryoablation Needle – Components

**NOTE:** The Active Zone Indicator is a marked band 20 mm from low-friction coated portion of the needle shaft. The location of the Active Zone Indicator guides needle withdrawal following Track Ablation (Cautery function).

| # | Features              | Description  |
|---|-----------------------|--|
| 1 | Heating Portion       | Located within the distal needle shaft below the coated area   |
| 2 | Low-friction Coating  | (A) Starts 4.5 mm from 1.5 CX tip; ends 50 mm from 1.5 CX tip<br>(B) Starts 4.5 mm from 1.5 CX tip; ends 35 mm from 1.5 CX tip |
| 3 | Thin Marks            | (A) Starts 55 mm from 1.5 CX tip, in 5 mm intervals<br>(B) Starts 45 mm from 1.5 CX tip, in 5 mm intervals                     |
| 4 | Thick Marks           | (A) Starts 60 mm from 1.5 CX tip, in 10 mm intervals<br>(B) Starts 40 mm from 1.5 CX tip, in 10 mm intervals                   |
| 5 | Active Zone Indicator | (A) Starts 70 mm from 1.5 CX tip; 10 mm band<br>(B) Starts 55 mm from 1.5 CX tip; 10 mm band                                   |
| 6 | Double-thick Mark     | Located 100 mm from 1.5 CX tip   |
| 7 | Triple-thick Mark     | Located 150 mm from 1.5 CX tip   |

**Figure 2. 1.5 CX Cryoablation Needle - Needle Features / Shaft Marks**


| # | Component                        | # | Component                   |
|---|----------------------------------|---|-----------------------------|
| 1 | Cutting Tip                      | 5 | Needle Handle (color coded) |
| 2 | Low-friction (non-stick) Coating | 6 | Gas Tubing                  |
| 3 | Active Zone Indicator            | 7 | Electrical Connector        |
| 4 | Needle Shaft (diameter 2.1 mm)   | 8 | Gas Connector               |

**Figure 3. 2.1 CX Cryoablation Needle - Components**

**NOTE:** The Active Zone Indicator is a marked band 20 mm from low-friction coated portion of the needle shaft. The location of the Active Zone Indicator guides needle withdrawal following Track Ablation (Cautery function).

BSC (MB IFU Template 8 x 8 G60ba1\_92310048M)\_IFU\_MB\_CX-Needle\_OUS\_514471619-01B

| # | Features              | Description  |
|---|-----------------------|--|
| 1 | Heating Portion       | Located within the distal needle shaft below the coated area   |
| 2 | Low-friction Coating  | (C) Starts 5.5 mm from 2.1 CX tip; ends 35 mm from 2.1 CX tip<br>(D) Starts 5.5 mm from 2.1 CX tip; ends 50 mm from 2.1 CX tip   |
| 3 | Single-thick Mark     | Located 50 mm from tip, (C) only   |
| 4 | Active Zone Indicator | (C) Starts 55 mm from 2.1 CX tip; 10 mm band<br>(D) Starts 70 mm from 2.1 CX tip; 10 mm band   |
| 5 | Thick Marks           | (C) Starts 40 mm from 2.1 CX tip, in 10 mm intervals<br>Starts 70 mm from 2.1 CX L tip, in 10 mm intervals<br>(D) Starts 60 mm from 2.1 CX tip, in 10 mm intervals<br>Starts 90 mm from 2.1 CX L tip, in 10 mm intervals |
| 6 | Thin Marks            | (C) Starts 45 mm from 2.1 CX tip, in 10 mm intervals<br>Starts 75 mm from 2.1 CX L tip, in 10 mm intervals<br>(D) Start 55 mm from 2.1 CX tip, in 10 mm intervals<br>Start 85 mm from 2.1 CX L tip, in 10 mm intervals   |
| 7 | Double-thick Mark     | Located 100 mm from tip for both 2.1 CX (C) and (D)<br>Located 200 mm from tip for both 2.1 CX L (C) and (D)   |
| 8 | Triple-thick Mark     | Located 150 mm from tip  |

Figure 4. 2.1 CX Cryoablation Needle – Needle Features / Shaft Marks

Table 1. CX Cryoablation Needle - Handle Colors

| Needle Brand Name                     | Configuration | Needle Shaft Length | Handle Color |
|---------------------------------------|---------------|---------------------|--------------|
| IceRod Cryoablation Needle 1.5 CX     | Angled 90°    | 175 mm              | Red          |
| IceSphere Cryoablation Needle 1.5 CX  | Angled 90°    | 175 mm              | Yellow       |
| IcePearl Cryoablation Needle 2.1 CX   | Straight      | 175 mm              | White        |
| IcePearl Cryoablation Needle 2.1 CX   | Angled 90°    | 175 mm              | White        |
| IcePearl Cryoablation Needle 2.1 CX L | Angled 90°    | 230 mm              | White        |

| Needle Brand Name                     | Configuration | Needle Shaft Length | Handle Color |
|---------------------------------------|---------------|---------------------|--------------|
| IceForce Cryoablation Needle 2.1 CX   | Straight      | 175 mm              | Gray         |
| IceForce Cryoablation Needle 2.1 CX   | Angled 90°    | 175 mm              | Gray         |
| IceForce Cryoablation Needle 2.1 CX L | Angled 90°    | 230 mm              | Gray         |

#### Materials

|                      |                            |  |
|----------------------|----------------------------|--|
| Materials            | Needle Shaft               | Stainless steel                        |
|                      | Low-friction Shaft Coating | Teflon type fluoropolymer              |
|                      | Needle Handle              | Brass (coated with heat shrink tubing) |
|                      | Gas Tubing                 | Polyurethane                           |
|                      | Connector                  | Polyoxymethylene                       |
| Sterilization Method | Ethylene Oxide (EO)        |  |

#### Operating Principle

The Boston Scientific ICEfx Cryoablation System or Visual-ICE Cryoablation System are intended for cryoablative tissue destruction using a minimally invasive procedure. The cryoablation system is controlled with a touch screen user interface that allows the user to control and monitor the procedure. The therapy delivered by the system is based on the Joule-Thomson effect displayed by compressed gases. The Joule-Thomson effect is a change in the temperature of a compressed gas as it flows through a narrow orifice and expands to a lower pressure. The cryoablation system uses high-pressure argon gas that circulates through closed-tip cryoablation needles to induce tissue freezing. Freezing is due to the Joule-Thomson effect augmented by a recuperative heat exchanger within the cryoablation needle.

Tissue thawing is achieved by controlling the heating element inside the CX Needles to provide active helium-free thawing (i-Thaw Function or FastThaw Function) and track ablation (Cautery function).

Cryoablative treatments are performed using a small-diameter cryoablation needle(s) inserted into the target lesion, usually under CT image guidance. When cryoablation needle(s) are placed into or near the target tissue and freezing is initiated, an iceball grows around the distal end of the needle shafts. In time, the iceball completely engulfs and destroys the target tissue. The shape of the iceball (ellipsoidal or spherical) depends on the type of needle used.

#### Needle Performance - Cryoablation

In-vivo iceball dimensions and resulting ablation zone are determined by the selected cryoablation needle (type and diameter), the number of needles placed, tissue and tumor characteristics, thermal heat sink from surrounding vasculature, and treatment duration. Monitoring iceball formation provides direct control throughout the procedure and is key to cryoablation success.

**NOTE:** Using ultrasound or CT visualization, monitor iceball formation throughout the cryoablation procedure.

Iceball shapes and dimensions are specific to the needle brand. The following laboratory models of iceball dimensions are provided to assist users in selecting the cryoablation needle(s) and needle placement to appropriately ablate the target

area. Typically, in-vivo dimensions are smaller than the dimensions generated in the laboratory. Laboratory testing was performed in room temperature gel; measurements were made after two (10 minute) Freeze cycles separated by a 5 minute passive Thaw cycle. Accuracy is  $\pm 3$  mm width,  $\pm 4$  mm length.

## 2.1 CX L 90° Porcine Testing/Lung

An in-vivo study was conducted in the porcine model to characterize the performance of the 2.1 CX L 90° cryoablation needles in lung tissue. The ablation zones produced by these needles were demarcated by histopathology and the dimensions of the ablated tissue were measured. The measurements were taken at the largest point of the ablated volume. The largest diameter of the ablated zone of lung tissue was measured to be 28 mm for the IcePearl Cryoablation Needle and 34 mm for the IceForce Cryoablation Needle. These results are consistent with laboratory measurements in iceballs produced by these needles in gel. The study concluded that the 2.1 CX L 90° IcePearl Cryoablation Needle and IceForce Cryoablation Needle are appropriate for use in cryoablation therapy in lung tissue.


Figure 5. 1.5 CX Cryoablation Needle Isotherm Data


Figure 6. IcePearl Cryoablation Needle 2.1 CX Isotherm Data


Figure 7. IceForce Cryoablation Needle 2.1 CX Isotherm Data

## Needle Performance - Active Thaw

CX Needles support active thawing using helium, or active thawing without helium using an i-Thaw Function internal heater. In addition, the i-Thaw Function heater can be activated to conduct the optional FastThaw Function. The FastThaw Function generates a temperature that is higher than that required for i-Thaw Function, resulting in a faster thaw time.

Laboratory testing was performed in 37 °C gel to evaluate the thawing performance of CX Needles and provide users with comparative thawing characteristics. A single (10 minute) Freeze cycle was conducted to create an iceball on each test needle. Activation of CX Needles FastThaw Function for a 5 minute thaw resulted in faster iceball thawing than CX Needles operated in i-Thaw Function for 5 minutes by approximately 15% faster for 1.5 CX Needles and 18% faster for 2.1 CX Needles.

**NOTE:** In-vivo thawing performance will be determined by numerous factors, including brand (type) and number of cryoablation needles used, needle location, tissue characteristics, and duration of thaw activation.

## Needle Performance - Track Ablation (Cautery function)

When CX Needles are connected to a Boston Scientific cryoablation system, activation of track ablation is an option. During track ablation, the needle shaft surface heats to a temperature of 85 °C to 100 °C, transferring heat along the distal shaft at a 30 mm segment of the IceRod Cryoablation Needle, a 14 mm segment of the IceSphere Cryoablation Needle, a 13 mm segment of the IcePearl Cryoablation Needle, and a 29 mm segment of the IceForce Cryoablation Needle.

Table 2. 1.5 CX Cryoablation Needle - Track Ablation

| Product Name                         | Track Ablation (30 seconds) |
|--------------------------------------|-----------------------------|
| IceRod Cryoablation Needle 1.5 CX    | ~2.3 mm $\pm$ 0.7 mm*       |
| IceSphere Cryoablation Needle 1.5 CX | ~1.7 mm $\pm$ 0.1 mm**      |
| *In-vivo testing                     |                             |
| **Ex-vivo testing                    |                             |

**Table 3. 2.1 CX Cryoablation Needle – Track Ablation**

| Product Name   | Track Ablation (30 seconds) | Track Ablation (3 minutes) |
|--|-----------------------------|----------------------------|
| IcePearl Cryoablation Needle 2.1 CX<br>IcePearl Cryoablation Needle 2.1 CX L | ~2.2 mm ± 0.6 mm            | ~2.1 mm ± 0.6 mm           |
| IceForce Cryoablation Needle 2.1 CX<br>IceForce Cryoablation Needle 2.1 CX L | ~2.5 mm ± 0.4 mm            | ~2.5 mm ± 0.4 mm           |

### CX Needles Porcine Testing/Liver

An in-vivo study was conducted in porcine liver tissue to evaluate the track ablation (Cautery function) performance of the IceRod Cryoablation Needle, IcePearl Cryoablation Needle, and IceForce Cryoablation Needle and to measure the depth of tissue necrosis. Cryoablation, followed by track ablation, was performed in the porcine tissue. Histological assessment demonstrated a fairly uniform zone of tissue necrosis surrounding the path of track ablation. The average width of the zone of histologically-confirmed tissue damage, measured from the edge of the needle track to the normal tissue interface. Refer to Table 2 and Table 3.

An ex-vivo study was conducted in freshly excised porcine liver tissue to evaluate track ablation performance of the IceSphere Cryoablation Needle and to measure the depth of tissue necrosis. Track ablation was performed in the porcine tissue. Triphenyltetrazolium chloride (TTC) staining demonstrated a contiguous zone of tissue necrosis surrounding the path of track ablation. The average width of the zone of TTC negative confirmed tissue necrosis, measured from the edge of the needle track to the normal tissue interface is described in the table. The exact size of the track ablation zone may vary in clinical settings due to perfusion or vasculature not present in ex-vivo tissue. Refer to Table 2.

### User Information

The intended users are medical professionals who have a thorough understanding of the technical principles, clinical applications, and risk associated with cryoablation procedures. Refer to the ICEfx Cryoablation System or Visual-ICE Cryoablation System User Manuals for optional education.

### INTENDED USE

The Boston Scientific disposable IceRod Cryoablation Needle, IceSphere Cryoablation Needle, IcePearl Cryoablation Needle, and IceForce Cryoablation Needle (CX Needles) are meant to be connected to a Boston Scientific ICEfx Cryoablation System or Visual-ICE Cryoablation System when performing cryoablative tissue destruction through application of extremely cold temperatures. The needles are intended to convert high-pressure gas to either a very cold Freezing application or to a warm Thawing application.

### Patient Groups

The intended population includes patients intended for cryoablative destruction of tissue during surgical procedures.

### INDICATIONS FOR USE

The Boston Scientific disposable IceRod Cryoablation Needle, IceSphere Cryoablation Needle, IcePearl Cryoablation Needle, and IceForce Cryoablation Needle (CX Needles), when used with a Boston Scientific ICEfx Cryoablation System or Visual-ICE Cryoablation System, are designed to destroy tissue by the application of extremely cold temperatures. A full list of specific indications can be found in the respective Boston Scientific cryoablation system user manual.

### Clinical Benefit Statement

The ICEfx Cryoablation System and Visual-ICE Cryoablation System (used in conjunction with Boston Scientific accessory products) are intended to destroy tissue (including prostate and kidney tissue, liver metastases, tumors, and skin lesions) by the application of extremely cold temperatures during minimally invasive procedures.

The clinical benefit is measured by overall clinical outcomes, including acceptable safety and performance outcomes specific to the target anatomy.

### CONTRAINDICATIONS

There are no known contraindications specific to use of CX Needles.

### WARNINGS

- A thorough understanding of the technical principles, clinical applications, and risk associated with cryoablation is necessary before using this product. Refer to the ICEfx Cryoablation System or Visual-ICE Cryoablation System User Manuals for optional education.
- Do not use this device for any purpose other than the stated intended use.
- Observe the expiration date of this product. Do not use past the listed expiration date.
- Select CX Needles appropriate for the application and tumor size. The iceball shape and size for CX Needles are described in the **Device Description** Section.
- Do not use CX Needles within a Magnetic Resonance Imaging (MRI) suite or environment.
- Use special care to ensure that a cryoablation needle does not come into contact with an implanted device.
- Cryoablation needles have sharp tips. Use care to ensure safe handling of a cryoablation needle to eliminate the risk of injury or possible exposure to blood-borne pathogens.
- The sterile field and sterility of the cryoablation needle should be maintained at all times. Do not contaminate the distal end of the sterile cryoablation needle.
- CX Needles are designed and indicated for Freezing and Thawing applications. These needles are not designed, indicated or tested for thermal protection. Serious injury to patient tissue may result if used for thermal protection.
- Inspect the packaging for damage. Do not use the needle if packaging appears opened or damaged; in the event of such occurrence, contact Boston Scientific Technical Assistance Center to arrange return of the complete package with the product.
- Do not use the needle if it is bent or damaged while attempting to unpack or use it. Never use a defective needle for a cryoablation procedure. A defective cryoablation needle that has a gas leak can cause a gas embolism or pneumatoxis in the patient.
- Avoid bending the needle shaft. Do not grasp needles with auxiliary instruments as this may cause damage to the needle shaft.
- During use, avoid damage to the needle from other surgical instruments.
- Do not kink, pinch, cut or pull excessively on the needle tubing. Damage to needle handle or tubing may cause the needle to become unusable.
- Ensure appropriate stability of tubing to avoid inadvertent movement of tubing or needle shifting during a procedure.
- Each needle must be locked into a needle channel before initiating a cryoablation procedure to avoid the risk of forceful ejection of the needles while under gas pressure.
- Prior to starting a cryoablation procedure, set up the cryoablation system (reference the SYSTEM SETUP Section within the cryoablation system User Manual) and then perform the Needle Integrity and Functionality Test. This test must be successfully completed in order to begin the procedure.

- Do not use the needle if there is no ice formation during the Freeze phase while performing the Needle Integrity and Functionality Test. Refer to cryoablation system user manual for troubleshooting. If the issue does not resolve, obtain a new needle and repeat the testing procedure.
- Do not use the needle if bubbles are seen escaping from the needle during Needle Integrity and Functionality Test.
- Ensure adequate measures are taken to protect organs and structures adjacent to the targeted tissue.
- Continuously monitor needle insertion, needle positioning, iceball formation and removal using image guidance (such as direct visualization, ultrasound, or Computed Tomography (CT)) to ensure adequate tissue coverage and to avoid damage to adjacent structures.
- Use Boston Scientific's Multi-Point 1.5 Thermal Sensor Device (MTS) to monitor the freeze / thaw temperatures for the intended treatment protocol or to monitor tissue temperature near critical structures.
- In the rare event that a needle breaks while inserted in the tissue, act immediately to remove needle parts from the patient's body and report such event to Boston Scientific.
- If a needle unintentionally strikes bone, do not start or continue the Freezing process as needle integrity may be compromised. Replace the needle prior to continuing the procedure.
- Needle handles and the needle shaft may frost during freezing. Avoid prolonged contact with frosted portions of the needle handle to avoid unintended thermal tissue damage to the patient or clinician. The patient's skin surface should be protected by warm saline irrigation or other means as determined by the physician.
- Needle tubing may become extremely cold when conducting freeze cycles during a cryoablation procedure. It is important that a patient's skin is protected from direct contact with needle tubing to avoid the potential for thermal injury to the patient. Ensure an appropriate insulating barrier is placed as needed (such as towels) or other method is employed to prevent needle tubing from touching a patient's skin.
- Ensure the Active Zone Indicator is not positioned outside the patient's skin when track ablation (Cautery function) is activated.
- The needle handle may become warm during active thawing. Prolonged contact with warm portions of the needle handle could cause unintended thermal tissue damage/burn to the patient or clinician.
- Active thawing produces heat along the distal needle shaft. Use care to avoid thermal injury/burn to nontargeted tissues.
- Ensure adequate thawing or cooling before attempting to remove needles from the patient
- Discontinue all needle operation prior to needle removal to minimize risk of thermal injury and/or tissue injury.
- When conducting FastThaw Function or track ablation (Cautery function), be alert for the Active Zone Indicator as when track ablation is activated and when the needle is withdrawn to prevent unintended tissue damage from the hot needle.
- Remove needles from the patient prior to disconnecting needles from the Boston Scientific cryoablation system.
- Dispose of needle(s) in accordance with the **Disposal** Section.
- No data regarding cryoablation in combination with other therapies is available from Boston Scientific.

## PRECAUTIONS

### General

- Boston Scientific CX Needles operate only with a Boston Scientific ICEFx Cryoablation System or Visual-ICE Cryoablation System.
- Boston Scientific recommends that only needles of the same type be placed together in a single channel. Using needles of differing types in a channel may affect the accuracy of the Gas Indicator.
- Confirm availability of sufficient gas to conduct the planned procedure: the number of needles, needle operations activated, gas cylinder size, pressure and gas flow affect the required gas volume.

### Handling

- Use of multiple needles is recommended to fully cover a target site and provide a suitable margin.
- Multiple needles placed in an adjacent configuration will typically create a large, coalesced iceball. Iceball formation must be monitored using image guidance to optimize a successful ablation procedure.
- Availability of a back-up needle is recommended should a replacement or additional needle be required during a procedure.
- If the Boston Scientific cryoablation system contains pressurized helium, track ablation (Cautery function), i-Thaw Function and FastThaw Function cannot be activated.

### Procedural

- In a safe and controlled fashion, twist the protective sleeve in a rotational motion while simultaneously pulling it away from the device handle. Take extra care when removing the protective sleeve from the device to avoid contact with the distal end of the needle.

### ADVERSE EVENTS

The potential adverse events associated with the device and/or cryoablation procedure include, but are not limited to:

- Allergic reaction (contrast, device, other)
- Angina
- Arrhythmia
- Atelectasis
- Bladder spasms
- Bleeding/hemorrhage
- Burn/frostbite
- Cerebrovascular accident (CVA)/stroke
- Cryoshock phenomenon (e.g. multi organ failure, severe coagulopathy, disseminated intravascular coagulation (DIC))
- Death
- Distension
- Edema/swelling
- Ejaculatory dysfunction
- Embolism (air, device, thrombus)
- Erectile dysfunction
- Fever

- Fistula
- Fracture
- Gastrointestinal symptoms (e.g. nausea, vomiting, diarrhea, constipation)
- Healing, impaired
- Hematoma
- Hematuria
- Hemothorax
- Hepatic dysfunction/failure
- Hernia
- Hypertension
- Hypotension
- Hypothermia
- Ileus
- Impotence
- Infection/abscess/sepsis
- Inflammation
- Muscle spasm
- Myocardial infarction
- Necrosis
- Need for additional intervention or surgery
- Nerve injury
- Neuropathy
- Obstruction
- Pain/discomfort
- Perforation (including organ and adjacent structures)
- Pericardial effusion
- Perirenal fluid collection
- Pleural effusion
- Pneumatisis (air or gas in an abnormal quantity and/or place in the body)
- Pneumothorax
- Post ablation syndrome (e.g. fever, pain, nausea, vomiting, malaise, myalgia)
- Renal insufficiency/failure
- Renal parenchymal or capsule fracture
- Respiratory distress/insufficiency/failure
- Scrotal edema
- Stenosis/stricture
- Subcutaneous emphysema
- Thrombosis/thrombus
- Tissue damage
- Transient ischemic attack (TIA)
- Tumor cell seeding
- Urethral sloughing
- Urinary frequency/urgency
- Urinary incontinence

- Urinary retention
- Urinary tract infection
- Vasovagal response
- Vessel trauma (e.g. dissection, injury, perforation, pseudoaneurysm, rupture, or other)
- Wound infection

#### CONFORMANCE TO STANDARDS

EMF output conforms with IEC 60601-1-2 "General requirements for safety – Collateral standard: Electromagnetic compatibility – Requirements and tests."

#### HOW SUPPLIED

##### Device Details

Boston Scientific disposable CX Needles are individually packaged in a sealed film-Tyvek Pouch. Each pouch is labeled Sterile, for SINGLE USE only.

---

**NOTE:** Boston Scientific cryoablation systems are supplied separately.

---

Needle connectors are protected by a rubber cap. This cap should be removed before connecting a needle to a cryoablation system.

Do not use if package is damaged or unintentionally opened before use.

Do not use if labeling is incomplete or illegible.

##### Handling and Storage

This product has no special handling or storage requirements.

#### OPERATIONAL INSTRUCTIONS

##### Needle Preparation

1. Using aseptic technique, carefully remove the cryoablation needle from the package and place in a sterile work area. In a safe and controlled fashion, remove the protective sleeve by twisting the protective sleeve in a rotational motion while simultaneously pulling it away from the device handle.
2. To conduct the Needle Integrity and Functionality Test, prepare a large basin containing sterile water or saline that is sufficiently able to accommodate the full needle shaft length.
3. Remove the connector cap, then connect the needle to the selected port.
4. Place the needles, individually or in groups, in the basin such that the full length of the needle is submerged in the sterile water or saline.

The needle is now ready for the Needle Integrity and Functionality Test.

---

**NOTE:** For detailed instructions on connecting needles to the system's Needle Connection Panel and performing the Needle Integrity and Functionality Test, refer to the respective Boston Scientific cryoablation system User Manual.

---

##### Needle Use

##### Needle Handling and Insertion

- Correct insertion of cryoablation needles into the target tissue is the responsibility of the physician.

---

**NOTE:** Although the needle has a sharp tip, a small skin incision may be made at the point of needle insertion.

---

- Always use two hands and support the needle mid-shaft with two fingers to eliminate the risk of bending. Do not insert the needle into tissue while holding the handle with one hand only.
- Insertion depth may be estimated using the needle marks on the shaft. Use image guidance to guide needle insertion and placement.
- Use image guidance as necessary to verify that the cryoablation needle is placed at the desired location prior to activating the needle.

**CAUTION:**

- After needle insertion, route the needle tubing in a manner that ensures it is not in contact with the patient's skin.
- Place an appropriate insulating barrier (such as towels) or employ other methods to prevent needle tubing from touching a patient's skin.
- The needle handle becomes cold during cryoablation. If the handle is in contact with skin, the skin surface should be protected by warm saline irrigation or other means as determined by the physician.

**NOTE:** See the respective Boston Scientific cryoablation system User Manual for instructions on system controls available to manage each Freeze cycle.

**CAUTION:**

- Continuously monitor iceball formation using direct visualization or image guidance such as Ultrasound or Computed Tomography (CT) to ensure adequate tissue coverage and to avoid damage to adjacent structures.
- A portion of the handle will remain frost-free to facilitate intraprocedural handling.
- The needle shaft becomes cold during cryoablation. Protect the patient's skin by using warm saline irrigation or other means as determined by the physician.
- Pay attention to needle handle position. Prolonged contact with frosted portions of the needle handle could cause unintended thermal tissue damage to the patient or clinician.
- Avoid allowing needle tubing to be in direct contact with a patient's skin. Direct contact with needle tubing could cause unintended thermal damage to the patient's skin. Place an appropriate insulating barrier (such as towels) or employ other methods to prevent needle tubing from touching a patient's skin.

**NOTE:** See the respective Boston Scientific cryoablation system User Manual for instructions on system controls available for i-Thaw Function and FastThaw Function options.

**CAUTION:**

- The distal portion of the needle handle may become warm during thawing.
- Pay attention to needle handle position. Prolonged contact with warm portions of the needle handle could cause unintended thermal tissue damage to the patient or clinician.
- Thaw thoroughly and stop all needle operation prior to removing needles to minimize the risk of tissue injury.
- If needle sticking is experienced, use a slight, gentle twist of the needle followed by slow withdrawal.

Notes for Performing Track Ablation (Cautery function)

- Track ablation can be activated at any time during a cryoablation procedure.
- Hold the needle in a stationary position during track ablation.
- Track ablation may be repeated as required. Prior to each reactivation, slowly withdraw per the table below, then activate track ablation.

**Table 4. 1.5 CX Cryoablation Needle – Withdrawal Distance**

| Product Name                         | Withdrawal Distance |
|--------------------------------------|---------------------|
| IceRod Cryoablation Needle 1.5 CX    | 30 mm               |
| IceSphere Cryoablation Needle 1.5 CX | 10 mm               |

**Table 5. 2.1 CX Cryoablation Needle – Withdrawal Distance**

| Product Name                          | Withdrawal Distance |
|---------------------------------------|---------------------|
| IcePearl Cryoablation Needle 2.1 CX   | 10 mm               |
| IcePearl Cryoablation Needle 2.1 CX L | 10 mm               |
| IceForce Cryoablation Needle 2.1 CX   | 25 mm               |
| IceForce Cryoablation Needle 2.1 CX L | 25 mm               |

- For detailed instructions on using and controlling the track ablation option, refer to the respective Boston Scientific cryoablation system User Manual.

**WARNING:**

- Ensure the Active Zone Indicator is not positioned outside the patient's skin when track ablation (Cautery function) is activated.
- Active thawing produces heat along the distal needle shaft. Use care to avoid thermal injury/burn to nontargeted tissues.
- The needle handle may become warm during active thawing. Prolonged contact with warm portions of the needle handle could cause unintended thermal tissue damage/burn to the patient or clinician.

**Needle Removal**

If track ablation (Cautery function) is performed:

- Do not remove the needle until needle cooling is complete.
- If needle sticking is experienced, use a slight, gentle twist of the needle followed by slow withdrawal.

If track ablation (Cautery function) is not performed:

- Thaw thoroughly and stop all needle operation prior to removing needles to minimize the risk of tissue injury.

**NOTE:** Boston Scientific's CX Needles are specially designed with a three-facet, trocar-like tip to minimize bleeding. However, some bleeding may occur. In the event of bleeding, apply treatment in accordance with good clinical practice and the hospital's treatment protocol. For example, following needle removal, hold compression until hemostasis is achieved; if necessary, place an appropriate dressing on the needle insertion site.

**CAUTION:** Removing the needle while it is still hot presents a risk of injury to adjacent tissue and/or organs.

## Disposal

This device contains a sharps hazard. Take precautions to ensure that sharps are handled properly. Dispose of all sharps directly into a sharps disposal container labeled with a biological hazard symbol. Sharps waste should be safely disposed of using available sharps waste channels in accordance with hospital, administrative, and/or local government policy.

## Post-Procedure

Any serious incident that occurs in relation to this device should be reported to the manufacturer and relevant local regulatory authority.

For customers in Australia, report any serious incident that occurs in relation to this device to Boston Scientific and to the Therapeutic Goods Administration (<https://www.tga.gov.au>).

## INFORMATION TO BRIEF THE PATIENT

The physician should consider the following points while counseling the patient on the use of the cryoablation system in association with the interventional procedure:

- Discuss the risks and benefits, including review of potential adverse events listed in this IFU and the cryoablation system User Manuals, of other interventional treatments likely to be employed.
- Discuss post-procedure instructions, including any lifestyle changes, medications, and home care and rehabilitation guidelines.

## WARRANTY

For device warranty information, visit ([www.bostonscientific.com/warranty](http://www.bostonscientific.com/warranty)).

EU Importer: Boston Scientific International B.V., Vestastraat 6, 6468 EX Kerkrade, The Netherlands

FastThaw, i-Thaw, IceForce, ICEfx, IcePearl, IceRod, IceSphere, Multi-Point 1.5 Thermal Sensor, and Visual-ICE are trademarks of Boston Scientific Corporation or its affiliates.

All other trademarks are the property of their respective owners.

## SYMBOL DEFINITIONS

Commonly used medical device symbols that appear on the labeling are defined at [www.bostonscientific.com/SymbolsGlossary](http://www.bostonscientific.com/SymbolsGlossary).

Additional symbols are defined at the end of this document.

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## INSTRUCTIONS FOR USE

The Instructions for Use (IFU) for this product are supplied in electronic form over the Internet.

Visit [www.IFU-BSCI.com](http://www.IFU-BSCI.com) to access the IFU in Adobe® Portable Document Format (PDF), and be sure to have the product label available for reference.

If you have difficulty accessing the IFU online, or would prefer to receive a paper copy, please contact Boston Scientific Customer Service or your local country contact. A copy will be sent to you at no charge and should arrive within seven days.

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All rights reserved.
Boston Scientific (Australia) Pty Ltd
© 2022 Boston Scientific Corporation or its affiliates.