ICEfx™ Cryoablation System

Cryoablation System

en	User's Manual	2
es	Manual del usuario	83
fr	Manuel d'utilisation	171
de	Bedienungsanleitung	257
it	Manuale dell'utente	343
nl	Gebruikershandleiding	428
pt-eu	Manual do Utilizador	514

TABLE OF CONTENTS

REUSE WARNING	6
DEVICE DESCRIPTION	6
System Description	6
Contents	6
ICEfx Cryoablation Console	7
Figure 1. Console Front View	7
Figure 2. Console Rear View	8
Touch Screen Monitor	8
ICEfx Cryoablation Cart	10
Figure 3. ICEfx Cryoablation Cart	10
Figure 4. ICEfx Cryoablation Cart - Wheel Lock Levers	10
Operating Principle	11
Materials	11
Non-pyrogenic	11
User Information	11
INTENDED USE	11
Patient Groups	11
INDICATIONS FOR USE	12
Clinical Benefit Statement	12
CONTRAINDICATIONS	12
WARNINGS	12
PRECAUTIONS	15
ADVERSE EVENTS	17
CONFORMANCE TO STANDARDS	19
Electrical Specifications	19
Electromagnetic Compatibility and Immunity (EMC and EMI)	19
Table 1. Cable Lengths	19
Table 2. Electromagnetic Emissions	20
Table 3. Electromagnetic Immunity	21
Table 4. Electromagnetic Immunity for Systems that are Not Life-supporting	22
Table 5. Recommended separation distances between portable and mobile RF communications equipment and the ICEfx Cryoablation System	23
HOW SUPPLIED	24
Device Details	24
Handling and Storage	24
Removing the Console from the Cart	24

OPERATIONAL INSTRUCTIONS	25
Additional Required Items	25
Installation, Calibration, and Servicing.....	25
Preparation	26
System Setup	26
Table 6. System Setup Flow	26
Figure 5. Position ICEfx Cryoablation Console on Cart.....	27
Figure 6. Tighten Retaining Screw on Cart	27
Figure 7. Argon Shutoff Valve - ON Position	28
Connecting the Argon Gas.....	30
Figure 9. Attach Safety Cable to Console.....	30
Figure 10. Connect Argon Gas Supply Line.....	30
Figure 11. Gas Cylinder Setup.....	31
Figure 12. Open Argon Gas Cylinder.....	31
Figure 13. EZ-Connect2 Dual Cylinder Adapter	32
Screen 1: Gas Indicator.....	32
Table 7. Working Argon Gas Pressure	32
Navigating the User Interface.....	33
Title Bar.....	33
Screen 2: Title Bar	33
Table 8. Title Bar.....	33
Options Menu	34
Screen 3: Options Menu	34
Table 9. Options Menu.....	34
Screen 4: Start Case Screen.....	35
Table 10. Start Case Screen Controls	35
Screen 5: Select Tissue Type	35
Table 11. Select Tissue Type	35
Screen 6: Edit Case Information Screen.....	36
Table 12. Edit Case Information Screen Controls.....	36
Screen 7: Initial Procedure Screen.....	37
Table 13. Common Procedure Controls	37
Screen 8: Timer During FastThaw Function	38
Table 14. Timer Information.....	38
MTS	39
Screen 9: Procedure Screen – MTS	39
Table 15. MTS Controls	39

Screen 10: Procedure Screen – Change MTS Identifier Name	40
Table 16. Change MTS Name.....	40
Screen 11: Procedure Screen during Testing	41
Table 17. Test Controls.....	41
Initiating a Freeze Cycle	42
Screen 12: Procedure Screen during Freezing.....	42
Table 18. Freeze Controls.....	42
Screen 13: Expanded Channel History.....	43
Initiating a Thaw Cycle (CX needles only)	43
Screen 14: Procedure Screen during Thawing (FastThaw Function).....	44
Table 19. Thaw Controls.....	44
Screen 15: Cautery Function Screen.....	45
Table 20. Cautery Function Controls.....	46
Advanced Channel Controls.....	46
Screen 16: Needle Type List.....	47
Program Cycles Control.....	47
Screen 17: Select Programmed Cycle.....	47
Screen 18: Program Cycles Dialog.....	48
Screen 19: Programming Cycles Screen.....	48
Screen 20: Save Cycle Sequence Screen.....	49
Reports.....	49
Screen 21: All Reports Screen.....	49
Screen 22: Procedure Report.....	50
Configuring Settings.....	51
Table 21. Settings Screen – Configurable Settings.....	51
Table 22. Settings Screen – Control Buttons.....	52
Administrative Options.....	52
Procedure.....	53
Performing The Cryoablation Procedure.....	53
Table 23. Cryoablation Procedure Flow.....	53
Pre-Procedure Needle/MTS Testing.....	54
Screen 23: Procedure Screen.....	54
Figure 14. Locking Needle into Channel.....	54
Initiating a Cryoablation Procedure	57
Screen 24. Argon Remaining	57
Adding a CX Needle during a Cryoablation Procedure	59
Changing Argon Cylinders during a Procedure.....	59

SYSTEM SHUTDOWN	61
Shutting Down the System	61
Figure 15. Vent Valve – Fully Open	61
Figure 16. (same as Figure 8) Vent Valve - Fully Closed	61
Post-Procedure	62
Cleaning the ICEfx Cryoablation System	62
Disposal	62
TROUBLESHOOTING	63
Software Recovery	63
Replacing Fuses	65
Displayed Messages	71
SYSTEM SPECIFICATIONS	81
Mechanical Specifications	81
Argon Gas Supply	81
Gas Cylinder Specifications	81
Accuracy of Displayed Values	81
Essential Performance	81
Boston Scientific Technical Assistance Center	82
PATIENT COUNSELING INFORMATION	82
WARRANTY	82
SYMBOL DEFINITIONS	82

⚠ ONLY

Caution: Federal Law (USA) restricts this device to sale by or on the order of a physician.

REUSE WARNING

The single-use only devices used with the ICEfx Cryoablation System are sterile. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

DEVICE DESCRIPTION

System Description

The ICEfx Cryoablation System consists of an ICEfx Cryoablation Console and an optional ICEfx Cryoablation Cart.

Contents

One (1) ICEfx Cryoablation Console

One (1) Optional ICEfx Cryoablation Cart: The cart provides storage for up to four cryoablation needle boxes and provides a hook to hold an accessory bag.

One (1) ICEfx Cryoablation System User Manual: The User Manual may be physical copy or available online at www.IFU-BSCI.com. The User Manual Describes the system and provides instructions for system operation and maintenance.

One (1) ICEfx Cryoablation System Quick Reference Guide: The Quick Reference Guide provides a visual overview of system setup and operation.

One (1) Power Cord

One (1) Console Cover: The Console Cover is used to protect the ICEfx Cryoablation Console during storage.

One (1) USB flash drive (4 GB) within an attached pouch: The USB flash drive is used to transfer procedure reports to a customer's computer for saving or printing.

One (1) Wrench: The Wrench is used when connecting or disconnecting the gas supply line.

One (1) Accessory Bag: The Accessory Bag is used for storing the system accessories and can be attached to the hook on the side of the ICEfx Cryoablation Cart.

One (1) EZ-Connect2 Dual Cylinder Adapter: The Dual Cylinder Adapter is a component used to connect two argon cylinders in tandem to the ICEfx Cryoablation System. The dual cylinder adapter consists of a four-way adapter assembly with argon pressure gauge and auxiliary gas supply connection attached to a long gas supply line and a short auxiliary gas supply line with gas cylinder connection. Reference the **Connecting the Argon Gas** Section for instructions on using the EZ-Connect2 Dual Cylinder Adapter.

- To accommodate procedure room variations, the gas supply line is available in alternative lengths to connect the argon cylinder to the Console. Reference Table 1.

ICEfx Cryoablation Console

The ICEfx Cryoablation Console is a portable helium-free desktop cryoablation unit. The front of the Console contains four needle connection channels (each supporting two needle ports), two Multi-Point 1.5 Thermal Sensor Device (MTS) ports, and a recessed 10.1 inch touch screen monitor.

Figure 1. Console Front View

- | | |
|-------------------------|----------------------------|
| 1. Touch Screen Monitor | 4. Locking Bars |
| 2. Visual Indicator | 5. Needle Connection Ports |
| 3. MTS Connection Ports | 6. Handles |

The rear of the Console (Figure 2) contains the power switch, one USB Port, one Ethernet Port (inactive), a software reset button, an argon gas inlet, a Vent Valve, and an Argon Shutoff Valve.


Figure 2. Console Rear View

- | | |
|-----------------------------|---------------------------------------|
| 1. Power Switch | 7. Vent Valve |
| 2. USB Port | 8. Argon Shutoff Valve |
| 3. Ethernet Port (inactive) | 9. Argon Gas Inlet (quick-connect) |
| 4. Software Reset Button | 10. Humidity Cover |
| 5. Service Door | 11. Gas Safety Clip |
| 6. Cooling Fan | 12. Grounding Post (select countries) |

Touch Screen Monitor

The cryoablation procedure is controlled by using the touch screen monitor. The monitor can be tilted and rotated to provide the optimal viewing and operating angle for the user. The touch screen monitor includes a virtual on-screen English QWERTY keyboard for entering procedure related data. The monitor folds flat into the monitor storage basin on the top of the unit for system storage.

Needle Connection Panel

The Needle Connection Panel contains two MTS Connection Ports and four numbered needle channels; each channel contains two ports to support connection of up to two cryoablation needles.

MTS Connection Ports

Two **MTS Connection Ports** are located above the Needle Connection Channels to support connection of the Boston Scientific MTS used to monitor temperature within and near the cryoablation target site.

Needle Channels

The Needle Connection Panel contains four numbered needle channels; each channel contains two ports to support connection of up to two cryoablation needles. Each channel operates independently of all other channels in either freezing or thawing mode.

Locking Bars

The locking bar on each needle channel locks the needles into the ports to secure them during the procedure.

Handles

Two handles, one on each side are available for lifting the Console. The handles are recessed cutouts in the underside of the base.

CAUTION: If the Console is attached to the ICEfx Cryoablation Cart, use the cart handles to lift the system.

Power Switch

The **Power Switch** powers ON the Console in preparation for a cryoablation procedure.

Communication Port

A **USB 2.0** communication port is located on the rear of the Console to allow the user to save reports to a Boston Scientific USB flash drive for download to another computer or for printing.

Software Reset

The **Software Reset** button is used to boot the ICEfx Cryoablation Console into a recovery mode if the software becomes corrupted (reference the **Software Recovery** Section).

Argon Shutoff Valve

The **Argon Shutoff Valve** is used to turn the gas supply to the Console ON or OFF. It is recommended to leave the valve in the GAS ON position. Only turn OFF argon gas in an emergency.

Vent Valve

The **Vent Valve** is used to vent the high-pressure argon gas from the ICEfx Cryoablation Console.

Argon Gas Inlet

The gas supply line connects the argon gas supply from the argon gas cylinder to the **Argon Gas Inlet** via a quick-connect connector.

ICEfx Cryoablation Cart

The optional ICEfx Cryoablation Cart provides a secure platform for the Console. The ICEfx Cryoablation Cart has four locking wheels to ensure the cart does not move during a cryoablation procedure. An opening inside the cart provides storage for up to four cryoablation needle boxes.

The cart has handles for easy movement and a hook for attaching an accessory bag.


Figure 3. ICEfx Cryoablation Cart

1. Storage Area
2. Locking Wheel
3. Handle
4. Accessory Bag Hook

Wheel Lock Levers

Use the Wheel Lock Levers to lock the ICEfx Cryoablation Cart in place during a procedure (Figure 4). The wheel lock levers can be locked from either side. Pressing the wheel lock lever DOWN locks the wheel and prevents movement during the procedure. Lifting the wheel lock lever UP unlocks the wheel.


Figure 4. ICEfx Cryoablation Cart - Wheel Lock Levers

Operating Principle

The ICEfx Cryoablation System is a portable system intended for cryoablative tissue destruction using a minimally invasive procedure. The system is computer-controlled with a touch screen user interface that allows the user to control and monitor the procedure. The therapy delivered by the system is based on the Joule-Thomson effect displayed by compressed gases. The Joule-Thomson effect is a change in the temperature of a compressed gas as it flows through a narrow orifice and expands to a lower pressure. Certain gases, such as argon, decrease in temperature due to the Joule-Thomson effect, while other gases, such as helium, increase in temperature.

The ICEfx Cryoablation System uses high-pressure argon gas that circulates through closed-tip cryoablation needles to induce tissue freezing. Freezing is due to the Joule-Thomson effect augmented by a recuperative heat exchanger within the cryoablation needle. Active tissue thawing is achieved by controlling the heating element inside cryoablation CX needles to provide active helium-free thawing (i-Thaw Function or FastThaw Function) and track ablation (Cautery function).

Cryoablative treatments are performed using a small-diameter cryoablation needle(s) inserted into the target lesion, usually under CT image guidance. When cryoablation needle(s) are placed into or near the target tissue and freezing is initiated, an iceball grows around the distal end of the needle shafts. In time, the iceball completely engulfs and destroys the target tissue.

The cryoablation needle design and the temperature of the argon gas as well as the freeze duration determine the size and shape of the iceball. Tissue ablation is achieved by repeated freeze and thaw cycles with both freezing and thawing contributing to cell death. Each freeze-thaw cycle consists of a freeze cycle followed by a thaw cycle. Generally, multiple freeze-thaw cycles are used to achieve complete destruction of the target tissue.

An important benefit of cryoablation is that imaging procedures, such as ultrasound and CT, are able to display the location and size of the iceball. This benefit of cryoablation is used for proper control of the therapy. During use, the procedure must be monitored using image guidance to ensure adequate tissue coverage and to avoid damage to adjacent structures.

In addition to image guidance, Boston Scientific provides temperature sensors to aid in monitoring tissue temperature near the target site and adjacent critical structures. These temperature sensors can provide quantitative data to supplement the qualitative information provided by the imaging modality. Needle temperature display for CX-type needles provides a visual means to monitor needle performance.

Materials

See Boston Scientific cryoablation needle and accessory product IFUs for specific information on materials.

Non-pyrogenic

See Boston Scientific cryoablation needle and accessory product IFUs for specific information on pyrogenicity.

User Information

The ICEfx Cryoablation System is designed to be operated by medical professionals who have a thorough understanding of the technical principles, clinical applications, and risks associated with cryoablation procedures. Optional education is available from your Boston Scientific representative.

INTENDED USE

The ICEfx Cryoablation System is intended for cryoablative destruction of tissue during minimally invasive procedures; various Boston Scientific accessory products are required to perform these procedures. The ICEfx Cryoablation System is intended for use as a cryosurgical tool in the fields of general surgery, dermatology, neurology (including cryoanalgesia), thoracic surgery (with the exception of cardiac tissue), gynecology, oncology, and urology. This system is designed to destroy tissue (including prostate and kidney tissue, liver metastases, tumors, and skin lesions) by the application of extremely cold temperatures.

Patient Groups

The intended population includes patients intended for cryoablative destruction of tissue during surgical procedures.

INDICATIONS FOR USE

The ICEfx Cryoablation System is indicated for use as a cryosurgical tool in the fields of general surgery, dermatology, neurology (including cryoanalgesia), thoracic surgery (with the exception of cardiac tissue), gynecology, oncology, and urology. This system is designed to destroy tissue (including prostate and kidney tissue, liver metastases, tumors, and skin lesions) by the application of extremely cold temperatures.

The ICEfx Cryoablation System has the following specific indications:

- Urology - Ablation of prostate tissue in cases of prostate cancer
- Oncology - Ablation of cancerous or malignant tissue and benign tumors, and palliative intervention
- Dermatology - Ablation or freezing of skin cancers and other cutaneous disorders
- Gynecology - Ablation of malignant neoplasia or benign dysplasia of the female genitalia
- General surgery - Palliation of tumors, recurrent cancerous lesions, and ablation of breast fibroadenomas
- Thoracic surgery - (with the exception of cardiac tissue)

Clinical Benefit Statement

The ICEfx Cryoablation System, when used with various Boston Scientific accessory products, is intended to destroy tissue (including prostate and kidney tissue, liver metastases, tumors, and skin lesions) by the application of extremely cold temperatures during minimally invasive procedures.

The clinical benefit is measured by overall clinical outcomes, including acceptable safety and performance outcomes specific to the target anatomy.

CONTRAINDICATIONS

There are no known contraindications specific to the use of the ICEfx Cryoablation System.

WARNINGS

General

- The ICEfx Cryoablation System is designed to be operated by medical professionals who have a thorough understanding of the technical principles, clinical applications, and risks associated with cryoablation procedures.
- See Boston Scientific cryoablation needle and accessory product IFUs for Warnings specific for these products.
- Do not use this device for any purpose other than the stated intended use and indications for use.
- Do not modify the ICEfx Cryoablation System in any way. Only authorized Boston Scientific personnel or Boston Scientific-trained authorized personnel are to service an ICEfx Cryoablation System.
- The ICEfx Cryoablation System should be inspected and serviced periodically according to the system specifications. Service must be performed by authorized service engineers. Reference the **Installation, Calibration, and Servicing** Section for detailed information.
- Do not use the ICEfx Cryoablation System if the system is visibly damaged, exposing any internal components or sharp edges.
- Do not use an ICEfx Cryoablation System near magnetic resonance imaging (MRI) equipment.
- The ICEfx Cryoablation System should not be used adjacent to or stacked with other equipment.
- If the ICEfx Cryoablation Console is mounted on a cart, lock the wheels on the cart prior to using the system to avoid inadvertent movement of the system during a procedure.
- To avoid risk of electric shock, this equipment must only be connected to a hospital grade electrical power outlet with a protective earth.
- Do not start a cryoablation procedure before verifying that the ICEfx Cryoablation System and all ancillary equipment are fully operational.

- Use of cables other than those specified, with the exception of those sold by Boston Scientific for use as replacement parts for internal components, may result in increased emissions or decreased immunity of the ICEfx Cryoablation System.
- Only use non-MRI needles with the ICEfx Cryoablation System.
- Boston Scientific recommends using the Boston Scientific MTS to monitor the freeze / thaw temperatures for the intended treatment protocol and to monitor temperatures in the adjacent organs and structures.
- Do not use the needle if it is bent or damaged while attempting to unpack or use it. Never use a defective needle for a cryoablation procedure. A defective cryoablation needle that has a gas leak can cause a gas embolism in the patient.
- Do not kink, pinch, cut or pull excessively on the needle tubing. Damage to needle handle or tubing may cause the needle to become unusable.
- Have sufficient argon gas available to conduct the planned cryoablation procedure: the number and type of needles, gas cylinder size, pressure and rate of gas flow affect the required gas volume (reference the **SYSTEM SPECIFICATIONS** Section for gas purity requirements). At least one full spare cylinder should be available for each treatment.
- High-pressure gas is dangerous if handled improperly. Local laws and safety rules regarding pressurized gas systems, reservoirs and components should always be observed.
- Ensure that the argon gas cylinder is chained to a wall or an approved cart to prevent inadvertent tipping of the cylinder.
- Do not connect the ICEfx Cryoablation System to a gas supply exceeding 6000 psi (414 bar, 41.4 MPa) to avoid damage to internal system components.
- The ICEfx Cryoablation System should not be operated in the presence of flammable fumes, e.g., flammable anesthetics or volatile substances.
- Do not bend or kink the gas supply line. Sharp bends or kinks may compromise the integrity of the gas supply line.
- Do not roll the ICEfx Cryoablation System over the gas supply line; such activity may damage the line.

Procedural

- Prior to starting a cryoablation procedure, set up the ICEfx Cryoablation System (reference the **System Setup** Section) and then perform Needle Integrity and Functionality Tests. Tests must be successfully completed in order to begin the procedure.
- Do not use the needle if there is no ice formation during the Freeze phase. Obtain a new needle and repeat the testing procedure.
- Do not use the needle if bubbles are seen escaping from the needle during Needle Integrity and Functionality Test.
- Ensure adequate measures are taken to protect organs and structures adjacent to the targeted tissue.
- Ensure that the MTSs are functioning properly before inserting into the patient by verifying that they are reading a reasonable room temperature.
- The sterile field and sterility of cryoablation needles should be maintained at all times. Do not contaminate the distal end of the sterile cryoablation needle.
- Avoid contact with the distal portion of the cryoablation needle to maintain sterility during testing.
- Continuously monitor needle insertion, needle positioning, iceball formation and removal using image guidance (such as direct visualization, ultrasound, or Computed Tomography (CT)) to ensure adequate tissue coverage and to avoid damage to adjacent structures.

- Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 in) to any part of the ICEfx Cryoablation System, including cables specified for use with the system. Otherwise, degradation of the performance of this equipment could result.
- Confirm the high-pressure argon gas supply is securely connected before opening the gas cylinder.
- Secure the safety cable at the end of the gas supply line to the Console before connecting the argon gas supply line to the argon gas inlet. The safety cable provides backup protection if the gas supply line becomes inadvertently disconnected from the system. Do not use a gas supply line with a missing safety cable. Doing so could compromise the safety of personnel in the room. Contact Boston Scientific Technical Assistance Center for further instructions.
- Each needle must be locked into a needle channel before initiating a cryoablation procedure to avoid the risk of forceful ejection of the needles while under gas pressure.
- If the needles are still connected, do not unlock the channels or disconnect the needles from the needle connection panel until all operations in the channel have completed.
- Use **Freeze** and **Thaw** operations only when the needle is placed in the target tissue.
- If a needle needs to be added during a procedure, connect the needle to an open channel. If there is no open channel available, reference the **Adding a CX Needle During a Cryoablation Procedure** Section for the instructions for adding a CX needle to a channel with an open port.
- Needle handles and the gas line may frost during freezing. Avoid prolonged contact with frosted portions of a needle handle to avoid unintended thermal tissue damage to the patient or clinician.
- Needle tubing may become extremely cold when conducting freeze cycles during a cryoablation procedure. It is important that a patient's skin is protected from direct contact with needle tubing to avoid the potential for thermal injury to the patient. Ensure an appropriate insulating barrier is placed as needed (such as towels) or other method is employed to prevent needle tubing from touching a patient's skin.
- The needle handle may become warm during active thawing. Pay attention to the needle handle position. Prolonged contact with warm portions of the needle handle could cause unintended thermal tissue damage/burn to the patient or clinician.
- Active thawing produces heat along the distal needle shaft. Use care to avoid thermal injury/burn to nontargeted tissues.
- Ensure adequate thawing or cooling before attempting to remove needles from the patient.
- Discontinue all needle operation prior to needle removal to minimize risk of thermal injury and/or tissue injury.
- When conducting **FastThaw** Function or when conducting **Cautery** function for Track Ablation, be alert for the Active Zone Indicator as the needle is withdrawn to prevent unintended tissue damage from the hot needle.
- If a needle without CX functionality is connected to the ICEfx Cryoablation System, active thawing is not possible. Needle release from an iceball will only be achievable with passive thaw. Use image guidance to ensure adequate thawing.
- Do not touch the ICEfx Cryoablation Console while touching the patient to avoid the risk of shocking the patient if an inadvertent electrical fault exists.
- Do not touch the screen if the touch screen monitor goes blank for more than five (5) seconds during a procedure. Immediately turn off power to the system and end the procedure to avoid inadvertent activation of needles.
- Warn procedure personnel prior to venting the ICEfx Cryoablation System to avoid startling them.
- If it is difficult to loosen the pressure gauge connected to the cylinder, or the argon gas supply line cannot be disconnected from the argon inlet connection, do not apply excessive force to release the gas supply line or to loosen the pressure gauge. The gas line may still be under pressure.
- Do not pull on the power cord. Grasp the plug, not the power cord, to disconnect the device from the wall socket.
- Dispose of device and accessories in accordance with the **Disposal** Section.

PRECAUTIONS

General

- Carefully read all instructions prior to use. Failure to observe all warnings and precautions may result in complications.
- Do not use the ICEfx Cryoablation System if any moisture or condensation is present on the surfaces of the system. Allow the system to dry completely for 12 hours prior to powering up the system. Powering up the system containing moisture or condensation could result in permanent damage to the electrical boards, causing the system to be inoperable.
- Take precautions to avoid potential electrostatic discharge. If an electrostatic discharge occurs after touching the monitor, the screen may flicker and MTS temperature readings may be inaccurate for a few seconds. The system will remain functional and the monitor will refresh momentarily.
- No data regarding cryoablation in combination with other therapies is available from Boston Scientific.
- Use the Boston Scientific supplied USB flash drive only to export reports or update software. Other data or software may corrupt the ICEfx Cryoablation System.
- Do not connect any other USB equipment to the ICEfx Cryoablation System USB port.
- Do not use a USB extension cable to connect the USB flash drive to the USB port. Connect the USB flash drive directly to the USB port provided on the ICEfx Cryoablation System. Use of a USB extension cable may result in electromagnetic emissions exceeding regulatory limits.
- Select a unique Patient ID that does not reveal the patient's identity to other system users.

Handling

- Handle the ICEfx Cryoablation System with care. Rough handling may damage the system and cause it to become inoperable. The Console should never be tilted. If the Console is not attached to a cart, ensure that the Console is placed on a flat, stable surface.
- Do not set food, drinks, or other objects on top of the Console. Doing so may damage the system.
- Do not set heavy objects on the monitor when in the down position or on the monitor storage basin when the monitor is in the up position.
- Prior to lowering the monitor, ensure that no objects are located in the monitor storage basin. Use care when lowering the monitor into the monitor storage basin; do not use excessive force to avoid damaging the monitor.
- Use care when lowering or swiveling the touch screen monitor to avoid the potential for pinching fingers.
- If the Console is attached to the ICEfx Cryoablation Cart, use the cart handles to lift the system to clear any threshold that is higher than 2 cm. Two people, one on each side, should use the handles to lift the system.
- Clean the ICEfx Cryoablation System by following the instructions in the **Cleaning the ICEfx Cryoablation System** Section. Do not use cleaning agents such as Betadine Antiseptic Solution or bleach solution, which may damage the touch screen.
- When placing the ICEfx Cryoablation Console on or removing the ICEfx Cryoablation Console from the ICEfx Cryoablation Cart, lift the Console using the handles on the underside of the base. Refer to the **Removing the Console from the Cart** Section for instructions on how to place or remove the Console from the ICEfx Cryoablation Cart.
- Position the argon cylinder close enough to the Console to ensure that the gas supply line is not stretched and does not present a trip hazard.
- Direct the high-pressure argon gas supply line towards the floor and secure the line with the clip located on the rear of the Console to minimize the potential for tripping.

Procedural

- The ICEfx Cryoablation System should be placed in close proximity for accessible needle connection and use.
- Turn the ICEfx Cryoablation Console ON before connecting the argon cylinder to ensure the proper diagnostic tests are conducted.
- Confirm that the Vent Valve is closed and Argon Shutoff Valve is in the ON position before connecting the gas line to the Console.
- If the Console produces a continuous hissing sound, verify that the Vent Valve is fully closed. If the Vent Valve is fully closed and the hissing sounds persists, shut OFF the system using the power switch located on the rear of the Console (Figure 2). Close the gas supply using the cylinder valve. Contact Boston Scientific Technical Assistance Center.
- Failure to operate the ICEfx Cryoablation System within the working pressure limits indicated on the user interface (Table 7) may affect iceball formation.
- Boston Scientific recommends that only needles of the same type be placed together in a single channel. Using needles of differing types in a channel may affect the accuracy of the **Gas Indicator**.
- During use, avoid damage to the needle from other surgical instruments.
- Do not allow the MTS to touch the cryoablation needle during use; erratic temperature measurements may be displayed.
- Needles without CX functionality do not support i-Thaw Function or FastThaw Function. Needle release from an iceball will only be achievable with passive thaw.
- Needles without CX functionality do not support cautery function for the purposes of track ablation.
- The Cautery function should not be initiated if the Active Zone Indicator is visible outside the patient's skin.
- If a needle appears to be blocked, press the **Thaw** button to thaw the needle for at least one minute to clear the blockage.
- If a needle without CX functionality appears to be blocked, wait for any formed ice to passively thaw. Do not use the needle. Use a new needle to continue the procedure.
- When the argon cylinder pressure falls below the lower working pressure limit, the system displays an alert message. To ensure optimal performance, replace an argon cylinder if the pressure falls below the lower working pressure limit.
- Depressurize the system after the cryoablation procedure is completed (reference the **Shutting Down the System** Section).

ADVERSE EVENTS

The potential adverse events associated with the device and/or cryoablation procedure include, but are not limited to:

- Angina
- Arrhythmia
- Atelectasis
- Bladder spasms
- Bleeding/hemorrhage
- Burn/frostbite
- Cerebrovascular Accident (CVA)/stroke
- Cryoshock Phenomenon (e.g. multi organ failure, severe coagulopathy, disseminated intravascular coagulation (DIC))
- Death
- Distension
- Edema/swelling
- Ejaculatory dysfunction
- Embolism (air, device, thrombus)
- Erectile dysfunction
- Fever
- Fistula
- Fracture
- Gastrointestinal symptoms (e.g. nausea, vomiting, diarrhea, constipation)
- Healing, impaired
- Hematoma
- Hematuria
- Hemothorax
- Hepatic dysfunction/failure
- Hernia
- Hypertension
- Hypotension
- Hypothermia
- Ileus
- Impotence
- Infection/abscess/sepsis
- Inflammation
- Muscle spasm
- Myocardial infarction
- Necrosis
- Need for additional intervention or surgery
- Nerve Injury

- Neuropathy
- Obstruction
- Pain/discomfort
- Perforation (including organ and adjacent structures)
- Pericardial Effusion
- Perirenal fluid collection
- Pleural effusion
- Pneumatosis (air or gas in an abnormal quantity and/or place in the body)
- Pneumothorax
- Post ablation syndrome (e.g. fever, pain, nausea, vomiting, malaise, myalgia)
- Renal insufficiency/failure
- Renal parenchymal or capsule fracture
- Respiratory distress/insufficiency/failure
- Scrotal edema
- Stenosis/stricture
- Subcutaneous emphysema
- Thrombosis/thrombus
- Tissue damage
- Transient ischemic attack (TIA)
- Tumor cell seeding
- Urethral sloughing
- Urinary frequency/urgency
- Urinary incontinence
- Urinary retention
- Urinary tract infection
- Vasovagal response
- Vessel trauma (e.g. dissection, injury, perforation, pseudoaneurysm, rupture, or other)
- Wound infection

CONFORMANCE TO STANDARDS

Electrical Specifications

- Input Voltage: 100 VAC to 240 VAC, single phase
- Input Frequency: 50 Hz to 60 Hz
- VA Rating: 250 VA
- IP Rating: IP10
- Fuse Rating: T 2AL
- Electrical Protection: Class I, Type BF protection against shock
- Signal Input/Output Ports: one (1) USB 2.0 full-speed port

EM disturbances beyond the compliance limits may cause loss of ability to control the touch screen.

Electromagnetic Compatibility and Immunity (EMC and EMI)

The ICEfx Cryoablation System needs special precautions regarding electromagnetic compatibility (EMC) and needs to be installed and placed into service according to the EMC information provided below.

The ICEfx Cryoablation System has been tested in the operating room environment for electromagnetic compatibility (EMC) and electromagnetic interference (EMI) compliance. The ICEfx Cryoablation System has been tested to be in compliance with IEC 60601-1-2 and EN 55011.

Portable and mobile radiofrequency (RF) communications equipment can affect the ICEfx Cryoablation System, causing it to operate improperly.

Table 1. Cable Lengths

Cable	Length
Power Cable	4.6 m (15 ft)
Gas Tubing (connected to needles)	2.5 m (8 ft)
Gas Supply Line (connected to argon cylinder)	Available lengths: 5 m (16 ft), 8 m (26 ft)

NOTE: The Gas Supply Line is available in more than one length to accommodate procedure room variations.

WARNING: Use of cables other than those specified, with the exception of those sold by Boston Scientific for use as replacement parts for internal components, may result in increased emissions or decreased immunity of the ICEfx Cryoablation System.

WARNING: The ICEfx Cryoablation System should not be used adjacent to or stacked with other equipment.

WARNING: Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the ICEfx Cryoablation System, including cables specified for use with the system. Otherwise, degradation of the performance of this equipment could result.

Table 2. Electromagnetic Emissions

Guidance and manufacturer's declaration – electromagnetic emissions		
The ICEfx Cryoablation System is intended for use in a professional healthcare facility electromagnetic environment with the compliance levels specified below. The customer or the user of the ICEfx Cryoablation System should assure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment – guidance
RF emissions CISPR 11	Group 1	The ICEfx Cryoablation System uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class A	
Harmonic emissions, IEC 61000-3-2	Class A	
Emissions test	Compliance	Electromagnetic environment – guidance
Voltage fluctuations/ flicker emissions, IEC 61000-3-3	Complies	
<p>NOTE: The emissions characteristics of this equipment make it suitable for use in industrial areas and hospitals (CISPR 11 class A). If it is used in a residential environment (for which CISPR 11 class B is normally required) this equipment might not offer adequate protection to radio-frequency communication services. The user might need to take mitigation measures, such as relocating or reorienting the equipment.</p>		

Zastarano uporabnik. Ne uporabljajte. Outdated version. No utilize. Versión obsoleta. Nemojte uporabljati. Version périmée. Ne pas utiliser. Παλιά έκδοση. Μην χρησιμοποιείτε. Outdated version. No utilize. Versión obsoleta. Nemojte uporabljati. Zastarajela verzija. Ne koristite. Úrelet útgáfa. Notið ekki. Versione obsoleta. Nemojte uporabljati. Novecojsi versija. Nemojte uporabljati. Pasenusi versija. Neizmantot. Elavult versio. Ne használjate. Dit is een verouderde versie. Niet gebruiken. Utdatert versjon. Skal ikke brukes. Wersja przeterminowana. Nie używać. Versão obsoleta. Não utilize. Versiune expirată. A nu se utiliza. Zastaraná verzia. Nepoužívať. Vanhentunut versio. Älä käytä. Föråldrad version. Använd ej. Güncel olmayan sürüm. Kullanmayın.

Table 3. Electromagnetic Immunity

Guidance and manufacturer's declaration – electromagnetic immunity			
The ICEfx Cryoablation System is intended for use in a professional healthcare facility electromagnetic environment with the immunity compliance levels specified below. The customer or the user of the ICEfx Cryoablation System should assure that it is used in such an environment.			
Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 8 kV contact ± 15 kV air	± 8 kV contact ± 15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	± 2 kV for power supply lines ± 1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 0.5 kV, ± 1 kV line(s) to line(s) ± 0.5 kV, ± 1 kV, ± 2 kV line(s) to earth	± 0.5 kV, ± 1 kV line(s) to line(s) ± 0.5 kV, ± 1 kV, ± 2 kV line(s) to earth	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	0 % UT; 0.5 cycle at 0°, 45°, 90°, 135°, 180°, 225°, 270°, and 315°. 0 % UT; 1 cycle 70 % UT; 25 cycles/30 cycles at 0° and 50Hz/60Hz. 0 % UT; 250 cycles/300 cycles at 50Hz/60Hz.	0 % UT; 0.5 cycle at 0°, 45°, 90°, 135°, 180°, 225°, 270°, and 315°. 0 % UT; 1 cycle 70 % UT; 25 cycles/30 cycles at 0° and 50Hz/60Hz. 0 % UT; 250 cycles/300 cycles at 50Hz/60Hz.	Mains power quality should be that of a typical commercial or hospital environment. If the user of the ICEfx Cryoablation System requires continued operation during power mains interruptions, it is recommended that the ICEfx Cryoablation System be powered from an uninterruptible power supply or a battery.
Power frequency (50 Hz/60 Hz) magnetic field IEC 61000-4-8	30 A/m	30 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
NOTE: UT is the a.c. mains voltage prior to application of the test level.			

Table 4. Electromagnetic Immunity for Systems that are Not Life-supporting

Guidance and manufacturer's declaration – electromagnetic immunity			
The ICEfx Cryoablation System is intended for use in a professional healthcare facility electromagnetic environment with the immunity compliance levels specified below. The customer or the user of the ICEfx Cryoablation System should assure that it is used in such an environment.			
Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic environment - guidance
Conducted RF IEC 61000-4-6	3 Vrms, 6 Vrms in ISM bands over 150 kHz to 80 MHz 80 % AM at 1 kHz	3 Vrms, 6 Vrms in ISM bands over 150 kHz to 80 MHz 80 % AM at 1 kHz	Portable and mobile RF communications equipment should be used no closer to any part of the ICEfx Cryoablation System, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance: $d = \left[\frac{3,5}{3} \right] \sqrt{P}$
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.7 GHz 80 % AM at 1 kHz	3 V/m 80 MHz to 2.7 GHz 80 % AM at 1 kHz	$d = \left[\frac{3,5}{3} \right] \sqrt{P}$ 80 MHz to 800 MHz $d = \left[\frac{7}{3} \right] \sqrt{P}$ 80 MHz to 2.5 GHz
Radiated RF Proximity fields IEC 61000-4-3 (per IEC 60601-1-2 Ed. 4)	9 V/m to 28 V/m per IEC 60601-1-2 Ed. 4, Table 9	9 V/m to 28 V/m per IEC 60601-1-2 Ed. 4, Table 9	where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey ^a , should be less than the compliance level in each frequency range ^b . Interference may occur in the vicinity of equipment marked with the following symbol: 
NOTE 1: At 80 MHz and 800 MHz, the higher frequency range applies.			
NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			
^a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the ICEfx Cryoablation System is used exceeds the applicable RF compliance level above, the ICEfx Cryoablation System should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the ICEfx Cryoablation System.			
^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.			

Table 5. Recommended separation distances between portable and mobile RF communications equipment and the ICEfx Cryoablation System

Recommended separation distances between portable and mobile RF communications equipment and the ICEfx Cryoablation System			
The ICEfx Cryoablation System is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the ICEfx Cryoablation System can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the ICEfx Cryoablation System as recommended below, according to the maximum output power of the communications equipment.			
Rated maximum output power of transmitter Watts (W)	Separation distance according to frequency of transmitter Meters (m)		
	150 kHz to 80 MHz	80 MHz to 800 MHz	800 MHz to 2.5 GHz
	$d = \left[\frac{3,5}{3} \right] \sqrt{P}$	$d = \left[\frac{3,5}{3} \right] \sqrt{P}$	$d = \left[\frac{7}{3} \right] \sqrt{P}$
0.01	0.117	0.117	0.233
0.1	0.369	0.369	0.738
1	1.167	1.167	2.333
10	3.689	3.689	7.379
100	11.667	11.667	23.333
For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer. NOTE 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies. NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			

HOW SUPPLIED

Device Details

The ICEfx Cryoablation System is supplied non-sterile. Boston Scientific accessory products required to perform the cryoablation procedure are supplied separately.

Do not use if package is damaged or unintentionally opened before use.

Do not use if labeling is incomplete or illegible.

Handling and Storage

Operating Conditions

- Temperature: 10 °C to 40 °C
- Relative Humidity: 30 % to 75 %

Storage Conditions

- Temperature: -15 °C to 50 °C
- Relative Humidity: 10 % to 90 %

Transportation Conditions

When shipping the ICEfx Cryoablation System, use the original shipping container to prevent damage from occurring during transport. If the original shipping container is not available, the customer takes responsibility to ensure proper transport conditions are satisfied or contacts Boston Scientific Technical Assistance Center to obtain the appropriate shipping container.

Removing the Console from the Cart

The following procedure describes how to detach the Console from the cart if required for storage, or if the Console is being returned to Boston Scientific for service.

1. Loosen the retaining screw that secures the cart to the Console.
2. Lift the Console off the cart and transfer to a stable surface.

CAUTION: Do not set food, drinks, or other objects on top of the Console. Doing so may damage the system.

OPERATIONAL INSTRUCTIONS

WARNING: The ICEfx Cryoablation System is designed to be operated by medical professionals who have a thorough understanding of the technical principles, clinical applications, and risks associated with cryoablation procedures. Optional education is available from your Boston Scientific representative.

CAUTION: Carefully read all instructions prior to use. Failure to observe all warnings and precautions may result in complications.

Additional Required Items

Accessory Products Used to Conduct Cryoablation Procedures

NOTE: See product specific Instructions for Use.

The following needles must be used with the ICEfx Cryoablation System:

- **Boston Scientific Cryoablation Needles:** The cryoablation needles are specifically designed for use with the Boston Scientific Cryoablation Systems and are available in a range of configurations that produce various iceball sizes and shapes, allowing the clinician to match the needles to the desired ablation zone. The cryoablation needles are supplied sterile.
-

WARNING: Only use non-MRI needles with the system.

Optional Accessory Items:

- **Cryoablation Needle Channel ID Stickers:** Cryoablation Needle Channel ID Stickers are affixed to the needle tubing to provide easier identification of the needles during a cryoablation procedure. Contact Boston Scientific Technical Assistance Center to order Cryoablation Needle Channel ID Stickers.
- **Boston Scientific MTS:** The MTS contain four sensor locations along the distal shaft of the needle to monitor temperature near the target site and adjacent critical structures.
- **Boston Scientific Urethral Warming Set:** The Urethral Warming Set is the conduit that circulates warm saline through the urethra during a prostate cryoablation procedure.

The following items used with the ICEfx Cryoablation System are re-usable and should be cleaned and/or sterilized in accordance with the Instructions for Use that accompanies each product.

The following items are needed to conduct cryoablation procedures and are not available from Boston Scientific.

- **Argon gas cylinder(s)**
-

NOTE: The argon gas must meet the purity requirements specified in the **SYSTEM SPECIFICATIONS** Section.

Installation, Calibration, and Servicing

Only Boston Scientific or authorized personnel are permitted to conduct service and system preventive maintenance on the system. Preventive maintenance on the ICEfx Cryoablation System is required every two years. Scheduled Preventative Maintenance must be completed to maintain system performance and safety.

WARNING: Do not modify the ICEfx Cryoablation System in any way. Only authorized Boston Scientific personnel or Boston Scientific-trained authorized personnel are to service an ICEfx Cryoablation System. If service is required, contact Boston Scientific Technical Assistance Center.

The ICEfx Cryoablation System will display a reminder on the screen approximately one month prior to when the preventive maintenance is due. If the reminder message appears and the preventive maintenance has not yet been scheduled, contact Boston Scientific Technical Assistance Center to schedule the service.

Preparation

System Setup

Table 6 depicts the order and steps of ICEfx Cryoablation System Setup. This section describes each step in-detail.

Table 6. System Setup Flow

1	Set Up Console	<ul style="list-style-type: none"> - Confirm the availability of gas, needles, and accessories - Ensure the Console is placed on a flat, stable surface (if not attached to the ICEfx Cryoablation Cart) - Lock the brakes on the cart (if using the ICEfx Cryoablation Cart) - Plug the power cord into a hospital grade power outlet - Ensure the Vent Valve is closed and the Argon Shut off Valve is ON - Power the ICEfx Cryoablation Console ON - Login
2	Connect Argon	<ul style="list-style-type: none"> - Attach the safety cable and connect the argon gas supply line between the Console and argon cylinder - OPEN the valve on the argon cylinder - Verify gas pressure is within working limits for the cryoablation procedure

Preparation for Use

Prior to using the ICEfx Cryoablation Console, inspect the Console, power cord, brakes, safety cable, gas supply line, and the monitor touch screen to ensure they are not damaged.

If any of the components are damaged, contact Boston Scientific Customer Service.

CAUTION: Do not use the ICEfx Cryoablation System if any moisture or condensation is present on the surfaces of the system. Allow the system to dry completely for 12 hours prior to powering up the system. Powering up the system containing moisture or condensation could result in permanent damage to the electrical boards, causing the system to be inoperable.

Attaching the Console to the ICEfx Cryoablation Cart

1. Position the Console on the ICEfx Cryoablation Cart. The depression on the bottom of the Console rests on the dome on top of the ICEfx Cryoablation Cart.

Figure 5. Position ICEfx Cryoablation Console on Cart

2. Tighten the retaining screw on the cart.

Figure 6. Tighten Retaining Screw on Cart

Setting up the Console

CAUTION: The Console should never be tilted. If the Console is not attached to a cart, ensure that the Console is placed on a flat, stable surface.

1. Verify that the power switch, Argon Shutoff Valve, and Vent Valve on the rear of the Console are easily accessible.

NOTE: Allocate space with adequate ventilation and free airflow. In order to guarantee proper ventilation, always keep the sides of the Console at least 0.5 m (20 in) from walls or other obstructions to airflow.

2. Lock the wheels on the cart.
 3. Plug the power cord into a hospital grade power outlet (electrical mains outlet) with a ground connection. Boston Scientific recommends using a stable and uninterruptible power outlet.
-

NOTE: If the power supply to the ICEfx Cryoablation Console is not stable or is noisy, the MTS temperature readings may be inaccurate.

WARNING: To avoid risk of electric shock, this equipment must only be connected to a hospital grade electrical power outlet with a protective earth.

OPTIONAL: When performing a prostate cryoablation procedure, set up a urethral warmer system by following the directions provided in the Urethral Warming Set Instructions for Use.

4. Confirm that the Argon Shutoff Valve on the rear of the Console is in the ON position. If needed, turn the knob counterclockwise to the ON position.


Figure 7. Argon Shutoff Valve - ON Position

5. Confirm the Vent Valve is fully closed. If needed, turn the knob clockwise to the CLOSED position.


Figure 8. Vent Valve - Fully Closed

6. Raise the monitor to the UP position and adjust it to a comfortable viewing angle.
-

CAUTION: Use care when swiveling the touch screen monitor to avoid the potential for pinching fingers.

7. Turn the Console on using the power switch located on the rear of the Console.

During boot-up, the system performs several diagnostic tests to verify that the hardware and software are operating properly.

NOTE: If the system was improperly shutdown after the previous procedure, the boot-up process could take up to 2 minutes.

NOTE: It is important to power up the system prior to connecting the gas to the system. If the system is not powered up prior to connecting the gas, the diagnostic tests will not be performed by the software.

The diagnostic tests check:

- The correct version of firmware is running on the system.
- Critical system components including the solenoid valves, internal power supplies, cooling fans, pressure transducers, and temperature measuring circuits.

A failure that prevents use of the system displays a message directing you to contact Boston Scientific Customer Service. Reference the **Displayed Messages** Section for information on displayed messages.

If the ICEfx Cryoablation System software detects pressurized gas in the system and the gas supply line is not connected, a message will display requesting that you manually vent the gas from the system.

Login

The *Login Screen* appears after the boot-up process is complete.

1. Enter your assigned **Username** and **Password** using the on-screen virtual keyboard.
2. Press **Login**.

NOTE:

- The username and password are not case-sensitive. Numbers are displayed when the keyboard is in upper case. To change text case, use the shift button on the virtual keyboard.
- If you leave the user interface idle for a pre-set time with no activity, the ICEfx Cryoablation System software requires you to re-enter your password to unlock the user interface.

Additional Login Options:

Forgot Password	If you forgot your password, contact your System Administrator and request that the administrator log in, go to the <i>Manage Users Screen</i> and change your password. Alternatively, press the Forgot Password button on the <i>Login Screen</i> . A displayed message provides a challenge to be relayed to Boston Scientific Technical Assistance Center. The Boston Scientific Technical Assistance Center will provide a response to enter on the screen using the virtual keyboard. Once a correct response is received, your password is reset and you are given the opportunity to change your password in the <i>Settings Screen</i> .
Emergency Login	In case of emergency, press the Emergency Login on the <i>Login Screen</i> . A message is displayed with a challenge. Call Boston Scientific Technical Assistance Center to receive the proper response to enter, then press the Login button. NOTE: This action does not reset your password. After successful login, the <i>Start Case Screen</i> appears (Screen 4).

Connecting the Argon Gas

WARNING: Do not connect the ICEfx Cryoablation System to a gas supply exceeding 6000 psi (414 bar, 41.4 MPa) to avoid damage to internal system components.

WARNING: Ensure that the argon gas cylinder is chained to a wall or an approved cart to prevent inadvertent tipping of the cylinder.

WARNING: Have sufficient argon gas available to conduct the planned cryoablation procedure: the number and type of needles, gas cylinder size, pressure and rate of gas flow affect the required gas volume (reference the **SYSTEM SPECIFICATIONS** Section for gas purity requirements). At least one full spare cylinder should be available for each treatment.

CAUTION: Turn the ICEfx Cryoablation Console ON before connecting the argon cylinder to ensure the proper diagnostic tests are conducted.

CAUTION: Position the argon cylinder close enough to the Console to ensure that the gas supply line is not stretched and does not present a trip hazard.

1. On the rear of the Console, ensure that the Vent Valve is in the CLOSED position (see Figure 8).
2. Remove the humidity cover from the argon inlet on the Console.
3. Secure the safety cable at the end of the gas supply line to safety clip on the Console (see Figure 9).

WARNING: Ensure that the safety cable is properly secured to the Console in case the gas supply line becomes accidentally disconnected.


Figure 9. Attach Safety Cable to Console

4. Connect the argon high-pressure gas supply line to the argon inlet of the Console using the quick-connect connector located at the rear of the Console (see Figure 10).


Figure 10. Connect Argon Gas Supply Line

- Connect the argon high-pressure gas supply line to the argon cylinder by securing the gauge assembly adapter onto the cylinder connection (see Figure 11).


Figure 11. Gas Cylinder Setup

- Carefully turn the cylinder valve on the argon gas cylinder counterclockwise for one quarter turn. Ensure that the pressure reading on the gauge responds immediately. Turn the cylinder valve further counterclockwise (approximately one full turn) to open the gas cylinder so that sufficient gas flow is present.


Figure 12. Open Argon Gas Cylinder

If no argon pressure is displayed on the system pressure gauge, ensure the Argon Shutoff Valve is in the ON position.

OPTIONAL:

The EZ-Connect2 Dual Cylinder Adapter connects two cylinders of argon gas to the Console to support a cryoablation procedure. A four-way adapter assembly with argon pressure gauge connects a gas supply line, the primary gas cylinder, and an auxiliary gas supply line.

If using the optional EZ-Connect2 Dual Cylinder Adapter, connect the gas supply line with four-way pressure gauge adapter assembly to the primary argon cylinder by securing the gauge assembly adapter onto the cylinder connection.

- Connect the end of the gas supply line to the argon inlet of the Console using the quick-connect connector.
- Connect the auxiliary gas supply line to the four-way adapter assembly using the quick-connect connector located on the end of the auxiliary gas supply line.
- Connect the opposite end of the auxiliary gas supply line to the second argon cylinder by securing the end of the auxiliary line to the cylinder connection.
- Open the primary cylinder valve first and use this cylinder until depleted. Do not open the cylinder valve on the second cylinder until the initial cylinder has been depleted.
- Reference the **Changing Argon Cylinders during a Procedure** Section for instructions on changing the gas cylinder during a procedure, if the second cylinder is also depleted during the procedure.



Figure 13. EZ-Connect2 Dual Cylinder Adapter

- 1 Auxiliary Gas Supply Line
- 2 Four-way Adapter Assembly with Pressure Gauge
- 3 Gas Supply Line

WARNING: Have sufficient argon gas available to conduct the planned cryoablation procedure: the number and type of needles, gas cylinder size, pressure and rate of gas flow affect the required gas volume (reference the **SYSTEM SPECIFICATIONS** Section for gas purity requirements). At least one full spare cylinder should be available for each treatment.

7. Ensure the Gas Indicator (Screen 1) displays the minimum working pressure before initiating a procedure (Table 7). The Gas Indicator should indicate the pressure is within the green range.

Screen 1: Gas Indicator

If the software detects that the pressure reading for the argon cylinder is less than 3.45 bar (50 psi), a message displays on the user interface indicating that no gas is connected. Connect the argon cylinder to the Console.

Table 7. Working Argon Gas Pressure

Nominal Working Pressure	Working Pressure Limits
3500 psi	3200 psi to 3800 psi
241 bar	221 bar to 262 bar
24.1 MPa	22.1 MPa to 26.2 MPa

CAUTION:

- When the argon cylinder pressure falls below the lower working pressure limit, the system displays an alert message. To ensure optimal performance, replace an argon cylinder if the pressure falls below the lower working pressure limit.
- Failure to operate the ICEfx Cryoablation Console within the working pressure limits may affect iceball formation.
- If the Console produces a continuous hissing sound, verify that the Vent Valve is fully closed. If the Vent Valve is fully closed and the hissing sounds persists, shut OFF the system using the power switch located on the rear of the Console (Figure 2). Close the gas supply using the cylinder valve. Contact Boston Scientific Technical Assistance Center.

Getting Started With The User Interface

The ICEfx Cryoablation System provides a graphical user interface that facilitates rapid communication between the user and the system through a touch screen interface.

This section provides information to help familiarize the user with the software user interface.

Refer to the **Performing The Cryoablation Procedure** Section for detailed guidance on how to perform Needle Integrity and Functionality Testing and conduct a cryoablation procedure.

Document Conventions

Throughout the User Manual typeset conventions represent various user interface sections, software buttons, positions, and steps.

- **Software Screen** section
- **Control** button
- ON position
- *OPTIONAL* = optional or alternate step

Logging In

When the system is powered on, the *Login Screen* appears after the boot-up process is complete.

To login, enter the **Username** and **Password**, and then press **Login**.

Navigating the User Interface

Setup and Gas Venting Wizards


The software provides wizards to guide the user through each of the steps involved with starting or ending a cryoablation procedure.

The Setup wizard appears when **START CASE** is pressed and guides the user through system setup including needle connection and testing. To skip the wizard, press **SKIP SETUP WIZARD**.

The Gas Venting wizard appears when **END CASE** is pressed and instructs the user on how to properly vent the system and disconnect the argon gas supply line. To skip the wizard, press **RETURN HOME**.

Title Bar

The Title Bar contains the Options Menu, the argon remaining, and the Pressure Gauge / Gas Indicator.

Screen 2: Title Bar


Table 8. Title Bar

#	Control	Description
1	Options Menu	Pressing the Options Menu icon displays an expanded menu of options. Reference the Options Menu Section for descriptions of the available options.
2	Argon Remaining	Displays an estimate of time remaining before depletion of the cylinder assuming that all connected needles are operated simultaneously.
3	Pressure Gauge / Gas Indicator	Displays the supplied pressure of the argon gas cylinder. NOTE: The ICEfx Cryoablation System contains an internal regulator that regulates the argon pressure to appropriate working limits. The Gas Indicator updates in real-time as needles are disconnected or additional needles are connected and as the freeze intensity is adjusted.

Options Menu

The Options Menu provides access to additional features, such as accessing reports and changing settings.

To view and select options, press the  icon in the upper left of the screen and select the desired option from the expanded menu.


Screen 3: Options Menu

Table 9. Options Menu

#	Button	Description
1	Settings	Configure various system settings (reference the Configuring Settings Section). NOTE: Some configuration parameters are limited to administrative and/or service users only.
2	Reports	When there is no active case, select Reports to view the contents of saved procedure reports and export reports to a USB flash drive. During a cryoablation procedure, select Reports to display details for the current case. NOTE: Administrative users may also delete reports.
3	User Manual	View information on how to access an electronic version of the User Manual.
4	Log Out	Log out of the system.
5	Shut Down	Shut down the system.

Starting the Procedure


After logging into the system, the *Start Case Screen* appears. Press **START CASE** to start a cryoablation procedure. Pressing **START CASE** displays the *Select Tissue Type Screen*. Select an organ and the Start Case wizard is activated.

Screen 4: Start Case Screen

Table 10. Start Case Screen Controls

#	Button	Description
1	START CASE	Press START CASE to initiate the Start Case wizard. The wizard guides the user through the steps for setting up the system. When setup completes, the <i>Procedure Screen</i> is displayed.
2	Shut Down	Press Shut Down to initiate shut down of the system.

Screen 5: Select Tissue Type


Table 11. Select Tissue Type

#	Button	Description
1	Tissue Type	Select kidney, lung, liver, prostate, other.
2	CONFIRM	Press CONFIRM to proceed.
3	CANCEL CASE	To select a different tissue type, press CANCEL CASE to return to the <i>Start Case Screen</i> .

Entering Case Information

Use the *Edit Case Information Screen* to provide additional information on the cryoablation procedure. This is an optional step.

The *Edit Case Information Screen* is displayed at the beginning and end of a cryoablation procedure. At the end of the procedure, an option is provided to save the case information.


Screen 6: Edit Case Information Screen

Table 12. Edit Case Information Screen Controls

#	Control	Description
1	Patient ID	Type the ID for the patient in this field. CAUTION: Select a unique Patient ID that does not reveal the patient's identity to other system users.
2	Physician	Type the name of the physician performing the procedure in this field.
3	Report Name	Specify a name for the report. The report name will be displayed on the <i>Reports Screen</i> .
4	Notes	Enter notes about the procedure.
5	Hospital Name	Type the name of the hospital in this field.
6	Hospital Address	Type the address of the hospital in this field.
7	CANCEL CASE	Press CANCEL CASE to discard the entered registration information and return to the <i>Start Case Screen</i> .
8	DONE	Press DONE to save the registration information.

Using the Procedure Screen



Use the *Procedure Screen* to control and monitor a cryoablation procedure. Initially, only the controls related to needle testing are displayed as shown in Screen 7. After successful completion of Needle Integrity and Functionality Testing, the freeze (Screen 12) and thaw (Screen 13) controls are enabled for the tested needles.

Screen 7: Initial Procedure Screen

Table 13 describes the controls common to the test, freeze, and thaw functions. The sections which follow describe the controls relevant to each of these functions.


Table 13. Common Procedure Controls

#	Control	Description
1		Use the ALL channel to perform a freeze, thaw, or test function on ALL active channels simultaneously. Press the button corresponding to the function you want to perform and then press the Start button. NOTE: Pressing the ALL button provides an option to program freeze-thaw cycles for all active channels.
2		Channel button – Identifies the active channel(s). Channels 1 through 4 are labeled individually and contain independent controls for the currently selected function (test, freeze, thaw, or cautery). NOTE: Channels without connected needles appear dimmed. <i>Advanced Channel Controls:</i> Pressing the Channel button provides options to change the needle type selected for that channel and program freeze thaw cycles for the channel.
3		Channel Port button – Identifies the channel ports. NOTE: Ports without connected needles appear dimmed.
4	Timer	The Timer displays the elapsed time for the currently selected function and provides additional data specific to the function. NOTE: During testing, the timer counts down from two minutes instead of showing the elapsed time.
5		END CASE — Press END CASE to end the cryoablation procedure and return to the <i>Start Case Screen</i> . Pressing this button generates a confirmation request and a request to save the procedure report.

Viewing the Timer

The timer displays the elapsed time during freezing or thawing. The timer counts down from two minutes during testing. The timer counts down from the selected duration during the cautery function. The upper left and upper right portions of the timer provide additional data for the freeze, thaw, and cautery functions.

Screen 8 provides an example of a timer displayed during FastThaw Function.

Screen 8: Timer During FastThaw Function

Table 14 provides additional information on the timer.

Table 14. Timer Information

#	Control	Description
1	Mode	The upper left portion of the timer displays information relevant to the current function: <ul style="list-style-type: none"> Freezing – Displays the freeze intensity. Thawing – Displays the thaw mode (i-Thaw Function or FastThaw Function). Cautery – Displays the cautery function duration.
2		Thaw button – Press the Thaw button to change from FastThaw Function to i-Thaw Function. The background color of the button is orange if FastThaw Function is selected and is yellow if i-Thaw Function is selected.
3	Progress Indicator Bar	A color coded progress indicator beneath the elapsed time indicates the type and duration of the functions performed on the respective channel. The segment color indicates the function: <ul style="list-style-type: none"> Blue = Freezing. The shade of blue indicates freeze intensity. Stick is a light blue color. Orange = FastThaw Function Yellow = i-Thaw Function Pink = Cautery function The length of a segment is relative to the elapsed time for the respective function. Clicking the progress indicator bar opens a more detailed view.
4	Timer	Displays the elapsed time during freezing or thawing. During testing, the timer counts down from two minutes. During the cautery function, the timer counts down from the selected duration. The color of the timer digits corresponds to the selected function: <ul style="list-style-type: none"> Green = Testing Blue = Freezing. Stick is a unique blue. Orange = FastThaw Function Yellow = i-Thaw Function Pink = Cautery function

#	Control	Description
5	Temperature Indicator	<p>For CX needles, the upper right portion of the timer displays temperature data for the current function:</p> <ul style="list-style-type: none"> • Testing – Displays the needle temperature. • Freezing – Displays the internal gas temperature at the needle tip. • Thawing – Displays the estimated temperature range of the needle shaft. <p>NOTE: During the needle warming phase, the timer displays a rotating warming indicator.</p> <ul style="list-style-type: none"> • Cautery function – Displays the estimated temperature range of the needle shaft. <p>NOTE: During the needle warming phase, the timer displays a rotating warming indicator.</p>

MTS

The *Procedure Screen* displays the MTS identifiers, sensor locations, and sensor temperatures for the connected MTS. The coldest sensor for each connected MTS is indicated by the blue indicator on the corresponding button (if the MTS Display Coldest setting is turned *On in the Settings Screen*).


Screen 9: Procedure Screen – MTS

Table 15. MTS Controls

#	Control	Description
1	MTS Identifier (MTS 1, MTS 2)	Identifies the needle.
2	Sensor Location	Distance of sensor from the needle tip: Orange - 35 mm, Pink - 25 mm, Green - 15 mm, Yellow - 5 mm
3	Sensor Temperature	The coldest temperature for each sensor is indicated by the blue on the corresponding button (if the MTS Display Coldest setting is turned On in the <i>Settings Screen</i>).

Changing the MTS Identifier

Click the MTS identifier (MTS 1, MTS 2) to open the *Procedure Screen - Change MTS Identifier Name*. This screen makes it possible to change the MTS name. In this example, the name is changed to A4.

Screen 10: Procedure Screen – Change MTS Identifier Name

Table 16. Change MTS Name


#	Control	Description
1	Change MTS name	In the field provided, enter the new name and click OK , CANCEL , or DONE .
2	OK	Changes the name and closes the screen.
3	CANCEL	Returns to the previous screen.
4	DONE	Changes the name and closes the screen.

Initiating Testing

WARNING: If a needle needs to be added during a procedure, connect the needle to an open channel. If there is no open channel available, reference the **Adding a CX Needle during a Cryoablation Procedure** Section for the instructions for adding a CX needle to a channel with an open port.

Press the **Start** button to initiate Needle Integrity and Functionality Testing. Testing can be initiated for a single active channel or for **ALL** active channels simultaneously.

Needle Integrity and Functionality Testing is required prior to using any cryoablation needle. Reference the **Pre-Procedure Needle/MTS Testing** Section for detailed information on how to conduct Needle Integrity and Functionality Testing.


Screen 11: Procedure Screen during Testing

Table 17. Test Controls

#	Control	Description
1		Test button – Indicates the channel is ready for testing. No other controls on a channel are enabled until needle testing is completed. Upon completion of Needle Integrity and Functionality Testing, the procedure controls for the channel become active.
2		Press Start to initiate testing. The timer counts down from two minutes when testing is initiated.


Initiating a Freeze Cycle

To start a freeze cycle, press the **Freeze** button to select the desired freeze intensity (default is 100 %), and then press the **Start** button. A freeze cycle can be initiated for a single active channel or for ALL active channels simultaneously. During freezing, the timer displays the elapsed time and the selected freeze intensity. For CX needles, the needle tip temperature during freezing is displayed in the upper right corner of the timer.


Screen 12: Procedure Screen during Freezing

Table 18. Freeze Controls

#	Control	Description
1		Freeze button – Press to select the freeze function and to select the freeze intensity. Or select Stick.
2		Thaw button – Press the Thaw button to change from FastThaw Function to i-Thaw Function. The background color of the button is orange if FastThaw Function is selected and is yellow if i-Thaw Function is selected.
3		Press Start to initiate freezing or select Stick.
4	Temperature Indicator	For CX needles, displays the needle tip temperature during freezing.
5	Freeze Intensity	Pressing the Freeze button displays the selectable freeze intensities. Select a freeze intensity from the menu or select “Stick” intensity.  NOTE: The ICEfx Cryoablation System controls freeze intensity by adjusting the duration of argon flow over each 10 second block of time (e.g., 30% freeze intensity freezes for 3 seconds and idles for 7 seconds). NOTE: When Stick is selected, argon gas is flowed through the cryoablation needle at a low duty cycle to create a very thin layer of ice around the needle shaft. The thin layer of ice secures the needle to prevent accidental displacement while the clinician places other needles.

Screen 13 shows the result of clicking the progress indicator bar.

Screen 13: Expanded Channel History

Initiating a Thaw Cycle (CX needles only)

CAUTION: Needles without CX functionality do not support i-Thaw Function or FastThaw Function. Needle release from an iceball will only be achievable with passive thaw.

WARNING: The needle handle and distal needle shaft may become warm during active thawing. Pay attention to the needle handle position. Prolonged contact with warm portions of the needle handle or distal needle shaft could cause unintended thermal tissue damage to the patient or clinician.

To initiate a thaw cycle:

NOTE: A thaw cycle can be initiated for a single active channel or for ALL active channels simultaneously.

1. Press the **Thaw** button to select the thaw mode (i-Thaw Function or FastThaw Function).

NOTE:

- When the **i-Thaw** Function is operational, selecting the **FastThaw** Function will change the status.
- Similarly, when the FastThaw Function is selected, selecting the **i-Thaw** Function will change the status.

2. Press the **Start** button to initiate thawing.


3. Press **Stop** button to stop the thaw cycle.

WARNING: Ensure adequate thawing or cooling before attempting to remove needles from the patient.

WARNING: Discontinue all needle operation prior to needle removal to minimize risk of thermal injury and/or tissue injury.

NOTE: If needle sticking is experienced, use a slight, gentle twist of the needle followed by slow withdrawal.

Thaw Controls

Screen 14: Procedure Screen during Thawing (FastThaw Function)

Table 19. Thaw Controls

#	Control	Description
1		Displays the currently selected thaw mode (i-Thaw Function or FastThaw Function). NOTE: The FastThaw Function generates a temperature that is higher than that generated for i-Thaw Function, resulting in a faster thaw time.
2		Displays the estimated temperature range of the needle shaft during thawing. NOTE: During the needle warming phase, the timer displays a rotating warming indicator.
3		Thaw button — Press the Thaw button to change from FastThaw Function to i-Thaw Function. The background color of the button is orange if FastThaw Function is selected and is yellow if i-Thaw Function is selected.
4		Press Stop to stop thawing.

Using Cautery Function (CX needles only)

Needle track ablation is achieved using the Cautery function. During cautery function, the needle track is ablated through the use of thermal energy. Each needle is operated independently using the **Cautery** button for the channel and port in which the needle is connected.

WARNING: Continuously monitor needle insertion, needle positioning, iceball formation and removal using image guidance (such as direct visualization, ultrasound, or Computed Tomography (CT)) to ensure adequate tissue coverage and to avoid damage to adjacent structures.

WARNING: When conducting **FastThaw** Function or when conducting **Cautery** function for Track Ablation, be alert for the Active Zone Indicator as the needle is withdrawn to prevent unintended tissue damage from the hot needle

CAUTION:

- Needles without CX functionality do not support cautery function for the purposes of track ablation.
- The Cautery function should not be initiated if the Active Zone Indicator is visible outside the patient's skin.

To use cautery function:

1. From the *Procedure Screen*, press the **CAUTERY** button to display the cautery function controls.
2. Press the desired needle port (A or B) on each channel containing a needle you wish to track ablate. Only one needle per channel may activate the cautery function at a given time.
3. Press the **CAUTERY** button (Screen 14) for the channel to select the cautery function duration.

NOTE: The duration for the cautery function phase may be selected from 30 seconds to 3 minutes, adjustable in 30 second increments.

- During the needle warming phase, the timer displays a rotating warming indicator.
- When the cautery function threshold is reached, Cauterization automatically begins.
 - The timer will display Cauterizing in the upper left.
 - The temperature indicator in the upper right of the timer will display the estimated temperature range for the needle shaft.

NOTE: The shaft temperature range is displayed as a temperature range because tissue and procedure values affect the temperature.


- When the cautery function has terminated, the timer indicates the status as Stopped and the temperature display indicates that the needle is Cooling.

4. Press the **Start** button for the selected channel to start cautery function.
5. If desired, press the **Start** button to re-initiate track ablation on additional sections of the needle track.
6. After completing track ablation, carefully remove the needle.

WARNING: Ensure adequate thawing or cooling before attempting to remove needles from the patient.



- If needle sticking is experienced, use a slight, gentle twist of the needle followed by slow withdrawal.
- During needle removal, watch for the 10 mm wide Active Zone Indicator on the needle to guide needle withdrawal. The Active Zone Indicator is a marked band along the needle shaft, located 20 mm distal to the heating portion of the needle.

Cautery Function Controls

Screen 15: Cautery Function Screen

Table 20. Cautery Function Controls


#	Control	Description
1		PROCEDURE button – Press to exit Cautery function mode and return to the <i>Procedure</i> Screen.
2	 	Cautery button – Press the Cautery button to select the cautery function duration.
3		Press Start to initiate cautery function.
4	Cautery Function Status	Displays the Cautery Function Status.
5	Temperature Indicator	Displays the estimated temperature range of the needle shaft during cautery function. NOTE: During the needle warming phase, the timer displays a rotating warming indicator.
6		Press Stop to stop cautery function.

Advanced Channel Controls

The Advanced Channel Controls provide options to change the needle type for the selected channel, retest a needle, and program multiple freeze-thaw cycles.

Select Needle Type Control

1. To select or change the needle type for a channel, press the **Channel** button to display the needle type list.

Screen 16: Needle Type List

2. Select the proper needle type from the list.

Retest a Needle

To retest a needle, press the **Channel** button and then press  to initiate a needle retest.


Program Cycles Control

Each freeze-thaw cycle consists of a freeze cycle followed by a thaw cycle. Using the *Program Cycles Control*, the user can program multiple freeze-thaw cycles for an individual channel or for all active channels simultaneously.

Each programmed cycle consists of one or more steps. Each step represents a freeze or thaw cycle that is set for a specific running time.

To run a saved cycle sequence:

1. To run a saved cycle sequence, press the individual **channel** button or the **ALL** channel button, and then press **PROGRAM CYCLES**.
2. The Program Cycles dialog displays. Select the desired cycle sequence from the list and press **LOAD**.


Screen 17: Select Programmed Cycle

3. Press **START** to initiate the cycle.


To program freeze-thaw cycles:

1. To program freeze-thaw cycles for a specific channel, press the individual channel button or the **ALL** channel button, and then press **PROGRAM CYCLES**.
2. Press **NEW** on the *Program Cycles* dialog.


Screen 18: Program Cycles Dialog

The *Programming Cycles* Screen displays for the selected channel.

Screen 19: Programming Cycles Screen

3. Press the plus sign  at the bottom left of the grid to create a new step.
4. For each step:
 - Choose between freeze or thaw
 - If freeze is selected, choose the freeze intensity
 - Select the time to run for the step
5. Repeat Steps 3 and 4 to add additional steps.
6. Use the UP and DOWN arrows to change the order of steps. To delete a step, press the delete icon.
7. Press **SAVE** to save the cycle sequence.

Screen 20: Save Cycle Sequence Screen

8. Provide a name for the cycle sequence and press SAVE.


Reports

Procedure reports provide a summary of a cryoablation procedure. The reports contain the information provided in the *Edit Case Information Screen*, details of the freeze-thaw cycles, number of MTS needles, a graphical history of MTS readings, and a graphical history of freeze, thaw, and cautery phases.

Viewing All Reports

Press  and select **Reports** from the expanded menu to view the saved procedure reports (when there is no active case).

The *All Reports Screen* displays a list of all saved reports in the ICEfx Cryoablation System. Select the report you want to view and press **VIEW REPORT**.

Screen 21: All Reports Screen

When a report is selected from the *All Reports Screen*, the report that displays provides options to delete or export the report.

NOTE: Users with an administrative login ID may delete any report.

Sorting a Report

To sort the list by Report Name, Hospital Name, Doctor Name, or Organ Type, press the appropriate header section of the report list.


Exporting a Report

From the *All Reports Screen*, press **VIEW REPORT** to display the report you would like to export.

From the *View Report Screen*, click **EXPORT** to export the report details and select the export file type and name. Reports can be exported in HTML or CSV formats. Wait for confirmation before removing the USB flash drive from the Console.

Viewing a Report during a Cryoablation Procedure

Press  and select **Report** from the expanded menu to view the current patient data during an active case. When a report is viewed during a cryoablation procedure, the report that displays provides an option to refresh the details. Options to save, delete, and export are not available during a cryoablation procedure.

Screen 22: Procedure Report

Configuring Settings

Press and select **Settings** from the expanded menu to display the *Settings Screen*. The *Settings Screen* allows the user to configure settings used during a cryoablation procedure.

Table 21 lists and describes the configurable settings. After changing settings, press **Confirm** to save the setting changes.

Table 21. Settings Screen – Configurable Settings

Setting	Description
Date	Set the values for the month, day, and year. The date can only be changed by service personnel.
Language	Select the language in which the software is displayed.
Custom Fields	Enter custom names to label the two custom fields (Custom Field 1 and Custom Field 2) available for entering supplementary information on the <i>Edit Case Information Screen</i> .
Time	Set the values for the hour, minute, and time period (AM or PM). The time can only be changed by service personnel.
Argon Cylinder Volume	Select the argon cylinder volume according to the standard in the geographical region. The cylinder volume can only be changed by administrative or service personnel.
Time Zone	The time zone can be changed by administrative or service personnel. The ICEfx Cryoablation System automatically adjusts for Daylight Saving Time.
Volume Units	Select the argon cylinder volume units according to the standard in the geographical region.
Pressure Units	Select the argon cylinder pressure units according to the standard in the geographical region.
Temperature Units	Select the temperature units displayed for the MTS and CX needles.
Clock Type	Select the clock type (12 hour or 24 hour). The clock type can be changed by administrative or service personnel.
Inactivity Timeout	Select the desired duration from 30 minutes to 180 minutes that the system can be inactive before requiring you to re-enter your password.
Idle Timer	Turn On or Off the timer that tracks system inactivity. Turn On this option if you want the system to timeout when it is inactive for a specified Inactivity Timeout value.
Upload Registration	Use the drop-down menu to enable or disable the option to upload registration data with the procedure reports. By default, registration data is not uploaded. This function is only available to administrative or service personnel.
Button Sound Click	Turn On or Off the clicking sound when a button is pressed.
MTS Display Coldest	Turn On or Off display of the coldest MTS location. If turned On, the coldest MTS is highlighted in blue.
Number MTS	Set the number of MTS locations to display.
Low Cylinder Alert	Select the desired reminder interval (0 minutes to 15 minutes) for the Gas Indicator to display an alert that the estimated remaining volume in the argon cylinder is low.
Displays Needle Temperatures	Turn On or Off the temperature display for CX needles.
Cycles	Select the option to program cycles.

Control Buttons

Table 22 describes the control buttons on the *Settings Screen*.

Table 22. Settings Screen – Control Buttons

Button	Description
ADMIN	Press ADMIN to access the administrative features. Reference the Administrative Options Section for more information. NOTE: The ADMIN options are only available to administrative and service users.
EXPORT LOGS	Press EXPORT LOGS to export the logs to a USB drive.
CLEAR HOSPITAL DATA	Press CLEAR HOSPITAL DATA to clear the hospital name, address, and physician name from the system history file.
CANCEL	Press CANCEL to discard any changes and to exit the screen.
CONFIRM	Press CONFIRM to confirm all setting changes.

Administrative Options

To display the administrative options:

1. Press  and select Settings from the expanded menu to display the *Settings Screen*.
2. Press the **ADMIN** button.

Software Update

Manual Software Update

Administrative and service users can manually update the ICEfx Cryoablation System software via the Boston Scientific supplied USB flash drive.

1. Press the **ADMIN** button on the *Settings Screen*.
2. Press **Software Update**.
3. Select the **USB** option to perform a manual software update.
4. Insert the USB flash drive containing the software update into the USB port on the rear of the Console.
5. Browse to select the update file and press **Update**. A confirmation message appears when the software update is complete.

Manage Users

Select **Manage Users** to add users, remove users, or change the password or access level of any user.

Demo Mode

Demo Mode is for training only.

Select **Demo Mode** and click under the Argon gas gauge to simulate argon pressure. No gas will flow and no ice ball forms. To remove simulated pressure, click again. Caution and Thaw functionality will not heat the needles.

Demo mode is reset to off when power is turned off or when the **Demo Mode** button is clicked again.

Procedure

Performing The Cryoablation Procedure

Table 23 outlines the steps for testing the cryoablation needles and MTS and initiating the cryoablation procedure. This section describes each step in-detail.


Table 23. Cryoablation Procedure Flow

1	Needle/MTS Testing	<ul style="list-style-type: none"> • Press the START CASE button • (Optional) Use the <i>Edit Case Information Screen</i> to enter patient treatment information • Select and prepare sterile needles and MTS for testing • Connect needles to the Needle Connection Panel and lock channels <p>NOTE: Only place needles of the same type in a single channel</p> <ul style="list-style-type: none"> • Connect MTS to the MTS ports, if applicable • Perform Needle Integrity and Functionality Testing
2	Perform Cryoablation Procedure	<ul style="list-style-type: none"> • Insert needles and MTS into target tissue • Initiate freeze activity by pressing the Freeze button to select the freeze intensity, and then pressing the Start button to start freezing. • Using image guidance, continuously monitor iceball formation throughout the procedure. • Initiate thaw activity: <ul style="list-style-type: none"> – For CX needles, press the Thaw button to select the thaw mode, and then press the Start button to begin thawing. – For non-CX needles, needle release from an iceball will only be achievable with passive thaw. • Press the Channel button and select Program Cycles to program freeze-thaw cycles
3	End Procedure	<ul style="list-style-type: none"> • Remove needles and MTS • End Procedure (Press the END CASE button on <i>Procedure Screen</i>) • View and save report if desired <ul style="list-style-type: none"> – Export report to a USB flash drive

Pre-Procedure Needle/MTS Testing

WARNING: Prior to starting a cryoablation procedure, set up the ICEfx Cryoablation System (reference the **System Setup** Section) and then perform Needle Integrity and Functionality Tests. Tests must be successfully completed in order to begin the to proceed with procedure.

1. On the touch screen monitor, press **START CASE**. The Setup wizard displays to demonstrate the setup steps.

Screen 23: Procedure Screen

2. **OPTIONAL:** On the *Edit Case Information* Screen that displays, enter the optional patient treatment information and any notes using the virtual keyboard. Available data entry fields include Patient ID, Hospital Name, Hospital Address, Physician Name, and Organ Type. If you need to enter any other information, additional custom fields can be labeled in the *Settings* Screen (reference the **System Setup** Section).

CAUTION: Select a unique Patient ID that does not reveal the patient's identity to other system users.

3. Using aseptic technique, carefully remove the cryoablation needle from the package and place in a sterile work area.
4. Remove the connector cap, and then connect the needle to the Needle Connection Port.

WARNING: Do not kink, pinch, cut or pull excessively on the needle tubing. Damage to needle handle or tubing may cause the needle to become unusable.

5. After a needle is inserted into the desired channel, lock the channel by sliding the locking bar away from the center of the system (Figure 14).


Figure 14. Locking Needle into Channel

- To allow for easier identification of the needle when multiple cryoablation needles are used during a cryoablation procedure, it is recommended to place a Cryoablation Needle Channel ID Sticker on the needle tubing.

NOTE: Contact Boston Scientific Technical Assistance Center to order Cryoablation Needle Channel ID Stickers for the cryoablation needles.

- Repeat Steps 3 through 6 for each cryoablation needle that will be tested.

CAUTION: Boston Scientific recommends that only needles of the same type be placed together in a single channel. Using needles of differing types in a channel may affect the accuracy of the **Gas Indicator**.

When a channel is locked, the software detects connection of the needle and the channel opens for testing. A dark gray channel button indicates a channel with needles connected.

For CX needles, the ICEfx Cryoablation System automatically detects the type of needle being used, the lot number, and the expiration date of the product. If the system detects that the expiration date has passed, a message appears indicating that the needle is expired and gas flow to that channel is disabled.

For non-CX needles, a Passive Thaw Required message box will appear.

- Press the **Channel** button to open **Advanced Channel Controls**, which allows you to change the needle type for a channel, as needed. Reference the **Advanced Channel Controls** Section for more information on Advanced Channel Controls.

- Prepare to conduct Needle Integrity and Functionality Testing.

WARNING: Avoid contact with the distal portion of the cryoablation needle to maintain sterility during testing.

- Secure the needle tubing to the sterile table prior to beginning the needle testing process.
- Fill a large basin (at least 30 cm in diameter) halfway with sterile water or saline.
- Place the needles, individually or in groups, in the basin such that the full length of the needle shaft is submerged in the sterile water or saline.

- Perform the Needle Integrity and Functionality Test on each needle by pressing the **Test** button on the channel containing the needle(s) and then pressing the **Start** button. The 120-second test automatically conducts a series of flush, freeze, and thaw phases. The durations of these phases are: 50 second flush with low-pressure argon, 15 second freeze with high-pressure argon, and 55 second thaw with low-pressure argon.

OPTIONAL: Alternatively, all needles may be tested simultaneously by pressing the **Start** button on the ALL channel. A message requests confirmation to test all needles. If appropriate, choose YES.

OPTIONAL: If additional testing is required, press the **Test** button again to repeat the test.

NOTE: For CX needles, if a previously tested needle is moved to a new channel during a procedure, the ICEfx Cryoablation System will recognize the needle has successfully completed the Needle Integrity and Functionality Test.

During testing, carefully watch each needle for the following:

Flush:

Ensure no bubbles are formed along the needle shaft and tip.

WARNING: A defective cryoablation needle that has a gas leak can cause a gas embolism in the patient. Never use a defective needle for a cryoablation procedure. Return defective needles to Boston Scientific for evaluation.

Freeze:

Ensure that ice begins to form around the tip of the needle.

WARNING: A needle is defective if there is no ice formation during the Freeze phase. Do not use a defective needle. Obtain a new needle and repeat the testing procedure.

Thaw:

For CX needles, ensure that the iceball detaches from the tip of the needle and that no bubbles escape from the needle tip.

For non-CX needles, ensure that no bubbles escape from the needle tip. A non-CX needle will not actively thaw the needle. Wait for the iceball to passively thaw and detach from the needle tip.

During the Needle Integrity and Functionality Test, the Gas Indicator for the argon cylinder provides an estimate of the time remaining before depletion of the cylinder assuming that all connected needles are operated simultaneously (reference the **Title Bar** Section).

Upon successful completion of the Needle Integrity and Functionality Test, the remaining control buttons on the channel become active. The needle(s) are ready for use.

10. When using MTS, prepare the desired number of sensors for testing. The ICEfx Cryoablation System supports use of two thermal sensors.
 11. Connect each thermal sensor to an MTS connection port by pushing the connector into the port. Ensure that the rectangular tongue on the MTS connector is aligned with the groove at the top of the connection port for proper MTS insertion.
 12. Confirm MTS functionality by verifying that the connected sensor produces a reasonable on-screen temperature reading (e.g., close to room temperature) on the *Procedure Screen*.
-

WARNING:

- Do not use an MTS that has failed the MTS Integrity and Functionality Test; erroneous temperature measurements may be displayed.
-

CAUTION:

- Do not allow the MTS to touch the cryoablation needle during use; erratic temperature measurements may be displayed.
-

Initiating a Cryoablation Procedure

WARNING: Continuously monitor needle insertion, needle positioning, iceball formation and removal using image guidance (such as direct visualization, ultrasound, or Computed Tomography (CT)) to ensure adequate tissue coverage and to avoid damage to adjacent structures.

WARNING: Do not touch the screen if the touch screen monitor goes blank for more than five (5) seconds during a procedure. Immediately turn off power to the system and end the procedure to avoid inadvertent activation of needles.

WARNING: Needle tubing may become extremely cold when conducting freeze cycles during a cryoablation procedure. It is important that a patient's skin is protected from direct contact with needle tubing to avoid the potential for thermal injury to the patient. Ensure an appropriate insulating barrier is placed as needed (such as towels) or other method is employed to prevent needle tubing from touching a patient's skin.

WARNING: If a needle without CX functionality is connected to the ICEfx Cryoablation System, active thawing is not possible. Needle release from an iceball will only be achievable with passive thaw. Use image guidance to ensure adequate thawing.

WARNING: If a needle needs to be added during a procedure, connect the needle to an open channel. If there is no open channel available, reference the **Adding a CX Needle during a Cryoablation Procedure** Section for the instructions for adding a CX needle to a channel with an open port.

1. Position the cryoablation needles and thermal sensors in the target tissue.

CAUTION: During use, avoid damage to the needle from other surgical instruments.

NOTE: Throughout the procedure, monitor the amount of argon remaining on the Gas Indicator (Screen 24). If it is necessary to change the argon cylinder during a procedure, follow the instructions provided in the **Changing Argon Cylinders during a Procedure** Section.

Screen 24. Argon Remaining

2. Press the **Freeze** button on selected channels containing needles to select the freeze intensity, and then press the **Start** button to begin the initial freeze phase of the procedure. Reference the **Initiating a Freeze Cycle** Section for more information on freezing.

To adjust the freeze intensity, press the **Freeze** button and select the desired intensity from the menu. The Freeze phase will continue at the selected freeze intensity until the freeze intensity is changed or the operation is stopped.

OPTIONAL: To initiate a freeze phase on all needles simultaneously, press the **Freeze** button on the channel labeled **ALL** and then press the **Start** button for that channel. Pressing any function button on the channel labeled **ALL** displays a message asking you to confirm simultaneous operation of all needles.

NOTE: Selecting **ALL** will initiate a Freeze phase at the intensity selected for each channel. To conduct freezing in all active channels at the same intensity, press the **Freeze** button for the **ALL** channel and select the freeze intensity to apply to all active channels. After selecting the freeze intensity, press the **Start** button for the **ALL** channel to initiate freezing on all active channels simultaneously.

3. Watch the timer to monitor the freeze phase elapsed time. When the desired freeze duration has elapsed, press the **Stop** button to enter an idle phase.
 4. If thermal sensors are used, monitor the tissue temperature using the **Temperature Sensors** section of the *Procedure Screen* (reference the **MTS** Section).
 5. **For CX needles**, to actively thaw the iceball, press the **Thaw** button on the channels containing needles and select the thaw mode (i-Thaw Function or FastThaw Function), and then press the **Start** button to start the thaw phase. Reference the **Initiating a Thaw Cycle (CX needles only)** Section for more information on thawing.
OPTIONAL: To initiate a thaw phase on all needles simultaneously, press the **Thaw** button on the channel labeled **ALL** and then press **Start**. Pressing any function button on the channel labeled **ALL** displays a message asking you to confirm simultaneous operation of all needles.
-

WARNING: The needle handle and distal needle shaft may become warm during active thawing. Pay attention to the needle handle position. Prolonged contact with warm portions of the needle handle or distal needle shaft could cause unintended thermal tissue damage to the patient or clinician.

6. **For CX needles**, watch the timer to monitor the thaw phase elapsed time (reference the **Initiating a Thaw Cycle (CX needles only)** Section for instructions on performing a timed thaw phase). When the desired thaw duration has elapsed, press the **Stop** button to enter an idle phase.
 7. **For non-CX needles**, active thawing is not possible. Allow time for passive thawing.
 8. Repeat Steps 2 through 6 until the desired number of freeze-thaw cycles has been completed.
 9. *OPTIONAL:* Press the **Channel** button and select **Program Cycles** to program freeze-thaw cycles for an individual channel. Reference the **Program Cycles Control** Section for information on programming freeze-thaw cycles.
-

NOTE: To program cycles for all active channels, press **ALL** and then press **Program Cycles** to set the programming options for all active channels.

WARNING: Ensure adequate thawing or cooling before attempting to remove needles from the patient.

WARNING: Discontinue all needle operation prior to needle removal to minimize risk of thermal injury and/or tissue injury.

10. Remove all needles and MTS from the patient.
11. Unlock the locking bar(s) and remove all needles and MTS from the needle connection panel.
12. See instructions in **Disposal** Section on how to dispose of used needles and MTS.
13. When the procedure is complete, press the **END CASE** button on the *Procedure Screen*. The Gas Venting wizard displays and details the shutdown process.

Follow the on-screen prompts:

- Provide confirmation to end the procedure – Press **Yes, End Case**.
- Request to save a report – Select the appropriate option for saving the information entered on the *Edit Case Information Screen*. Reference the **Reports** Section for more information on reports.
- Option to view the Gas Venting wizard which describes how to depressurize the system. Press **GAS VENTING WIZARD** to view the steps describing how to depressurize the system.

14. Depressurize the system:

- a. Turn the shutoff valve on the argon cylinder clockwise to close the cylinder.
- b. Turn the Vent Valve on the rear of the Console to the VENT position to release the high-pressure gas from the system.
- c. Turn the Vent Valve to the CLOSED position after gas venting completes.

WARNING: If the needles are still connected, do not unlock the channels or disconnect the needles from the needle connection panel while venting is in process.

15. If you are ready to shutdown the system, reference the **Shutting Down the System** Section for the system shutdown procedure.

Adding a CX Needle during a Cryoablation Procedure

If an additional CX needle needs to be added during a cryoablation procedure, add the CX needle to an open channel, lock the needle in the channel, and perform testing.

If there are no open channels, use the following procedure to add a new CX needle.

CAUTION: Boston Scientific recommends that only needles of the same type be placed together in a single channel. Using needles of differing types in a channel may affect the accuracy of the Gas Indicator.

1. Select the channel for the new needle and stop all operations for the existing needle in the channel.
2. Unlock the channel, disconnect the needle from the Console, and set the needle aside.
3. Insert the new needle into the channel and lock the channel.
4. Test the new needle.
5. When the test completes, unlock the channel and reconnect the needle removed from the Console in Step 2.
6. Lock the channel and continue with the cryoablation procedure.

Adding a non-CX Needle during a Cryoablation Procedure

If an additional needle needs to be added during a cryoablation procedure, add the needle to an open channel, lock the needle in the channel, and perform testing.

Changing Argon Cylinders during a Procedure

If it becomes necessary to replace an argon cylinder during a procedure, stop all freezing and thawing operations.

Standard Argon Cylinder Setup

1. Plan the appropriate time to change a cylinder by estimating the argon required to complete the procedure. The **Gas Indicator** indicates how many minutes of argon is remaining based on the selected gas flow intensity, type and number of needles being used. Also, consider the number of planned freeze-thaw cycles for the procedure.
2. Safely position a full argon cylinder, with the required gas type and purity, near the empty cylinder.
3. Close and tighten the argon cylinder valve.
4. Slowly open the Vent Valve to discharge gas from the system and the high-pressure gas supply line. Wait until all the pressure is released and the pressure gauge on the gas supply line indicates zero pressure.
5. Use the appropriate wrench to remove the gauge assembly from the empty argon cylinder.
6. Connect the gauge assembly to the full argon cylinder.
7. Close and tighten the Vent Valve.
8. Carefully turn the cylinder valve on the argon cylinder counterclockwise for one quarter turn. Ensure that the pressure reading on the gauge responds immediately. Turn the cylinder valve further counterclockwise to open the gas cylinder so that sufficient gas flow is present.
If no gas pressure is displayed on the **Gas Indicator**, ensure the Argon Valve is OPEN.
9. Continue the cryoablation procedure at the next planned freeze or thaw phase.

Dual Argon Cylinder Connection

1. Safely position a full argon cylinder, with the required purity, near the empty cylinder.
2. Close and tighten the cylinder valve of the empty argon cylinder.
3. Open the Vent Valve to discharge gas from the system and the high-pressure gas supply line. Wait until all the pressure is released and the pressure gauge indicates zero pressure.
4. Close the Vent Valve on the Console.
5. Connect the auxiliary gas supply line to the EZ-Connect2 Dual Cylinder Adapter using the quick-connect connectors.
6. Connect the opposite end of the auxiliary gas supply line to the new argon cylinder.
7. Carefully turn the cylinder valve on the new argon cylinder counterclockwise for one quarter turn. Ensure that the pressure reading on the gauge responds immediately. Turn the cylinder valve further counterclockwise to open the argon cylinder so that sufficient gas flow is present.

SYSTEM SHUTDOWN

Shutting Down the System

WARNING: Warn procedure personnel prior to venting the ICEfx Cryoablation System to avoid startling them.

1. If the Console was not vented at the end of the procedure, perform these steps to vent the system.
 - a. Turn the shutoff valve on the argon cylinder clockwise to close the cylinder.
 - b. Turn the Vent Valve on the Console to the VENT position to release the high-pressure gas from the system.


Figure 15. Vent Valve – Fully Open

- c. Turn the Vent Valve to the CLOSED position after venting completes.


Figure 16. (same as Figure 8) Vent Valve - Fully Closed

2. Disconnect the high-pressure gas supply line from the Console and from the argon cylinder. Store the gas supply line and gauge assembly in the accessory bag.

WARNING: If it is difficult to loosen the pressure gauge connected to the cylinder, or the argon gas supply line cannot be disconnected from the argon inlet connection, do not apply excessive force to release the gas supply line or to loosen the pressure gauge. The gas line may still be under pressure.

3. Cover the argon inlet with the humidity plug.
4. Press  and select **Logout** from the expanded menu to log out of the system.
5. Press the **Shutdown** button on the *Login Screen* to power down the system. A message appears requesting confirmation to shut the system down.

6. Wait until the screen goes black. Turn the Power Switch on the rear of the Console to the OFF position.
7. Unplug the Console.

WARNING: Do not pull on the power cord. Grasp the plug, not the power cord, to disconnect the device from the wall socket.

8. Clean the Console after each use using the instructions provided in the **Cleaning the ICEfx Cryoablation System** Section. Ensure the system is dry prior to storage.
9. Lower the Touch Screen Monitor into the monitor storage basin prior to system storage.

CAUTION: Prior to lowering the monitor, ensure that no objects are located in the monitor storage basin. Use care when lowering the monitor into the monitor storage basin; do not use excessive force to avoid damaging the monitor.

CAUTION: Use care when lowering the touch screen monitor to avoid the potential for pinching fingers.

Post-Procedure

Any serious incident that occurs in relation to this device should be reported to the manufacturer and to the relevant local regulatory authority.

For customers in Australia, report any serious incident that occurs in relation to this device to Boston Scientific and to the Therapeutic Goods Administration (<https://www.tga.gov.au>).

Cleaning the ICEfx Cryoablation System

Clean the ICEfx Cryoablation System after each use by following the steps below.

1. Clean the Touch Screen Monitor when the ICEfx Cryoablation Console is powered OFF.
 - Wipe the screen gently with a damp gauze wipe.
 - Use water or isopropyl alcohol cleaning solutions.
 - Do not use cleaning agents such as Betadine Antiseptic Solution or bleach solution.
2. Clean the Console by wiping with a damp gauze wipe.
 - Use soap and water or isopropyl alcohol cleaning solutions.
 - Do not use cleaning agents such as Betadine Antiseptic Solution or bleach solution.
 - Do not allow water or other liquid to drip or seep inside the needle connection ports. The needle connection ports must remain completely dry at all times. Avoid system exposure to liquids.
3. Ensure that cleaned surfaces are dry prior to closing or turning on the system.

Disposal

All external and accessible surfaces of this device should be cleaned per the Cleaning the ICEfx Cryoablation System instructions included in the user manual. Include any common detachable cables (power cord, video cables, patch cords, etc.). Review the user manual to identify if any hazardous materials are present.

If placing the unit into an electronics recycling stream, notify the receiver of the presence of such materials. Use of recycling service suppliers familiar with medical electrical equipment is recommended, but not required. Do not dispose of by incineration, burial or placement into municipal waste stream.

Device should be safely disposed of in accordance with hospital, administrative, and/or local government policy or returned to Boston Scientific. For a Returned Product Kit, contact Boston Scientific Technical Assistance Center.

Dispose of all sharps directly into a sharps disposal container labeled with a biological hazard symbol. Sharps waste should be safely disposed of using available sharps waste channels in accordance with hospital, administrative, and/or local government policy.

TROUBLESHOOTING

Boston Scientific suggests the following options for troubleshooting the ICEfx Cryoablation System. If the suggested approaches do not resolve the issue or you experience an issue not referenced below, contact Boston Scientific Technical Assistance Center.

Software Recovery

In the event of software corruption or failure, the software may be restored to the previous software version.

Administrative and service users may update the software with an appropriate USB flash drive.

1. Shutdown the Console.
2. Hold down the **Software Reset** button by inserting a straightened paper clip into the Software Reset hole while simultaneously powering up the system. The system displays the *Software Recovery Screen*.
3. Press the **Rollback** button to restore the software to the previous software version.
4. *OPTIONAL*: Press the **Load** button to update the software if a message displayed on the *Login Screen* indicating that the software configuration is invalid.
5. If updating the software to a newer version available on a USB flash drive.
 - Login as an administrative user.
 - Select **Settings** from the Options Menu.
 - Press the **ADMIN** button on the *Settings Screen*.
 - Press the **Software Update** button.
 - Select the manual software update option.
 - Insert the USB flash drive.

NOTE: Wait 20 seconds for the system to recognize the flash drive.

- Press the **Browse** button.
 - Select the file to conduct the upgrade.
 - Press the **Update** button.
-

NOTE:



- Wait until a displayed message confirms completion of the update.
 - The upgrade may take half an hour to complete.
-


Electronics, Electrical and User Error Related Issues

Symptom	Potential Causes / Solutions
System does not turn ON (i.e., the fan is not running) or power is lost during procedure	<ol style="list-style-type: none"> 1. The power switch on the rear panel is OFF (Figure 2). Turn the power ON. 2. The power cable to the ICEfx Cryoablation Console is disconnected from the power outlet or from the rear panel of the Console. Connect the power cable to the Console, checking that the power cable is fully seated. Connect the power cable to the power outlet. 3. There is no power at the wall outlet. Ensure the power outlet is switched ON. Call the hospital Biomedical Engineer for assistance, if needed. 4. A fuse may have blown. Reference Replacing Fuses Section for information on how to replace fuses in the system
Channel or needle is not recognized by the system and not available for use	<ol style="list-style-type: none"> 1. Check the corresponding Channel Locking Bar and ensure it is fully in the locked position. 2. At least one needle must be inserted into a channel for the channel to be available for use. 3. If steps 1 and 2 are unsuccessful, select a new needle for use. Contact Boston Scientific Technical Assistance Center to arrange return of the non-operational needle. 4. The channel may be defective. Do not use this channel. Move the needle(s) to a different channel. Conduct the Needle Integrity and Functionality Test again.
USB flash drive is not functioning Or The USB flash drive is not recognized by the system	<ol style="list-style-type: none"> 1. The USB flash drive is not properly connected to the USB port. Remove the USB flash drive from the USB port on the Console. Wait a few seconds and reconnect the USB flash drive to the USB port. 2. If the problem persists, try using a different Boston Scientific supplied USB flash drive. 3. The USB flash drive is faulty. Replace the USB flash drive with a new Boston Scientific supplied USB flash drive.
The <i>Login Screen</i> appeared after the system was left idle for more than 2 hours while on the <i>Procedure Screen</i>	Enter the appropriate password to return to the <i>Procedure Screen</i> .
Touch screen does not respond	Turn OFF and restart the system using the power switch on the rear of the system (Figure 2).
Touch Screen goes blank during a procedure	<ol style="list-style-type: none"> 1. The video cable may be disconnected. WARNING: Do not touch the screen if the touch screen monitor goes blank for greater than five (5) seconds during a procedure. Immediately turn off power to the system and end the procedure to avoid inadvertent activation of needles. 2. Contact Boston Scientific Technical Assistance Center.

Replacing Fuses

Instruction	Photo
Contact Boston Scientific Technical Assistance Center to obtain the replacement fuses required for this procedure.	
Turn the Power Switch OFF on the rear of the ICEfx Cryoablation Console.	
Pull on the red cable retention mechanism to disconnect the power cable from the power inlet.	

Instruction	Photo
<p>Insert a small screwdriver into the slot on the top of the fuse holder and flip down the cover of the fuse holder.</p>	
<p>Insert a small screwdriver underneath the edge of the fuse holder to remove the fuse holder from the power outlet.</p>	
<p>While holding your hand under the fuse holder, carefully slide the fuse holder out of the power inlet. NOTE: There are two fuses in the fuse holder.</p>	

Instruction	Photo
<p>Catch the fuse holder as it is removed from the power inlet. The two fuses that are in the fuse holder are the fuses that are part of the system circuit.</p>	
<p>Replace the fuses in the fuse holder with the two replacement fuses supplied from Boston Scientific Technical Assistance Center. NOTE: Use only Boston Scientific specified fuses in the ICEfx Cryoablation Console.</p>	
<p>Plug in the power cord.</p>	
<p>Call Boston Scientific Technical Assistance Center to schedule a service call to determine the cause of the blown fuses, to determine if service is required, and to order replacement fuses.</p>	

Gas Issues

Symptom	Potential Causes / Solutions
The ICEfx Cryoablation System does not allow testing of a needle in a locked channel	<ol style="list-style-type: none"> 1. Verify that the system is at operating pressure. Reference Table 7 for information on working gas pressures 2. The Argon Shutoff Valve may be in the OFF position. Verify that the Argon Shutoff Valve (Figure 2) is in the ON position to allow sufficient gas flow. 3. If necessary, further open the argon cylinder valve by turning the valve counterclockwise to allow sufficient gas flow. Verify that appropriate pressure is displayed on the Gas Indicator.
Needle is not freezing during the Needle Integrity and Functionality Test	<ol style="list-style-type: none"> 1. The argon cylinder valve may be closed. Open the gas cylinder valve by turning the gas cylinder valve counterclockwise to allow sufficient gas flow. Verify that appropriate pressure is displayed on the Gas Indicator. 2. Confirm that the argon cylinder is connected to the argon inlet. 3. The needle may be clogged (by dust or ice). Try to test it again. 4. If the needle still does not freeze, <ul style="list-style-type: none"> • Press the Stop button to stop all activity on the channel. • Hold the needle connector firmly with one hand and unlock the channel to disconnect the needle. • Move the needle to a different channel and run the test again. <p>NOTE: If only one needle is connected to the channel, there may be some residual pressure behind the needle connector.</p> <ol style="list-style-type: none"> 5. If the problem persists, replace the needle with a new needle and contact Boston Scientific Technical Assistance Center.
The Vent Valve is leaking gas	The Vent Valve may be open. Fully close the Vent Valve.
The pressure displayed on the Gas Indicator indicates the gas pressure is too low (Table 9)	<ol style="list-style-type: none"> 1. Verify the Argon Shutoff Valve is open. 2. Verify that the argon cylinder valve is sufficiently opened to allow gas flow. If necessary, open the valve approximately another half turn. 3. Verify that the cylinder contains sufficient pressure using the gauge on the cylinder. 4. Replace the cylinder, if needed.
It is difficult to loosen the pressure gauge connected to the cylinder, or the argon gas supply line cannot be disconnected from the argon inlet connection	<p>The gas line was not vented and the gas line is still under pressure.</p> <p>WARNING: If it is difficult to loosen the pressure gauge connected to the cylinder, or the argon gas supply line cannot be disconnected from the argon inlet connection, do not apply excessive force to release the gas supply line or to loosen the pressure gauge. The gas line may still be under pressure.</p> <ol style="list-style-type: none"> 1. Verify that the argon cylinder is CLOSED. 2. Verify the pressure gauge at the gas cylinder is 0. 3. If in the <i>Procedure Screen</i>, verify the gas pressure display shows no gas connected. 4. If still unable to disconnect the gas supply line, or if the system is OFF, open the Vent Valve on the rear of the system to completely vent the system. 5. Upon completion, close the Vent Valve.
Gas begins leaking through a needle port after gas flow was started using the Test , Freeze or Thaw buttons	<p>The channel may have a loose or broken socket.</p> <ol style="list-style-type: none"> 1. Disconnect the needle and move the needle to a different channel. 2. Conduct the Needle Integrity and Functionality Test again on the needle. 3. Contact Boston Scientific Technical Assistance Center.

Mechanical Issues

Symptom	Potential Causes / Solutions
Needle cannot be firmly connected in the needle port	<ol style="list-style-type: none"> 1. Verify that the locking bar is in the UNLOCKED position. 2. Needle connector may be defective. Use a different needle. 3. Residual gas pressure may be in the needle port. Use a different channel. 4. Check the Gas Indicator display. If there is pressure in the system, end the procedure and open the Vent Valve on the rear of the system to completely vent the system. 5. Upon completion, close the Vent Valve.
Locking bar on the Needle Connection Panel cannot be positioned in the LOCKED position	<ol style="list-style-type: none"> 1. Verify that all needles in the channel are fully inserted in the needle connection ports. 2. The Locking Bar may be defective. Move the needle to a different channel. Contact Boston Scientific Technical Assistance Center to schedule a service call. 3. Check the Gas Indicator display and verify the system is not pressurized. If there is pressure in the system, end the procedure and open the Vent Valve on the rear of the system to completely vent the system. Upon completion, close the Vent Valve.
Cart does not move freely	Release the wheel brake levers.
The Vent Valve cannot be used to depressurize the system.	<p>When viewing the Gas Venting Wizard, if a dialog displays indicating the pressure is not decreasing and the argon cylinder is closed, the Vent Valve may be plugged.</p> <p>If the Vent Valve is plugged:</p> <ol style="list-style-type: none"> 1. Press CANCEL. 2. Press START CASE. 3. Select SKIP SETUP WIZARD. 4. Perform a freeze on one or more channels. 5. Monitor the pressure on the Gas Indicator. 6. When the pressure goes below 3.45 bar (50 psi), which could take more than 10 minutes, stop the freeze and then press END CASE. 7. Follow the instructions presented on the wizard steps.

Gas Cylinder and Gas Supply Line

Symptom	Potential Causes / Solutions
Safety cable is missing at either the cylinder side or the Console side of the gas supply line	Do not use a gas supply line with a missing safety cable. Doing so could compromise the safety of personnel in the room. Contact Boston Scientific Technical Assistance Center for further instructions.
A gauge or a gas supply line is damaged	Do not use any damaged product. Contact Boston Scientific Technical Assistance Center for new accessories.
A gas leak was detected between the gauge adapter and the cylinder valve	<ul style="list-style-type: none"> • Tighten the connection using the wrench provided with the ICEfx Cryoablation System. • Close the cylinder valve and vent the gas from the ICEfx Cryoablation System and the gas supply line using the Vent Valve (Figure 3). Verify the system is depressurized. Loosen and remove the gauge assembly adapter. Verify no debris is on the gas cylinder connection point; clean the sealing surface as needed to remove any debris. Reposition and tighten the gauge assembly adapter onto the cylinder valve, using the wrench provided with the ICEfx Cryoablation System.

Needles

Symptom	Potential Causes / Solutions
During the freeze phase or following a thaw phase, no iceball or small iceball forms on the needle(s) in a given channel	<p>For CX needles:</p> <ol style="list-style-type: none"> 1. Perform the following steps in the order in which they are described: <ul style="list-style-type: none"> • Stop freezing/thawing actions for all channels. • Thaw the problematic needle(s) for at least one minute. • Freeze the needle(s) to verify proper function. 2. If the problem persists, connect a new needle into a different channel and test the needle. Continue the procedure with the newly tested needle. 3. Contact Boston Scientific Technical Assistance Center for further instructions upon completion of the procedure. <p>For non-CX needles:</p> <ol style="list-style-type: none"> 1. Perform the following steps in the order in which they are described: <ul style="list-style-type: none"> • Stop freezing/thawing actions for all channels. • Wait for any formed ice to passively thaw. Do not use the needle. 2. Use a new needle to continue the procedure. 3. Contact Boston Scientific Technical Assistance Center for further instructions upon completion of the procedure.
During the Needle Integrity and Functionality test bubbles are seen escaping from the needle	<p>WARNING: Do not use the needle.</p> <ul style="list-style-type: none"> • Disconnect and segregate the needle. • Return the needle to Boston Scientific for evaluation. • Use a new needle to continue the procedure. • Test the new needle to confirm the integrity and functionality of the new needle. • Contact Boston Scientific Technical Assistance Center for further instructions upon completion of the procedure.
Needle is bent or damaged while attempting to unpack or use it	<p>WARNING: Do not use the needle.</p> <ul style="list-style-type: none"> • Segregate the needle. • Use a different needle to continue the procedure. • Contact Boston Scientific Technical Assistance Center for further instructions upon completion of the procedure.

CX Needles

Symptom	Potential Causes / Solutions
A needle is not recognized by the system.	<ol style="list-style-type: none"> 1. Move the needle to a different port. 2. If the system still does not recognize the needle, press the Channel button to display the needle type drop-down list. Select the needle type from the list.

MTS

Symptom	Potential Causes / Solutions
Not all of the expected temperature sensor locations on the MTS needle are registered on the <i>Procedure Screen</i> .	<ol style="list-style-type: none"> 1. On the <i>Settings Screen</i>, select the desired number of sensors to display. 2. Move the MTS needle to a different port. 3. Replace the MTS needle. 4. Contact Boston Scientific Technical Assistance Center.
MTS needle is connected but temperature is not registered on the <i>Procedure Screen</i> .	<ol style="list-style-type: none"> 1. Move the MTS needle to a different port. 2. Replace the MTS needle. 3. Contact Boston Scientific Technical Assistance Center.

Displayed Messages

The ICEfx Cryoablation System displays a message on the user interface when a user requests assistance or when user, needle or system errors are detected.

NOTE: Record and report the message number (e.g., 10-01, 80-02) if assistance is required from Boston Scientific Technical Assistance Center.

Login Messages

Message	Reason for Appearance / Solutions
10-01 Login <ul style="list-style-type: none"> You have not entered the correct Login Name Reenter your Login Name Contact your System Administrator if assistance is required Contact Boston Scientific Technical Assistance Center if further assistance is required 	Reason for Appearance: <ul style="list-style-type: none"> No name was entered. OR <ul style="list-style-type: none"> The name entered did not match assigned names in the system. Solutions: <ol style="list-style-type: none"> Reenter your login credentials. Contact the System Administrator for further assistance. Contact Boston Scientific Technical Assistance Center.
10-02 Login <ul style="list-style-type: none"> You have not entered the correct Password Reenter your Password Contact your System Administrator if assistance is required Contact Boston Scientific Technical Assistance Center if further assistance is required 	Reason for Appearance: <ul style="list-style-type: none"> No password was entered. OR <ul style="list-style-type: none"> The password entered did not match the password associated with the login name. Solutions: <ol style="list-style-type: none"> Reenter the Password. Contact the System Administrator for further assistance. Contact Boston Scientific Technical Assistance Center.
10-03 Reset Password Challenge <ul style="list-style-type: none"> To reset your password, contact Boston Scientific Technical Assistance Center Relay the Challenge on the screen below Enter the Response provided by Customer Service Press the Reset button 	Reason for Appearance: <p>User forgot their password, pressed the Forgot Password button and received a challenge to be relayed to Customer Service.</p> Solutions: <ol style="list-style-type: none"> Contact Boston Scientific Technical Assistance Center and relay the Challenge displayed on the screen. Enter the Response provided by Customer Service. Press Reset button.
10-04 Password Reset <ul style="list-style-type: none"> Enter your new password below 	Reason for Appearance: <p>User correctly entered the reset password challenge and should now set the new password</p> Solutions: <ol style="list-style-type: none"> Enter the new password in the New Password field. Enter the same password in the Confirm Password field.

Message	Reason for Appearance / Solutions
10-05 Emergency Login <ul style="list-style-type: none"> To obtain an emergency login, contact Boston Scientific Technical Assistance Center Relay the Challenge on the screen below Enter the Response provided by Customer Service Press the Login button 	Reason for Appearance: User requested an emergency login and received a challenge to be relayed to Customer Service Solutions: <ol style="list-style-type: none"> Contact Boston Scientific Technical Assistance Center and relay the Challenge displayed on the screen. Enter the Response provided by Customer Service. Press the Login button.
10-06 Incorrect Response <ul style="list-style-type: none"> You have not entered the correct response Contact Boston Scientific Technical Assistance Center for a Response to the on-screen Challenge 	Reason for Appearance: User tried to access the Emergency Login and did not respond correctly to the challenge. Solutions: <ol style="list-style-type: none"> Contact Boston Scientific Technical Assistance Center and relay the Challenge displayed on the screen. Enter the Response provided by Customer Service. Press the Login button. NOTE: This action does not reset a password.
10-07 Incorrect Response <ul style="list-style-type: none"> You have not entered the correct response Contact Boston Scientific Technical Assistance Center for a Response to the on-screen Challenge 	Reason for Appearance: User tried to reset their password and did not respond correctly to the challenge Solutions: <ol style="list-style-type: none"> Contact Boston Scientific Technical Assistance Center and relay the Challenge displayed on the screen. Enter the Response provided by Customer Service. Press the Login button.

Procedure Messages

Message	Reason for Appearance / Solutions
20-01 Cannot Start Test <ul style="list-style-type: none"> The gas pressure is too low/high to begin a procedure Check that the gas cylinders have appropriate pressure to start the procedure 	Reason for Appearance: User pressed the Test button when either the gas cylinder was not yet connected or the gas pressure was below operating pressure (reference Table 7). Solutions: <ol style="list-style-type: none"> Verify that the gas cylinder is connected (reference the Connecting the Argon Gas Section). Check that the gas cylinder valve is sufficiently opened to allow gas flow (reference the Connecting the Argon Gas Section). If necessary, open the valve approximately another half turn. Verify that the gas cylinder contains sufficient pressure using the gauge on the cylinder. Replace the gas cylinder, if needed.
20-02 Argon Shut Off Valve <ul style="list-style-type: none"> The argon shut off valve may be closed Check and open if necessary 	Reason for Appearance: At system start-up, the system detected gas was connected, but no gas is entering the system. The argon shut off valve may be closed. Solutions: Verify that the argon shutoff valve is in the ON position.

Message	Reason for Appearance / Solutions
<p>20-03 Test All</p> <ul style="list-style-type: none"> You selected to control all needles simultaneously Do you want to initiate test in all needles now? <p style="text-align: center;">YES NO</p>	<p>Reason for Appearance: User selected Test button in ALL channel to control all active channels.</p> <p>Solutions: Confirm testing of all connected needles.</p>
<p>20-04 Freeze All Channels</p> <ul style="list-style-type: none"> You selected to control all needles simultaneously Do you want to initiate [x%] freeze for all active needles now? <p style="text-align: center;">YES NO</p>	<p>Reason for Appearance: User selected Freeze button in ALL channel to control all active channels.</p> <p>Solutions: Confirm initiation of freeze phase on all active needles.</p>
<p>20-06 Thaw All Channels</p> <ul style="list-style-type: none"> The argon shut off valve may be closed Check and open if necessary 	<p>Reason for Appearance: User selected Thaw button in ALL channel to control all active channels.</p> <p>Solutions: Confirm initiation of thawing on all active needles.</p>
<p>20-07 Stop All Channels</p> <ul style="list-style-type: none"> You selected to control all needles simultaneously Do you want to stop activity in all needles now? <p style="text-align: center;">YES NO</p>	<p>Reason for Appearance: User selected Stop button in ALL channel to control all active channels.</p> <p>Solutions: Select whether to stop activity for all active needles.</p>
<p>20-18 System Shutdown</p> <ul style="list-style-type: none"> You selected to control all needles simultaneously Do you want to initiate [x%] freeze for all active needles now? <p style="text-align: center;">YES NO</p>	<p>Reason for Appearance: User selected Shutdown in the <i>Login Screen</i> to power down the system.</p> <p>Solutions: Confirm shutdown of the system. Reference the Shutting Down the System Section for the shutdown procedure.</p>
<p>20-19 Procedure Timeout</p> <ul style="list-style-type: none"> The procedure has exceeded the allowable time Procedure will be terminated <p style="text-align: center;">YES NO</p>	<p>Reason for Appearance: The procedure exceeded the allowable 8 hour duration.</p> <p>Solutions:</p> <ul style="list-style-type: none"> If the message appeared because the system was left unattended, shutdown the system. <p>OR</p> <ul style="list-style-type: none"> If more time is required for the procedure, enter the <i>Procedure Screen</i> again and continue with the procedure.

Gas Messages

Message	Reason for Appearance / Solutions
30-02 Vent System <ul style="list-style-type: none"> The system was not vented at the end of the last case Vent the system 	Reason for Appearance: The user did not vent the system completely before the last shutdown. Solutions: To clear the remaining gas pressure, manually vent the gas using the Vent Valve.
30-03 Low Gas Level <ul style="list-style-type: none"> Low level of argon gas remains in the cylinder Replace the gas cylinder as soon as feasible with a new cylinder 	Reason for Appearance: Remaining gas in the cylinder was calculated to last only X more minutes. Where X represents the Low Cylinder Alert set on the <i>Settings Screen</i> (reference the Configuring Settings Section). Solutions: Replace the gas cylinder as soon as possible..

Needle Messages (CX Needles)

Message	Reason for Appearance / Solutions
40-01 Unsupported Needle <ul style="list-style-type: none"> The system was not vented at the end of the last case Vent the system 	Reason for Appearance: A needle type not supported by the software configuration was connected in the channel. The channel will be disabled until an appropriate needle is connected. OR Software identified needles not regulatory cleared or supported in the country of operation. The channel will be disabled until an appropriate needle is connected. Solutions: <ol style="list-style-type: none"> 1. Disconnect and segregate the needle. 2. Connect a supported needle type and select the needle type from the drop down menu. 3. Contact Boston Scientific Technical Assistance Center if further guidance is needed.
40-03 Recalled Needle <ul style="list-style-type: none"> This needle lot number is identified as part of a recall and is unavailable for use Return the needle to Boston Scientific Connect a needle from a different lot number to continue the procedure 	Reason for Appearance: The needle lot number was identified as a lot number in a regulatory recall. The channel will be disabled until a needle from a different lot number is connected. Solutions: <ol style="list-style-type: none"> 1. Disconnect and segregate the needle. 2. Return the needle to Boston Scientific. 3. Connect a needle from a different lot number to continue the procedure

Message	Reason for Appearance / Solutions
<p>40-04 Expired Needle</p> <ul style="list-style-type: none"> The needle 'Use by' date is in the past Replace with a new needle 	<p>Reason for Appearance: The needle expiration date was identified as an expired date. The channel will be disabled until an appropriate needle is connected.</p> <p>Solutions:</p> <ol style="list-style-type: none"> Disconnect and segregate the needle. Replace with a new needle.
<p>40-05 Used Needle</p> <ul style="list-style-type: none"> This needle has been previously used Replace with a new needle 	<p>Reason for Appearance: The software prevents needle reuse and identified a needle previously used. The channel will be disabled until a new needle is connected.</p> <p>Solutions:</p> <ol style="list-style-type: none"> Disconnect and segregate the needle. Replace with a new needle.
<p>40-06 Corrupt Memory</p> <ul style="list-style-type: none"> Needle Memory chip is corrupt on Channel X Manually select the needle type 	<p>Reason for Appearance: User connected a needle containing a bad memory chip.</p> <p>Solutions: Manually select the needle type using the drop-down menu.</p>
<p>40-07 Passive Thaw Required</p> <ul style="list-style-type: none"> The currently selected needle on channel [x] does not have i-Thaw Function capability Passive thaw is required on this channel 	<p>Reason for Appearance:</p> <ul style="list-style-type: none"> Combinations of needles that are i-Thaw Function capable and that are not i-Thaw Function capable were connected to the needle connection panel. <p>OR</p> <ul style="list-style-type: none"> Passive thaw is required. <p>Solutions:</p> <ul style="list-style-type: none"> Wait for ice to passively thaw.
<p>40-08 i-Thaw Error</p> <ul style="list-style-type: none"> Channel [Xa or Xb] i-Thaw Function needle is defective for electrical thawing Use passive thaw 	<p>Reason for Appearance: An i-Thaw Function needle was identified as not i-Thaw Function operational.</p> <p>Solutions:</p> <ul style="list-style-type: none"> Use passive thawing.

Software Messages

Message	Reason for Appearance / Solutions
60-01 Update Failed <ul style="list-style-type: none"> An error occurred during the update process Retry update 	Reason for Appearance: During user update of the system, an error occurred that prevented completion of the update. Solutions: <ol style="list-style-type: none"> Retry the update. Contact Boston Scientific Technical Assistance Center if the update fails.
60-02 Update Failed <ul style="list-style-type: none"> An error occurred during the upload process Retry upload 	Reason for Appearance: During log data upload of the system, an error occurred that prevented completion of the upload. Solutions: <ol style="list-style-type: none"> Retry the update. Contact Boston Scientific Technical Assistance Center if the update fails.
60-03 Incompatible Software <ul style="list-style-type: none"> The ICEfx Cryoablation System software is not compatible with regulatory approvals If problem persists, contact Boston Scientific Technical Assistance Center 	Reason for Appearance: The software was checked against approved software versions in the regulatory files for each respective market. An incompatibility with regulatory approval was identified. Solutions: <ol style="list-style-type: none"> Contact Boston Scientific Technical Assistance Center if further assistance is required.
60-04 Invalid Configuration <ul style="list-style-type: none"> A software file is nonoperational If problem persists, contact Boston Scientific Technical Assistance Center 	Reason for Appearance: There was a problem with the software's configuration files. Solutions: Contact Boston Scientific Technical Assistance Center.
60-05 Software Recovery <ul style="list-style-type: none"> An error occurred during the upload process Retry upload 	Reason for Appearance: User pressed the Software Recovery button then selected Rollback. Activation will rollback the software to the previous software version. Solutions Confirm software rollback. Reference the Software Recovery Section for information on software recovery.
60-06 Software Recovery <ul style="list-style-type: none"> Are you sure you want to restore the default settings for all configurations? <p style="text-align: center;">YES NO</p>	Reason for Appearance: User pressed the Software Recovery button then selected Load. Activation will restore the system default settings for all configurations. Solutions: Confirm restoration of default settings.
60-07 Incompatible Hardware <ul style="list-style-type: none"> The hardware is not compatible with the current software Try downloading a software update from the USB flash drive 	Reason for Appearance: When the system performed self-checks at start up, an incompatibility between the hardware and software was detected. Solutions: <ul style="list-style-type: none"> Manually update the software using the provided USB flash drive.
60-08 Unable to export logs to the USB flash drive <ul style="list-style-type: none"> The USB flash drive may be disconnected or full. 	Reason for Appearance: During log export, an error occurred that prevented completion of the export. Solutions: <ol style="list-style-type: none"> Verify that the USB flash drive is connected and has enough space. Contact Boston Scientific Technical Assistance Center if further assistance is required.

Reports Messages

Message	Reason for Appearance / Solutions
70-03 Report Error <ul style="list-style-type: none"> Errors occurred while compiling the report The report may be incomplete 	Reason for Appearance: User chose to access a report during a procedure or to save data to a report at the end of a procedure. Errors occurred which may affect the completeness of the report. Solutions: Contact Boston Scientific Technical Assistance Center.
70-04 Report Saved <ul style="list-style-type: none"> Report saved successfully 	Reason for Appearance: Report successfully saved to a USB flash drive. Solutions: No user action is required.
70-05 Duplicate Filename <ul style="list-style-type: none"> The chosen filename already exists on the USB flash drive Choose a different filename 	Reason for Appearance: User tried to export a report to a USB flash drive using a filename duplicate to a filename on the flash drive. Solutions: To export the report, choose a different filename.
70-06 Report Error <ul style="list-style-type: none"> Unable to export the report to the USB flash drive The USB flash drive may be disconnected or full 	Reason for Appearance: User selected Save Reports to Flash Drive. No flash drive was detected or there was insufficient space on the flash drive. Solutions: <ol style="list-style-type: none"> Remove the USB flash drive from the USB port on the Console. Wait a few seconds and reconnect the USB flash drive to the USB port. If the problem persists, contact Boston Scientific Technical Assistance Center for a replacement USB flash drive.

System Messages

NOTE: System Error messages are displayed in the center of the Title Bar. Press on the error message to view further details.

Message	Reason for Appearance / Solutions
80-01 Communication Failure <ul style="list-style-type: none"> Internal communication failed Reconnection attempt failed Restarting ICEfx Cryoablation System If problem persists, contact Boston Scientific Technical Assistance Center 	Reason for Appearance: Software was unable to connect with hardware after an attempt to reinitialize communication. If the restart fails, the system is not usable. Solutions: <ol style="list-style-type: none"> Restart the ICEfx Cryoablation Console. Contact Boston Scientific Technical Assistance Center if the problem persists.
80-02 Startup Failure <ul style="list-style-type: none"> System self-checks failed Restart the system If problem persists, contact Boston Scientific Technical Assistance Center 	Reason for Appearance: Software self-checks located a failure that required the system to be restarted. Solutions: <ol style="list-style-type: none"> Restart the ICEfx Cryoablation Console. Contact Boston Scientific Technical Assistance Center if the problem persists.
80-03 Pressure Alert <ul style="list-style-type: none"> Pressure exceeds safe operating limits Close the gas cylinders and vent gas from the system Procedure will be terminated 	Reason for Appearance: System detected that the internal pressure was above safe limits. The system will terminate the procedure. The user must vent gas from the system. Solutions: <ol style="list-style-type: none"> Close the argon cylinder. Vent the system using the Vent Valve.
80-04 Temperature Warning <ul style="list-style-type: none"> The internal temperature of the ICEfx Cryoablation System exceeds appropriate operating limits Discontinue the cryoablation procedure as soon as safe to do so Contact Boston Scientific Technical Assistance Center 	Reason for Appearance: The internal temperature of the system exceeded appropriate operating limits. Solutions: <ol style="list-style-type: none"> End the cryoablation procedure as soon as it is safe to do so. Contact Boston Scientific Technical Assistance Center.
80-05 Service Due <ul style="list-style-type: none"> Low battery detected Contact Boston Scientific Technical Assistance Center to schedule service 	Reason for Appearance: System detected a low battery. System operation may be impacted if operated with a low battery. Solutions: Contact Boston Scientific Technical Assistance Center.
80-30 System Error <ul style="list-style-type: none"> Gas pressure display may be inaccurate. Monitor procedure carefully using image guidance. Use passive thaw. At the conclusion of the procedure, Boston Scientific Technical Assistance Center. 	Reason for Appearance: Internal pressure checks were inconsistent and may produce an inaccurate display on the Pressure Gauge. Solutions: <ol style="list-style-type: none"> Continuously monitor the procedure using image guidance. After completing the procedure, contact Boston Scientific Technical Assistance Center.

Message	Reason for Appearance / Solutions
<p>80-31 System Error</p> <ul style="list-style-type: none"> Gas cylinder shut off valve is not open enough to provide adequate flow. If necessary, open the valve approximately another half turn. 	<p>Reason for Appearance: Gas flow from the argon cylinder was insufficient. To improve gas flow, the cylinder valve should be opened further.</p> <p>Solutions: Check that the argon cylinder valve is sufficiently opened to allow gas flow (reference the Reports Section). If necessary, open the valve approximately another half turn.</p>
<p>80-32 System Error</p> <ul style="list-style-type: none"> Channel X is defective. Choose another channel. At the conclusion of the procedure, contact Boston Scientific Technical Assistance Center. 	<p>Reason for Appearance: A defective solenoid on Channel X was detected; another channel should be selected.</p> <p>Solutions:</p> <ol style="list-style-type: none"> Select a different channel. After completing the procedure, contact Boston Scientific Technical Assistance Center.
<p>80-33 System Error</p> <ul style="list-style-type: none"> Defective MTS in channel [X]. Replace with a new MTS. Contact Boston Scientific Technical Assistance Center 	<p>Reason for Appearance: The MTS temperature display was out of the expected range when initially connected to the needle connection panel.</p> <p>Solutions:</p> <ol style="list-style-type: none"> Disconnect and segregate the MTS Connect a new MTS
<p>80-34 System Error</p> <ul style="list-style-type: none"> Fan X failure. At the conclusion of the procedure, contact Boston Scientific Technical Assistance Center. 	<p>Reason for Appearance: Fan X was identified as not operating.</p> <p>Solutions: After completing the procedure, contact Boston Scientific Technical Assistance Center</p>
<p>80-35 System Error</p> <ul style="list-style-type: none"> The gas flow rate in channel XX exceeds the system operation requirements and may impact performance. Reduce the number of active needles. 	<p>Reason for Appearance: The calculated flow rate for a given channel exceeded the recommended flow rate for a channel. The number of active needles should be reduced.</p> <p>Solutions:</p> <ol style="list-style-type: none"> Reduce the number of active needles. Contact Boston Scientific Technical Assistance Center for further guidance.
<p>80-36 System Error</p> <ul style="list-style-type: none"> The gas flow rate exceeds the system operation requirements and may impact performance. Reduce the number of active needles. 	<p>Reason for Appearance: The cumulative flow rate for all channels exceeded the recommended flow rate for the system. The number of active needles should be reduced.</p> <p>Solutions:</p> <ol style="list-style-type: none"> Reduce the number of active needles. Contact Boston Scientific Technical Assistance Center for further guidance.

Service Messages

Message	Reason for Appearance / Solutions
90-01 Service Due <ul style="list-style-type: none"> • Service the ICEfx Cryoablation System soon • Service must be completed by [DATE]. • Contact Boston Scientific Technical Assistance Center to schedule service. 	Reason for Appearance: User was reminded of the deadline to complete scheduled system service. Reminder begins four weeks prior to the service due date. Solutions: Contact Boston Scientific Technical Assistance Center to schedule service.
90-02 Service Due <ul style="list-style-type: none"> • Service for the ICEfx Cryoablation System is past due • Contact Boston Scientific Technical Assistance Center to schedule service. 	Reason for Appearance: The ICEfx Cryoablation System was not serviced on schedule. During subsequent start-ups, the user is reminded the service is past due. Solutions: Contact Boston Scientific Technical Assistance Center to schedule service.
90-03 System End of Life <ul style="list-style-type: none"> • The ICEfx Cryoablation System is at the end of operational life • Contact Boston Scientific Technical Assistance Center to arrange return of this system to be refurbished, replaced or disposed. 	Reason for Appearance: The ICEfx Cryoablation System reached the end of its operating life. Solutions: Contact Boston Scientific Technical Assistance Center to schedule service.

SYSTEM SPECIFICATIONS

Mechanical Specifications

ICEfx Cryoablation Console

- Weight: 20 kg (44 lb)
- Height: 28 cm (11 in), monitor down 48 cm (19 in), monitor up
- Foot Print: 36 cm x 53 cm (14 in x 21 in)

ICEfx Cryoablation Cart

- Weight: 23 kg (50 lb)
- Height: 91 cm (36 in)
- Foot Print: 45 cm x 53 cm (18 in x 21 in)

Argon Gas Supply

- Purity Level: 99.998 % or higher
- Solid particle size: < 5 μm

Gas Cylinder Specifications

- Maximum Pressure: 6000 psi (414 bar, 41.4 MPa)
- Recommended Volume of Gas Cylinder: 42 L to 50 L

Accuracy of Displayed Values

• Temperature Accuracy:

- o $\pm 3^{\circ}\text{C}$ over range of -60°C to 40°C

• Supplied Gas Pressure Accuracy:

- o ± 50 psi, over range of 1000 psi to 6000 psi
- o ± 3.4 bar over range of 69 bar to 414 bar
- o ± 0.314 MPa over range of 6.9 MPa to 41.4 MPa

• Built-in Regulator Gas Pressure:

- o ± 50 psi, over range of 1000 psi to 4000 psi
- o ± 3.4 bar over range of 69 bar to 276 bar
- o ± 0.314 MPa over range of 6.9 MPa to 27.6 MPa

• Time Intervals:

- o ± 5 seconds over any 10 minute interval

Essential Performance

The ICEfx Cryoablation System essential performance is defined as:

- The ability to turn on and turn off the flow or argon to the connected cryoablation needles via command inputs from the user interface
- The ability to maintain uninterrupted gas flow to user-selected needles when gas flow is turned on by the user
- The ability to uninterruptedly prevent gas flow to user-selected needles when gas flow is turned off by the user
- The ability to turn on and turn off electrical power to connected CX cryoablation needles via command inputs from the user interface
- The ability to maintain uninterrupted electrical power to the needle heaters of user-selected CX needles when heaters are activated via command inputs from the user interface
- The ability to reliably prevent electrical power flow to user-selected CX needles when electrical heating is discontinued via command inputs from the user interface
- The ability to continuously monitor and display temperatures measured by user-selected thermal sensors

Boston Scientific Technical Assistance Center

Region	Contact Number	email
United States	+1.800.949-6708	CEtechsupportUSA@bsci.com
APAC (Asia, Middle East)	+65.64.18.8878	CETechSupportAPAC@bsci.com
Japan	+81.120.177.779	JapanCESTAC@bsci.com
China	+86 400.801.6908	ChinaCapitalService@bsci.com
Australia/New Zealand	+61 1800.676133 - option 5	CapitalequipmentANZ@bsci.com
Brazil	+55 1155459063	CEBrazilTeam@bsci.com
Mexico	+52 15559924100	CESupportMXC@bsci.com
Europe (See Below for Individual countries)	+31 45.546.7707	CEtechsupportEMEA@bsci.com
Austria	+43 16 08 10 37	CEtechsupportEMEA@bsci.com
Denmark	+45 80 25 34 29	CEtechsupportEMEA@bsci.com
Czech Republic	0800 142942	CEtechsupportEMEA@bsci.com
Finland	+358 800770055	CEtechsupportEMEA@bsci.com
France	+33 139 304 971	CEtechsupportEMEA@bsci.com
Germany	+49 815 126 86118	CEtechsupportEMEA@bsci.com
Italy	+39 022 698 3218	CEtechsupportEMEA@bsci.com
Netherlands	+31 45 5467707	CEtechsupportEMEA@bsci.com
Norway	0800 14236	CEtechsupportEMEA@bsci.com
Spain	+34 917 619 999	CEtechsupportEMEA@bsci.com
Sweden	+46 020790077	CEtechsupportEMEA@bsci.com
UK	+44 1442 411 686	CEtechsupportEMEA@bsci.com

PATIENT COUNSELING INFORMATION

The physician should consider the following points while counseling the patient on the use of the ICEfx Cryoablation System in association with the interventional procedure:

- Discuss the risks and benefits, including review of potential adverse events listed in the IFUs for the ICEfx Cryoablation System and accessory products used to conduct Cryoablation procedures and for other interventional treatments likely to be employed.
- Discuss post-procedure instructions, including any lifestyle changes, medications, and home care and rehabilitation guidelines.

WARRANTY

For device warranty information, visit (www.bostonscientific.com/warranty).

SYMBOL DEFINITIONS

Commonly used medical device symbols that appear on the labeling are defined at www.bostonscientific.com/SymbolsGlossary.

Additional symbols are defined at the end of this document.

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INSTRUCTIONS FOR USE

en

The Instructions for Use (IFU) for this product are supplied in electronic form over the Internet.

Visit www.IFU-BSCI.com to access the IFU using your internet browser. Information from the product labeling such as Product Code, UDI-DI or Product Name may be used to search for the specific IFU.

If you have difficulty accessing the IFU online, or would prefer to receive a paper copy, please contact Boston Scientific Customer Service or your local country contact. A copy will be sent to you at no charge and should arrive within seven days.

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