Visual-ICE™

Cryoablation System

en	User's Manual	2
es	Manual del usuario	89
fr	Manuel d'utilisation	189
de	Bedienungsanleitung	285
it	Manuale dell'utente	384
nl	Gebruikershandleiding	478
pt-eu	Manual do Utilizador	574

TABLE OF CONTENTS

REUSE WARNING	6
DEVICE DESCRIPTION	6
System Description	6
Contents	6
Figure 1. Visual-ICE Cryoablation System Front View	7
Figure 2. Visual-ICE Cryoablation System Rear View	7
Touch Screen Monitor	8
Communication Ports	8
Figure 3. Monitor Storage Basin	8
Mouse Track Pad	8
Storage Compartment	8
Brake Pedal	8
Argon Shutoff Valve	8
Gas Inlets	8
Manual Vent Valve	9
Needle Connection Panel	9
Figure 4. Visual-ICE Cryoablation System Needle Connection Panel	9
Figure 5. Visual-ICE Cryoablation System Needle Channel	9
MTS Connection Ports	9
Software Reset	9
Power Control Knob	10
Needle Channels	10
Operating Principle	10
Materials	10
Non-pyrogenic	10
User Information	10
INTENDED USE	11
INDICATIONS FOR USE	11
Clinical Benefit Statement	11
CONTRAINDICATIONS	11
WARNINGS	11
PRECAUTIONS	14
ADVERSE EVENTS	16
CONFORMANCE TO STANDARDS	18
Table 1. Cable Lengths	18
Table 2. Electromagnetic Emissions	19
Table 3. Electromagnetic Immunity	20

Table 4. Electromagnetic Immunity for Systems that are Not Life-supporting	21
Table 5. Recommended separation distances between portable and mobile RF communications equipment and the Visual-ICE Cryoablation System	22
HOW SUPPLIED	22
Device Details.....	22
Handling and Storage	22
Storage Conditions.....	22
OPERATIONAL INSTRUCTIONS	23
Additional Required Items.....	23
Installation, Calibration, and Servicing.....	24
PREPARATION	24
System Operation	24
Table 6. Cryoablation Procedure Flow	24
Preparation for Use	24
System Set-up.....	25
Screen 1. Disabled Channel.....	26
Screen 2. Vent Gas Message.....	26
Screen 3. Login Screen.....	27
Screen 4. Incorrect Login.....	27
Screen 5. Reset Password Challenge.....	28
Screen 6. Password Reset.....	28
Screen 7. Emergency Login.....	29
Screen 8. Startup Screen.....	29
Figure 6. Visual-ICE Cryoablation System Gas Connections.....	30
Figure 7. Gas Cylinder Set-up	31
Figure 8. EZ-Connect2 Dual Cylinder Adapter.....	32
Screen 9. No Gas Connected message.....	32
Table 7. Working Gas Pressures.....	33
Screen 10. Optimizing Lines message.....	33
Pre-Procedure Testing.....	34
Screen 11. Procedure Screen.....	34
Figure 9. Locking Needle into Channel	35
Screen 12. Thawing Mode message.....	35
Screen 13. Needle Past Expiration Date message.....	36
Screen 14. Select Needle Type Menu	36
Figure 10. MTS Connection	38
Navigating The User Interface.....	38
Screen 15. Login Screen	39

Startup Screen.....	39
Screen 16. Startup Screen	39
Table 8. Startup Screen Buttons.....	39
Screen 17. Procedure Screen.....	40
Navigation Tool Bar.....	40
Screen 18. Navigation Tool Bar	40
Table 9. Navigation Tool Bar	41
Table 10. Channel Controls.....	42
Channel Status	43
Screen 19. Channel Controls and Channel Status Section	43
Screen 20. Active Thaw Channel Status Section	43
Screen 21. Needle Tip Temperature Display.....	43
Screen 22. Enlarged Timer	44
Screen 23. Repositioned Enlarged Timers	44
Temperature Sensors	45
Screen 24. Temperature Sensors Section.....	45
Organ Map	46
Screen 25. Organ Map Section	46
Table 11. Organ Map Controls	47
Screen 26. Procedure Report Example.....	48
Screen 27. View Reports Screen.....	49
Screen 28. Export Report Screen.....	50
Configuring Settings.....	51
Screen 29. Configure Settings.....	51
Table 12. Configure Settings Options.....	51
PROCEDURE	52
Performing a Cryoablation Procedure.....	52
Screen 30. Gas Time Remaining	52
Reports.....	54
Screen 31. Export Report Screen	54
Screen 32. Report Exported message	55
System Shutdown	55
Changing Gas Cylinders during a Procedure.....	56
Standard Gas Cylinder Setup	56
Dual Gas Cylinder Connection.....	56
i-Thaw Function and FastThaw Function Control for CX-type Needles.....	57
Screen 33. Advanced Thaw Controls FastThaw Function	57
Screen 34. Advanced Thaw Controls i-Thaw Function	58

Screen 35. Advanced Thaw - Needle Warming.....	58
Screen 36. Advanced Thaw - Temperature Display During Thawing.....	59
Screen 37. Cautery Confirmation message.....	59
Screen 38. 1.5 CX Needle Warming.....	60
Screen 39. 1.5 CX Needle Cautery In Progress.....	60
Screen 40. Advanced Channel Controls.....	62
Screen 41. Linked Channels.....	62
Cycle Programming Control.....	63
Screen 42. Advanced Cycle Controls.....	63
Screen 43. Saved Sequence Controls.....	64
Screen 44. Saved Sequence Controls.....	65
Advanced Temperature Sensor Controls.....	66
Screen 45. Advanced Temperature Sensor Controls.....	66
ADMINISTRATIVE FUNCTIONS.....	67
Configure Settings.....	67
Screen 46. Configure Settings.....	67
Table 13. Configure Settings Controls.....	68
Manual Software Update.....	69
Screen 47. Software Update Confirmation.....	69
POST-PROCEDURE.....	69
Cleaning the Visual-ICE Cryoablation System.....	69
Disposal.....	69
TROUBLESHOOTING.....	70
Software Recovery.....	70
Screen 48. Software Recovery Screen.....	70
Screen 49. Invalid Configuration Message.....	71
Electronics, Electrical and User Error Related Issues.....	72
Replacing Fuses.....	73
Gas Issues.....	74
Mechanical Issues.....	76
Gas Cylinder and Gas Supply Line.....	76
Needles.....	77
Displayed Messages.....	78
SYSTEM SPECIFICATIONS.....	87
External Gas Supply.....	87
PATIENT COUNSELING INFORMATION.....	88
WARRANTY.....	88
SYMBOL DEFINITIONS.....	88

Rx ONLY

Caution: Federal Law (USA) restricts this device to sale by or on the order of a physician.

REUSE WARNING

The single-use only devices used with the Visual-ICE Cryoablation System are sterile. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

DEVICE DESCRIPTION

System Description

The Visual-ICE Cryoablation System is a mobile cryoablation system that contains ten needle connection channels (each supporting two needle ports), four Multi-Point 1.5 Thermal Sensor Device (MTS) ports, one argon inlet connector, one helium inlet connector, a recessed 19 inch touch screen monitor, a mouse track pad, a USB port, and an Ethernet port (inactive). The system hardware and operating software are housed within the system.

The system is mounted on four swivel caster wheels for system mobility. The system contains a three-way brake pedal mechanism in front to immobilize the two front wheels of the system during a procedure and to provide directional control during movement of the system. The rear wheels have individual brake pedals. Gas supply line clips at the rear of the system (Figure 2) are used to direct the gas supply lines toward the floor, minimizing tripping hazards. A compartment in the lower section provides storage for system accessories.

Contents

One (1) Visual-ICE Cryoablation System

One (1) Visual-ICE Cryoablation System User Manual: The User Manual may be physical copy or available online at www.IFU-BSCI.com. The User Manual Describes the system and provides instructions for system operation and maintenance.

One (1) Visual-ICE Cryoablation System Quick Reference Guide: The Quick Reference Guide summarizes key steps in system operation.

One (1) Wrench

One (1) Console Cover: The Console Cover is used to protect the Visual-ICE Cryoablation System during storage.

One (1) USB flash drive (4 GB) within an attached pouch: The USB flash drive is used to transfer procedure reports to a customer's computer for saving or printing.

TWO (2) Single Cylinder Adapters: Single Cylinder Adapters are comprised of a high-pressure gas supply line attached to a pressure gauge.

One (1) Argon Single Cylinder Adapter: The Argon Single Cylinder Adapter is a high-pressure argon gas supply line with attached pressure gauge.

- To accommodate procedure room variations, the gas supply line is available in alternative lengths to connect the argon cylinder to the Visual-ICE Cryoablation System. Reference Table 1.

One (1) Helium Single Cylinder Adapter: The Helium Single Cylinder Adapter is a high-pressure helium gas supply line with attached pressure gauge.

- To accommodate procedure room variations, the gas supply line is available in alternative lengths to connect the helium cylinder to the Visual-ICE Cryoablation System. Reference Table 1.

Optional

One (1) EZ-Connect2 Dual Cylinder Adapter: The EZ-Connect2 Dual Cylinder Adapter is an optional component used to connect two gas cylinders in tandem to the Visual-ICE Cryoablation System. The dual cylinder adapter consists of a four-way adapter assembly with argon pressure gauge and long gas supply line with system connection, short gas supply line with cylinder connection. Reference the **Dual Gas Cylinder Connection** Section for instructions on using the EZ-Connect2 Dual Cylinder Adapter.


Figure 1. Visual-ICE Cryoablation System Front View

- | | | |
|-------------------------|---------------------------|--------------------|
| 1 Touch Screen Monitor | 4 Needle Connection Panel | 7 Brake Pedal |
| 2 Monitor Storage Basin | 5 Power Control Knob | 8 Wheels (casters) |
| 3 Mouse Track Pad | 6 Storage Compartment | |


Figure 2. Visual-ICE Cryoablation System Rear View

- | | | |
|--------------------------|----------------------------------|---------------------------|
| 1 Argon Shutoff Valve | 5 Cord Wrap Bracket | 9 Manual Vent Valve |
| 2 Cooling Fan | 6 Power Switch | 10 Rear Wheel Brake Pedal |
| 3 Argon Inlet Connector | 7 Ground Post (select countries) | |
| 4 Helium Inlet Connector | 8 Gas Supply Line Clip | |

Touch Screen Monitor

The cryoablation procedure is controlled through the use of the touch screen monitor. The monitor can be tilted and rotated to provide the optimal viewing and operating angle for the user. The touch screen monitor includes a virtual on-screen English QWERTY keyboard for entering procedure related data and may be operated by finger touch. The monitor folds flat into the monitor storage basin on the top of the unit for system storage.

Communication Ports

2 communication ports are located on the back panel of the monitor storage basin (Figure 3).

- The Ethernet port is inactive
- The USB 2.0 port allows you to save reports to a USB flash drive for download to another computer or for printing.


Figure 3. Monitor Storage Basin

- 1 Ethernet Port (Inactive)
- 2 USB 2.0 Port

Mouse Track Pad

The mouse track pad is located within the recessed monitor storage basin. The mouse track pad provides an alternative to the touch screen as a means of interacting with the system. Use the track pad to move and position the cursor on the monitor. To press a button on the screen, position the cursor over the button and press the left button on the track pad.

Storage Compartment

The storage compartment may be used to store Visual-ICE Cryoablation System accessories, such as the gas supply lines and tools. Do not place very heavy objects in the storage compartment. The weight limit is 23 kg (50 lb). Do not store liquids in the compartment. Liquids spilled in the storage compartment may drip into the system; the compartment is not watertight.

Brake Pedal

The brake pedal operates on the front two wheels of the Visual-ICE Cryoablation System. Position the brake in the UP position to prevent the front two wheels from swiveling during transportation. Position the brake in the DOWN position to lock the front two wheels into place. When the brake pedal is in the middle position, the front two wheels are free to roll and swivel. If the flooring is uneven, it may be necessary to lock the two rear wheels in addition to locking the front wheels. Lock each rear wheel by using the individual locking pedal on each wheel.

Argon Shutoff Valve

The Argon Shutoff Valve is used to turn the gas supply to the Visual-ICE Cryoablation System ON or OFF. It should be left in the **Argon ON** position and used to turn **OFF** argon gas only in an emergency.

Gas Inlets

The gas supply lines connect the argon and helium gas supply from the respective gas cylinders to the argon and helium gas inlets. The argon inlet is a plug connector; the helium inlet is a socket connector.

Manual Vent Valve

The Manual Vent Valve is used to vent the high-pressure gas from the Visual-ICE Cryoablation System if the automatic venting feature is not used.

Needle Connection Panel


Figure 4. Visual-ICE Cryoablation System Needle Connection Panel


Figure 5. Visual-ICE Cryoablation System Needle Channel

- | | | | |
|-----------------------|----------------------|--|------------------|
| 1 MTS Connection Port | 3 Power Control Knob | 5 Locking Bar | 7 Channel Number |
| 2 Software Reset | 4 Needle Channel | 6 Needle Port with Electrical Connection | |

MTS Connection Ports

Four MTS Connection Ports are located above the Needle Connection Channels. The MTS Connection Ports support connection of the Boston Scientific Multi-Point 1.5 Thermal Sensor Devices that are used to monitor temperature within and near the cryoablation target site.

Software Reset

The **Software Reset** button is used to boot the Visual-ICE Cryoablation System into a recovery mode if the software becomes corrupted (reference the **Software Recovery** Section).

Power Control Knob

The Power Control Knob powers ON the Visual-ICE Cryoablation System in preparation for a procedure.

Needle Channels

The Needle Connection Panel contains ten numbered needle channels; each channel contains two ports to support connection of up to two cryoablation needles. Each channel operates independently of all other channels in either freezing or thawing mode. The electrical connection is used for needles with a memory chip, i-Thaw Function, FastThaw Function or cautery function capabilities. The locking bar on each channel locks the needles into the ports to secure them during the procedure.

Operating Principle

The Visual-ICE Cryoablation System is a mobile system intended for cryoablative tissue destruction using a minimally invasive procedure. The system is computer-controlled with a touch screen user interface that allows the user to control and monitor the procedure. Innovative gas dryers produce consistent iceballs and boost freezing performance for all needles.

The therapy delivered by the system is based on the Joule-Thomson effect displayed by compressed gases. The Joule-Thomson effect is a change in the temperature of a compressed gas as it flows through a narrow orifice and expands to a lower pressure. Certain gases, such as argon, decrease in temperature due to the Joule-Thomson effect, while other gases, such as helium, increase in temperature.

The Visual-ICE Cryoablation System uses high-pressure argon gas that circulates through closed-tip cryoablation needles to induce tissue freezing. Active tissue thawing is achieved by circulating helium gas through the needles or, alternatively, using a heating element inside the cryoablation needle which can be energized to cause thawing (i-Thaw Function). The Visual-ICE Cryoablation System also controls the heating element inside cryoablation CX needles to provide active helium-free thawing (i-Thaw Function or FastThaw Function) and track ablation (cautery function).

Tissue ablation is achieved by repeated freeze and thaw cycles with both freezing and thawing contributing to cell death. Generally, multiple freeze-thaw cycles are used to achieve complete destruction of the target tissue.

When multiple cryoablation needles are placed into or near the target tissue and freezing is initiated, an iceball grows around the distal end of the needle shafts. In time, the iceballs coalesce and completely engulf the target tissue. An important benefit of cryoablation is that imaging procedures, such as ultrasound and CT, are able to display the location and size of the iceball. This benefit of cryoablation is used for proper control of the therapy. During use, the procedure must be monitored via image guidance to ensure adequate tissue coverage and to avoid damage to adjacent structures.

In addition to image guidance, Boston Scientific provides temperature sensors to aid in monitoring tissue temperature near the target site and adjacent critical structures. These temperature sensors can provide quantitative data to supplement the qualitative information provided by the imaging modality. Needle tip temperature display for CX-type needles provides a visual means to monitor needle performance.

Materials

See Boston Scientific cryoablation needle and accessory product IFUs for specific information on materials.

Non-pyrogenic

See Boston Scientific cryoablation needle and accessory product IFUs for specific information on pyrogenicity.

User Information

The Visual-ICE Cryoablation System is designed to be operated by medical professionals who have a thorough understanding of the technical principles, clinical applications, and risks associated with cryoablation procedures. Optional education is available from your Boston Scientific representative.

INTENDED USE

The Visual-ICE Cryoablation System is intended for cryoablative destruction of tissue during minimally invasive procedures; various Boston Scientific accessory products are required to perform these procedures. The Visual-ICE Cryoablation System is intended for use as a cryosurgical tool in the fields of general surgery, dermatology, neurology (including cryoanalgesia), thoracic surgery (with the exception of cardiac tissue), gynecology, oncology, and urology. This system is designed to destroy tissue (including prostate and kidney tissue, liver metastases, tumors, and skin lesions) by the application of extremely cold temperatures.

Patient Groups

The intended population includes patients intended for cryoablative destruction of tissue during surgical procedures.

INDICATIONS FOR USE

The Visual-ICE Cryoablation System is indicated for use as a cryosurgical tool in the fields of general surgery, dermatology, neurology (including cryoanalgesia), thoracic surgery (with the exception of cardiac tissue), gynecology, oncology, and urology. This system is designed to destroy tissue (including prostate and kidney tissue, liver metastases, tumors, and skin lesions) by the application of extremely cold temperatures.

The Visual-ICE Cryoablation System has the following specific indications:

- Urology - Ablation of prostate tissue in cases of prostate cancer
- Oncology - Ablation of cancerous or malignant tissue and benign tumors, and palliative intervention
- Dermatology - Ablation or freezing of skin cancers and other cutaneous disorders
- Gynecology - Ablation of malignant neoplasia or benign dysplasia of the female genitalia
- General surgery - Palliation of tumors, recurrent cancerous lesions, and ablation of breast fibroadenomas
- Thoracic surgery - (with the exception of cardiac tissue)

Clinical Benefit Statement

The Visual-ICE Cryoablation System, when used with various Boston Scientific accessory products, is intended to destroy tissue (including prostate and kidney tissue, liver metastases, tumors, and skin lesions) by the application of extremely cold temperatures during minimally invasive procedures.

The clinical benefit is measured by overall clinical outcomes with acceptable safety specific to the target anatomy and indication.

CONTRAINDICATIONS

There are no known contraindications specific to the use of the Visual-ICE Cryoablation System.

WARNINGS

General

- The Visual-ICE Cryoablation System is designed to be operated by medical professionals who have a thorough understanding of the technical principles, clinical applications, and risks associated with cryoablation procedures.
- See Boston Scientific cryoablation needle and accessory product IFUs for Warnings specific for these products.
- Do not use this device for any purpose other than the stated intended use and indications for use.
- Do not modify the Visual-ICE Cryoablation System in any way. Only authorized Boston Scientific personnel or Boston Scientific-trained authorized personnel are to service Visual-ICE Cryoablation Systems.
- The Visual-ICE Cryoablation System should be inspected and serviced periodically according to the system specifications. Service must be performed by authorized service engineers. Reference the **Installation, Calibration, and Servicing** Section for detailed information.
- Do not use the Visual-ICE Cryoablation System if the system is visibly damaged, exposing any internal components or sharp edges.

- Do not use the Visual-ICE Cryoablation System near magnetic resonance imaging (MRI) equipment.
- The Visual-ICE Cryoablation System should not be used adjacent to or stacked with other equipment.
- Lock the wheels on the Visual-ICE Cryoablation System prior to using the system to avoid inadvertent movement of the system during a procedure.
- To avoid risk of electric shock, this equipment must only be connected to a hospital grade electrical power outlet with a protective earth.
- Do not start a cryoablation procedure before verifying that the Visual-ICE Cryoablation System and all ancillary equipment are fully operational.
- Use of cables other than those specified, with the exception of those sold by Boston Scientific for use as replacement parts for internal components, may result in increased emissions or decreased immunity of the Visual-ICE Cryoablation System.
- Only use non-MRI needles with the Visual-ICE Cryoablation System.
- Boston Scientific recommends using the Boston Scientific MTS to monitor the freeze / thaw temperatures for the intended treatment protocol and to monitor temperatures in the adjacent organs and structures.
- Do not use the needle if it is bent or damaged while attempting to unpack or use it. Never use a defective needle for a cryoablation procedure. A defective cryoablation needle that has a gas leak can cause a gas embolism in the patient.
- Do not kink, pinch, cut or pull excessively on the needle tubing. Damage to needle handle or tubing may cause the needle to become unusable.
- Have sufficient argon gas available to conduct the planned cryoablation procedure: the number and type of needles, gas cylinder size, pressure and rate of gas flow affect the required gas volume (reference the **SYSTEM SPECIFICATIONS** Section for gas purity requirements). At least one full spare cylinder should be available for each treatment.
- High-pressure gas is dangerous if handled improperly. Local laws and safety rules regarding pressurized gas systems, reservoirs and components should always be observed.
- Ensure that the gas cylinders are chained to a wall or an approved cart to prevent inadvertent tipping of the cylinders.
- Do not connect the Visual-ICE Cryoablation System to a gas supply exceeding 6000 psi (414 bar, 41.4 MPa) to avoid damage to internal system components.
- The Visual-ICE Cryoablation System should not be operated in the presence of flammable fumes, e.g., flammable anesthetics or volatile substances.
- Do not bend or kink the gas supply line. Sharp bends or kinks may compromise the integrity of the gas supply line.
- Do not roll the Visual-ICE Cryoablation System over the gas supply line; such activity may damage the line.

Procedural

- Prior to starting a cryoablation procedure, set up the Visual-ICE Cryoablation System (reference the **System Set-up** Section) and then perform Needle Integrity and Functionality Tests. Tests must be successfully completed in order to begin the procedure.
- Do not use the needle if there is no ice formation during the Freeze phase. Obtain a new needle and repeat the testing procedure.
- Do not use the needle if bubbles are seen escaping from the needle during Needle Integrity and Functionality Test.
- Ensure adequate measures are taken to protect organs and structures adjacent to the targeted tissue.
- Ensure that the MTSs are functioning properly before inserting into the patient by verifying that they are reading a reasonable room temperature.
- The sterile field and sterility of cryoablation needles should be maintained at all times. Do not contaminate the distal end of the sterile cryoablation needle.

- Avoid contact with the distal portion of the cryoablation needle to maintain sterility during testing.
- Continuously monitor needle insertion, needle positioning, iceball formation and removal using image guidance (such as direct visualization, ultrasound, or Computed Tomography (CT)) to ensure adequate tissue coverage and to avoid damage to adjacent structures.
- Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the Visual-ICE Cryoablation System, including cables specified for use with the system. Otherwise, degradation of the performance of this equipment could result.
- Confirm the high-pressure gas supply line(s) is securely connected before opening the gas cylinder(s).
- Secure the safety cable at the end of the gas supply line to the System before connecting the argon gas supply line to the argon gas inlet. The safety cable provides backup protection if the gas supply line becomes inadvertently disconnected from the system. Do not use a gas supply line with a missing safety cable. Doing so could compromise the safety of personnel in the room. Contact Boston Scientific Technical Assistance Center for further instructions.
- Each needle must be locked into a needle channel before initiating a cryoablation procedure to avoid the risk of forceful ejection of the needles while under gas pressure.
- If the needles are still connected, do not unlock the channels or disconnect the needles from the needle connection panel until all operations in the channel have completed.
- Use **Freeze** and **Thaw** operations only when the needle is placed in the target tissue.
- Needle handles and the gas line may frost during freezing. Avoid prolonged contact with frosted portions of a needle handle to avoid unintended thermal tissue damage to the patient or clinician.
- Needle tubing may become extremely cold when conducting freeze cycles during a cryoablation procedure. It is important that a patient's skin is protected from direct contact with needle tubing to avoid the potential for thermal injury to the patient. Ensure an appropriate insulating barrier is placed as needed (such as towels) or other method is employed to prevent needle tubing from touching a patient's skin.
- The needle handle may become warm during active thawing. Pay attention to the needle handle position. Prolonged contact with warm portions of the needle handle could cause unintended thermal tissue damage/burn to the patient or clinician.
- Active thawing produces heat along the distal needle shaft. Use care to avoid thermal injury/burn to nontargeted tissues.
- Ensure adequate thawing or cooling before attempting to remove needles from the patient.
- Discontinue all needle operation prior to needle removal to minimize risk of thermal injury and/or tissue injury.
- When conducting **FastThaw** Function or when conducting **cautery** function for Track Ablation, be alert for the Active Zone Indicator as the needle is withdrawn to prevent unintended tissue damage from the hot needle.
- Do not touch the Visual-ICE Cryoablation System while touching the patient to avoid the risk of shocking the patient if an inadvertent electrical fault exists.
- Do not touch the screen if the touch screen monitor goes blank for more than five (5) seconds during a procedure. Immediately turn off power to the system and end the procedure to avoid inadvertent activation of needles.
- Warn procedure personnel prior to venting the Visual-ICE Cryoablation System to avoid startling them.
- If it is difficult to loosen the pressure gauge connected to the cylinder, or the high-pressure gas supply line(s) cannot be disconnected from the inlet connections, do not apply excessive force to release the gas supply line or to loosen the pressure gauge. The gas line may still be under pressure.
- Do not pull on the power cord. Grasp the plug, not the power cord, to disconnect the device from the wall socket.
- Dispose of device and accessories in accordance with the **Disposal** Section.

PRECAUTIONS

General

- Carefully read all instructions prior to use. Failure to observe all warnings and precautions may result in complications.
- Do not use the Visual-ICE Cryoablation System if any moisture or condensation is present on the surfaces of the system. Allow the system to dry completely for 12 hours prior to powering up the system. Powering up the system containing moisture or condensation could result in permanent damage to the electrical boards, causing the system to be inoperable.
- Take precautions to avoid potential electrostatic discharge. If an electrostatic discharge occurs after touching the monitor, the screen may flicker and MTS temperature readings may be inaccurate for a few seconds. The system will remain functional and the monitor will refresh momentarily.
- Use care to avoid electrostatic discharge (ESD) events when removing the cover from the Visual-ICE Cryoablation System. Boston Scientific recommends that the operator touch one or more metal parts on the rear of the system prior to touching anything on the Needle Connection Panel.
- No data regarding cryoablation in combination with other therapies is available from Boston Scientific.
- Use the Boston Scientific supplied USB flash drive only to export reports or update software. Other data or software may corrupt the Visual-ICE Cryoablation System.
- Do not connect any other USB equipment to the Visual-ICE Cryoablation System USB port.
- Do not use a USB extension cable to connect the USB flash drive to the USB port. Connect USB flash drive directly to the USB port provided on the Visual-ICE Cryoablation System. Use of a USB extension cable may result in electromagnetic emissions exceeding regulatory limits.
- Select a unique Patient ID that does not reveal the patient's identity to other system users.

Handling

- Handle the Visual-ICE Cryoablation System with care. Rough handling may damage the system and cause it to become inoperable. The system should never be tilted.
- Maneuver the Visual-ICE Cryoablation System by pulling on the system using the rear handle.
- Do not set food, drinks, or other objects on top of the system. Doing so may damage the system.
- Do not store liquids in the storage compartment. The storage compartment is not watertight.
- Do not set heavy objects on the monitor when in the down position or on the monitor storage basin when the monitor is in the up position. The weight limit is 9 kg (20 lb).
- Prior to lowering the monitor, ensure that no objects are located in the monitor storage basin. Use care when lowering the monitor into the monitor storage basin; do not use excessive force to avoid damaging the monitor.
- Use care when lowering or swiveling the touch screen monitor to avoid the potential for pinching fingers.
- Lift the Visual-ICE Cryoablation System to clear any threshold that is higher than 1 cm. Two people, one on each side, should use the handles to lift the system.
- Clean the Visual-ICE Cryoablation System by following the instructions in the **Handling and Storage** Section. Do not use cleaning agents such as Betadine Antiseptic Solution or bleach solution, which may damage the touch screen.
- Position the argon cylinder close enough to the system to ensure that the gas supply line is not stretched and does not present a trip hazard.
- Direct the high-pressure gas supply lines toward the floor and secure the lines with the clips located on the rear of the Visual-ICE Cryoablation System to minimize the potential for tripping.

Procedural

- The Visual-ICE Cryoablation System should be placed in close proximity for accessible needle connection and use.
- Turn the Visual-ICE Cryoablation System ON before connecting the gas cylinders to ensure the proper diagnostic tests are conducted.
- Confirm that the Manual Vent Valve is closed and Argon Shutoff Valve is in the ON position before connecting the gas line to the system.
- If the system produces a continuous hissing sound, verify that the Manual Vent Valve is fully closed. If the Manual Vent Valve is fully closed and the hissing sounds persists, shut OFF the system using the power control knob located on the front of the system (Figure 1). Close the gas supplies using the cylinder valves. Contact Boston Scientific Technical Assistance Center.
- Failure to operate the Visual-ICE Cryoablation System within the working pressure limits indicated on the user interface (Table 7) may affect the iceball formation.
- Boston Scientific recommends that only needles of the same type be placed together in a single channel. Do not combine i-Thaw Function capable needles with non-i-Thaw Function capable needles in the same channel. Additionally, using needles of differing types in a channel may affect the accuracy of the **Gas Indicator**.
- During use, avoid damage to the needle from other surgical instruments.
- Do not allow the MTS to touch the cryoablation needle during use; erratic temperature measurements may be displayed.
- The cautery function should not be initiated if the Active Zone Indicator is visible outside the patient's skin.
- If a needle appears to be blocked, press the **Thaw** button to thaw the needle for at least one minute to clear the blockage.
- When the argon cylinder pressure falls below the lower working pressure limit, the system displays an alert message. To ensure optimal performance, replace an argon cylinder if the pressure falls below the lower working pressure limit.
- Depressurize the system after the cryoablation procedure is completed (reference the **System Shutdown** Section).
- If the Visual-ICE Cryoablation System detects that helium gas is connected, the system defaults to helium thaw mode and i-Thaw Function, FastThaw Function, and cautery function capabilities are disabled.
- Ice formation during the Flush and Thaw phases indicates that argon gas is connected to the helium inlet. Before proceeding change the cylinders and ensure that each gas supply line is connected to the correct cylinder (reference the **Standard Gas Cylinder Setup** Section).
- Any interruption of a programmed phase immediately terminates that phase and the programmed cycle.
- Activation of thermal cautery function immediately ceases all other cryoablation operations for the duration of cautery function activity.

ADVERSE EVENTS

The potential adverse events associated with the device and/or cryoablation procedure include, but are not limited to:

- Angina
- Arrhythmia
- Atelectasis
- Bladder spasms
- Bleeding/hemorrhage
- Burn/frostbite
- Cerebrovascular Accident (CVA)/stroke
- Cryoshock Phenomenon (e.g. multi organ failure, severe coagulopathy, disseminated intravascular coagulation (DIC))
- Death
- Distension
- Edema/swelling
- Ejaculatory dysfunction
- Embolism (air, device, thrombus)
- Erectile dysfunction
- Fever
- Fistula
- Fracture
- Gastrointestinal symptoms (e.g. nausea, vomiting, diarrhea, constipation)
- Healing, impaired
- Hematoma
- Hematuria
- Hemothorax
- Hepatic dysfunction/failure
- Hernia
- Hypertension
- Hypotension
- Hypothermia
- Ileus
- Impotence
- Infection/abscess/sepsis
- Inflammation

- Muscle spasm
- Myocardial infarction
- Necrosis
- Need for additional intervention or surgery
- Nerve Injury
- Neuropathy
- Obstruction
- Pain/discomfort
- Perforation (including organ and adjacent structures)
- Pericardial Effusion
- Perirenal fluid collection
- Pleural effusion
- Pneumatosis (air or gas in an abnormal quantity and/or place in the body)
- Pneumothorax
- Post ablation syndrome (e.g. fever, pain, nausea, vomiting, malaise, myalgia)
- Renal insufficiency/failure
- Renal parenchymal or capsule fracture
- Respiratory distress/insufficiency/failure
- Scrotal edema
- Stenosis/stricture
- Subcutaneous emphysema
- Thrombosis/thrombus
- Tissue damage
- Transient ischemic attack (TIA)
- Tumor cell seeding
- Urethral sloughing
- Urinary frequency/urgency
- Urinary incontinence
- Urinary retention
- Urinary tract infection
- Vasovagal response
- Vessel trauma (e.g. dissection, injury, perforation, pseudoaneurysm, rupture, or other)
- Wound infection

CONFORMANCE TO STANDARDS

Electrical Specifications:

- Input Voltage: 100 VAC to 240 VAC, single phase
- Input Frequency: 50 Hz – 60 Hz
- VA Rating: 250 VA
- IP Rating: IP10
- Fuse Rating: T 3.15AL
- Electrical Protection: Class I, Type BF protection against shock
- Signal Input/Output Ports: one (1) Ethernet port (inactive), one (1) USB 2.0 port

Electromagnetic Compatibility and Immunity (EMC and EMI)

The Visual-ICE Cryoablation System needs special precautions regarding electromagnetic compatibility (EMC) and needs to be installed and placed into service according to the EMC information provided below.

The Visual-ICE Cryoablation System has been tested in the operating room environment for electromagnetic compatibility (EMC) and electromagnetic interference (EMI) compliance. The Visual-ICE Cryoablation System has been tested to be in compliance with IEC 60601-1-2 and EN 55011.

Portable and mobile radiofrequency (RF) communications equipment can affect the Visual-ICE Cryoablation System, causing it to operate improperly.

Table 1. Cable Lengths

Cable	Length
Power Cable	4.6 m (15 ft)
Gas Tubing (connected to needles)	2.5 m (8 ft)
Gas Supply Line (connected to argon cylinder)	Available lengths: 5 m (16 ft), 8 m (26 ft)

NOTE: The Gas Supply Line is available in more than one length to accommodate procedure room variations.

WARNING: Use of cables other than those specified, with the exception of those sold by Boston Scientific for use as replacement parts for internal components, may result in increased emissions or decreased immunity of the Visual-ICE Cryoablation System.

WARNING: The Visual-ICE Cryoablation System should not be used adjacent to or stacked with other equipment.

WARNING: Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 in) to any part of the Visual-ICE Cryoablation System, including cables specified for use with the system. Otherwise, degradation of the performance of this equipment could result.

Table 2. Electromagnetic Emissions

Guidance and manufacturer's declaration – electromagnetic emissions		
The Visual-ICE Cryoablation System is intended for use in a professional healthcare facility electromagnetic environment with the compliance levels specified below. The customer or the user of the Visual-ICE Cryoablation System should assure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment – guidance
RF emissions CISPR 11	Group 1	The Visual-ICE Cryoablation System uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class A	
Harmonic emissions, IEC 61000-3-2	Class A	
Voltage fluctuations/ flicker emissions, IEC 61000-3-3	Complies	
<p>NOTE: The emissions characteristics of this equipment make it suitable for use in industrial areas and hospitals (CISPR 11 class A). If it is used in a residential environment (for which CISPR 11 class B is normally required) this equipment might not offer adequate protection to radio-frequency communication services. The user might need to take mitigation measures, such as relocating or reorienting the equipment.</p>		

Table 3. Electromagnetic Immunity

Guidance and manufacturer's declaration – electromagnetic immunity			
The Visual-ICE Cryoablation System is intended for use in a professional healthcare facility electromagnetic environment with the immunity compliance levels specified below. The customer or the user of the Visual-ICE Cryoablation System should assure that it is used in such an environment.			
Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 8 kV contact ± 15 kV air	± 8 kV contact ± 15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	± 2 kV for power supply lines ± 1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 0.5 kV, ± 1 kV line(s) to line(s) ± 0.5 kV, ± 1 kV, ± 2 kV line(s) to earth	± 0.5 kV, ± 1 kV line(s) to line(s) ± 0.5 kV, ± 1 kV, ± 2 kV line(s) to earth	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	0 % U_T ; 0.5 cycle at 0°, 45°, 90°, 135°, 180°, 225°, 270°, and 315° 0 % U_T ; 1 cycle 70 % U_T ; 25 cycles/30 cycles at 0° and 50Hz/60Hz. 0 % U_T ; 250 cycles/300 cycles at 50Hz/60Hz.	0 % U_T ; 0.5 cycle at 0°, 45°, 90°, 135°, 180°, 225°, 270°, and 315°. 0 % U_T ; 1 cycle 70 % U_T ; 25 cycles/30 cycles at 0° and 50Hz/60Hz. 0 % U_T ; 250 cycles/300 cycles at 50Hz/60Hz.	Mains power quality should be that of a typical commercial or hospital environment. If the user of the Visual-ICE Cryoablation System requires continued operation during power mains interruptions, it is recommended that the Visual-ICE Cryoablation System be powered from an uninterruptible power supply or a battery.
Power frequency (50 Hz/60 Hz) magnetic field IEC 61000-4-8	30 A/m	30 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
NOTE: U_T is the a.c. mains voltage prior to application of the test level.			

Table 4. Electromagnetic Immunity for Systems that are Not Life-supporting

Guidance and manufacturer's declaration – electromagnetic immunity			
The Visual-ICE Cryoablation System is intended for use in a professional healthcare facility electromagnetic environment with the immunity compliance levels specified below. The customer or the user of the Visual-ICE Cryoablation System should assure that it is used in such an environment.			
Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic environment - guidance
Conducted RF IEC 61000-4-6	3 Vrms, 6 Vrms in ISM bands over 150 kHz to 80 MHz 80 % AM at 1 kHz	3 Vrms, 6 Vrms in ISM bands over 150 kHz to 80 MHz 80 % AM at 1 kHz	Portable and mobile RF communications equipment should be used no closer to any part of the Visual-ICE Cryoablation System, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance: $d = \left[\frac{3,5}{3} \right] \sqrt{P}$ $d = \left[\frac{3,5}{3} \right] \sqrt{P} \quad 80 \text{ MHz to } 8200 \text{ MHz}$ $d = \left[\frac{7}{3} \right] \sqrt{P} \quad 800 \text{ MHz to } 2.5 \text{ GHz}$ where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey ^a , should be less than the compliance level in each frequency range ^b . Interference may occur in the vicinity of equipment marked with the following symbol: 
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.7 GHz 80 % AM at 1 kHz	3 V/m 80 MHz to 2.7 GHz 80 % AM at 1 kHz	
Radiated RF Proximity fields IEC 61000-4-3 (per IEC 60601-1-2 Ed. 4)	9 V/m – 28 V/m per IEC 60601-1-2 Ed. 4, Table 9	9 V/m – 28 V/m per IEC 60601-1-2 Ed. 4, Table 9	
<p>NOTE 1: At 80 MHz and 800 MHz, the higher frequency range applies.</p> <p>NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.</p> <p>^a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the Visual-ICE Cryoablation System is used exceeds the applicable RF compliance level above, the Visual-ICE Cryoablation System should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the Visual-ICE Cryoablation System.</p> <p>^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.</p>			

Table 5. Recommended separation distances between portable and mobile RF communications equipment and the Visual-ICE Cryoablation System

Recommended separation distances between portable and mobile RF communications equipment and the Visual-ICE Cryoablation System			
The Visual-ICE Cryoablation System is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the Visual-ICE Cryoablation System can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Visual-ICE Cryoablation System as recommended below, according to the maximum output power of the communications equipment.			
Rated maximum output power of transmitter Watts (W)	Separation distance according to frequency of transmitter Meters (m)		
	150 kHz to 80 MHz $d = \left[\frac{3,5}{3} \right] \sqrt{P}$	80 MHz to 800 MHz $d = \left[\frac{3,5}{3} \right] \sqrt{P}$	800 MHz to 2.5 GHz $d = \left[\frac{7}{3} \right] \sqrt{P}$
0.01	0.117	0.117	0.233
0.1	0.369	0.369	0.738
1	1.167	1.167	2.333
10	3.689	3.689	7.379
100	11.667	11.667	23.333
For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.			
NOTE 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.			
NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			

HOW SUPPLIED

Device Details

The Visual-ICE Cryoablation System is supplied non-sterile and is intended for multiple use. Boston Scientific accessory products required to perform the cryoablation procedure are supplied separately.

Do not use if package is damaged or unintentionally opened before use.

Do not use if labeling is incomplete or illegible.

Handling and Storage

Operating Conditions

- Temperature: 10 °C to 40 °C
- Relative Humidity: 30 % to 75 %

Storage Conditions

- Temperature: -15 °C to 50 °C
- Relative Humidity: 10 % to 90 %

Transportation Conditions

When shipping the Visual-ICE Cryoablation System, use the original shipping container to prevent damage from occurring during transport. If the original shipping container is not available, the customer takes responsibility to ensure proper transport conditions are satisfied or contacts Boston Scientific Technical Assistance Center to obtain the appropriate shipping container.

CAUTION: Do not set food, drinks, or other objects on top of the System. Doing so may damage the system.

OPERATIONAL INSTRUCTIONS

WARNING: The Visual-ICE Cryoablation System is designed to be operated by medical professionals who have a thorough understanding of the technical principles, clinical applications, and risks associated with cryoablation procedures. Optional education is available from your Boston Scientific representative.

CAUTION: Carefully read all instructions prior to use. Failure to observe all warnings and precautions may result in complications.

Additional Required Items

Accessory Products Used to Conduct Cryoablation Procedures

NOTE: See product specific Instructions for Use.

The following needles must be used with the Visual-ICE Cryoablation System:

- **Boston Scientific Cryoablation Needles:** The cryoablation needles are specifically designed for use with the Boston Scientific Cryoablation Systems and are available in a range of configurations that produce various iceball sizes and shapes, allowing the clinician to match the needles to the desired ablation zone. The cryoablation needles are supplied sterile.
-

WARNING: Only use non-MRI needles with the system.

Optional Accessory Items:

- **Cryoablation Needle Channel ID Stickers:** Cryoablation Needle Channel ID Stickers are affixed to the needle tubing to provide easier identification of the needles during a cryoablation procedure. Contact Boston Scientific Technical Assistance Center to order Cryoablation Needle Channel ID Stickers.
- **Boston Scientific MTS:** The MTS contain four sensor locations along the distal shaft of the needle to monitor temperature near the target site and adjacent critical structures.
- **Boston Scientific Urethral Warming Set:** The Urethral Warming Set is the conduit that circulates warm saline through the urethra during a prostate cryoablation procedure.

The following items used with the Visual-ICE Cryoablation System are re-usable and should be cleaned and/or sterilized in accordance with the Instructions for Use that accompanies each product.

- **IV Pole and Bracket for use with a fluid warmer and pump system**

The following items are needed to conduct cryoablation procedures and are not available from Boston Scientific.

- **Argon gas cylinder(s)**
 - **Helium gas cylinder(s) if using helium for thawing**
-

NOTE: The argon gas must meet the purity requirements specified in the **SYSTEM SPECIFICATIONS** Section.

- **A fluid warmer and pump system** when using the Boston Scientific Urethral Warming Set

Boston Scientific recommends the use of a sterile drape (customer supplied) to cover the touch screen if the system will be operated by members of the sterile team.

Installation, Calibration, and Servicing

Only Boston Scientific or authorized personnel are permitted to conduct service and system preventive maintenance on the system. Preventive maintenance on the Visual-ICE Cryoablation System is required every two years. Scheduled Preventative Maintenance must be completed to maintain system performance and safety.

WARNING: Do not modify the Visual-ICE Cryoablation System in any way. Only authorized Boston Scientific personnel or Boston Scientific-trained authorized personnel are to service Visual-ICE Cryoablation Systems. If service is required, contact Boston Scientific Technical Assistance Center.

The Visual-ICE Cryoablation System will display a reminder on the screen approximately one month prior to when the preventive maintenance is due. If the reminder message appears and the preventive maintenance has not yet been scheduled, contact Boston Scientific Technical Assistance Center to schedule the service.

PREPARATION

System Operation

Table 6 depicts the order and steps of the Visual-ICE Cryoablation System set-up and treatment procedure. Each step is described in detail in this chapter.

Table 6. Cryoablation Procedure Flow

1	Set-up System	<ul style="list-style-type: none"> • Confirm availability of gas, needles, and accessories • Position Visual-ICE Cryoablation System and lock brake • Ensure Manual Vent Valve is <u>closed</u> and Argon Shutoff Valve is ON • Power Visual-ICE Cryoablation System ON • Login
2	Connect Gas Cylinders	<ul style="list-style-type: none"> • Connect helium gas, if applicable • Connect argon gas • Ensure safety cables are connected • OPEN valve on helium, then argon
3	Needle / MTS Testing	<ul style="list-style-type: none"> • Enter patient treatment information in Registration Button • Prepare sterile needles for testing • Connect needles to Visual-ICE Cryoablation System and lock channels • Perform Needle Integrity and Functionality Test • Test MTSs, if applicable
4	Conduct Cryoablation Procedure	<ul style="list-style-type: none"> • Insert needles into target tissue • Deliver the cryoablation treatment and remove needles • View and save reports, if desired • Depressurize the system and shutdown

Preparation for Use

Prior to using the Visual-ICE Cryoablation System, inspect the chassis, power cord, brakes, safety cables, gas supply lines, gas connections, and the monitor touch screen to ensure they are not damaged. If any of the components are damaged, contact Boston Scientific Technical Assistance Center.

CAUTION: Do not use the Visual-ICE Cryoablation System if any moisture or condensation is present on the surfaces of the system. Allow the system to dry completely for 12 hours prior to powering up the system. Powering up the system containing moisture or condensation could result in permanent damage to the electrical boards, causing the system to be inoperable.

Prior to starting a cryoablation procedure, set up the Visual-ICE Cryoablation System, connect the gas cylinders, and conduct functionality tests on each cryoablation needle and thermal sensor (reference the **Pre-Procedure Testing** Section).

System Set-up

1. Position the Visual-ICE Cryoablation System alongside the patient table. Ensure that the needle gas tubing has sufficient length to reach the patient. Ensure that the power switch and the power control knob (Figure 1 and Figure 2) are both easily accessible.

NOTE: Allocate space with adequate ventilation and free airflow. In order to guarantee proper ventilation, always keep the sides of the Console at least 0.5 m (20 in) from walls or other obstructions to airflow.

2. Lock the two front wheels using the Brake Pedal on the Visual-ICE Cryoablation System. If necessary, lock the two rear wheels using the individual brakes on each wheel.
3. Plug the power cord into a hospital grade power outlet (electrical mains outlet) with a ground connection. Boston Scientific recommends using a stable and uninterruptable power outlet.

NOTE: If the power supply to the Visual-ICE Cryoablation System is not stable or is noisy, the MTS temperature readings may be inaccurate.

WARNING: To avoid risk of electric shock, this equipment must only be connected to a hospital grade electrical power outlet with a protective earth.

WARNING: Do not touch the chassis of the Visual-ICE Cryoablation System while touching the patient to avoid the potential risk of shocking the patient if an inadvertent electrical fault exists.

OPTIONAL: When performing a prostate cryoablation procedure, set up a urethral warmer system by following directions provided in the Urethral Warming Set Instructions for Use.

4. Ensure that the power switch located on the rear of the system is in the ON position (Figure 2). This power switch should remain ON at all times. The Visual-ICE Cryoablation System will not turn ON if this power switch is in the OFF position.
5. Confirm that the Argon Shutoff Valve on the Visual-ICE Cryoablation System is in the Argon ON position. Turn it to the Argon ON position, if needed.
6. Confirm the Manual Vent Valve is fully closed (Figure 2). If needed, turn the knob clockwise until it is fully closed.
7. Raise the monitor to the UP position and adjust it to a comfortable viewing angle.

CAUTION: Use care when swiveling the touch screen monitor to avoid the potential for pinching fingers.

8. Turn the system ON using the power control knob located near the Needle Connection Panel (Figure 4). During boot-up, the system performs several diagnostic tests to verify that the hardware and software are operating properly. The system may produce a series of clicking sounds when the system performs self-diagnostic tests. This boot-up process is completed in approximately 45 seconds.


NOTE: If the system was improperly shutdown after the previous procedure, the boot-up process could take up to 2 minutes.

NOTE: It is important to power up the system prior to connecting the gas to the system. If the system is not powered up prior to connecting the gas, the diagnostic tests will not be performed by the software.

The diagnostic tests check:

- The correct version of firmware is running on the system
- Critical system components including the solenoid valves, internal power supplies, cooling fans, pressure transducers, and temperature measuring circuits.


If the system detects a failure of an individual channel, the channel is disabled and the needle type display window, bordered in red, notes the channel is defective. A message indicating the failure is displayed in the upper right-hand corner of the Navigation Tool Bar (see Screen 1).

Screen 1. Disabled Channel


A failure that prevents use of the system displays a message directing you to contact Boston Scientific Technical Assistance Center (reference the **Displayed Messages** Section).

If the Visual-ICE Cryoablation System software detects pressurized gas in the system and the gas supply is not connected, a message appears requesting you to vent the gas from the system.

Screen 2. Vent Gas Message

The *Login Screen* appears after the boot-up process is complete.


Screen 3. Login Screen

9. Enter your assigned login name and password using the on-screen virtual keyboard.

NOTE: The login name and password are not case-sensitive. Numbers are displayed when the keyboard is in upper case. To change text case, use the shift button on the virtual keyboard.

NOTE: If you leave the user interface idle for a pre-set time with no activity, the Visual-ICE Cryoablation System software requires you to re-enter your password to unlock the user interface (reference the **Configure Settings** Section).


Screen 4. Incorrect Login

Additional Login Options:


If you forgot your password, contact your System Administrator and request that the administrator log in, go to the *Manage Users Screen* and change your password.

Alternatively, press the **Forgot Password** button at the top of the *Login Screen* (Screen 3). A displayed message provides a challenge to be relayed to Boston Scientific Technical Assistance Center (Screen 5).

Screen 5. Reset Password Challenge


The Technical Assistance Center will provide a response to be entered on the screen using the virtual keyboard. Your password will be reset (Screen 6) and you will have an opportunity to change your password in the *Configuration Screen*.

Screen 6: Password Reset

In case of emergency, press the **Emergency Login** button at the top of the screen (Screen 3). A message is displayed with a challenge. Call Boston Scientific Technical Assistance Center to receive the proper response to enter, then press the **Login** button (Screen 7).

NOTE: This action does not reset your password.

Screen 7. Emergency Login

After successful login, the *Startup Screen* appears (Screen 8).

Screen 8. Startup Screen

Connecting Gas Cylinders

WARNING: Do not connect the Visual-ICE Cryoablation System to a gas supply exceeding 6000 psi (414 bar, 41.4 MPa) to avoid damage to internal system components.

WARNING: Ensure that the gas cylinders are chained to a wall or an approved cart to prevent inadvertent tipping of the cylinders.

WARNING: Have sufficient argon gas available to conduct the planned cryoablation procedure: the number and type of needles, gas cylinder size, pressure and rate of gas flow affect the required gas volume (reference the **External Gas Supply** Section for gas purity requirements). At least one full spare cylinder should be available for each treatment.

CAUTION: Turn the Visual-ICE Cryoablation System ON before connecting the gas cylinders to ensure the proper diagnostic tests are conducted.

NOTE: If cryoablation needles with i-Thaw Function are being used, only argon gas is required to be connected. If helium gas is connected, i-Thaw Function is disabled and only thawing with helium gas is available.

1. Position the gas cylinder(s) close enough to the Visual-ICE Cryoablation System to ensure that the gas supply line is not stretched and does not present a trip hazard.
 2. On the rear of the Visual-ICE Cryoablation System, ensure that the Manual Vent Valve is in the CLOSED position.
 3. Remove the humidity covers from the helium and argon inlets on the Visual-ICE Cryoablation System.
 4. Secure the safety cable at the end of the gas supply line to the system.
-

WARNING: Ensure that the safety cable is properly secured to the system in case the gas supply line becomes accidentally disconnected.

5. Connect the helium high-pressure gas supply line to the helium inlet of the Visual-ICE Cryoablation System using the quick-connect connector located at the rear of the system.


Figure 6. Visual-ICE Cryoablation System Gas Connections

- 1 Safety cable
 - 2 Quick-connect connectors
 - 3 Humidity cover
6. Route the helium gas supply line through the supply line clip on the System.
 7. Connect the helium high-pressure gas supply line to the helium cylinder by securing the gauge assembly adapter onto the cylinder connection (Figure 7).
-

NOTE: Gas cylinder connections have left-handed threads.


Figure 7. Gas Cylinder Set-up

- 1 Gauge assembly adapter 2 Cylinder valve

8. Carefully turn the cylinder valve on the helium gas cylinder counter-clockwise for one quarter turn. Ensure that the pressure reading on the gauge responds immediately. Turn the cylinder valve further counter-clockwise (approximately one full turn) to fully open the gas cylinder so that sufficient gas flow is present.
9. Repeat the procedure outlined in Steps 4 through 8 to connect the argon gas cylinder to the Visual-ICE Cryoablation System, using the argon gas supply line.

If no argon pressure is displayed on the system pressure gauge, ensure the Argon Shutoff Valve is in the Argon ON position.

OPTIONAL:

The EZ-Connect2 Dual Cylinder Adapter connects two cylinders of argon gas to the Visual ICE Cryoablation System to support a cryoablation procedure. A four-way adapter assembly with argon pressure gauge connects a gas supply line, the primary gas cylinder, and an auxiliary gas supply line.

If using the optional EZ-Connect2 Dual Cylinder Adapter, connect the gas supply line with four-way pressure gauge adapter assembly to the primary argon cylinder by securing the gauge assembly adapter onto the cylinder connection.

- Connect the end of the gas supply line to the argon inlet of the Visual-ICE Cryoablation System using the quick-connect connector.
- Connect the auxiliary gas supply line to the four-way adapter assembly using the quick-connect connector located on the end of the auxiliary gas supply line.
- Connect the opposite end of the auxiliary gas supply line to the second argon cylinder by securing the end of the auxiliary line to the cylinder connection.
- Open the primary cylinder valve first and use this cylinder until depleted. Do not open the cylinder valve on the second cylinder until the initial cylinder has been depleted.
- Reference the **Changing Gas Cylinders during a Procedure** Section for instructions on changing the gas cylinder during a procedure, if the second cylinder is also depleted during the procedure.



Figure 8. EZ-Connect2 Dual Cylinder Adapter

- 1 Auxiliary Gas Supply Line 2 Four-way Adapter Assembly with Pressure Gauge 3 Gas Supply Line

WARNING: Have sufficient argon gas available to conduct the planned cryoablation procedure: the number and type of needles, gas cylinder size, pressure and rate of gas flow affect the required gas volume (reference the **External Gas Supply** Section for gas purity requirements). At least one full spare cylinder should be available for each treatment.

10. Ensure the **Gas Indicator** (Screen 9) displays the minimum working pressure before initiating a procedure (Table 7). The **Gas Indicator** should indicate the pressure to be in the green range. If the system detects that the pressure reading for either gas cylinder is less than 50 psi (3.4 bar, 0.344 MPa), a message displays on the *Navigation Tool Bar* (Screen 9). Connect the gas cylinders to the Visual-ICE Cryoablation System.

Screen 9. No Gas Connected message

Table 7. Working Gas Pressures

Gas	Nominal Working Pressure	Working Pressure Limits
Argon	3500 psi 241 bar 24.1 MPa	3200 psi to 3800 psi 221 bar to 262 bar 22.1 MPa to 26.2 MPa
Helium	2200 psi 152 bar 15.2 MPa	1800 psi to 2500 psi 124 bar to 172 bar 12.4 MPa to 17.2 MPa

NOTE: When helium is not connected, active thawing can be conducted with an i-Thaw Function capable cryoablation needle. The system will flow argon gas at 500 psi (34 bar, 3.4 MPa) through the needle and will energize the internal heating element of the i-Thaw Function capable needle when the **Thaw** button is pressed.

CAUTION:

- When the gas cylinder pressure falls below the lower working pressure limit, the system displays an alert message in the *Navigation Tool Bar*. To ensure optimal performance, replace a gas cylinder if the pressure falls below the lower working pressure limit.
- Failure to operate the Visual-ICE Cryoablation System within the working pressure limits may affect the cryoablation procedure.
- If the system produces a continuous hissing sound, verify that the Manual Vent Valve is fully closed. If the Manual Vent Valve is fully closed and the hissing sounds persists, shut OFF the system using the power control knob located on the front of the system (Figure 1). Close the gas supplies using the cylinder valves. Contact Boston Scientific Technical Assistance Center.

Gas Line Flushing


After argon is connected and the gas lines pressurized, the Visual-ICE Cryoablation System automatically conducts gas line flushing to reduce potential gas line contaminants.

- When argon is connected to the system and the gas lines pressurized prior to pressing **Start Procedure**, the automated flushing will occur before the *Procedure Screen* is displayed.
- When the *Procedure Screen* is displayed before the argon is connected, automated flushing will occur after the connected gas pressurizes the lines.

During the flushing procedure, a sequence of gas flow, followed by an idle period, will occur three times.

- Gas will flow for 3 seconds.
- An idle period of 30 seconds will follow the gas flow.
- The automated flushing process will be completed in ~ 90 seconds.

As the system conducts gas line flushing, the screen displays *Optimizing Lines*.


Screen 10. Optimizing Lines message

Pre-Procedure Testing

WARNING: Prior to starting a cryoablation procedure, set up the Visual-ICE Cryoablation System and perform Needle Integrity and Functionality Tests on each cryoablation needle and thermal sensor.

1. On the touch screen monitor, press **Start Procedure**. The *Procedure Screen* appears (Screen 11).

Screen 11. Procedure Screen

2. Using aseptic technique, carefully remove the cryoablation needle from the package and place in a sterile work area.
3. Remove the connector cap, and then connect the needle to the Visual-ICE Cryoablation System needle connection panel (Figure 4).

WARNING: Do not kink, pinch, cut or pull excessively on the needle tubing. Damage to needle handle or tubing may cause the needle to become unusable.

4. After a needle(s) is inserted into the desired channel, lock the channel by sliding the locking bar away from the center of the system (Figure 9).


Figure 9. Locking Needle into Channel

- To allow for easier identification of the needle when multiple cryoablation needles are used during a cryoablation procedure, it is recommended to place a Needle Channel ID sticker on the needle tubing.

NOTE: Contact Boston Scientific Technical Assistance Center to order Needle Channel ID sticker for the cryoablation needles.

- Repeat Steps 2 through 5 for each cryoablation needle that will be tested.

CAUTION: Boston Scientific recommends that only needles of the same type be placed together in a single channel. Do not combine i-Thaw Function capable needles with non-i-Thaw Function capable needles in the same channel. Additionally, using needles of differing types in a channel may affect the accuracy of the **Gas Indicator**.


CAUTION: If the Visual-ICE Cryoablation System detects that helium gas is connected, the system defaults to helium thaw mode and i-Thaw Function, FastThaw Function, and cautery function capabilities are disabled.

NOTE: The Visual-ICE Cryoablation System displays a message if one or more needles without i-Thaw Function capability are connected and helium gas is not connected to the system. Ensure helium gas is connected to the system to use active thaw without i-Thaw Function.


Screen 12. Thawing Mode message

When the channel has been locked, the software detects that a needle has been connected and the channel opens for testing. A dark gray channel button indicates a channel with needles connected. If the needle contains a memory chip, the Visual-ICE Cryoablation System automatically detects the type of needle being used, the lot number, and the expiration date of the product. If the system detects that the expiration date has passed, a message appears and gas flow to that channel is disabled.

Screen 13. Needle Past Expiration Date message

- If the needle does not contain a memory chip, a menu appears with a selection of needle types (Screen 14). Select the proper needle type from the drop-down menu.

Screen 14. Select Needle Type Menu

- After selecting the first needle, the selection for subsequent needles defaults to the initial selection. Confirm needle type displayed in each channel matches the needle type connected.
- Press and hold the **Channel** button to open Advanced Channel Controls, which allows you to change the needle type for a channel, as needed.
- Prepare to conduct Needle Integrity and Functionality Testing.

WARNING: The sterile field and sterility of cryoablation needles should be maintained at all times. Do not contaminate the distal end of the sterile cryoablation needle. Avoid contact with the distal portion of the cryoablation needle to maintain sterility during testing.

- Secure the needle tubing to the sterile table prior to beginning the needle testing process.
 - Fill a large basin (at least 30 cm in diameter) halfway with sterile water or saline.
 - Place the needles, individually or in groups, in the basin such that the full length of the needle shaft is submerged in the sterile water or saline.
11. Perform the Needle Integrity and Functionality Test on each needle by pressing the **Test** button on the channel containing the needle(s). The 90 second test automatically conducts a series of flush, freeze, and thaw phases. The durations of these phases are: 45 seconds flush with helium, 15 seconds freeze with argon, and 30 seconds thaw with helium.

OPTIONAL: Alternatively, all needles may be tested simultaneously by pressing the **Test** button on the channel labeled **ALL**. A message requests confirmation to test all needles. If appropriate, choose YES.

OPTIONAL: If additional testing is required, press the **Test** button again to repeat the test.

NOTE: When helium is not connected, the two-minute test consists of 50 seconds low-pressure argon flow, 15 seconds freeze with high-pressure argon, and 55 seconds of low-pressure argon flow.

NOTE: If a previously tested needle is moved to a new channel at any time during a procedure, the Needle Integrity and Functionality Test must be conducted again on that needle.

NOTE: If a previously tested needle containing a memory chip is moved to a new channel during a procedure, the Visual-ICE Cryoablation System will recognize the needle has successfully completed the Needle Integrity and Functionality Test.

During testing, carefully watch each needle for the following:

Flush: Ensure no bubbles are formed along the needle shaft and tip. Ensure that no iceball is formed during the Flush phase.

WARNING: A defective cryoablation needle that has a gas leak can cause a gas embolism in the patient. Never use a defective needle for a cryoablation procedure. Return defective needles to Boston Scientific for evaluation.

CAUTION: Ice formation during the Flush phase indicates that argon gas is connected to the helium inlet. Before proceeding change the cylinders and ensure that each gas supply line is connected to the correct cylinder (reference the **Standard Gas Cylinder Setup** Section).

Freeze: Ensure that ice begins to form around the tip of the needle.

WARNING: A needle is defective if there is no ice formation during the Freeze phase. Do not use a defective needle. Obtain a new needle and repeat the testing procedure.

Thaw: Ensure that the iceball detaches from the tip of the needle and that no bubbles escape from the needle tip.

CAUTION: Ice formation during the Thaw phase indicates that the argon gas is connected to the helium inlet. Before proceeding change the cylinders and ensure that each cylinder is connected to the correct inlet (reference the **Standard Gas Cylinder Setup** Section).

During the Needle Integrity and Functionality Test, the Gas Indicators for both gases provide estimates of the time remaining before depletion of the cylinders assuming that all connected needles are operated simultaneously (reference the **Navigation Tool Bar** Section).

Upon successful completion of the Needle Integrity and Functionality Test, the **Test** button displays a green check and remaining control buttons on the channel become active. The needle(s) are ready for use.

12. When using MTSs, prepare the desired number of sensors for testing. The Visual-ICE Cryoablation System supports use of four thermal sensors. Connect each thermal sensor to an MTS connection port by pushing the connector into the port. Ensure that the rectangular tongue on the MTS connector is aligned with the groove at the top of the connection port for proper MTS insertion (Figure 10).

Figure 10. MTS Connection

13. Confirm MTS functionality by verifying that the connected sensor produces a reasonable onscreen temperature reading (e.g., close to room temperature) in the **Temperature Sensor** Section of the software screen.

WARNING: Do not use an MTS that has failed the MTS Integrity and Functionality Test; erroneous temperature measurements may be displayed.

CAUTION: Do not allow the MTS to touch the cryoablation needle during use; erratic temperature measurements may be displayed.

Navigating The User Interface


Throughout the User Manual typeset conventions represent various user interface sections, software buttons, positions, and steps.

- *Software Screen* section
- **Control** button
- ON position
- *OPTIONAL* = optional or alternate step

The Visual-ICE Cryoablation System provides a graphical user interface that facilitates rapid communication between the user and the system through a touch screen interface.

Login Screen

When the system is powered on, the *Login Screen* appears after the boot-up process is complete (reference the **System Set-up** Section).

Screen 15. Login Screen

Startup Screen

After logging into the system, the *Startup Screen* displays several options.

Screen 16. Startup Screen

Table 8. Startup Screen Buttons

Button	Description
Start Procedure	Transition to the <i>Procedure Screen</i> to start a cryoablation procedure.
Logout	Log out of the system.
View Reports	View the contents of a report and export reports to a USB flash drive. NOTE: Administrative users may also delete reports.
Configure Settings	Configure various system settings (reference the Configure Settings Section). NOTE: Some configuration parameters are limited to administrative and/or service users only.
User Manual	View information on how to access an electronic version of the User Manual.
Service	Field Service personnel login to modify configuration settings, perform and record preventive maintenance activities. NOTE: This option is only available to authorized field service personnel.

Procedure Screen

The Visual-ICE Cryoablation System *Procedure Screen* provides a single-screen view to control and monitor a cryoablation procedure. The *Procedure Screen* is divided into sections for the Navigation Tool Bar, Context Sensitive Help, Channel Controls, Channel Status, Temperature Sensors, and Organ Map. The title bar of each section on the *Procedure Screen* provides user selected self-help for that section.

Screen 17. Procedure Screen

- | | | | | | |
|---|------------------------|---|-----------------|---|---------------------|
| 1 | Navigation Tool Bar | 3 | Channel Control | 5 | Temperature Sensors |
| 2 | Context Sensitive Help | 4 | Channel Status | 6 | Organ Map |

Navigation Tool Bar

Screen 18. Navigation Tool Bar

The Navigation Tool Bar contains the Pressure Gauge / Gas Indicator and procedure buttons that you can select to enter registration information, configure procedure settings, enter procedure notes, view and export reports, and end the procedure. On occasion, error messages may be displayed in lieu of the logo.

Table 9. Navigation Tool Bar

Button	Description
Pressure Gauge / Gas Indicator	Displays the working pressure of the argon and helium gases within the system. NOTE: The Visual-ICE Cryoablation System contains internal regulators that regulate the gas pressure to appropriate working limits. The pressure displayed on the Gas Indicator is the internal, regulated pressure and not the cylinder gas pressure. Pressing the Pressure Gauge toggles the pressure gauge to display the estimated procedure time remaining before the gas cylinders are depleted. The estimated times are displayed in hours:minutes:seconds. During needle testing, both pressure gauges display the estimated time remaining. The initial estimates during needle testing are based on the assumption that all connected needles are operating simultaneously at 100 % freeze intensity. The Gas Indicator updates in real time as needles are disconnected or additional needles are connected and as the freeze intensity is adjusted. Pressing the Gas Indicator toggles the display back to Pressure Gauge.
Registration	Provides optional data entry fields to record Patient ID, Hospital Name, Hospital Address, Physician Name, and Organ Type. Two custom fields are provided for additional information. The names of the custom fields can be specified in the <i>Configure Settings Screen</i> (reference the Configure Settings Section)
Notes	A location to input text. Selecting this button displays the on-screen keyboard for data entry. Procedure notes entered in this location are included in the procedure report (reference <i>Configure Settings Screen</i> (reference the Startup Screen Section).
Settings and Procedure	Displays the <i>Procedure Screen</i> to begin a cryoablation procedure.
Report	Displays a report of all procedure data that has been entered and captured for the current procedure. The report can be saved to the USB flash drive. Pressing the Report button during a procedure displays all procedure information saved up to that point in time.
End Procedure	Ends the current procedure and returns to the <i>Startup Screen</i> . Pressing this button generates a confirmation request, a request to save the report, and an option to automatically vent the system.

Context Sensitive Help

The *Context Sensitive Help* displays a summary of procedural steps to guide you through a cryoablation procedure. It is intended to provide only a general overview; step-by-step instructions are provided in the **Preparation for Use** Section.

User-Selected Self Help

The title bar of each section provides access to additional help information. Press the title bar to access an explanation of the buttons and fields available in each section of the *Procedure Screen*.

Channel Controls


Channels 1 through 10 are labeled individually and contain independent controls for **Test, Freeze, Freeze Intensity, Thaw, and Stop**. Each individual channel displays the needle type of the needles connected adjacent to the channel controls (Screen 19). The Channel labeled **ALL** operates all active channels simultaneously.

Table 10. Channel Controls

Button	Description
	Channel button – Identifies the active channel(s). > <i>Advanced Channel Controls</i> : Pressing and holding the Channel button provides options to change the needle type selected for that channel, link two adjacent channels together for simultaneous operation, and program freeze – thaw cycles.
	Channel labeled ALL – Allows testing, freezing, and thawing on ALL active channels simultaneously. Press the button for the desired function (Test , Freeze , or Thaw) on this channel to activate that function on all needles simultaneously.
	Test button – Begins the Needle Integrity and Functionality Testing that is required prior to using any cryoablation needle. No other controls are enabled until needle testing is complete.
	Tested button - Following completion of Needle Integrity and Functionality Testing, the button displays a check mark and the remaining control buttons on the channel become active.
	Freeze button – Begins a Freeze phase at the selected freeze intensity.
	Freeze Intensity drop-down menu – Provides an option to adjust the freeze intensity from 100 % to 5 % or to select “Stick” intensity. NOTE: The Visual-ICE Cryoablation System controls freeze intensity by adjusting the duration of argon flow over each 10 second block of time (e.g., 30 % freeze intensity freezes for 3 seconds and idles for 7 seconds).
	
	Thaw button – Begins a Thaw phase. > <i>Advanced Thaw Controls</i> : Pressing and holding the Thaw button provides an option for enabling FastThaw and an option for enabling track ablation via the Cautery Function button. The FastThaw Function and cautery function are available only on CX-type cryoablation needles.
	FastThaw Function button – Begins a FastThaw phase. > <i>Advanced Thaw Controls</i> : Pressing and holding the FastThaw Function button provides an option to change from FastThaw Function to i-Thaw Function. The FastThaw Function is available only on CX-type cryoablation needles.
	Stop button – Stops any activity.

Channel Status

Channel Status shows the status of each freeze, thaw, and idle phase with numeric and color-coded displays on the progress indicator. Variations in blue shading visually represent the selected freeze intensity. The **Timer** button to the right of the progress indicator displays the elapsed time of the current phase.

Screen 19. Channel Controls and Channel Status Section

Variations in orange shading visually distinguish active thaw (helium or i-Thaw Function), FastThaw Function and cautery function phases (Screen 20).

Screen 20. Active Thaw Channel Status Section

For i-Thaw Function and CX-type needles, the *Channel Status* also displays the internal gas temperature at the needle tip during freeze phase and the estimated temperature range for the needle shaft during active thaw phase. The temperature display is updated every 2 seconds.

NOTE: During the needle warming phase for CX-type needles, the channel status will display a rotating warming indicator.


Screen 21. Needle Tip Temperature Display

Enlarge and Reposition Timers

During a needle test, freeze, thaw, or idle phase, press the **Timer** button to enlarge the timer display (Screen 22). The enlarged Timer displays the channel number in the upper left corner of the timer window, elapsed time, and, when freezing, the selected freeze intensity. When CX-type needles are connected, the enlarged timer also displays the internal needle tip temperature during the freeze phase and the estimated temperature range of the needle shaft during the thaw phase.

NOTE: During the needle warming phase for CX-type needles, the timer will display a rotating warming indicator.

Timers for three selected channels can be simultaneously enlarged. Press the timer to return it to its original size.

Screen 22. Enlarged Timer

Change the location of an enlarged timer by dragging the timer to a new screen location.

Screen 23. Repositioned Enlarged Timers


To briefly display the time associated with a completed cycle, press the section of the status bar for the selected operation. Press the **Scale** button to adjust the graphical display of the channel status such that all operations are visible. Press the **Scroll** button to adjust the graphical display to 5-minute increments; the display scrolls throughout the procedure. Press the **Maximize** button (+) to enlarge the graphical display. Press the **Minimize** button (-) to reduce the display to the original size

Temperature Sensors

The **Temperature Sensors** Section (Screen 24) displays the temperature measured for each sensor location on the connected MTS needles and a black triangle highlights the coldest temperature on each MTS.

Temperature vs. Time can be plotted for any four selected sensor locations from the connected MTS needles. Press the button corresponding to the selected sensor location on an MTS to select that sensor location to be graphed. Press the respective **MTS Channel** button to graphically display all four sensor locations of an individual MTS. The color of the selected temperature sensor location corresponds to the color of the temperature line on the graph. Press the **Thermometer** button to graph the coldest temperature on a selected MTS. Press the **Thermometer** button again to return to the original display.

NOTE: A maximum of four sensor locations can be displayed on the graph at any one time. The four selected sensor locations can be from any of the connected MTS needles.

Screen 24. Temperature Sensors Section

- | | | | | | |
|---|---------------------------|---|---------------------------|---|--|
| 1 | MTS sensor locations | 3 | Thermometer button | 5 | Button corresponding to 5 mm sensor location |
| 2 | MTS Channel button | 4 | Maximize button | 6 | Triangle indicates coldest sensor |


Select the **Scale** button to adjust the graphical display of the temperature measurements such that the temperature throughout the entire procedure is visible. Select the **Scroll** button to display a scroll bar to scroll through the graphical data during the procedure.

To enlarge the *Temperature Sensors* graphical display, press the **Maximize** button (+). Press the **Minimize** button (-) to return the display to the original size.

Advanced Temperature Sensors Controls: Pressing and holding an **MTS Channel** button provides an option to set visual alarms when a chosen temperature measurement falls below a desired value or when the rate of temperature drop is greater than a desired value. This advanced control also provides an option to label the MTS location (reference the **Advanced Temperature Sensor Controls** Section).

Organ Map

The *Organ Map* may be used to provide a pictorial representation of tumor location, needle placement and cycle activity for a selected image (kidney, prostate, liver, or lung) or for a drawn image (other). Tumor shapes may be drawn within the displayed organ, and needles placed, moved, or removed within the organ map. A blank drawing canvas is available for customized display. Needles may be controlled from the organ map by using the **Control Needle** button. Needle activity status (Freeze, Thaw, or Idle) is represented by color coding.

Screen 25. Organ Map Section

Table 11. Organ Map Controls

Button	Description
	Select Organ button – Select an organ map (kidney, prostate, liver, lung, or other) to be displayed on the drawing area
	Draw Line button – Draw free-form on the organ map or drawing area. Line thickness and line color are selectable.
	Draw Circle button – Draw a circle of various sizes on the organ map or drawing area. Line thickness and line color are selectable.
	Move Circle button – Move a circle to another area on the organ map. This button appears only when a circle has been selected.
	Erase button – Erase lines previously added to the drawing space. The size of the eraser is selectable.
	Place Needle button – Select and place needles and MTS on the organ map. Needles are selectable and identified by channel and needle port (e.g., 1a, 1b, 2a, etc.). MTS are selectable by MTS port. Drag the needle or MTS to the desired location on the Organ Map to place it. Place each needle on the Organ Map as it is inserted into the target tissue to avoid potential mix-up of needle locations. The status of the needle(s) placed on the organ map displays the state of the needle (Freeze, Thaw, Idle) by color coding, matching that identified in the Channel Status Section.
	Control Needle button – Initiate freeze and thaw cycles of the selected needles.
	Needle Lock button – Lock a needle in place on the organ map so that it may not be inadvertently moved.
	Clear All button – Clear all drawings and needle placements from the drawing canvas.
	Trash button – Delete a needle on the organ map.
	Previous Image button – Change the view of the selected organ.
	Next Image button – Change the view of the selected organ.

To view a report that has been saved to the Visual-ICE Cryoablation System, press the **View Reports** button from the *Startup Screen* (Screen 16).

The *View Reports Screen* displays a list of all saved procedure reports in the Visual-ICE Cryoablation System (Screen 27). You may select a report to be viewed or exported, or you may delete your own reports. Users with an administrative login ID may delete any report.

The screenshot shows the 'View Reports' screen of the Visual-ICE Cryoablation System. At the top right, the time is 6:40 am and the date is 22 Aug 2021. A 'Back' button is located at the top left. Below the title bar, the text 'Startup >> View Reports' is displayed. The main area contains a table with the following data:

	Report Name	Hospital Name	Doctor Name	Procedure Type
1	2021_Aug_04_5_16_am	No information entered	No information entered	Other
2	2021_Jul_16_9_30_am	No information entered	No information entered	Other
3	2021_Jul_18_8_15_am_Tim	ST. Mary's	Dr. Thomas	Lung
4	2021_Jul_18_9_09_am_tim	ST. Mary's	Dr. Jones	Prostate
5	2021_Jul_18_9_24_am_Tim	ST. Mary's	Dr. Thomas	Kidney

Below the table, the selected report '2021_Jul_18_8_15_am_Tim' is displayed in detail. The report title is '2021_Jul_18_8_15_am_Tim'. On the left, there is a 'Report Sections' menu with options: Registration (selected), Procedure, Temperature Sensors, Organ Map, and Notes. Below the menu are 'Export Report' and 'Delete Report' buttons. The main content area shows the 'Lung Procedure' title, the time '8:15 am 18 Jul 2021', and the 'Registration' section with the following information:


- Patient ID: JohnKI-9001
- Hospital Name: ST. Mary's
- Hospital Address: 500 Peach St
- Doctor Name: Dr. Thomas

The 'Procedure' section is partially visible at the bottom.

Screen 27. View Reports Screen

To sort the list by Report Name, Hospital Name, Physician Name, or Procedure Type, press the appropriate header section of the report list.

The **Export Report** button displays a window to choose the Export File Type and File Name for exporting the report. Reports can be exported in HTML, PDF, or CSV formats. An additional option is available to export raw MTS data for further analysis.

Screen 28. Export Report Screen

Elä...
Dit is e...
Utdatert v...
Wersja przete...
Versão obsoleta.
Versiune expirată. A...
Zastaraná verzia. Nepou...
Zastarela različica. Ne upou...
Vanhentunut versio. Älä käytä.
Föråldrad version. Använd ej.
Güncel olmayan sürüm. Kullanmayın.

Configuring Settings

The *Configure Settings Screen* allows selection of settings used during a cryoablation procedure. Settings that may be changed include MTS Needle Locations, System, Procedure and Registration Settings and Units (reference the **Configure Settings** Section).

Control buttons have options to Manage Users and Manual Software Update (reference the **Configure Settings** Section). The Manual Software Update buttons are available only to system administrators and service personnel. Only Service personnel have the option to adjust the system time and date.

Screen 29. Configure Settings

Table 12. Configure Settings Options

Button	Description
Manage Users	Change your password. Administrative users may add users, remove users, or change the password of any user.
Manual Software Update	Install a software update via a USB flash drive. NOTE: This function is available only to administrative and service users.

Service Screen

The *Service Screen* is available only to Boston Scientific trained, authorized service personnel with a service login ID. The *Service Screen* provides service users with the ability to run system diagnostics, enable or disable system features, adjust minimum and maximum gas pressures, view event logs, and perform manual system configuration.

PROCEDURE

Performing a Cryoablation Procedure

WARNING: Do not touch the screen if the touch-screen monitor goes blank for more than five (5) seconds during a procedure. Immediately turn off power to the system and end the procedure to avoid inadvertent activation of needles.

1. **OPTIONAL:** On the *Procedure Screen*, select the **Registration** button to enter optional patient treatment information. Use your finger to enter information on the virtual keyboard. Available data entry fields are Patient ID, Hospital Name, Hospital Address, Physician Name, and Organ Type. If you need to enter any other registration information, two custom fields can be labeled in the *Configure Settings Screen* (reference the **Configure Settings** Section).

NOTE: The Visual-ICE Cryoablation System automatically loads the appropriate Organ Map based on the Organ Type selected on the **Registration** button.

CAUTION: Select a unique Patient ID that does not reveal the patient's identity to other system users.

2. **OPTIONAL:** Select the **Notes** button on the *Procedure Screen* to enter additional procedure notes. Notes can be entered at any time during the cryoablation procedure.
3. Position the cryoablation needles and thermal sensors in the target tissue.

CAUTION: During use, avoid damage to the needle from other surgical instruments.

OPTIONAL: Reference the **Organ Map** Section for instructions on placing needles and thermal sensors on the Organ Map.

WARNING: Use image guidance to verify that the cryoablation needles are placed at the desired location prior to activating a needle.

4. Select the desired Freeze Intensity using the drop-down menu.

NOTE: Throughout the procedure, monitor the amount of gas time remaining in the cylinders using the **Gas Indicator** on the Navigation Tool Bar (Screen 30). If it is necessary to change the gas cylinders during a procedure, follow the instructions provided in the **Changing Gas Cylinders during a Procedure** Section.

Screen 30. Gas Time Remaining

5. Press the **Freeze** button on selected channels containing needles to begin the initial freeze phase of the procedure. To adjust the freeze intensity, press the **Freeze Intensity** button and select the desired intensity from the drop-down menu. The Freeze cycle will continue at the selected freeze level until this operation is changed or stopped.

WARNING: Continuously monitor iceball formation using image guidance such as direct visualization, ultrasound, or Computed Tomography (CT) to ensure adequate tissue coverage and to avoid damage to adjacent structures.

OPTIONAL: To initiate a freeze phase on all needles simultaneously, press the **Freeze** button on the channel labeled ALL. Pressing any function button on the channel labeled ALL displays a message asking you to confirm simultaneous operation of all needles.

NOTE: Selecting **ALL** will initiate a Freeze phase at the intensity selected for each channel. To conduct freezing in all active channels at the same intensity, select the intensity in the **ALL** channel before pressing the **Freeze** button.

-
6. Watch the timer to monitor the freeze phase elapsed time (reference the **Channel Status** Section for instructions on enlarging the timer display). When the desired freeze duration has elapsed, press the **Stop** button to enter an idle phase.
 7. If thermal sensors are used, monitor the tissue temperature using the **Temperature Sensors** Section of the *Procedure Screen* (reference the **Temperature Sensors** Section).
 8. To actively thaw the iceball, press the **Thaw** button on the channels containing needles to start the thaw phase. If the Visual-ICE Cryoablation System detects that helium gas is connected, the system defaults to helium thaw mode. If needles of mixed type (some helium thaw and some i-Thaw Function) are connected, the system also defaults to helium thaw mode.

NOTE: When thawing with CX-type needles, thawing is limited to a maximum of 7 needles simultaneously activated. Thawing with the FastThaw Function is limited to a maximum of 4 needles simultaneously (reference the **i-Thaw Function and FastThaw Function Control for CX-type Needles** Section for instructions on using i-Thaw Function and FastThaw Function).

OPTIONAL: To initiate a thaw phase on all needles simultaneously, press the **Thaw** button on the channel labeled ALL. Pressing any function button on channel labeled **ALL** displays a message asking you to confirm simultaneous operation of all needles.

-
9. Watch the timer to monitor the thaw phase elapsed time (reference the **Cycle Programming Control** Section for instructions on performing a timed thaw phase). When the desired thaw duration has elapsed, press the **Stop** button to enter an idle phase.
 10. Repeat steps 4 through 9 until the desired number of freeze-thaw cycles has been completed.

WARNING: Ensure adequate thawing or cooling before attempting to remove needles from the patient.

-
11. Remove all needles and MTS from the patient.
 12. Unlock the locking bar(s) and remove all needles and MTS from the needle connection panel.
 13. Dispose of used needles and MTS in a biohazard container in accordance with hospital and safety regulations.
 14. When the procedure is complete, press the **End Procedure** button on the *Procedure Screen*. Three messages appear requesting action:
 - Confirmation to end the procedure – Press the **Yes** button to end the procedure.
 - Request to save a report – Press the **Yes** button to save a report.
 - Request to automatically vent the high-pressure gas – Press the **Yes** button to automatically vent the system. The system prompts you to close the gas supplies before venting. The automatic venting takes approximately 1.5 minutes. Before initiating automatic venting, warn others in the vicinity to expect a venting noise.

WARNING: If the needles are still connected, do not unlock the channels or disconnect the needles from the needle connection panel until all operations in the channel have completed.


-
15. If you are ready to shut down the system, reference the **Shutting Down the System** Section for the system shutdown procedure.

Reports

At any time during a procedure, press the **Report** button on the *Procedure Screen* to view a summary of the report information that has been saved up to that point.

At the end of a cryoablation procedure, a report summarizing the complete procedure can be saved to the system and exported for use on a personal computer.


1. Press the **Report** button on the *Procedure Screen*.
2. When viewing a report, you may scroll through the report using the scroll bar on the right-hand side of the screen or you may select a section to view by pressing the report section name on the left-hand side of the screen.
3. Press the **Export Report** button to save the report to the USB flash drive. A window appears for you to select the file format and file name. Enter the file name using the onscreen virtual keyboard. You may also choose to export the temperature sensor data in a file to be used for further analysis.

Screen 31. Export Report Screen

CAUTION: Use only a Boston Scientific supplied USB flash drive with the Visual-ICE Cryoablation System. Do not use that flash drive for purposes unrelated to Visual-ICE Cryoablation System data and reports.

4. Press the **Export** button to begin exporting the file. Wait for confirmation before removing the USB flash drive from the system.

Screen 32. Report Exported message

System Shutdown

WARNING: Warn operating room personnel prior to venting the Visual-ICE Cryoablation System to avoid startling them.

1. If you did not choose to automatically vent the Visual-ICE Cryoablation System, turn the shutoff valve on the gas cylinders clockwise to close the gas cylinders.
2. Turn the Manual Vent Valve on the Visual-ICE Cryoablation System to the OPEN position to release the high-pressure gas from the system.
3. Disconnect the high-pressure gas supply lines from the Visual-ICE Cryoablation System and from the gas cylinders. Store the gas supply lines and gauge assemblies in the storage compartment provided on the system (Figure 1).

WARNING: If it is difficult to loosen the pressure gauge connected to the cylinder, or the high-pressure gas supply line(s) cannot be disconnected from the inlet connections, do not apply excessive force to release the gas supply line or to loosen the pressure gauge. The gas line may still be under pressure.

4. Turn the Manual Vent Valve to the CLOSED position after the gas has been vented.
5. Press the **Logout** button on the *Startup Screen* to log out of the system.
6. Press the **Shutdown** button on the *Login Screen* to power down the system. A message appears requesting confirmation to shut the system down.
7. Wait until the screen goes black. Turn the Power Control Knob to the OFF position.
8. Unplug the Visual-ICE Cryoablation System and wrap the power cord around the cord wrap on the rear of the system.

WARNING: Do not pull on the power cord. Grasp the plug, not the power cord, to disconnect the device from the wall socket.

9. Cover the helium and argon inlets with the humidity plugs.
10. Clean the system after each use using the instructions provided in the **Cleaning** Section. Ensure the system is dry prior to storage.
11. Lower the Touch Screen Monitor into the monitor storage basin prior to system storage.

CAUTION: Prior to lowering the monitor, ensure that no objects, such as the USB flash drive, are located in the monitor storage basin. Use care when lowering the monitor into the monitor storage basin; do not use excessive force to avoid damaging the monitor.

CAUTION: Use care when lowering the touch screen monitor to avoid the potential for pinching fingers.

12. Cover the Visual-ICE Cryoablation System with the Console cover.

Changing Gas Cylinders during a Procedure

If it becomes necessary to replace a gas cylinder during a procedure, stop all freezing and thawing operations.

Standard Gas Cylinder Setup

1. Plan the appropriate time to make a cylinder change by estimating the gas required to complete the procedure. The **Gas Indicator** on the Navigation Tool Bar indicates how much time is remaining in each gas cylinder based on the selected gas flow intensity, type and number of needles being used. Also, consider the number of planned freeze-thaw cycles for the procedure.
2. Safely position a full gas cylinder, with the required gas type and purity, near the empty cylinder.
3. Close and tighten the cylinder valves of both gas cylinders.
4. Slowly open the Manual Vent Valve to discharge gas from the system and the high-pressure gas supply line. Wait until all the pressure is released and both gauges on the gas supply lines indicate zero pressure.
5. Use the wrench to remove the gauge assembly from the empty cylinder.
6. Connect the gauge assembly to the full cylinder.
7. Close and tighten the Manual Vent Valve.
8. Carefully turn the cylinder valve on the helium gas cylinder counterclockwise for one quarter turn. Ensure that the pressure reading on the gauge responds immediately. Turn the cylinder valve further counterclockwise to open the gas cylinder so that sufficient gas flow is present.
9. Carefully turn the cylinder valve on the argon gas cylinder counterclockwise for one quarter turn. Ensure that the pressure reading on the gauge responds immediately. Turn the cylinder valve further counterclockwise to open the gas cylinder so that sufficient gas flow is present. If no argon pressure is displayed on the Gas Indicator, ensure the Argon Shutoff Valve is OPEN.
10. Continue the cryoablation procedure at the next planned freeze or thaw phase.

Disconnecting a Helium Cylinder

If helium is connected, i-Thaw Function, FastThaw Function and cautery functions are disabled. To use i-Thaw Function, FastThaw Function, and cautery functions, helium must be disconnected and the lines cleared of helium.

1. Close both the argon and helium cylinder valves.
2. Open the Manual Vent Valve to discharge gas from the system and the high-pressure gas supply lines. Wait until all the pressure is released and both pressure gauges on the Navigation Tool Bar indicate zero pressure.
3. Close the Manual Vent Valve.
4. Slowly open the argon cylinder valve one quarter to one half of a turn. Allow pressure to build on the argon pressure gauge. Open the argon cylinder valve fully for sufficient argon flow.

Dual Gas Cylinder Connection

1. Safely position a full argon gas cylinder, with the required purity, near the empty cylinder.
2. Close and tighten the cylinder valve of the empty gas cylinder.
3. Open the Manual Vent Valve to discharge gas from the system and the high-pressure gas supply lines. Wait until all the pressure is released and the pressure gauge(s) on the Navigation Tool Bar indicate zero pressure.

4. Close the Manual Vent Valve.
5. Connect the auxiliary gas supply line to the EZ-Connect2 Dual Cylinder Adapter using the quick-connect connectors.
6. Connect the opposite end of the auxiliary gas supply line to the new cylinder.
7. Carefully turn the cylinder valve on the new gas cylinder counterclockwise for one quarter turn. Ensure that the pressure reading on the gauge responds immediately. Turn the cylinder valve further counterclockwise to open the gas cylinder so that sufficient gas flow is present.

Advanced Thaw Controls

The Visual-ICE Cryoablation System provides options to select helium-free thawing (i-Thaw Function and FastThaw Function) and to conduct track ablation (**cautery**).

NOTE: i-Thaw Function, FastThaw Function and cautery functions are only available when needles with that capability are connected.

NOTE: Argon gas only must be connected to access these functions. Connecting helium gas to the Visual-ICE Cryoablation System disables these functions.

i-Thaw Function and FastThaw Function Control for CX-type Needles

CAUTION: Active thawing produces heat along the distal needle shaft. Use care to avoid thermal injury to non-targeted tissues.

WARNING: The needle handle may become warm during active thawing. Pay attention to the needle handle position. Prolonged contact with warm portions of the needle handle could cause unintended thermal tissue damage/burn to the patient or clinician.

1. Press and hold the **Thaw** button to access *Advanced Thaw Controls* (Screen 33). The *Advanced Thaw Controls* menu provides access to i-Thaw Function / FastThaw Function and cautery functions.
 - When i-Thaw Function is operational, pressing the **Change to FastThaw** Function button in the *Advanced Thaw Controls* will change the thaw from i-Thaw Function to FastThaw Function (Screen 33).

Screen 33. Advanced Thaw Controls FastThaw Function

- Similarly, when **FastThaw** Function is selected, pressing the **Change to i-Thaw** Function button will change the thaw from FastThaw Function to i-Thaw Function (Screen 34).


Screen 34. Advanced Thaw Controls i-Thaw Function

- After selecting the type of thaw, press the **Thaw** or **FastThaw** Function button in the Channel Control to initiate thawing.

NOTE: The FastThaw Function generates a temperature that is higher than that generated for i-Thaw Function, resulting in a faster thaw time.

NOTE: A maximum of 7 needles can thaw simultaneously using the i-Thaw Function; a maximum of 4 needles can thaw simultaneously using the FastThaw Function.


- During the needle warming phase for CX-type needles, the channel status will display a rotating warming indicator (Screen 35).

Screen 35. Advanced Thaw - Needle Warming

- When the threshold temperature for i-Thaw Function or FastThaw Function is reached, the channel status will display the estimated temperature range for the needle shaft (Screen 36).

NOTE: The shaft temperature is displayed as a temperature range because tissue and procedural variables will affect the temperature.

Screen 36. Advanced Thaw - Temperature Display During Thawing

- Thaw thoroughly and stop all needle operation prior to removing needles to minimize the risk of tissue injury.
 - If needle sticking is experienced, use a slight, gentle twist of the needle followed by slow withdrawal.

Cautery Control for Track Ablation

Needle track ablation is achieved using the cautery function in the *Advanced Thaw Controls*. The needle track is ablated through the use of thermal energy. Each cautery function enabled needle (1.5 CX type needles and 2.1 CX-type needles) is operated independently using the **Cautery** function button in the *Advanced Thaw Controls* for the channel and port in which the needle is connected.

WARNING: Use image guidance to monitor needle location and position to avoid the risk of thermal injury/burn to adjacent tissue/organs.

WARNING: Track ablation produces heat along the distal needle shaft. Use care to avoid thermal injury/burn to adjacent tissue/organs.

WARNING: The needle handle may become warm during track ablation. Pay attention to the needle handle position. Prolonged contact with warm portions of the needle handle could cause unintended thermal injury/burn to the patient or clinician.


CAUTION: The cautery function should not be initiated if the Active Zone Indicator is visible outside the patient's skin.

NOTE: While the Visual-ICE Cryoablation System is in cautery function mode, freezing and thawing are not allowed on any other channels.

NOTE: The cautery function is not available on the channel **ALL**.

Cautery Control for 1.5 CX Needles

- Press and hold the **Thaw** button to access Advanced Thaw Controls.
- Press the **Cautery** function button on the *Advanced Thaw Controls* menu to open the *Cautery Control Screen* (Screen 33). A confirmation message appears (Screen 37).


Screen 37. Cautery Confirmation message

- Press the desired needle port (A or B) on each channel containing a needle you wish to track ablate. Only one needle per channel may activate the cautery function at a given time.

NOTE: The Visual-ICE Cryoablation System uses a pre-set duration of 30 seconds for each cautery phase for 1.5 CX-type needles. This value is not adjustable, but the duration can be stopped in advance of 30 seconds by pressing the **Stop** button.

NOTE: A maximum of four (4) needles may use the cautery function simultaneously.


- Press the **Start** button to initiate track ablation (cautery).
 - During the needle warming phase, the *Cautery Control Screen* will display a rotating warming indicator (Screen 38).

Screen 38. 1.5 CX Needle Warming

- When the cautery threshold temperature is reached, Cauterization automatically begins.
 - The timer will display Cauterizing (Screen 39).
 - The *Cautery Control Screen* will display the estimated temperature range for the needle shaft.

NOTE: The shaft temperature is displayed as a temperature range because tissue and procedural variables affect the temperature.

Screen 39. 1.5 CX Needle Cautery In Progress

- When the cautery function has terminated, the timer indicates the status as Stopped and the temperature display indicates that the needle is Cooling.
- If desired, press the **Start** button to re-initiate track ablation on additional sections of the needle track.
 - After completing track ablation, carefully remove the needle.

WARNING: Ensure adequate thawing or cooling before attempting to remove needles from the patient.

- If needle sticking is experienced, use a slight, gentle twist of the needle followed by slow withdrawal.
- During needle removal, watch for the 10 mm wide Active Zone Indicator on the needle to guide needle withdrawal. The Active Zone Indicator is a marked band along the needle shaft, located 20 mm distal to the heating portion of the needle.

Cautery Control for 2.1 CX Needles

1. Press and hold the **Thaw** button to access **Advanced Thaw Controls**.
2. Press the **Cautery** function button on the **Advanced Thaw Controls** menu to open the *Cautery Control Screen* (Screen 33). A confirmation message appears (Screen 37).
3. Press the desired needle port (A or B) on each channel containing a needle you wish to track ablate. Only one needle per channel may activate the cautery function at a given time.

NOTE: A maximum of four (4) needles may use the cautery function simultaneously.

NOTE: The duration for the cautery function phase for 2.1 CX-type needles may be selected from 30 seconds to 3 minutes, adjustable in 30 second increments. If a 2.1 CX-type is connected, a drop down menu will provide cautery function duration choices.

4. Select the cautery function duration from the drop-down menu.
5. Press the **Start** button to initiate track ablation (cautery).
 - During the needle warming phase, the *Cautery Control Screen* will display a rotating warming indicator (Screen 38).
 - When the cautery function threshold temperature is reached, Cauterization automatically begins.
 - o The timer will display Cauterizing (Screen 39).
 - o The *Cautery Control Screen* will display the estimated temperature range for the needle shaft.

NOTE: The shaft temperature is displayed as a temperature range because tissue and procedural variables affect the temperature.

- When the cautery function has terminated, the timer indicates the status as Stopped and the temperature display indicates that the needle is Cooling.
6. If desired, press the **Start** button to re-initiate track ablation on additional sections of needle track.
 7. After completing track ablation, carefully remove the needle.

WARNING: Ensure adequate thawing or cooling before attempting to remove needles from the patient.

- If needle sticking is experienced, use a slight, gentle twist of the needle followed by slow withdrawal.
- During needle removal, watch for the 10 mm wide Active Zone Indicator on the needle to guide needle removal. The Active Zone Indicator is a marked band along the needle shaft, located 20 mm distal to the heating portion of the needle.

Advanced Channel Controls

Advanced Channel Controls for each channel provide options to change the needle type for the selected channel, link two channels together, and program multiple freeze-thaw cycles.

Select Needle Type Control

1. To change the needle type for a channel, press and hold the **Channel** button to enter Advanced Channel Controls for that channel (Screen 40).
2. Select the proper needle type from the drop-down menu.
3. Press the **OK** button.

Screen 40. Advanced Channel Controls

Link Channels Control

1. Press and hold the **Channel** button to enter *Advanced Channel Controls* for that channel (Advanced Channel Controls).
2. Press the **Link** button to link two channels together for simultaneous operation. When two channels are linked, the **Channel** button shows both channels (Screen 41).

NOTE: This function is not available on the channel labeled **ALL**. You can only link channels that are on the same horizontal plane on the Needle Connection Panel (e.g., 1 and 2, 3 and 4, 5 and 6).


Screen 41. Linked Channels

3. Press the **Unlink** button (accessed by pressing and holding the **Channel** button) to remove the link between two channels so that each operates independently.

Cycle Programming Control

1. Press and hold the **Channel** button to enter *Advanced Channel Controls* for that channel.
2. Press the **Cycles** button (Advanced Channel Controls) to enter *Advanced Cycle Controls* to program freeze-thaw cycle(s) (Screen 42).


Screen 42. Advanced Cycle Controls

3. Select the desired freeze intensity from the *Freeze* controls using the drop-down menu and the freeze phase duration using the appropriate arrow or the scroll bar.
4. Add the programmed freeze cycle to the *Cycle Sequence* menu using the left **Arrow** button next to the *Freeze* controls.
5. Select the desired thaw by clicking on the available options in the *Thaw* controls. Select the thaw duration using the appropriate arrows or the scroll bar.
6. Add the programmed thaw cycle to the *Cycle Sequence* menu using the left **Arrow** button next to the *Thaw* controls.
7. Program additional cycles by repeating steps 3 – 6 as appropriate.
8. Arrange the cycle sequence by highlighting a programmed cycle on the *Cycle Sequence* controls. Using the **Up** or **Down** button, move the cycle into the desired sequence.
9. Remove a cycle from the *Cycle Sequence* by highlighting the cycle then pressing the **Stop** button.
10. Press the **Go** button to initiate the cryoablation procedure with the programmed cycles.

CAUTION: Any interruption of a programmed phase immediately terminates that phase and the programmed cycle.


11. Repeat steps 1 to 10 to program additional channels.

NOTE: Programmed sequences may be saved by selecting the **Save Cycle** button. Name the sequence, then press the **Save** button (Screen 43). To run a saved sequence, enter *Advanced Channel Controls* for the selected channel, press the **Cycles** button and press the **Load Cycle** button. Using the drop-down menu, select the saved sequence, press the **Load** button, then press the **Go** button (Screen 44).

Screen 43. Saved Sequence Controls

To run a saved sequence, enter *Advanced Channel Controls* for the selected channel, press the **Cycles** button and press the **Load Cycle** button. Using the drop-down menu, select the saved sequence, press the **Load** button, then press the **Go** button (Screen 44).


Screen 44. Saved Sequence Controls

Advanced Temperature Sensor Controls

The advanced temperature sensor controls provide an option to set visual alarm notifications when a chosen temperature measurement on an MTS falls below a desired value or when the rate of temperature drop is greater than a desired value for a select MTS sensor location.

1. Press and hold the **MTS Channel** button (Screen 24) for a selected MTS to enter *Advanced Temperature Sensor Settings*.

Screen 45. Advanced Temperature Sensor Controls

2. Select the desired alarm type and desired sensor location using the drop-down menus.
3. Press the **Up** or **Down** button and adjust the temperature to define the desired alarm limit.
4. Press the **Set** button to set the alarm.

OPTIONAL: Select a name to be displayed above the MTS channel location using the Needle Location drop-down menu. The list of needle names provided is derived from the list available in the *Configure Settings Screen* (reference the **Configure Settings** Section) and is associated with the selected organ map.

ADMINISTRATIVE FUNCTIONS

Configure Settings

The *Configure Settings Screen* allows you to change system settings used during a cryoablation procedure. A maximum of five (5) user accounts can be configured for each Visual-ICE Cryoablation System.

Settings that may be changed include MTS Needle Locations, System, Procedure and Registration Settings and Units (reference Table 13). Only Service personnel have the option to adjust the system time and date; System Administrators can change the Time Zone.

When settings have been changed, press the **Back** button to return to the *Startup Screen*. A message is displayed summarizing the changes made to the settings and asking for confirmation to save the settings. Press **Yes** to save the settings, **No** to exit the screen without saving your changes, or **Cancel** to return to the *Configure Settings Screen* and continue making changes.

6:42 am 22 Aug 2021

Back Manage Users Manual Software Update Configure Ethernet

Startup >> Configuration

Visual ICE™
Cryoablation System

MTS Needle Locations

Kidney

MTS Locations

Add Location

Remove Location

System Settings

Argon Cylinder Volume: 42.0000

Helium Cylinder Volume: 42.0000

Units: Liters (checked), Cubic Meters

Inactivity Timeout (minutes): 120

Language: English (English)

Procedure Settings

Maximize by Default

Channel Status

Temperature Sensors

Organ Map: None (checked)

Low Cylinder Alert (minutes): 10

Link all channels

Passive thaw timer count up

Active Flush

Display Sensor Temperatures

FastThaw in Channel Controls

Automatic Flush

Registration Settings

Custom Fields

Upload Registration: Disabled

Clear hospital name, address and doctor name history

Clear Hospital Information

Units

Pressure Units: bar, psi (checked), MPa

Temperature Units: Celsius (checked), Fahrenheit

Next maintenance due on: 04 Oct 2022

Export Logs

Time

Date: 22 Aug 2021 (The time and date can only be changed by service personnel)

Time: 6:41 am

Time Zone: (UTC-6:00) Central Time

Use 24 hour clock

Screen 46. Configure Settings

Table 13. Configure Settings Controls

Setting	Description
MTS Needle Locations	Customize the list of names to be provided for the MTS needles for each associated organ map. Press Add Location to add a new name to the list. Press Remove Location to delete a name from the list.
Cylinder Volume	Select the gas cylinder volume and units of measure according to the standard in the geographical region. The gas cylinder volume and units can only be changed by administrative or service personnel.
Inactivity Timeout	Select the desired duration from 30 minutes to 180 minutes that the system can be inactive before requiring you to re-enter your password. The default inactivity timeout is two hours.
Language	Select the language in which the software is displayed.
Maximize by Default	Select a section of the <i>Procedure Screen</i> to be maximized by default at login.
Low Cylinder Alert	Select the desired reminder interval (0 minutes to 15 minutes) for the Gas Indicator to display an alert that the estimated remaining volume in the gas cylinder is low.
Link all channels	Check this box to automatically link all adjacent channels for simultaneous operation (e.g., 1 and 2, 3 and 4, 5 and 6, ...)
Passive thaw timer count up	Check this box to automatically display elapsed time during passive thaw. The digital timer reads Stopped and the time elapsed in passive thaw is indicated.
Active Flush	Unless disabled by removing the check from this box, an automated active 30 second flush will follow a freeze cycle conducted at 50 % or higher intensity and lasting greater than 3 minutes.
Display Sensor Temperatures	Check this box for a Channel Status display of the internal gas temperature at the needle tip during freeze phase and the estimated temperature range of the needle shaft during active thaw phases for i-Thaw Function and CX-type needles.
FastThaw in Channel Controls	Check this box for display of the FastThaw Function icon in the <i>Channel Controls</i> when the system is in i-Thaw Function mode and a FastThaw Function capable needle is connected.
Automatic Flush	Unless disabled by removing the check from this box, an automated gas line flushing will occur after argon is connected and the gas lines pressurized.
Custom Fields	Enter custom names to label the two custom fields available for entering information in the <i>Registration Screen</i> .
Upload Registration	Use the drop-down menu to enable or disable the option to upload registration data with the procedure reports. By default, registration data is not uploaded. This function is only available to administrative or service personnel.
Clear Hospital Information	Clear the hospital name, address, and physician name from the system history file.
Time Zone	The time zone can be changed by administrative or service personnel. The Visual-ICE Cryoablation System automatically adjusts for Daylight Saving Time.
Pressure Units	Select the units of pressure that the Gas Indicator displays.
Temperature Units	Select the temperature units that the Temperature Sensors Section displays and graphs.

Control buttons at the top of the screen provide options to Manage Users and Manual Software Update.


Manage Users: Select your username and press the **Change Password** button to change your password. Administrative users may add users, remove users, or change the password of any user.

Manual Software Update: Press the **Manual Software Update** button to install a software update via a USB flash drive. This function is available only to administrative and service users.

Manual Software Update

Administrative and service users can manually update the Visual-ICE Cryoablation System software via a USB flash drive.

1. Press the **Manual Software Update** button on the *Configure Settings Screen* (Screen 46).
2. Press the **Browse** button to select the update file and press **Update**. A confirmation message appears when the software update is complete (Screen 47).

Screen 47. Software Update Confirmation

POST-PROCEDURE

Any serious incident that occurs in relation to this device should be reported to the manufacturer and to the relevant local regulatory authority.

For customers in Australia, report any serious incident that occurs in relation to this device to Boston Scientific and to the Therapeutic Goods Administration (<https://www.tga.gov.au>).

Cleaning the Visual-ICE Cryoablation System

Clean the Visual-ICE Cryoablation System after each use by following the steps below.

1. Clean the Touch Screen Monitor when the Visual-ICE Cryoablation System is powered OFF.
 - Wipe the screen gently with a damp gauze wipe.
 - Use water or isopropyl alcohol cleaning solutions.
 - Do not use cleaning agents such as Betadine Antiseptic Solution or bleach solution.
2. Clean the System by wiping with a damp gauze wipe.
 - Use soap and water or isopropyl alcohol cleaning solutions.
 - Do not use cleaning agents such as Betadine Antiseptic Solution or bleach solution.
 - Do not allow water or other liquid to drip or seep inside the needle connection ports. The needle connection ports must remain completely dry at all times.
3. Ensure that cleaned surfaces are dry prior to closing or turning on the system.

Disposal

All external and accessible surfaces of this device should be cleaned per the Cleaning the Visual-ICE Cryoablation System instructions included in the user manual. Include any common detachable cables (power cord, video cables, patch cords, etc.). Review the user manual to identify if any hazardous materials are present.

If placing the unit into an electronics recycling stream, notify the receiver of the presence of such materials. Use of recycling service suppliers familiar with medical electrical equipment is recommended, but not required. Do not dispose of by incineration, burial or placement into municipal waste stream.

Device should be safely disposed of in accordance with hospital, administrative, and/or local government policy or returned to Boston Scientific. For a Returned Product Kit, contact Boston Scientific Technical Assistance Center.

Dispose of all sharps directly into a sharps disposal container labeled with a biological hazard symbol. Sharps waste should be safely disposed of using available sharps waste channels in accordance with hospital, administrative, and/or local government policy.

TROUBLESHOOTING


Boston Scientific suggests the following options for troubleshooting the Visual-ICE Cryoablation System. If the suggested approaches do not resolve the issue or you experience an issue not referenced below, contact Boston Scientific Technical Assistance Center.

Software Recovery

In the event of software corruption or failure, the software may be restored to the previous software version.

Administrative and service users may update the software with an appropriate USB flash drive.


1. Shutdown the Visual-ICE Cryoablation System.
2. Hold down the **Software Reset** button by inserting a straightened paper clip into the Software Reset hole while simultaneously powering up the system. The system displays the *Software Recovery Screen*.

Screen 48. Software Recovery Screen

3. Press the **Rollback** button to restore the software to the previous software version.

4. *OPTIONAL*: Press the **Load** button to update the software if the screen displayed a message indicating that the software configuration is invalid on the *Login Screen* (Screen 49).

Screen 49. Invalid Configuration Message

5. If updating the software to a newer version available on a USB flash drive.
- Login as an administrative user.
 - Press the **Configure Settings** button on *Startup Screen* (Screen 16).
 - Press the **Manual Software Update** button on the *Configure Settings Screen* (Screen 29).
 - Insert the USB flash drive.

NOTE: Wait 20 seconds for the system to recognize the flash drive.

- Press the **Browse** button.
 - Select the file to conduct the upgrade.
 - Press the **Update** button.
-

NOTE:


- Wait until a displayed message confirms completion of the update.
 - The upgrade may take half an hour to complete.
-

Electronics, Electrical and User Error Related Issues

Symptom	Potential Causes / Solutions
System does not turn ON (i.e., the fan is not running) or power is lost during procedure	<ol style="list-style-type: none"> 1. The power control knob on the front panel of the system or power switch on the rear panel is OFF (Figure 1 and Figure 2). Turn the power ON. 2. The power cable to the Visual-ICE Cryoablation System is disconnected from the power outlet or from the rear panel of the system. Connect the power cable to the Visual-ICE Cryoablation System, checking that the power cable is fully seated. Connect the power cable to the power outlet. 3. There is no power at the wall outlet. Ensure the power outlet is switched ON. Call the hospital Biomedical Engineer for assistance, if needed. 4. A fuse may have blown. Spare fuses are located in the power cord inlet on the system (Figure 2). Reference the Replacing Fuses Section for instructions on how to replace fuses in the system.
Channel or needle is not recognized by the system and not available for use	<ol style="list-style-type: none"> 1. Check the corresponding Channel Locking Bar and ensure it is fully in the locked position. 2. At least one needle must be inserted into a channel for the channel to be available for use. 3. If a needle with a memory chip is being used and steps 1 - 2 are unsuccessful, select a new needle for use. Contact Boston Scientific Technical Assistance Center to arrange return of the non-operational needle. 4. The channel may be defective. Do not use this channel. Move the needle(s) to a different channel. Conduct the Needle Integrity and Functionality Test again.
Touch screen does not respond	<ol style="list-style-type: none"> 1. The track pad may be used to control the system. 2. Turn OFF and restart the system using the power control knob on the front of the system (Figure 1).
USB flash drive is not functioning Or The USB flash drive is not recognized by the system	<ol style="list-style-type: none"> 1. The USB flash drive is not connected to the USB port. Connect the USB flash drive to the port with the USB port icon (Figure 3). 2. The USB flash drive is not properly connected to USB port. Remove the USB flash drive from the designated USB port on the system. Wait a few seconds and reconnect the USB flash drive to the designated USB port. 3. If the problem persists, try using a different USB flash drive. 4. The USB flash drive is faulty. Replace the USB flash drive with a new USB flash drive.
The <i>Login Screen</i> appeared after the system was left idle for more than 2 hours while on the <i>Procedure Screen</i>	Enter the appropriate password to return to the <i>Procedure Screen</i> .
Touch screen goes blank during a procedure	<p>The video cable may be disconnected.</p> <p>WARNING: Do not touch the screen if the touch screen monitor goes blank for greater than five (5) seconds during a procedure. Immediately turn off power to the system and end the procedure to avoid inadvertent activation of needles.</p> <p>Contact Boston Scientific Technical Assistance Center.</p>

Replacing Fuses

Instruction	Photo
<p>1. Turn the Power Switch OFF on the rear of the Visual-ICE Cryoablation System. Loosen the thumb screw on the cord retention clip.</p>	
<p>2. Remove the power cord from the retention clip. Remove the two screws that secure the retention clip and remove the retention clip from the power inlet.</p>	
<p>3. Insert a small screwdriver into the slot on the bottom of the fuse holder to begin to snap the fuse holder out of the power inlet.</p>	
<p>4. While holding your hand under the fuse holder, carefully slide the fuse holder out of the power inlet. NOTE: There are four fuses in the fuse holder.</p>	
<p>5. Catch the fuse holder and fuses as the fuse holder is removed from the power inlet. The two fuses that remain in the fuse holder are the fuses that are part of the system circuit.</p>	
<p>6. Replace the fuses in the fuse holder with the two loose fuses. NOTE: Use only Boston Scientific specified fuses in the Visual-ICE Cryoablation System.</p>	

Instruction	Photo
7. Snap the fuse holder back into the power inlet. Replace the retention clip, plug in the power cord, and tighten the thumb screw on the retention clip.	
8. Call Boston Scientific Technical Assistance Center to schedule a service call to determine the cause of the blown fuses, to determine if service is required and to provide replacement fuses.	

Gas Issues

Symptom	Potential Causes / Solutions
The Visual-ICE Cryoablation System does not allow testing of a needle in a locked channel	The Argon Shutoff Valve may be in the Argon OFF position. Verify that the Argon Shutoff Valve (Figure 2) is in the Argon ON position to allow sufficient gas flow. Verify that the gas cylinder valve is fully open.
Needle is not freezing during the Needle Integrity and Functionality Test	<ol style="list-style-type: none"> 1. The argon cylinder valve may be closed. Open the gas cylinder valve by turning the gas cylinder valve counterclockwise to allow sufficient gas flow. Verify that appropriate pressure is displayed on the Gas Indicator. 2. Confirm that the argon cylinder is connected to the argon inlet. 3. The needle may be clogged (by dust or ice). Try to test it again. 4. If the needle still does not freeze <ul style="list-style-type: none"> • Press the Stop button to stop all activity on the channel. • Hold the needle connector firmly with one hand and unlock the channel to disconnect the needle. • Move the needle to a different channel and run the test again. <p>NOTE: If only one needle is connected to the channel, there may be some residual pressure behind the needle connector.</p> <ol style="list-style-type: none"> 5. If the problem persists, replace the needle with a new needle and contact Boston Scientific Technical Assistance Center.
No helium gas is delivered to the needle	<ol style="list-style-type: none"> 1. The helium gas cylinder valve may be closed. Open the gas cylinder valve by turning the valve counterclockwise to allow sufficient gas flow. Verify that appropriate pressure is displayed on the Gas Indicator. 2. The needle may be clogged. <ul style="list-style-type: none"> • Press the Stop button to stop all activity on the channel. • Hold the needle connector firmly with one hand and unlock the channel to disconnect the needle. • Move the needle to a different channel and run the test again. <p>NOTE: If only one needle is connected to the system, there may be some residual pressure behind the needle connector.</p> <ol style="list-style-type: none"> 3. If the problem persists, replace the needle with a new needle.

Symptom	Potential Causes / Solutions
The Manual Vent Valve is leaking gas	The Manual Vent Valve may be open. Fully close the Manual Vent Valve.
There is a hissing sound before the needles are connected	<ol style="list-style-type: none"> 1. Check that the regulated pressures are within the working pressure limits (green range on the Gas Indicator display). The system may be venting to reduce the pressure to below 4200 psi (289.6 bar, 28.96 MPa) to prevent damage to the system. If the pressure drops to within working pressure range, the system will operate normally. 2. The Automatic Vent Valve may be stuck open. If the manual vent valve is fully closed and the hissing persists, shut down the system and contact Boston Scientific Technical Assistance Center.
The pressure displayed on the Gas Indicator indicates the gas pressure is too low (Table 7)	<ol style="list-style-type: none"> 1. Verify the Argon Shutoff Valve is open. 2. Verify that the argon cylinder valve is sufficiently opened to allow gas flow. If necessary, open the valve approximately another half turn. 3. Verify that the cylinder contains sufficient pressure using the gauge on the cylinder. 4. Replace the cylinder, if needed.
During the Needle Integrity and Functionality Test, the needle is freezing in the first 45 seconds of the test instead of thawing, and then starts thawing for 15 seconds instead of freezing	<p>The gases are connected incorrectly (e.g., the helium gas supply line is connected to the argon cylinder and vice versa).</p> <ul style="list-style-type: none"> • End the procedure. • Vent the high pressure gas from the system. • Disconnect the gas supply lines and re-connect them to the correct cylinder. • Start a new procedure. • Run the needle test again.
It is difficult to loosen the pressure gauge connected to the cylinder, or the argon gas supply line cannot be disconnected from the argon inlet connection	<p>The gas lines were not vented and the gas lines are still under pressure.</p> <p>WARNING: If it is difficult to loosen the pressure gauge connected to the cylinder, or the argon gas supply line cannot be disconnected from the argon inlet connection, do not apply excessive force to release the gas supply line or to loosen the pressure gauge. The gas line may still be under pressure.</p> <ul style="list-style-type: none"> • Verify that the gas cylinders are CLOSED. • Verify the pressure gauge at the gas cylinder is 0 psi (0 bar, 0 MPa). • If in the <i>Procedure Screen</i>, verify the gas pressure display shows no gas connected. • If the Visual-ICE Cryoablation System is ON, end the procedure and vent the system using the automatic vent function. • If still unable to disconnect the gas supply lines, or if the system is OFF, open the Manual Vent Valve on the rear of the system to completely vent the system. • Upon completion, close the Manual Vent Valve.
Gas begins leaking through a needle port after gas flow was started using the Test , Freeze or Thaw buttons	<p>The channel may have a loose or broken socket.</p> <ul style="list-style-type: none"> • Disconnect the needle and move the needle to a different channel. • Conduct the Needle Integrity and Functionality Test again on the needle. • Contact Boston Scientific Technical Assistance Center.

Mechanical Issues

Symptom	Potential Causes / Solutions
Needle cannot be firmly connected in the needle port	<ol style="list-style-type: none"> 1. Verify that the locking bar is in the UNLOCKED position. 2. Needle connector may be defective. Use a different needle. 3. Residual gas pressure may be in the needle port. Use a different channel. 4. Check the Gas Indicator display. If there is pressure in the system, end the procedure and vent the system using the automatic vent function.
Locking bar on the Needle Connection Panel cannot be positioned in the LOCKED position	<ol style="list-style-type: none"> 1. Verify that all needles in the channel are fully inserted in the needle connection ports. 2. The Locking Bar may be defective. Move the needle to a different channel. Contact Boston Scientific Technical Assistance Center to schedule a service call. 3. Check the Gas Indicator display and verify the system is not pressurized. If there is pressure in the system, end the procedure and vent the system using the automatic vent function.
System does not move freely	<ol style="list-style-type: none"> 1. Release the brake to unlock the front wheels. 2. Check the individual rear wheel brakes and ensure the brakes are released.

Gas Cylinder and Gas Supply Line

Symptom	Potential Causes / Solutions
Safety cable is missing at either the cylinder side or the system side of the gas supply line	Do not use a gas supply line with a missing safety cable. Doing so could compromise the safety of personnel in the room. Contact Boston Scientific Technical Assistance Center for further instructions.
A gauge or a gas supply line is damaged	Do not use any damaged product. Contact Boston Scientific Technical Assistance Center for new accessories.
A gas leak was detected between the gauge adapter and the cylinder valve	<ul style="list-style-type: none"> • Tighten the connection using the wrench provided with the Visual-ICE Cryoablation System. • Close the cylinder valve and vent the gas from the Visual-ICE Cryoablation System and the gas supply lines using the Manual Vent Valve (Figure 2). Verify the system is depressurized. Loosen and remove the gauge assembly adapter. Verify no debris is on the gas cylinder connection point; clean the sealing surface as needed to remove any debris. Reposition and tighten the gauge assembly adapter onto the cylinder valve, using the wrench provided with the Visual-ICE Cryoablation System.

Needles

Symptom	Potential Causes / Solutions
A needle with a memory chip is not recognized by the system	<ol style="list-style-type: none"> 1. Move the needle to a different port. 2. If the system still does not recognize the needle, select the needle type from the drop-down menu displayed by the system. 3. If the needle is still not recognized, check the system's software version. (The software version is displayed on the <i>Startup Screen</i>). Contact Boston Scientific Technical Assistance Center for further instructions upon completion of the procedure.
During the freeze phase or following a thaw phase, no iceball or small iceball forms on the needle(s) in a given channel	<ol style="list-style-type: none"> 1. Perform the following steps in the order in which they are described: <ul style="list-style-type: none"> • Stop freezing/thawing actions for all channels. • Thaw the problematic needle(s) for at least one minute. • Freeze the needle(s) to verify proper function. 2. If the problem persists, connect a new needle into a different channel and test the needle. Continue the procedure with the newly tested needle. 3. Contact Boston Scientific Technical Assistance Center for further instructions upon completion of the procedure.
Only some of the temperature sensor locations on the MTS needle are registered on the Temperature Sensors display	<ol style="list-style-type: none"> 1. Move the MTS needle to a different port. 2. Replace the MTS needle. 3. Contact Boston Scientific Technical Assistance Center for further instructions upon completion of the procedure.
MTS needle is connected but temperature is not registered on the Temperature Sensors display	<ol style="list-style-type: none"> 1. Move the MTS needle to a different port. 2. Replace the MTS needle. 3. Contact Boston Scientific Technical Assistance Center for further instructions upon completion of the procedure.
During the Needle Integrity and Functionality test bubbles are seen escaping from the needle	<p>WARNING: Do not use the needle.</p> <ul style="list-style-type: none"> • Disconnect the needle from the Visual-ICE Cryoablation System and segregate the needle. • Return the needle to Boston Scientific for evaluation. • Use a new needle to continue the procedure. • Test the new needle to confirm the integrity and functionality of the new needle. • Contact Boston Scientific Technical Assistance Center for further instructions upon completion of the procedure.
Needle is bent or damaged while attempting to unpack or use it	<p>WARNING: Do not use the needle.</p> <ul style="list-style-type: none"> • Segregate the needle. • Use a different needle to continue the procedure. • Contact Boston Scientific Technical Assistance Center for further instructions upon completion of the procedure.

Displayed Messages

The Visual-ICE System displays a message on the user interface when a user requests assistance or when user, needle or system errors are detected.

NOTE: Record and report the message number (e.g., 10-01, 80-02) if assistance is required from Boston Scientific Technical Assistance Center.

LOGIN

Message	Reason for Appearance / Solutions
10-01 Login <ul style="list-style-type: none"> You have not entered the correct Login Name Reenter your Login Name Contact your System Administrator if assistance is required Contact Boston Scientific Technical Assistance Center if further assistance is required 	No name was entered. OR The name entered did not match assigned names in the system.
10-02 Login <ul style="list-style-type: none"> You have not entered the correct Password Reenter your Password Contact your System Administrator if assistance is required Contact Boston Scientific Technical Assistance Center if further assistance is required 	No password was entered. OR The password entered did not match the password associated with the login name.
10-03 Reset Password Challenge <ul style="list-style-type: none"> To reset your password, contact Boston Scientific Technical Assistance Center Relay the Challenge on the screen below Enter the Response provided by Technical Assistance Center Press the Reset button 	User forgot their password, pressed the Forgot Password button and received a challenge to be relayed to the Technical Assistance Center.
10-04 Password Reset <ul style="list-style-type: none"> Your password has been reset to XXX Change your password in the configuration screen when convenient 	User correctly entered the reset password challenge and should now set the new password.
10-05 Emergency Login <ul style="list-style-type: none"> To obtain an emergency login, contact Boston Scientific Technical Assistance Center Relay the Challenge on the screen below Enter the Response provided by Technical Assistance Center Press the Login button 	User requested an emergency login and received a challenge to be relayed to the Technical Assistance Center.
10-06 Incorrect Response <ul style="list-style-type: none"> You have not entered the correct response Contact Boston Scientific Technical Assistance Center for a Response to the on-screen Challenge 	User tried to access the Emergency Login and did not respond correctly to the challenge. User must contact to the Technical Assistance Center in order to receive an emergency login response. This action does not reset a password.
10-07 Incorrect Response <ul style="list-style-type: none"> You have not entered the correct response Contact Boston Scientific Technical Assistance Center for a Response to the on-screen Challenge 	User tried to reset their password and did not respond correctly to the challenge. User must contact to the Technical Assistance Center in order to reset their password.

PROCEDURE

Message	Reason for Appearance / Solutions
20-01 Cannot Start Test <ul style="list-style-type: none"> The gas pressure is too low/high to begin a procedure Check that the gas cylinders have enough pressure to start the procedure 	User pressed the Test button when either gas cylinders were not yet connected or the gas pressure is below operating pressure (Reference Table 7). To proceed, gas cylinders with sufficient pressure must be connected.
20-02 Argon Shut Off Valve <ul style="list-style-type: none"> The argon shut off valve may be closed Check and open if necessary 	At system start-up, the system detected gas connected, but no gas was entering the system. The argon shut off valve may be closed. To proceed, the argon shut off valve must be open.
20-03 Test All <ul style="list-style-type: none"> You selected to control all needles simultaneously Place the needles, individually or in groups, in the basin such that the full length of the needle shaft is submerged in the water or saline. As the test cycles run, look at the needles to ensure that there are no air bubbles, and that a small ice ball forms during the freeze portion of the sequence. Do you want to initiate test in all needles now? <p style="text-align: center;">YES NO</p>	User selected Test button in ALL channel to control all active channels. To proceed, the user must confirm test for all connected needles.
20-03 Test All <ul style="list-style-type: none"> You selected to control all needles simultaneously Place the needles, individually or in groups, in the basin such that the full length of the needle shaft is submerged in the water or saline. As the test cycles run, look at the needles to ensure that there are no air bubbles, and that a small ice ball forms during the freeze portion of the sequence. Do you want to initiate test on all channels or just the untested channels? <p style="text-align: center;">YES NO</p>	User selected Test button in ALL channel to control all active channels. To proceed, the user must confirm test for all connected needles or just the untested channels.
20-04 Freeze All Channels <ul style="list-style-type: none"> You selected to control all needles simultaneously Do you want to initiate freeze for all active needles now? <p style="text-align: center;">YES NO</p>	User selected Freeze button in ALL channel to control all active channels. To proceed, the user must confirm freezing of all active needles.
20-05 Freeze All Intensity <ul style="list-style-type: none"> You selected to control all needles simultaneously. Do you want [x%] freeze intensity applied to all needles now? <p style="text-align: center;">YES NO</p>	User selected a freeze intensity in ALL channel to control all active channels. To proceed, the user must confirm the selected intensity is to be applied to all active needles.
20-06 Thaw All Channels <ul style="list-style-type: none"> You selected to control all needles simultaneously Do you want to initiate thaw for all active needles now? <p style="text-align: center;">YES NO</p>	User selected Thaw button in ALL channel to control all active channels. To proceed, the user must confirm thawing of all active needles.

Message	Reason for Appearance / Solutions
20-07 Stop All Channels <ul style="list-style-type: none"> You selected to control all needles simultaneously Do you want to stop activity in all needles now? <p style="text-align: center;">YES NO</p>	User selected Stop button in ALL channel to control all active channels. To proceed, the user must confirm stop for all active needles.
20-10 Cautery <ul style="list-style-type: none"> CAUTION: Activation of thermal cautery immediately ceases all other cryoablation operations for the duration of cautery activity. Do you want to proceed with cautery? <p style="text-align: center;">YES NO</p>	User pressed Thaw button to open the Advanced Thaw Controls and the thermal cautery option. The CAUTION is displayed upon user selecting advanced option for thermal cautery.
20-11 End Procedure <ul style="list-style-type: none"> Are you sure you want to end the procedure? <p style="text-align: center;">YES NO</p>	User selected End Procedure and must confirm desire to end the procedure.
20-12 Automatic Vent <ul style="list-style-type: none"> Do you want to automatically vent high pressure gas from the system? <p style="text-align: center;">YES NO</p>	User was provided an option for automatically venting the high pressure gas in the system.
20-13 Automatic Vent <ul style="list-style-type: none"> Is the gas supply closed? <p style="text-align: center;">YES CANCEL</p>	If user selected the option for automatic gas venting, the user will need to confirm the gas supply is closed before automatic venting is activated.
20-14 Automatic Vent <ul style="list-style-type: none"> The gas pressure is not dropping Check that the gas cylinder shut off valve is closed 	User selected to automatically vent high pressure gas at the end of a procedure, but the pressure was not dropping. User must ensure the shut off valve is closed.
20-15 Automatic Vent <ul style="list-style-type: none"> Venting is in progress If needles are still connected, do not unlock channels or disconnect needles until venting is complete 	User selected to automatically vent the high pressure gas at the end of the procedure.
20-16 Automatic Vent <ul style="list-style-type: none"> Automatic venting successfully completed 	User chose to automatically vent the high pressure gas at the end of the procedure.
20-17 Gas Vent <ul style="list-style-type: none"> Before disconnecting the gas hose, manually vent the Visual-ICE Cryoablation System using the Manual Vent Valve on the rear of the machine 	User selected to not use the feature to automatically vent high pressure gas from the System.
20-18 System Shutdown <ul style="list-style-type: none"> Do you want to shut down the system? <p style="text-align: center;">YES NO</p>	User selected Shutdown in the login screen to power down the system.

Message	Reason for Appearance / Solutions
20-19 Procedure Timeout <ul style="list-style-type: none"> The procedure has exceeded the allowable time Procedure will be terminated 	The procedure exceeded the allowable 8 hour duration.
20-20 Maximum Activity for FastThaw Function <ul style="list-style-type: none"> Maximum limit for FastThaw Function reached Wait until thawing is complete before activating additional FastThaw Function capable needles This channel will use passive thaw 	User attempted to activate more than four FastThaw Function capable needles. The Visual-ICE Cryoablation System supports simultaneous operation of no more than four FastThaw Function capable needles. After completing the thaw phase with four needles, additional needles may be activated for FastThaw Function.
20-21 Maximum Activity for i-Thaw Function <ul style="list-style-type: none"> Maximum limit for i-Thaw Function reached Wait until thawing is complete before activating additional i-Thaw Function capable needles This channel will use passive thaw 	Thawing with i-Thaw Function is limited to a maximum of 9 needles simultaneously activated in thaw mode. When thawing with CX-type needles, thawing is limited to a maximum of 7 needles simultaneously activated. After completing the thaw phase, additional needles may be activated.
20-22 Maximum Limit for Active Thaw <ul style="list-style-type: none"> Maximum limit for i-Thaw Function and FastThaw Function reached Wait until thawing is complete before activating thaw on additional needles 	User attempted to activate more needles than are supported for active thaw. After completing the thaw phase, additional needles may be activated.
20-23 Maximum Activity for Test <ul style="list-style-type: none"> Maximum limit for test reached Wait until testing is complete before activating additional needles for testing 	User attempted to initiate testing on more needles than are supported for simultaneous testing.
20-24 Advanced Thaw Unavailable <ul style="list-style-type: none"> Needle(s) in channel X are not capable of advanced thaw Advanced thaw activity is unavailable 	User attempted to initiate Advanced Thaw in a channel. The channel did not contain i-Thaw Function or FastThaw Function capable needles.
20-25 Advanced Thaw Unavailable <ul style="list-style-type: none"> Advanced thaw activity is unavailable on channel X while helium is connected Use helium thaw or disconnect helium 	User attempted to initiate Advanced Thaw while helium was connected. If helium is connected, the system defaults to helium thaw mode. If mixed needle types are connected (some helium thaw and some i-Thaw Function), the system also defaults to helium thaw mode.

GAS

Message	Reason for Appearance / Solutions
30-01 Automatic Vent <ul style="list-style-type: none"> Gas was not appropriately vented from the Visual-ICE Cryoablation System when the system was last used Vent the system, using either the automatic vent option or the Manual Vent Valve 	Upon startup, there was gas pressure remaining in the system which will make it difficult to connect needles.
30-02 Gas Vent <ul style="list-style-type: none"> Manually vent the Visual-ICE Cryoablation System using the valve on the rear of the machine 	User pressed Close button after receiving message that gas had not been appropriately vented when the system was last used. To clear the remaining gas pressure, user must manually vent the gas using the Manual Vent Valve.
30-03 Low Gas Level <ul style="list-style-type: none"> Low level of [helium / argon] gas remains in the cylinder Replace the gas cylinder as soon as feasible with a new cylinder 	The system displayed an alert that the estimated remaining volume in the gas cylinder is low. The alert can be configured by the user to be displayed at an interval of 0 minutes to 15 minutes. (Reference Table 13, Low Cylinder Alert).

NEEDLES

Message	Reason for Appearance / Solutions
40-01 Unsupported Needle <ul style="list-style-type: none"> The connected needle is not supported by the software Choose and connect a different needle type 	A needle type not supported by the software configuration was connected in the channel. The channel will be disabled until an appropriate needle is connected.
40-03 Recalled Needle <ul style="list-style-type: none"> This needle lot number is identified as part of a recall and is unavailable for use Return the needle to Boston Scientific Connect a needle from a different lot number to continue the procedure 	The needle lot number was identified as a lot number in a regulatory recall. The channel will be disabled until a needle from a different lot number is connected.
40-04 Expired Needle <ul style="list-style-type: none"> The needle 'Use by' date is in the past Replace with a new needle 	The needle expiration date was identified as an expired date. The channel will be disabled until an appropriate needle is connected.
40-05 Used Needle <ul style="list-style-type: none"> This needle has been previously used Replace with a new needle 	The software prevents needle reuse and identified a needle previously used. The channel will be disabled until a new needle is connected.
40-06 Corrupt Memory <ul style="list-style-type: none"> Needle Memory chip is corrupt on Channel X Manually select the needle type 	User connected a needle containing a bad memory chip. User may select needle type from the drop-down menu.
40-07 Passive Thaw Required <ul style="list-style-type: none"> The currently selected needle on channel [x] does not have i-Thaw Function capability Since helium is not connected, passive thaw is required on this channel 	Combinations of needles that are i-Thaw Function capable and that are not i-Thaw Function capable were connected to the needle connection panel; no helium gas was connected, so passive thaw was required.
40-08 i-Thaw Error <ul style="list-style-type: none"> Channel [Xa or Xb] i-Thaw Function capable needle is defective for electrical thawing Use passive thaw or connect helium gas 	An i-Thaw Function capable needle was identified as not i-Thaw Function operational. Passive thaw or connection of helium for active thaw is required.

TEMPERATURE SENSOR / ORGAN MAP

Message	Reason for Appearance / Solutions
50-01 Alarm <ul style="list-style-type: none"> Sensor point (5, 15, 25 or 35) on MTS (#X) met selected alarm limit 	The MTS detected the selected alarm temperature.
50-02 Alarm <ul style="list-style-type: none"> Sensor point (5, 15, 25 or 35) on MTS (#X) met selected alarm limit 	The MTS detected the selected alarm rate of temperature change within 30 seconds.
50-03 MTS Disconnected <ul style="list-style-type: none"> An MTS is disconnected Reconnect the MTS to continue to monitor temperature in that location 	An MTS became disconnected during a procedure.
50-05 Clear Drawing <ul style="list-style-type: none"> Are you sure you want to erase the entire drawing? YES NO 	User selected the Clear Drawing button. When selected, all information placed on the organ map, except needle placement, will be erased.
50-06 Organ Type <ul style="list-style-type: none"> All annotation and needles will be removed from the organ before proceeding Would you like to proceed with the organ type change? YES NO 	User changed the organ type.

SOFTWARE

Message	Reason for Appearance / Solutions
60-04 Update Failed <ul style="list-style-type: none"> An error occurred during the update process Retry update 	During User update of the system, an error occurred that prevented completion of the update. Another update should be attempted.
60-04 Update Failed <ul style="list-style-type: none"> An error occurred during the upload process Retry upload 	During User update of the system, an error occurred that prevented completion of the upload. Another upload should be attempted.
60-05 Incompatible Hardware <ul style="list-style-type: none"> The hardware is not compatible with the current software Contact Boston Scientific Technical Assistance Center to schedule service 	When the system performed self-checks at start up, an incompatibility between the hardware and software was detected. Contact Boston Scientific Technical Assistance Center.
60-06 Incompatible Software <ul style="list-style-type: none"> The Visual-ICE Cryoablation System software is not compatible with regulatory approvals Contact Boston Scientific Technical Assistance Center to schedule service 	The software was checked against approved software versions in the regulatory files for each respective market. Incompatibility with regulatory approval was identified. The system should be updated with the appropriate software. Contact Boston Scientific Technical Assistance Center.
60-07 Invalid Configuration <ul style="list-style-type: none"> A software file is nonoperational. If problem persists, contact Boston Scientific Technical Assistance Center 	There was a problem with the software's configuration files. Contact Boston Scientific Technical Assistance Center.
60-09 Software Recovery <ul style="list-style-type: none"> This will rollback the software on this system Are you sure you want to do this? YES NO 	User pressed the Software Recovery button then selected Rollback. Activation will rollback the software to the previous software version.

Message	Reason for Appearance / Solutions
60-10 Software Recovery <ul style="list-style-type: none"> Are you sure you want to restore the default settings for all configurations? <p style="text-align: center;">YES NO</p>	User pressed the Software Recovery button then selected Load. Activation will restore the system default settings for all configurations.

REPORTS

Message	Reason for Appearance / Solutions
70-01 Save Report <ul style="list-style-type: none"> Do you want to save the report to the Visual-ICE Cryoablation System? <p style="text-align: center;">YES NO</p>	User selected End Procedure and was provided an option to save the report before exiting from procedure.
70-02 System is Busy <ul style="list-style-type: none"> Finalizing Procedure 	System activity was displayed during process of report saving.
70-03 Report Error <ul style="list-style-type: none"> Errors occurred while compiling the report The report may be incomplete 	User chose to access a report during a procedure or to save data to a report at the end of a procedure. Errors occurred which may affect the completeness of the report.
70-04 Report Saved <ul style="list-style-type: none"> Report saved successfully 	Report successfully saved to a USB
70-05 Duplicate Filename <ul style="list-style-type: none"> The chosen filename already exists on the USB flash drive Choose a different filename 	User tried to export a report to a USB flash drive using a filename duplicate to a filename on the flash drive. To export the report, a different filename must be chosen.
70-06 Report Error <ul style="list-style-type: none"> Unable to export the report to the USB flash drive The USB flash drive may be disconnected or full 	User selected Save Reports to Flash Drive. No flash drive was detected or there was insufficient space on the flash drive.

SYSTEM

Message	Reason for Appearance / Solutions
80-01 Communication Failure <ul style="list-style-type: none"> Internal communication failed Reconnection attempt failed Restarting Visual-ICE Cryoablation System If problem persists, contact Boston Scientific Technical Assistance Center 	Software was unable to connect with hardware after attempt to reinitialize communication. If the restart fails, the system is not usable.
80-02 Startup Failure <ul style="list-style-type: none"> System self-checks failed Restart the system If problem persists, contact Boston Scientific Technical Assistance Center 	Software self-checks located a failure that required the system to be restarted.
80-03 Pressure Alert <ul style="list-style-type: none"> Pressure exceeds safe operating limits Close the gas cylinders Procedure will be terminated and gas vented from the system 	System detected the internal pressure was above safe limits. The system will terminate the procedure and vent gas from the system.
80-04 Temperature Warning <ul style="list-style-type: none"> The internal temperature of the Visual-ICE Cryoablation System exceeds appropriate operating limits Discontinue the cryoablation procedure as soon as safe to do so Contact Boston Scientific Technical Assistance Center 	The internal temperature of the system exceeded appropriate operating limits.
80-05 Service Due <ul style="list-style-type: none"> Low battery detected Contact Boston Scientific Technical Assistance Center to schedule service 	System detected a low battery. System operation may be impacted if operated with a low battery.
80-30 System Error* <ul style="list-style-type: none"> Gas pressure display may be inaccurate. Monitor procedure carefully using image guidance. Use passive thaw. At the conclusion of the procedure, contact Boston Scientific Technical Assistance Center. 	* System Error messages are displayed in the right corner of the navigation tool bar. The User is offered details of the error by "Press here" instruction. Internal pressure checks were inconsistent and may produce an inaccurate display on the Pressure Gauge .
80-31 System Error* <ul style="list-style-type: none"> Gas cylinder shut off valve is not open enough to provide adequate flow. If necessary, open the valve approximately another half turn. 	* System Error messages are displayed in the right corner of the navigation tool bar. The User is offered details of the error by "Press here" instruction. Gas flow from gas cylinder was insufficient. To improve gas flow, the cylinder valve should be opened further.
80-32 System Error* <ul style="list-style-type: none"> Channel X is defective. Choose another channel. At the conclusion of the procedure, contact Boston Scientific Technical Assistance Center. 	* System Error messages are displayed in the right corner of the navigation tool bar. The User is offered details of the error by "Press here" instruction. A defective solenoid on Channel X was detected; another channel should be selected.

Message	Reason for Appearance / Solutions
80-33 System Error* <ul style="list-style-type: none"> Defective MTS in channel [X]. Replace with a new MTS. 	* System Error messages are displayed in the right corner of the navigation tool bar. The User is offered details of the error by "Press here" instruction. The MTS temperature display was out of the expected range when initially connected to the needle connection panel.
80-34 System Error* <ul style="list-style-type: none"> Fan X failure. At the conclusion of the procedure, contact Boston Scientific Technical Assistance Center. 	* System Error messages are displayed in the right corner of the navigation tool bar. The User is offered details of the error by "Press here" instruction. Fan X was identified as not operating.
80-35 System Error* <ul style="list-style-type: none"> The gas flow rate in channel XX exceeds the system operation requirements and may impact performance. Reduce the number of active needles. 	* System Error messages are displayed in the right corner of the navigation tool bar. The User is offered details of the error by "Press here" instruction. The calculated flow rate for a given channel exceeded the optimum system requirement. The number of active needles should be reduced.
80-36 System Error* <ul style="list-style-type: none"> The gas flow rate exceeds the system operation requirements and may impact performance. Reduce the number of active needles. 	* System Error messages are displayed in the right corner of the navigation tool bar. The User is offered details of the error by "Press here" instruction. The cumulative flow rate for all channels exceeded the optimum system requirement. The number of active needles should be reduced.
80-37 System Error* <ul style="list-style-type: none"> Defective MTS in channel 1, 2, 3, 4 	* System Error messages are displayed in the right corner of the navigation tool bar. The User is offered details of the error by "Press here" instruction. Contact Boston Scientific Technical Assistance Center.
80-38 System Error* <ul style="list-style-type: none"> Passive thaw required 	* System Error messages are displayed in the right corner of the navigation tool bar. The User is offered details of the error by "Press here" instruction. Contact Boston Scientific Technical Assistance Center.

SERVICE

Message	Reason for Appearance / Solutions
90-01 Service Due <ul style="list-style-type: none"> Service the Visual-ICE Cryoablation System soon Service must be completed by [DATE]. Contact Boston Scientific Technical Assistance Center to schedule service. 	User was reminded of the deadline to complete scheduled system service. Reminder begins four weeks prior to service due date.
90-02 Service Due <ul style="list-style-type: none"> Service for the Visual-ICE Cryoablation System is past due Contact Boston Scientific Technical Assistance Center to schedule service. 	The Visual-ICE Cryoablation System was not serviced on schedule. During subsequent startups, the User is reminded the service is past due.
90-03 System End of Life <ul style="list-style-type: none"> The Visual-ICE Cryoablation System is at the end of operational life Contact Boston Scientific Technical Assistance Center to arrange return of this system to be refurbished, replaced or disposed. 	The Visual-ICE Cryoablation System reached the end of its operating life. Contact Boston Scientific Technical Assistance Center to schedule service.

SYSTEM SPECIFICATIONS

Mechanical Specifications

Visual-ICE Cryoablation System

- Weight: 100 kg (220 lb)
- Height: 107 cm (42 in), monitor down
157 cm (62 in), monitor up
- Foot Print: 56 cm x 66 cm (22 in x 26 in)
- Storage Compartment Weight Capacity: 22 kg (50 lb)
- Monitor Storage Basin Weight Capacity: 9 kg (20 lb)
- Closed Monitor Weight Capacity: 9 kg (20 lb)

External Gas Supply

- Argon Cylinder:
 - Purity Level: 99.998 % or higher
 - Solid particle size: < 5 µm
- Helium Cylinder:
 - Purity Level: 99.995 % or higher
 - Solid particle size: < 5 µm

Gas Cylinder Specifications

- Maximum Pressure: 6000 psi (414 bar, 41.4 MPa)
- Recommended Volume of Gas Cylinder: 42 L - 50 L

Accuracy of Displayed Values

- **Temperature Accuracy:**
 - ± 3 °C over range of -60 °C to 40 °C
- **Supplied Gas Pressure Accuracy:**
 - ± 50 psi, over range of 1000 psi to 6000 psi
 - ± 3.4 bar over range of 69 bar to 414 bar
 - ± 0.344 MPa over range of 6.9 MPa to 41.4 MPa
- **Built-in Regulator Gas Pressure:**
 - ± 50 psi, over range of 1000 psi to 4000 psi
 - ± 3.4 bar over range of 69 bar to 276 bar
 - ± 0.344 MPa over range of 6.9 MPa to 27.6 MPa
- **Time Intervals:**
 - ± 5 seconds over any 10 minute interval

Essential Performance

The Visual-ICE Cryoablation System essential performance is defined as:

- The ability to turn on and turn off the flow of argon or helium to the connected cryoablation needles via command inputs from the user interface
- The ability to maintain uninterrupted gas flow to user-selected needles when gas flow is turned on by the user
- The ability to uninterruptedly prevent gas flow to user-selected needles when gas flow is turned off by the user
- The ability to energize the electrical resistance heater on user-selected i-Thaw Function/cautery function capable needles
- The ability to continuously monitor and display temperatures measured by user-selected thermal sensors

Boston Scientific Technical Assistance Center:

Region	Contact Number	email
United States	+1.800.949-6708	CEtechsupportUSA@bsci.com
APAC (Asia, Middle East)	+65.64.18.8878	CETechSupportAPAC@bsci.com
Japan	+81.120.177.779	JapanCESTAC@bsci.com
China	+86 400.801.6908	ChinaCapitalService@bsci.com
Australia/New Zealand	+61 1800.676133 - option 5	CapitalequipmentANZ@bsci.com
Brazil	+55 1155459063	CEBrazilTeam@bsci.com
Mexico	+52 15559924100	CESupportMXC@bsci.com
Europe (See Below for Individual countries)	+31 45.546.7707	CEtechsupportEMEA@bsci.com
Austria	+43 16 08 10 37	CEtechsupportEMEA@bsci.com
Denmark	+45 80 25 34 29	CEtechsupportEMEA@bsci.com
Czech Republic	0800 142942	CEtechsupportEMEA@bsci.com
Finland	+358 800770055	CEtechsupportEMEA@bsci.com
France	+33 139 304 971	CEtechsupportEMEA@bsci.com
Germany	+49 815 126 86118	CEtechsupportEMEA@bsci.com
Italy	+39 022 698 3218	CEtechsupportEMEA@bsci.com
Netherlands	+31 45 5467707	CEtechsupportEMEA@bsci.com
Norway	0800 14236	CEtechsupportEMEA@bsci.com
Spain	+34 917 619 999	CEtechsupportEMEA@bsci.com
Sweden	+46 020790077	CEtechsupportEMEA@bsci.com
UK	+44 1442 411 686	CEtechsupportEMEA@bsci.com

PATIENT COUNSELING INFORMATION

The physician should consider the following points while counseling the patient on the use of the Visual-ICE Cryoablation System in association with the interventional procedure:

- Discuss the risks and benefits, including review of potential adverse events listed in the IFUs for the Visual-ICE Cryoablation System and accessory products used to conduct Cryoablation procedures and for other interventional treatments likely to be employed.
- Discuss post-procedure instructions, including any lifestyle changes, medications, and home care and rehabilitation guidelines.

WARRANTY

For device warranty information, visit (www.bostonscientific.com/warranty).

SYMBOL DEFINITIONS

Commonly used medical device symbols that appear on the labeling are defined at www.bostonscientific.com/SymbolsGlossary.

Additional symbols are defined at the end of this document.

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10-03 **Reset Password Challenge**

- To reset your password, contact Boston Scientific Technical Assistance Center
 - +1 800 949 6708
 - CEtechsupportUSA@bsci.com
- Relay the Challenge on the screen below
- Enter the Response provided by Technical Assistance Center
- Press the Reset button

User Name: A
Challenge: GFEECG
Response:

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10-04 **Password Reset**

- Your password has been reset to BSC
- Change your password in the configuration screen when convenient

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The screenshot displays the Visual-ICE Cryoablation System interface. At the top, there are navigation tabs for 'Registration', 'Settings & Procedures', 'Notes', 'Report', and 'End Procedure'. The main display area is divided into several sections: 'Connect Needles' with instructions, 'Channel Control' and 'Channel Status' panels with buttons for 10 channels, 'Temperature Sensors' with a graph showing temperature over time, and 'Organ Map' with a grid of 10 points. The interface also shows pressure gauges for 'Freeze' (24.2 MPa) and 'Thaw' (15.7 MPa), and a clock showing 6:54 on 01 Sep 2021.

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- | | | |
|------------------------------|------------------------------------|--|
| 1 Navigations-Werkzeugleiste | 3 Channel Control (Kanalsteuerung) | 5 Temperature Sensors (Temperatursensoren) |
| 2 Kontextsensitive Hilfe | 4 Channel Status (Kanalstatus) | 6 Organ Map (Organkarte) |

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Registration | **Settings & Procedures** | **Notes** | **Report** | **End Procedure**

Est. Remaining Time 00:05:32 (Freeze) | **Est. Remaining Time 00:05:34** (Thaw)

7:38 am | 26 Aug 2021

VISUAL ICE™
Cryoablation System

Report Sections

- Registration
- Procedure
- Temperature Sensors
- Organ Map
- Notes

Refresh | **Export Report**

Needle Types

Channel	A	B	A	B	Channel
1	IcePearl 2.1 CX				2
3	IcePearl 2.1 CX				4
5		IcePearl 2.1 CX			6
7					8
9					10

Channel 1 Timeline

0 sec | 2 min 33 sec

100% Freeze 1 min 54 sec | Stopped 39 sec

Number of Freeze Thaw cycles: 1

Activity Durations

Freeze Durations	Thaw Durations
1 min 54 sec	39 sec

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The screenshot shows the 'View Reports' interface. At the top right, the system name 'VISUAL ICE™ Cryoablation System' is displayed along with the time '6:40 am' and date '22 Aug 2021'. A 'Back' button is located in the top left. Below the navigation bar is a breadcrumb trail: 'Startup >> View Reports'.

The main content area features a table with the following columns: Report Name, Hospital Name, Doctor Name, and Procedure Type. The table contains five entries, with the third entry selected:

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Below the table, the selected report '2021_Jul_18_8_15_am_Tim' is displayed in a detailed view. On the left, a 'Report Sections' sidebar lists: Registration (selected), Procedure, Temperature Sensors, Organ Map, and Notes. At the bottom of this sidebar are 'Export Report' and 'Delete Report' buttons. The main report content area shows the 'VISUAL ICE™ CRYOABLATION SYSTEM' logo and 'Boston Scientific' branding. The title is 'Lung Procedure' with a timestamp of '8:15 am 18 Jul 2021'. Under the 'Registration' section, the following patient information is listed:

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The 'Procedure' section is currently empty.

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The screenshot displays the Visual ICE Cryoablation System interface. At the top right, the time is 6:41 am and the date is 22 Aug 2021. A 'Back' button is visible in the top left. The main area shows a table of reports with columns for Report Name, Hospital Name, Doctor Name, and Procedure Type. The third row is selected, showing a report for 'Lung' procedure at 'ST. Mary's' hospital by 'Dr. Thomas' on '2021 Jul 18 8 15 am Tim'. An 'Export Report' dialog box is open, allowing the user to choose a file type (PDF) and a file name (2021_Jul_18_8_15_am_Tim). Below the dialog, there are buttons for 'Export' and 'Cancel'. The background also shows a detailed view of the selected report, including patient information (Patient ID: JohnKI-9001, Hospital Name: ST. Mary's, Hospital Address: 500 Peach St, Doctor Name: Dr. Thomas) and a section for 'Procedure'.

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The screenshot displays the Visual ICE Cryoablation System interface. At the top, there's a navigation bar with 'Registration' and 'Settings & Procedures' tabs, and buttons for 'Notes', 'Report', and 'End Procedure'. Two pressure gauges are shown: 'Freeze' at 24.2 MPa and 'Thaw' at 15.7 MPa. The main interface is divided into several sections: 'Connect Needs' with instructions, 'Channel Control' with a 'Scale Scroll' slider, 'Channel Status' with a 'Scale Scroll' slider, 'Temperature Sensors' with a graph showing temperature over time, and 'Organ Map' with a grid of points labeled A and B. A numeric keypad is located at the bottom right.

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The screenshot displays the Visual ICE Cryoablation System interface, focusing on the navigation bar. It includes 'Registration' and 'Settings & Procedures' tabs, and buttons for 'Notes', 'Report', and 'End Procedure'. Two pressure gauges are shown: 'Freeze' at 24.2 MPa and 'Thaw' at 15.7 MPa. The top right corner shows the time '6:54' and date '01 Sep 2021'.

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INSTRUCTIONS FOR USE

en

The Instructions for Use (IFU) for this product are supplied in electronic form over the Internet.

Visit www.IFU-BSCI.com to access the IFU using your internet browser. Information from the product labeling such as Product Code, UDI-DI or Product Name may be used to search for the specific IFU.

If you have difficulty accessing the IFU online, or would prefer to receive a paper copy, please contact Boston Scientific Customer Service or your local country contact. A copy will be sent to you at no charge and should arrive within seven days.

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Visual-ICE, EZ-Connect2, FastThaw, i-Thaw, and Multi-Point 1.5 Thermal Sensor are trademarks of Boston Scientific Corporation or its affiliates.
All other trademarks are the property of their respective owners.
All rights reserved.
Boston Scientific Limited
Boston Scientific (Australia) Pty Ltd
Boston Scientific Corporation