Encore™ 26

Inflation Device

Directions for Use	2
Instrucciones de uso	4
Mode d'emploi	6
Gebrauchsanweisung	8
Istruzioni per l'uso	10
Gebruiksaanwijzing	12
Instruções de Utilização	14

90998329-01

2014-11

Encore™ 26

Inflation Device

Rx ONLY

Caution: Federal Law (USA) restricts this device to sale by or on the order of a physician.

WARNING

Contents supplied STERILE using an ethylene oxide (EO) process. Do not use if sterile barrier is damaged. If damage is found, call your Boston Scientific representative.

For single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

After use, dispose of product and packaging in accordance with hospital, administrative and/or local government policy.

DEVICE DESCRIPTION

The Encore 26 Inflation Device employs a threaded plunger and locking mechanism activated by a finger latch. When the finger latch is depressed, the threaded plunger is unlocked and may be advanced or withdrawn as required. When the finger latch is released, the device is in the locked position. This allows for the generation and monitoring of pressure in atmospheres (atm) and kilopascal (kPa). The device incorporates a pressure gauge with a 0-26 atm/bar (0-2634 kPa) dial, a 20 ml (20 cc) syringe, and a connecting tube. The gauge accuracy is 3-2-3, indicating an accuracy of $\pm 3\%$ in the upper and lower quarter of the scale and is $\pm 2\%$ in the middle half of the scale.

Contents

Inflation Device

INDICATIONS FOR USE/INTENDED USE

The Encore 26 Inflation Device is recommended for use with balloon dilatation catheters to create and monitor pressure in the balloon, and to deflate the balloon.

CONTRAINDICATIONS

None known.

WARNINGS

None known.

PRECAUTIONS

None known.

HOW SUPPLIED

Do not use if package is opened or damaged.


Do not use if labeling is incomplete or illegible.

Handling and Storage

Store in a cool, dry, dark place.

DIRECTIONS FOR USE

1. Hold the device in one hand and depress the finger latch to unlock or disengage the threaded plunger. **(Figure 1)**.
2. Pull back on the threaded plunger to fill the syringe with contrast and saline (see catheter manufacturer's recommendations).
3. Aspirate 5 ml (5 cc) to 8 ml (8 cc) into syringe, then hold upright and eliminate any air that may have collected in the syringe and tubing.
4. Place a meniscus of fluid on the balloon luer port of the prepared catheter and attach the connecting tube to the luer port.
5. A negative pressure may be applied to the catheter to maintain deflated balloon profile by pulling the plunger back to the desired volume and locking it in place by releasing the finger latch **(Figure 2)**.
6. With the finger latch in the locked position (out), pressure is increased by rotating the plunger handle in a clockwise direction. Pressure can be decreased instantaneously by depressing the finger latch to the unlocked position (in) while slightly rotating the plunger handle.

WARRANTY

Boston Scientific Corporation (BSC) warrants that reasonable care has been used in the design and manufacture of this instrument. **This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether express or implied by operation of law or otherwise, including, but not limited to, any implied warranties of merchantability or fitness for a particular purpose.** Handling, storage, cleaning and sterilization of this instrument as well as other factors relating to the patient, diagnosis, treatment, surgical procedures and other matters beyond BSC's control directly affect the instrument and the results obtained from its use. BSC's obligation under this warranty is limited to the repair or replacement of this instrument and BSC shall not be liable for any incidental or consequential loss, damage or expense directly or indirectly arising from the use of this instrument. BSC neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this instrument. **BSC assumes no liability with respect to instruments reused, reprocessed or resterilized and makes no warranties, express or implied, including but not limited to merchantability or fitness for a particular purpose, with respect to such instruments.**

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All rights reserved.
Boston Scientific Limited
Boston Scientific (Australia) Pty Ltd
© 2014 Boston Scientific Corporation or its affiliates.