# LithoVue™

Single-Use Digital Flexible  
Ureteroscope

en	Instructions for Use	2
es	Instrucciones de uso	13
fr	Mode d'emploi	24
de	Gebrauchsanweisung	35
it	Istruzioni per l'uso	46
nl	Gebruiksaanwijzing	57
pt	Instruções de Utilização	68

51116905-01

2023-04

## TABLE OF CONTENTS

<b>REUSE WARNING</b> .....	3
<b>PRODUCT DESCRIPTION</b> .....	3
Operating Principle .....	3
<b>DEVICE DESCRIPTION</b> .....	3
LithoVue Device Model Numbers .....	3
LithoVue Device Specifications .....	3
Precaution .....	3
Figure 1. LithoVue Single-Use Digital Flexible Ureteroscope .....	4
User Training, Knowledge, and Skills .....	4
Contents .....	4
Ureteroscope Features and Controls .....	4
Figure 2. Main features of the LithoVue Single-Use Digital Flexible Ureteroscope .....	4
<b>STANDARD DEFLECTION</b> .....	5
<b>REVERSE DEFLECTION</b> .....	5
Figure 3. Articulation control can be standard or reversed depending upon model used .....	5
Figure 4. Distal Tip – The distal tip includes 3 elements: a.) the opening of the working channel, b.) illumination optics, and c.) the digital imaging sensor .....	5
Conformance to Standards .....	5
<b>INTENDED USE/INDICATIONS FOR USE</b> .....	6
Clinical Benefit Statement .....	6
<b>CONTRAINDICATIONS</b> .....	6
<b>WARNINGS</b> .....	6
<b>PRECAUTIONS</b> .....	7
Figure 5. DO NOT BEND SHARPLY .....	7
<b>ADVERSE EVENTS</b> .....	8
<b>HOW SUPPLIED</b> .....	8
Device Details .....	8
Handling and Storage .....	8
<b>DEVICE PREPARATION</b> .....	8
Unpacking .....	8
Inspection and Operational Checks .....	8
<b>OPERATIONAL INSTRUCTIONS</b> .....	8
Ureteroscope Compatibility .....	8
Fluid Compatibility .....	9
Setup .....	9
How to Articulate the Distal Tip .....	9
How to Connect an Irrigation/Contrast Source .....	9
<b>ACCESS, VISUALIZATION AND APPLICATION OF THERAPY</b> .....	9
Access and Visualization .....	9
Figure 6. Articulation Section .....	10
Inserting and Removing an Accessory and Application of Therapy .....	10
Removing the Ureteroscope from the Patient .....	10
Post-Procedure .....	11
<b>DISPOSAL OF THE URETEROSCOPE AND PACKING MATERIALS</b> .....	11
<b>TROUBLESHOOTING</b> .....	11
Loss of Image .....	11
Loss of Articulation Control .....	11
<b>INFORMATION TO BRIEF THE PATIENT</b> .....	12
<b>WARRANTY</b> .....	12
<b>SYMBOL DEFINITIONS</b> .....	12

# LithoVue™

## Single-Use Digital Flexible Ureteroscope

### Rx ONLY

**Caution:** Federal Law (USA) restricts this device to sale by or on the order of a physician.

#### REUSE WARNING

For single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure, which in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

#### PRODUCT DESCRIPTION

The LithoVue Single-Use Digital Flexible Ureteroscope is a sterile, single-use device comprised of two main components: a handle with articulation controls and accessory access ports, and a flexible shaft portion.

The LithoVue Single-Use Digital Flexible Ureteroscope is also referred to as the LithoVue Flexscope or ureteroscope in these instructions.

#### Operating Principle

The LithoVue Flexscope is used by physicians to access, visualize, and perform procedures in the urinary tract. The ureteroscope enables delivery and use of accessories such as biopsy forceps, laser fibers, guidewires, graspers and retrieval baskets at a surgical site. The distal tip of the ureteroscope articulates to 270 degrees in two directions, and the distal tip can be rotated a total of 360 degrees by rotating the handle. In addition, the shaft of the ureteroscope features a secondary deflection that is passive.

#### DEVICE DESCRIPTION

The flexible shaft of the ureteroscope includes one working channel. The working channel enables the delivery of surgical accessories and irrigation/contrast solutions to the distal tip and surgical field.

The handle includes an articulation lever that controls articulation of the distal tip. The handle also includes two access ports. One port accepts surgical accessories and one port provides a connection point for irrigation/contrast solutions.

#### LithoVue Device Model Numbers

Standard Deflection Model	M0067913500
Reverse Deflection Model	M0067913600


#### LithoVue Device Specifications

Specification	Measure
Optical Working Distance	2 mm - 50 mm
Direction of View	0 degrees
Field of View	85 degrees in air
Insertion Portion Width (Distal Face)	7.7F (2.57 mm)
Maximum Insertion Portion Width (Overall Shaft Diameter)	9.5F (nominal) 3.23 mm (maximum)
Shaft Working Length	68 cm
Minimum Instrument Channel Width	3.6F (nominal) 1.15 mm (minimum)
Instrument Channel Length	82 cm
Maximum Angle of Deflection	270 degrees
Minimum Bend Radius	8 mm
Degree of Protection Against Ingress of Liquids, LithoVue Flexscope	IPX2

#### Precaution

- There is no guarantee that instruments selected solely using the minimum instrument channel width will be compatible with the LithoVue Flexscope.

- There is no guarantee that instruments selected solely using the maximum insertion portion width and working length will be compatible with the LithoVue Flexscope.

**Figure 1. LithoVue Single-Use Digital Flexible Ureteroscope**

### User Training, Knowledge, and Skills

Read this IFU and the LithoVue System Workstation User's Manual before using the LithoVue System. Study labeling thoroughly for safe handling and storage. Use the ureteroscope as intended.

Only physicians who have ureteroscopic diagnostic and therapeutic training should use the LithoVue Flexscope. A thorough understanding of the techniques, principles, clinical applications and risks associated with ureteroscopy is necessary before using the device. This includes, but is not limited to: devices used in laser lithotripsy, administration of irrigation/contrast, and use of procedural accessories such as baskets/forceps, access sheaths, renal sheaths, guidewires and antiretropulsion devices.

### Contents

LithoVue Single-Use Digital Flexible Ureteroscope

### Ureteroscope Features and Controls

The features and controls on the ureteroscope include the following:


**Handle** – The handle enables the user to insert, withdraw and rotate the flexible shaft (Figure 2). The handle includes an articulation lever, an accessory access port, an irrigation port, and a connector cable with plug. The accessory access and irrigation ports, respectively, enable the insertion of surgical accessories and the introduction of irrigation/contrast fluid through the flexible shaft.

**Articulation lever** – The articulation lever controls the up and down articulation of the distal tip. Depending on the model, the articulation control can be standard (lever up = tip up) or reverse (lever up = tip down) (Figure 3). Returning the articulation lever to the Neutral position causes the distal tip to return to its unarticulated position.

**Cable and plug** – The connector cable transmits power from the LithoVue System Workstation to the ureteroscope and transmits video signals from the ureteroscope to the LithoVue System Workstation.

**Accessory access port** – The accessory access port uses a 6 % taper luer-type connector and is the point of access for delivering surgical accessories into the working channel of the ureteroscope.

**Irrigation port** – The irrigation port uses a 6 % taper luer-type connector and is the point of access for delivering irrigation/contrast solutions to the working channel of the ureteroscope.


**Figure 2. Main features of the LithoVue Single-Use Digital Flexible Ureteroscope**


**Articulation section** – This is the portion of the flexible shaft that can be deflected using the articulation lever.

**STANDARD DEFLECTION**


**REVERSE DEFLECTION**

**Figure 3. Articulation control can be standard or reversed depending upon model used**

**Flexible shaft** – The flexible shaft is the insertion portion of the ureteroscope (See Figure 4).

It includes: (a) a working channel for delivery of surgical accessories (surgical tools) and irrigation/contrast solutions to the surgical field, (b) an optical conduit that delivers light to the surgical field, and (c) a digital imaging sensor at the distal tip to transmit real time images of the operative field to the LithoVue System Workstation. The working channel is 3.6F (1.2 mm) in diameter.

**Figure 4. Distal Tip – The distal tip includes 3 elements: a.) the opening of the working channel, b.) illumination optics, and c.) the digital imaging sensor.**

**Conformance to Standards**

The LithoVue Flexible Ureteroscope complies with IEC-60601, and CAN/CSA C22.2 No 60601-1.

## INTENDED USE/INDICATIONS FOR USE

The LithoVue System is intended to be used to visualize organs, cavities and canals in the urinary tract (urethra, bladder, ureter, calyces and renal papillae) via transurethral or percutaneous access routes. It can also be used in conjunction with endoscopic accessories to perform various diagnostic and therapeutic procedures in the urinary tract.

### Clinical Benefit Statement

LithoVue is a single-use ureteroscope, which is a supportive device to endoscopic diagnosis and treatments and contributes to procedural benefits by enabling navigation to and visualization of the organs, cavities, and canals in the urinary tract via transurethral or percutaneous access routes. It can also be used in conjunction with endoscopic accessories to perform various diagnostic and therapeutic procedures in the urinary tract.

### CONTRAINDICATIONS

Contraindications for this device are those specific to urinary tract endoscopy.

Diagnostic or therapeutic ureteroscopy is contraindicated in people with an untreated urinary tract infection.

Other contraindications to therapeutic ureteroscopy (e.g., lithotripsy, endopyelotomy, tumor therapy) are more numerous and can mirror those associated with the corresponding open surgical interventions. Patients on anticoagulants or with coagulopathies should be managed appropriately.

### WARNINGS

- Do not use electromedical energy sources in the presence of flammable detergents, anesthetics, nitrous oxide (N<sub>2</sub>O), or oxygen.
- Consult the operating manuals of all electromedical energy sources used with endoscopic instruments for appropriate instructions, warnings, and cautions prior to use. Such sources of energy include electrical, electrohydraulic, electrosurgical, heat, hydraulic, laser, light, pressure, sound, ultrasound, and vacuum.
- Do not insert or advance the ureteroscope unless there is a clear live endoscopic view of the lumen through which the scope is being advanced (or confirm with visualization by other imaging modalities). Doing so can cause patient injury such as perforation, avulsion, hemorrhage, or urothelial damage.
- In the event that the live endoscopic image is lost, do not advance or insert the ureteroscope and do not insert, advance or actuate accessories. Doing so can cause accessory damage or patient injury such as perforation, avulsion, hemorrhage or urothelial damage. See troubleshooting section for more information.
- Do not use excessive force while advancing or withdrawing the scope. Doing so can cause patient injury such as perforation, hemorrhage or urothelial damage or damage to the ureteroscope. If resistance is felt during advancement or withdrawal of the scope, investigate the source of resistance and take remedial action (e.g., fluoroscopy, contrast injection).
- Do not force the distal tip of the ureteroscope against the sidewall of the ureter or renal pelvis. Doing so can cause patient injury such as perforation, avulsion, hemorrhage or urothelial damage.
- Do not use excessive force when advancing or withdrawing an accessory within the ureteroscope. Doing so can cause patient injury such as perforation, avulsion, hemorrhage, urothelial damage or damage to the ureteroscope.
- When inserting or using accessories, maintain continuous visualization of the distal tip. Ensure that the distance between the distal tip of the ureteroscope and the object in view is greater than the ureteroscope's minimum visible distance. Failure to do so may result in accessory damage or patient injury such as perforation, hemorrhage or urothelial damage.
- Do not withdraw a laser fiber back into the ureteroscope while the laser is firing. Doing so may cause patient injury and/or scope damage.
- Do not look directly into the light emitted from the ureteroscope. Doing so can result in eye injury.
- Verify ground isolation when setting up and using accessories from different manufacturers. Failure to do so can result in shocks and accessory malfunction leading to patient injury.
- Do not open the handle of the ureteroscope. Doing so can compromise fluid ingress protection and result in risk of electric shock.
- The ureteroscope is a single-use device and there are no serviceable parts. Do not repair damaged or non-operating ureteroscopes. Do not use the ureteroscope if damage is discovered or suspected.
- Do not excessively bend the flexible shaft or the articulating section of the ureteroscope as this may break or kink the shaft.
- If damage to the ureteroscope occurs or it stops functioning during a procedure, stop using the ureteroscope immediately. See troubleshooting section for more information. Continue the procedure with a new ureteroscope, as appropriate.

- Prior to use of a cardiac defibrillator, remove the ureteroscope from the patient. Failure to remove the ureteroscope from a patient during use of a cardiac defibrillator could result in patient injury or damage to the system due to the discharge of the cardiac defibrillator.

#### PRECAUTIONS

- Only use the LithoVue Single-Use Digital Flexible Ureteroscope in conjunction with the LithoVue System Workstation. Connection to other devices may cause device or property damage or operator injury.
- Only physicians with adequate ureteroscopic training should perform procedures with the ureteroscope. Consult the medical literature regarding techniques, complications, and hazards prior to any procedure.
- Use the ureteroscope with caution in patients who have undergone previous urinary tract reconstructive surgery or with known strictures. These conditions may prevent passage of the flexible scope shaft.
- The LithoVue System has been tested and shown to be compatible with laser lithotripsy devices. The use of other energized procedural devices may cause loss of image, device damage, or patient injury.
- To ensure satisfactory performance, perform the prescribed inspections and operational checks described in the device preparation section before use.
- The distal tip of the ureteroscope should be straight when inserting and withdrawing accessories. Follow the accessory instructions for use regarding inserting the accessory into a flexible ureteroscope. Failure to follow the accessory instructions for use may result in patient injury, or accessory and ureteroscope damage or malfunction.
- The ureteroscope features a strain relief at the transition from the handle to the shaft. The strain relief protects the device during use. To prevent damage to the shaft and/or light fiber, do not bend the shaft sharply (Figure 5).


Figure 5. DO NOT BEND SHARPLY

- Use only those fluids/lubricants recommended in the ureteroscope compatibility section.
- Do not use accessories that fail to meet the compatibility requirements of the ureteroscope as stated on the ureteroscope labeling. Doing so may cause damage to the ureteroscope and/or accessory.
- Failure to thoroughly understand and follow all instructions, cautions and warnings provided in this IFU and the LithoVue System Workstation User's Manual may result in injury to the patient and/or user; and/or may result in damage to, or malfunction of, this equipment. Additionally, damage to other equipment or property may result. Follow all instructions, cautions and warnings provided with all products and equipment to be used in conjunction with the LithoVue System to avoid any possible hazards due to device incompatibility.
- Placing the LithoVue System Workstation near other medical electrical equipment may result in electromagnetic interference (EMI) which may degrade the video image. Additionally, EMI from the LithoVue System Workstation may interfere with other equipment in the operating room. Verify that all electrical equipment is working properly before starting the procedure. Failure to do so may result in equipment not working, possibly resulting in either a delay of the procedure or an adverse event.
- Do not insert a wet, contaminated, or damaged connector cable plug into the LithoVue System Workstation as poor video performance or damage to the system may result.

- Do not remove the connector cable from the LithoVue System Workstation by pulling on the cable as poor video performance or damage to the system may result. Slide the locking collar on the cable plug toward the cable and pull the plug out to remove the cable.

#### ADVERSE EVENTS

Possible complications include, but may not be limited to:

- Bleeding
- Avulsion
- Stenosis/Stricture
- Laceration
- Pain
- Discomfort
- Urinoma
- Extravasation
- Fever
- Sepsis
- Inflammation
- Perforation (ureter, renal pelvis or bladder)
- Ureteral Reflux
- Hematoma
- Urothelial damage
- Infection

#### HOW SUPPLIED

##### Device Details

Package contents include one LithoVue Single-Use Digital Flexible Ureteroscope. The ureteroscope is supplied STERILE using an ethylene oxide (EO) process in a sealed pouch and for single use only. One device in pouch and one IFU is packaged in a carton.

Do not use if sterile barrier is damaged.

Do not use if package is damaged or unintentionally opened before use.

Do not use if labeling is incomplete or illegible.

##### Handling and Storage

This product has no special handling or storage requirements.

#### DEVICE PREPARATION

##### Unpacking

Upon receipt, examine the shipping carton and its contents for signs of damage. Confirm that the sterile barrier is intact and check the expiration date. Do not use a damaged or expired ureteroscope. Do not attempt to repair the ureteroscope.

Inside the shipping carton, the ureteroscope is packaged in a tray which is sealed within a sterile barrier. Open the sterile barrier by peeling back the header on the handle side of the bag. Remove the ureteroscope and tray together. Remove the ureteroscope from the tray by gripping the handle and connector cable and lifting the ureteroscope out of the tray while gently sliding the ureteroscope shaft out of the enclosed portion of the tray. It is not necessary to open the enclosed portion of the tray to remove the device. After removing the ureteroscope from the packaging, perform the following operational checks.

##### Inspection and Operational Checks

**WARNING Failure to perform inspection and operational checks may result in patient injury and/or damage to the device and accessories.**

- Inspect the entire surface of the flexible shaft visually and by feeling the shaft with your fingertips. Examine the handle, articulation lever, and access ports. Ensure no components are loose, bent or broken.
- Visually inspect the distal tip for damage including dents, protrusions, tears and holes.
- Examine the ureteroscope in up, down, and straight articulated modes. Confirm the plane of articulation. Confirm that the scope tip articulates smoothly in both directions. DO NOT force the flexible tip into a straight or flexed position while holding the articulation lever. Doing so can damage the control mechanism.
- Visually examine the connector cable and plug for damage or contamination.

#### OPERATIONAL INSTRUCTIONS

##### Ureteroscope Compatibility

The recommended media for fluoroscopic visualization is radiopaque contrast diluted per normal procedure.

The LithoVue Single-Use Digital Flexible Ureteroscope is compatible with the following ancillary devices and accessories:

- LithoVue System Workstation
- Accessory devices with a minimum working length of 82 cm (32.3 in) and compatible with a 3.6F (1.2 mm) working channel diameter
- Irrigation/contrast pumps, syringe, or hanging bag with a maximum flow rate of 650 ml/min at a maximum head pressure of 40 psi (276 kPa) and a luer type fitting on the supply tubing
- Guidewires with a maximum outside diameter of 0.97 mm (0.038 in)

- Laser lithotripsy systems featuring aiming beams with adjustable intensity or pulse option

#### Fluid Compatibility

- The LithoVue Flexscope is known to withstand exposure to the following procedural fluids: blood, urine, saline, contrast media, petroleum-based and water-based lubricants. In the event of poor image quality due to debris covering the distal tip of the LithoVue Flexscope, a 15 to 70 percent alcohol in purified water solution may be used to wipe the distal tip with a cotton swab.

#### Setup

1. Setup and power up the LithoVue System Workstation following the LithoVue System Workstation User's Manual.
2. Insert the connector cable plug from the ureteroscope into the receptacle on the front of the LithoVue System Workstation until the cable plug is fully inserted.

---

**CAUTION:** The ureteroscope cable connector should connect to the cable receptacle easily. If the cable connector does not connect, verify that the arrows on the cable connector and cable receptacle line up. Forcing the cable connector into the cable receptacle may damage the scope and/or System Workstation.

---

3. Hold the distal tip of the ureteroscope close to an object (about the same distance as the width of the shaft) and ensure the video monitor displays a clear, sharply-focused, high-quality image.
4. Verify light output at the tip and adjust the brightness, as necessary. (See the LithoVue System Workstation User's Manual).

---

**Note:** White balancing is automatically done by the LithoVue System Workstation. Manual white balancing is not required.

---

5. Connect a luer-type cap or luer-type sealing device to the accessory access port on the handle to prevent fluid from leaking from the port during the procedure. If desired, connect the irrigation/contrast supply tubing from a compatible irrigation/contrast source to the irrigation port on the handle of the ureteroscope.

#### How to Articulate the Distal Tip

To articulate the distal tip, move the articulation lever on the handle with your thumb.

#### Note:

- For standard-type ureteroscopes, moving the articulation lever up moves the distal tip up. Moving the articulation lever down moves the distal tip down (Figure 3).
- For reverse-type ureteroscopes, moving the articulation lever up moves the distal tip down. Moving the articulation lever down moves the distal tip up (Figure 3).

#### How to Connect an Irrigation/Contrast Source

Connect the irrigation/contrast supply tubing of a compatible irrigation/contrast source (syringe, gravity-fed bag, or pump) to the irrigation port on the ureteroscope using a luer-type fitting (Figure 2).

---

**Note:** Connect a luer-type cap or luer-type sealing device to the accessory access port on the handle to prevent fluid from leaking from the port prior to beginning flow of irrigation into the ureteroscope.

---

**CAUTION:** The recommended media for fluoroscopic visualization is radiopaque contrast diluted per normal procedure.

---

### ACCESS, VISUALIZATION AND APPLICATION OF THERAPY

#### Access and Visualization

1. Gently advance the scope per standard practice to the desired treatment area.

---

**CAUTION:** If using an access sheath or renal sheath, do not articulate the distal tip inside the sheath.

---

**CAUTION:** Do not apply excessive force to the articulation lever.

---


**CAUTION:** Grasp the ureteroscope shaft close to the urethral meatus and advance with short strokes to prevent buckling of the shaft. If passing the scope through a nephrostomy cannula grasp the ureteroscope shaft close to the cannula and advance with short strokes.

---

- Utilize imaging to confirm position and articulation of the distal tip of the ureteroscope.

**CAUTION:** Confirm via the live endoscopic video image or other imaging modality that the distal tip of the ureteroscope is advancing when the ureteroscope shaft is being guided through the anatomy.

**Note:** The articulation section of the flexible shaft is the most radiopaque portion of the device. The distal tip extends approximately 5.5 mm beyond the end of the articulation section.

**Figure 6. Articulation Section**

#### Inserting and Removing an Accessory and Application of Therapy

**Note:** If using a laser lithotripsy system, select a system which features an aiming beam with adjustable intensity or pulse option. Adjust aiming beam intensity or pulsation to allow clear live video image of the lithotripsy target. Do not actuate any accessory device, including laser lithotripsy devices, without a clear endoscopic live video image. If the laser aiming beam cannot be adjusted to provide a satisfactory image, abort the procedure.

- Read, understand and follow all instructions for use supplied with accessory devices used in conjunction with the LithoVue Flexscope. Prior to advancing an accessory, ensure that the tip of the scope is not in contact with the sidewall of the collecting system. Failure to leave a space between the tip of the scope and the sidewall may result in an unintended perforation when the instrument is advanced beyond the tip of the scope.
- Follow the accessory instructions for use regarding the proper scope and accessory configuration for inserting and advancing that accessory through a flexible ureteroscope.
- Insert the desired accessory device into the ureteroscope using the accessory access port (Figure 2). While doing this, be sure that the scope is stabilized to prevent inadvertent movement which might cause injury to the patient or ureteroscope.
- Advance the accessory slowly while observing the live video image for the initial entry of the accessory into the field of view.
- Once the accessory is in the field of view, use the articulation lever to articulate the tip and accessory to achieve treatment objective.

**Note:** If you encounter resistance when advancing or withdrawing an accessory through the ureteroscope, visually confirm that the tip of the ureteroscope is not in contact with tissue. Then, ensure the articulation lever is in the Neutral position so that the distal tip of the ureteroscope is in its straight, non-articulated position. Also, ensure the accessory is in the proper configuration for passage through the ureteroscope. For example, confirm that the stone basket is in a closed configuration, the shaft is not kinked, etc. If you still encounter resistance, carefully withdraw the accessory, inspect it for any defects and if none are found, reintroduce it.

- Before withdrawing an accessory, visually confirm that it is in the proper configuration for passage through the ureteroscope before pulling it back.

#### Removing the Ureteroscope from the Patient

- Return the articulation lever to the Neutral position to allow the ureteroscope distal tip to relax into its straight, non-articulated position.

**CAUTION:** Ensure the articulation lever is in the Neutral position when removing the scope from the patient.

2. Under direct visualization, slowly withdraw the ureteroscope from the patient. If you feel resistance, investigate the source of the resistance using another imaging modality before continuing to withdraw the ureteroscope.

#### Post-Procedure

1. Unplug the ureteroscope from the LithoVue System Workstation by grasping the grey collar on the cable plug and pulling it backwards toward the user.
2. Dispose of the ureteroscope. (See "DISPOSAL OF THE URETEROSCOPE AND PACKING MATERIALS").

Any serious incident that occurs in relation to this device should be reported to the manufacturer and relevant local regulatory authority.

#### DISPOSAL OF THE URETEROSCOPE AND PACKING MATERIALS

To minimize the risk of infection or microbial hazards after use, dispose device and packaging as follows:

After use, device may contain biohazardous substances. The device and packaging should be treated and disposed of as biohazardous waste or have them treated and disposed of in accordance with any applicable hospital, administrative, and/or local government regulations. Use of a biohazardous container with biological hazard symbol is recommended. Untreated biohazardous waste should not be disposed of in the municipal waste system.

#### TROUBLESHOOTING

##### Loss of Image

---

**CAUTION:** In the event that the live video image becomes lost or unusable, during the procedure, **STOP - DO NOT ADVANCE THE URETEROSCOPE, DO NOT INSERT, ADVANCE, OR ACTUATE PROCEDURAL DEVICES IN THE WORKING CHANNEL.** Blind manipulation of the ureteroscope can result in patient injury and/or scope damage.

---

If a User Message is present on the LithoVue System Workstation, refer to the LithoVue System Workstation User's Manual for the potential cause and troubleshooting steps. If no user message is displayed, conduct the following steps to correct the issue:

1. Disconnect the Ureteroscope from the LithoVue System Workstation.
2. Ensure that the Ureteroscope cable and plug and LithoVue System Workstation receptacle are undamaged.
3. Ensure that the Ureteroscope plug and the LithoVue System Workstation receptacle are clean and dry.
4. Reconnect the Ureteroscope to the LithoVue System Workstation, firmly pushing in the connector until it is fully inserted.
5. Check for live video image. If video image is unusable, continue to Step 6.
6. Power-off the LithoVue System Workstation.
7. Power-on the LithoVue System Workstation.

If live image cannot be restored, discontinue use of the ureteroscope. Conduct the following steps to remove the ureteroscope:

1. Return the articulation lever to Neutral position – this will cause the distal tip to return to an unarticulated position. The articulation state of the ureteroscope can be confirmed with imaging.
2. Slowly withdraw the ureteroscope from the patient. If you feel resistance, investigate the source of resistance with imaging before continuing to withdraw the ureteroscope.
3. Dispose of the ureteroscope (See "DISPOSAL OF THE URETEROSCOPE AND PACKING MATERIALS") and use another scope to complete the procedure.

##### Loss of Articulation Control

---

**CAUTION:** In the event of a loss of articulation control, **STOP - DO NOT ADVANCE THE URETEROSCOPE, DO NOT INSERT, ADVANCE OR ACTUATE PROCEDURAL DEVICES IN THE WORKING CHANNEL.**

---

In the event that the articulation control becomes unresponsive, discontinue use of the Ureteroscope by conducting the following steps:

1. Return the articulation lever to the Neutral position.
2. Use imaging to confirm that the distal tip is in a non-articulated state.
3. If the distal tip remains articulated, utilize the articulation lever and imaging to straighten the tip.
4. If the articulation lever cannot be used to straighten the distal tip, an 0.035in (0.89 mm) or larger super stiff guidewire may be passed, floppy end first, through the working channel to aid in straightening the distal tip using the following steps:
  - a. Stabilize the handle of the ureteroscope and insert the guidewire into the ureteroscope using the accessory access port (Figure 2).

- b. Advance the guidewire slowly, observing the live video image for initial entry of the guidewire into the field of view.
  - c. Stop advancing the guidewire when the stiff portion of the guidewire emerges from the distal tip and is visible on the live video image.
5. Use fluoroscopy or other imaging modality to confirm that the ureteroscope shaft has straightened.
  6. Slowly withdraw the ureteroscope from the patient. Monitor the withdrawal using spot fluoroscopic imaging or other imaging modality. If you feel resistance, investigate the source of resistance and take remedial action before continuing to withdraw the ureteroscope.
  7. Dispose of the ureteroscope (See "DISPOSAL OF THE URETEROSCOPE AND PACKING MATERIALS").

#### INFORMATION TO BRIEF THE PATIENT

Prior to use of this device, brief the patient on any contraindications, warnings, precautions, adverse events, and post-procedure instructions outlined in this document that are of relevance to the patient. Potential adverse events associated with ureteroscopy procedures include avulsion, bleeding, fever, discomfort/pain, hematoma, infection, inflammation, laceration, perforation, sepsis, stenosis/stricture, ureteral reflux, urinoma, extravasation, and urothelial damage.

#### WARRANTY

For device warranty information, visit ([www.bostonscientific.com/warranty](http://www.bostonscientific.com/warranty)).

EU Importer: Boston Scientific International B.V., Vestastraat 6, 6468 EX Kerkrade, The Netherlands

LithoVue is a registered trademark of Boston Scientific Corporation or its affiliates.

All other trademarks are the property of their respective owners.

#### SYMBOL DEFINITIONS

Commonly used medical device symbols that appear on the labeling are defined at [www.bostonscientific.com/SymbolsGlossary](http://www.bostonscientific.com/SymbolsGlossary).

Additional symbols are defined at the end of this document.

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## INSTRUCTIONS FOR USE

en

The Instructions for Use (IFU) for this product are supplied in electronic form over the Internet.

Visit [www.IFU-BSCI.com](http://www.IFU-BSCI.com) to access the IFU using your internet browser. Information from the product labeling such as Product Code, UDI-DI or Product Name may be used to search for the specific IFU.

If you have difficulty accessing the IFU online, or would prefer to receive a paper copy, please contact Boston Scientific Customer Service or your local country contact. A copy will be sent to you at no charge and should arrive within seven days.

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All rights reserved.
Boston Scientific Limited
Boston Scientific (Australia) Pty Ltd
© 2023 Boston Scientific Corporation or its affiliates.