


# TrapEase™ Polyp Traps



Designed for safe specimen handling, easy and convenient removal of Polyps and to help prevent loss of Polyps during suction.

# TrapEase™ Polyp Traps

## Quad Trap



**Four Removable Polyp Chambers:** Designed to maximize efficiency by limiting the need for specimen transfer during a procedure




**Parabola Slot Design:** Chambers stay in place during procedure and remove easily for Polyp transfer



**Safety Screen:** Designed to prevent loss of specimen and allows collection of a fifth specimen



**Large Numbered Chambers**



Formalin Solution 10% Neutral Buffered  
CAUTION! Contains Formaldehyde, which is a  
Potential Cancer Hazard. Repeated or prolonged  
exposure increases the risk. Toxic by ingestion,  
and if swallowed. Irritating to eyes, respiratory  
system & skin. May cause sensitization by  
inhalation or skin contact. Risk of severe  
damage to eyes. May be FATAL if swallowed.  
Target Organs: respiratory system, eyes, skin.  
For more information see MSDS. For Laboratory  
Use Only! Formaldehyde CAS #50-00-0.

# TrapEase™ Polyp Traps

## Single Trap

**Single Chamber**

**Convenient Removal:**

Lift Polyps directly off the screen

**Safety Screen:** Designed to prevent loss of specimen



Quad Trap

Single Trap

Complementary Devices

Ordering Information




# TrapEase™ Polyp Traps


## Complementary Devices

Boston Scientific offers a wide portfolio of products for tissue resection and Polypectomy


**A Comprehensive Offering of Snares:**  
Including the Captivator™ II Single-Use Snare



**Retrieval Devices:**  
Including the RescueNet™ Retrieval Net



**Complication Management:**  
Including the Resolution™ and Resolution 360™ Clip





# TrapEase™ Polyp Traps

## Ordering Information

Order Number	Description	Packaging
STE-298-25	TrapEase Quad Trap	Box 25
STE-297-50	TrapEase Single Trap	Box 50

For more detailed information on devices, please visit [www.BostonScientific.com/gastro](http://www.BostonScientific.com/gastro)





**Boston  
Scientific**  
Advancing science for life™

All images taken by Boston Scientific.

Registered trademarks are the property of their respective owners.

Indications, contraindications, warnings and instructions for use can be found in the product labeling supplied with each device.

**CAUTION: Federal (USA) law restricts these devices to sale by or on the order of a physician.**

**CAUTION: The law restricts these devices to sale by or on the order of a physician. Information for use only in countries with applicable health authority registrations. Material not intended for use in France.**

Boston Scientific Corporation  
300 Boston Scientific Way  
Marlborough, MA 01752-1234  
[www.bostonscientific.com/gastro](http://www.bostonscientific.com/gastro)

©2017 Boston Scientific Corporation  
or its affiliates. All rights reserved.

ENDO-492405-AA August 2017