

Going PRO

**Boston
Scientific**
Advancing science for life™

What is the CRE™ PRO Wireguided Balloon Dilatation Catheter?

The CRE PRO Wireguided Balloon Dilatation Catheter is the CRE Wireguided Balloon Dilatation Catheter PLUS embedded platinum iridium RO markers and a streamlined profile for optimal passability.

- **Optimized Passability**—Made with a new balloon folding process for a lower profile leading to consistent ease of passage through the scope.
- **Simplified Scope Compatibility**—Compatible with 2.8 mm+ working channel endoscopes and 3.7 mm+ working channel duodenoscopes. Smaller sizes (6-8, 8-10, 10-15mm) are also compatible with 3.2 mm working channel duodenoscopes.
- **Fluoroscopic Visualization**—Embedded RO markers to facilitate balloon placement using fluoroscopy

CRE PRO Wireguided Balloon Dilatation Catheter was the first balloon dilator indicated for use in the DASE procedure for removal of difficult biliary stones.

To learn more about the CRE PRO Wireguided Balloon Dilatation Catheter visit www.bostonscientific.com/gastro

CRE™ PRO Wireguided Balloon Dilatation Catheter Ordering Information

Boston Scientific
Advancing science for life™

CRE PRO Wireguided Esophageal/Pyloric/Biliary Balloon Dilatation Catheters

Order Number	Balloon Inflated O.D.		Balloon Length (cm)	Inflation Pressures		Minimum Working Channel (mm)	Catheter Size (Fr)	Working Length (cm)
	(mm)	(Fr)		(ATM)	(kPa)			
M005 58600	6-7-8	18-21-24	5.5	3-6-10	304-608-1013	2.8/3.2*	7.5	180
M005 58610	8-9-10	24-27-30	5.5	3-5-9	304-557-912	2.8/3.2*	7.5	180
M005 58620	10-11-12	30-33-36	5.5	3-5-8	304-507-811	2.8/3.2*	7.5	180
M005 58630	12-13.5-15	36-40.5-45	5.5	3-4.5-8	304-456-811	2.8/3.2*	7.5	180
M005 58640	15-16.5-18	45-49.5-54	5.5	3-4.5-7	304-456-709	2.8/3.7*	7.5	180
M005 58650 **	18-19-20	54-57-60	5.5	3-4.5-6	304-456-608	2.8/3.7*	7.5	180

CRE PRO Wireguided Esophageal/Pyloric/Colonic/Biliary Balloon Dilatation Catheters

Order Number	Balloon Inflated O.D.		Balloon Length (cm)	Inflation Pressures		Minimum Working Channel (mm)	Catheter Size (Fr)	Working Length (cm)
	(mm)	(Fr)		(ATM)	(kPa)			
M005 58660	6-7-8	18-21-24	5.5	3-6-10	304-608-1013	2.8/3.2*	7.5	240
M005 58670	8-9-10	24-27-30	5.5	3-5-9	304-557-912	2.8/3.2*	7.5	240
M005 58680	10-11-12	30-33-36	5.5	3-5-8	304-507-811	2.8/3.2*	7.5	240
M005 58690	12-13.5-15	36-40.5-45	5.5	3-4.5-8	304-456-811	2.8/3.2*	7.5	240
M005 58700	15-16.5-18	45-49.5-54	5.5	3-4.5-7	304-456-709	2.8/3.7*	7.5	240
M005 58710 ***	18-19-20	54-57-60	5.5	3-4.5-6	304-456-608	2.8/3.7*	7.5	240

* Duodenoscope Channel Size

** Esophageal and Biliary application only.

*** Esophageal, Colonic and Biliary application only.

Recommended Inflation Device: Alliance™ II Integrated Inflation Handle Order #5062, Sterile Syringe Order #5060. Guidewire included and preloaded in balloon catheter.

All trademarks are the property of their respective owners. Indications, Contraindications, Warnings and Instructions for Use can be found in the product labeling supplied with each device.

CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician.

Boston Scientific Corporation
300 Boston Scientific Way
Marlborough, MA 01752
www.bostonscientific.com/gastro

©2015 Boston Scientific Corporation or its affiliates. All rights reserved.
END0-347803-AB December 2015