

| ImageReady™ MR-Conditional* Model List*

Therapy	Device Name	Device Model Numbers	MRI System Configuration
ICM	LUX-Dx™	M301	1.5T or 3T
	LUX-Dx II™	M302	
	LUX-Dx II+™	M312	
S-ICD	EMBLEM™	A209	1.5T
	EMBLEM™ MRI	A219	

The following S-ICD electrodes are labeled as MR-Conditional*
S-ICD ELECTRODES: 3010, 3400, 3401, 3501

Brady Pulse Generators			
Therapy	Device Name	Device Model Numbers	MRI System Configuration
Pacemakers	ACCOLADE™ MRI	L310, L311, L331	1.5T or 3T
	ESSENTIO™ MRI	L110, L111, L131	
	Proponent™ MRI	L210, L211, L231	
	VITALIO™ MRI ¹	K275, K277	1.5T
CRT-P	VISIONIST™ X4	U228	1.5T or 3T
	VALITUDE™ X4	U128	

The following Pacemaker and CRT-P leads and accessories are labeled as MR-Conditional*

Right Atrium and Right Ventricle

INGEVITY™ MRI: 7731, 7732, 7735, 7736, 7740, 7741, 7742
INGEVITY+™: 7840, 7841, 7842
Suture Sleeves for INGEVITY™ MRI/INGEVITY+™ leads: 6402
FINELINE™ II: 4469, 4470, 4471, 4472, 4473, 4474, 4479, 4480
 (Not valid with VITALIO MRI.)
FINELINE™ II: 4456, 4457, 4458, 4459 (Not valid with VITALIO MRI.)
Suture Sleeves for Fineline™ leads: 6220, 6221

Port Plugs: 7145

¹ VITALIO™ MRI is only MR-Conditional with INGEVITY MRI/INGEVITY+ leads only.

Left Ventricle

ACUITY™ X4: 4671, 4672, 4674, 4675, 4677, 4678
Suture Sleeves: 4603
Port Plugs: 7148

This document contains Boston Scientific CRM devices that are approved by the FDA as MR-Conditional as of January 2024.
 Visit <http://www.bostonscientific.com/imageready> for additional information cardiology/radiology checklists, patient resources, and the MRI technical manual.

All trademarks are the property of their respective owners. Rx only.

For a full list of Conditions of Use, and potential adverse events applicable when the Conditions of Use are met or not met, refer to the MRI Technical Guide. All items on the full list of Conditions of Use must be met in order for an MRI scan to be considered MRI Conditional.



ImageReady™ MR-Conditional* Model List*

Tachy Pulse Generators			
Therapy	Device Name	Device Model Numbers	MRI System Configuration
ICDs	RESONATE™ HF	D532, D533	1.5T or 3T
	RESONATE™	D432, D433	
	VIGILANT™	D232, D233	
	PERCIVA™ DF4	D412, D413	
	PERCIVA™ DF1	D400, D401	1.5T
	MOMENTUM™	D120, D121	
	AUTOGEN™	D160, D161, D162, D163	
	DYNAGEN™	D020, D021, D022, D023, D150, D151, D152, D153	
	INOGEN™	D010, D011, D012, D013, D140, D141, D142, D143	
CRT-Ds	RESONATE™ HF	G547	1.5T or 3T
	RESONATE™ X4	G447	
	VIGILANT™	G247	
	VIGILANT™	G248	1.5T
	MOMENTUM™	G124, G125, G138	1.5T
	AUTOGEN™	G160, G161, G166, G168	
	DYNAGEN™	G150, G151, G156, G158	
	INOGEN™	G140, G141, G146, G148	

The following ICD and CRT-D leads and accessories are labeled as MR-Conditional*

Right Atrium

INGEVITY™ MRI: 7735, 7736, 7740, 7741, 7742

INGEVITY+™: 7840, 7841, 7842

Suture Sleeves for INGEVITY™ MRI/INGEVITY+ leads: 6402

FINELINE™ II: 4469, 4470, 4471, 4472, 4473, 4474, 4479, 4480

Suture Sleeves for FINELINE leads: 6220, 6221

Port Plugs: 7145

Right Ventricle

ENDOTAK RELIANCE 4-SITE™ (DF4): 0262, 0263, 0265, 0266, 0272, 0273, 0275, 0276, 0282, 0283, 0285, 0286, 0292, 0293, 0295, 0296

ENDOTAK RELIANCE™ (DF1): 0127, 0128, 0129, 0137, 0138, 0139, 0143, 0147, 0148, 0149, 0153, 0157, 0158, 0159, 0170, 0171, 0172, 0173, 0174, 0175, 0176, 0177, 0180, 0181, 0182, 0183, 0184, 0185, 0186, 0187 (1.5 T only applies to all)

RELIANCE 4-FRONT™: 0636, 0650, 0651, 0652, 0653, 0654, 0655, 0657, 0658, 0662, 0663, 0665, 0672, 0673, 0675, 0676, 0682, 0683, 0685, 0686, 0692, 0693, 0695, 0696

Suture Sleeves for RELIANCE 4-FRONT™: 6403

DF-1 Lead Port Plug for ENDOTAK RELIANCE™ (IS-1) Leads- Single Coil: 6996 (1.5 T only)

*When conditions are met.

Left Ventricle

ACUITY™ X4: 4671, 4672, 4674, 4675, 4677, 4678

Suture Sleeves: 4603

IS4/DF4 Port Plugs: 7148

ACUITY™ Spiral: 4591, 4592, 4593 (1.5 T only)

Suture Sleeves for ACUITY Spiral Leads: 6100 (1.5 T only applies to all)

EASYTRAK™ 2: 4542, 4543, 4544 (1.5 T only)

Suture Sleeves for EASYTRAK™ 2 Leads: 6773 (1.5 T only applies to all)

IS-1 Port Plugs: 7145 (1.5 T only)



Advancing science for life™

Cardiology

300 Boston Scientific Way
Marlborough, MA 01752-1234
www.bostonscientific.com

Medical Professionals:
1.800.CARDIAC (227.3422)
Patients and Families:
1.866.484.3268

© 2024 Boston Scientific Corporation
or its affiliates. All rights reserved.

CRM-507108-AL

This document contains Boston Scientific CRM devices that are approved by the FDA as MR-Conditional as of January 2024. Visit <http://www.bostonscientific.com/imageready> for additional information cardiology/radiology checklists, patient resources, and the MRI technical manual.

All trademarks are the property of their respective owners. Rx only.