

Radial Jaw™ 4

Single-Use Biopsy Forceps

Available



with
needle



without
needle



Radial Jaw 4
Standard Capacity



Radial Jaw 4
Large Capacity



Radial Jaw 4
Jumbo



Radial Jaw 4
Gastropediatric

Samples Like Never Before™



Jaw Configuration Size

Designed to provide:

- ▶ Tissue specimens for excellent sample handling, preparation and accurate histological diagnosis
- ▶ Consistent sample retention for precise sample identification
- ▶ Enhanced passability through tortuous anatomy

Micro-Mesh Tooth Design

- ▶ Engineered for a clean, precise bite



Dual Pull Wire Design

Intended to:

- ▶ Increase bite force
- ▶ Enable the jaws to pivot onto the targeted biopsy site for accurate tissue acquisition



Positioning Markers

- ▶ Engineered to assist in the insertion and withdrawal of forceps to enhance procedural safety and efficiency

Streamlined Catheter

Designed to provide:

- ▶ The right balance of columnar strength and flexibility for excellent pushability and control during scope passage



Visible Distal End Tube

- ▶ Improved visibility
- ▶ Prevents inadvertent lodging of the tube in the scope working channel

Single-Use

- ▶ Designed to eliminate the risk of transmitting patient-to-patient disease



Multibite™
Single-Use
Biopsy Forceps

- ▶ The Multibite Multiple Sample Biopsy Forceps are designed specifically to collect tissue endoscopically for histologic examinations

ORDERING INFORMATION



Radial Jaw™ 4 Standard Capacity Biopsy Forceps



Order Number	Description	Jaw O.D. (mm)	Working Length (cm)	Minimum Working Channel (mm)	Color Code	Units
M00513383	Radial Jaw 4 Standard Capacity	2.2	160	2.8	Yellow	Box 5
M00513381	Radial Jaw 4 Standard Capacity	2.2	160	2.8	Yellow	Box 20
M00513382	Radial Jaw 4 Standard Capacity	2.2	160	2.8	Yellow	Box 40
M00513393	Radial Jaw 4 Standard Capacity w/Needle	2.2	160	2.8	Yellow	Box 5
M00513391	Radial Jaw 4 Standard Capacity w/Needle	2.2	160	2.8	Yellow	Box 20
M00513392	Radial Jaw 4 Standard Capacity w/Needle	2.2	160	2.8	Yellow	Box 40
M00513403	Radial Jaw 4 Standard Capacity	2.2	240	2.8	Orange	Box 5
M00513401	Radial Jaw 4 Standard Capacity	2.2	240	2.8	Orange	Box 20
M00513402	Radial Jaw 4 Standard Capacity	2.2	240	2.8	Orange	Box 40
M00513413	Radial Jaw 4 Standard Capacity w/Needle	2.2	240	2.8	Orange	Box 5
M00513411	Radial Jaw 4 Standard Capacity w/Needle	2.2	240	2.8	Orange	Box 20
M00513412	Radial Jaw 4 Standard Capacity w/Needle	2.2	240	2.8	Orange	Box 40



Radial Jaw 4 Large Capacity Biopsy Forceps



Order Number	Description	Jaw O.D. (mm)	Working Length (cm)	Minimum Working Channel (mm)	Color Code	Units
M00513303	Radial Jaw 4 Large Capacity	2.4	160	2.8	Yellow	Box 5
M00513301	Radial Jaw 4 Large Capacity	2.4	160	2.8	Yellow	Box 20
M00513302	Radial Jaw 4 Large Capacity	2.4	160	2.8	Yellow	Box 40
M00513313	Radial Jaw 4 Large Capacity w/Needle	2.4	160	2.8	Yellow	Box 5
M00513311	Radial Jaw 4 Large Capacity w/Needle	2.4	160	2.8	Yellow	Box 20
M00513312	Radial Jaw 4 Large Capacity w/Needle	2.4	160	2.8	Yellow	Box 40
M00513323	Radial Jaw 4 Large Capacity	2.4	240	2.8	Orange	Box 5
M00513321	Radial Jaw 4 Large Capacity	2.4	240	2.8	Orange	Box 20
M00513322	Radial Jaw 4 Large Capacity	2.4	240	2.8	Orange	Box 40
M00513333	Radial Jaw 4 Large Capacity w/Needle	2.4	240	2.8	Orange	Box 5
M00513331	Radial Jaw 4 Large Capacity w/Needle	2.4	240	2.8	Orange	Box 20
M00513332	Radial Jaw 4 Large Capacity w/Needle	2.4	240	2.8	Orange	Box 40



For single use only. Do not re-use.

All trademarks are the property of their respective owners.

CAUTION: The law restricts this device to sale by or on the order of a physician. Indications, Contraindications, Warnings and Instructions for Use can be found in the product labeling supplied with each device. Products shown for INFORMATION purposes only and may not be approved for sale in certain countries. This material not intended for use in France. RX Only.



Radial Jaw 4 Jumbo Biopsy Forceps



Order Number	Description	Jaw O.D. (mm)	Working Length (cm)	Minimum Working Channel (mm)	Color Code	Units
M00513343	Radial Jaw 4 Jumbo	2.8	160	3.2	Yellow	Box 5
M00513353	Radial Jaw 4 Jumbo w/Needle	2.8	160	3.2	Yellow	Box 5
M00513363	Radial Jaw 4 Jumbo	2.8	240	3.2	Orange	Box 5
M00513361	Radial Jaw 4 Jumbo	2.8	240	3.2	Orange	Box 20
M00513362	Radial Jaw 4 Jumbo	2.8	240	3.2	Orange	Box 40
M00513373	Radial Jaw 4 Jumbo w/Needle	2.8	240	3.2	Orange	Box 5
M00513371	Radial Jaw 4 Jumbo w/Needle	2.8	240	3.2	Orange	Box 20
M00513372	Radial Jaw 4 Jumbo w/Needle	2.8	240	3.2	Orange	Box 40



Radial Jaw 4 Gastropediatric Biopsy Forceps



Order Number	Description	Jaw O.D. (mm)	Working Length (cm)	Minimum Working Channel (mm)	Color Code	Units
M00513443	Radial Jaw 4 Gastropediatric	1.8	160	2.0	Yellow	Box 5
M00513453	Radial Jaw 4 Gastropediatric w/Needle	1.8	160	2.0	Yellow	Box 5



Multibite™ Multiple Sample Biopsy Forceps



Order Number	Description	Jaw O.D. (mm)	Working Length (cm)	Minimum Working Channel (mm)	Color Code	Units
M00510101	Multibite Multiple Sample Biopsy Forceps	2.4	160	2.8	Yellow	Box 5
M00510121	Multibite Multiple Sample Biopsy Forceps	2.4	240	2.8	Orange	Box 5

Boston Scientific

Advancing science for life™

Boston Scientific Corporation
300 Boston Scientific Way
Marlborough, MA 01752-1234
www.bostonscientific.com/gastro

©2016 Boston Scientific Corporation or its affiliates. All rights reserved.

ENDO-392404-AA June 2016

CE 2797