

Hemostasis

Hemostasis

Clipping Devices

Resolution™ Clip..... 41

Bipolar Electrohemostasis Catheters

Gold Probe™ Single-Use Electrohemostasis Catheters..... 41

Injection Gold Probe™ Single-Use Electrohemostasis /
Injection Catheters 41

Needles

Interject™ Clear Single-Use Injection Therapy Needle Catheters 42

Interject Contrast Single-Use Injection Therapy Needle Catheters 42

Optiflo™ Single-Use Hemostasis Catheters..... 43

Multiple Band Ligators

Speedband Superview Super 7™ Multiple Band Ligators..... 43


Generator

Endostat™ III Bipolar / Monopolar Electrosurgical Generator 44

Endostat Generator Accessories 44


Electrosurgical Bipolar Cable Adapter

Electrosurgical Bipolar Cable Adapter 45




Hemostasis

Clipping Devices





 Magnetic Resonance Conditional¹

Resolution™ Clip 				
Order Number	Working Length (cm)	Minimum Working Channel (mm)	Clip Opening (mm)	Packaging
M00522600	155	2.8	11	Each
M00522601	155	2.8	11	Box 10
M00522602	155	2.8	11	Box 20
M00522610	235	2.8	11	Each
M00522611	235	2.8	11	Box 10
M00522612	235	2.8	11	Box 20

¹Product with this symbol may be used in Magnetic Resonance environment according to the conditions described in Directions for Use.


Bipolar Electrohemostasis Catheters




Injection Gold Probe™ Single-Use Electrohemostasis / Injection Catheters 						
Order Number	Description	O.D. (Fr)	Needle Size Gauge	Working Length (cm)	Type Generator / Connector	Packaging
M00560150	Injection Gold Probe Catheter	7	25	210	Standard Plug*	Each
M00560160	Injection Gold Probe Catheter	10	25	210	Standard Plug*	Each

*Compatible with Endostat™ Bipolar-Monopolar Electrosurgical Generators.

For connections to other generators, please refer to the Electrosurgical Bipolar Cable Adapters, (see page 45).




Gold Probe™ Single-Use Electrohemostasis Catheters 						
Order Number	Description	O.D. (Fr)	Working Length (cm)	Type Generator / Connector	Packaging	
M00560070	Gold Probe Catheter	7	300	Standard Plug*	Each	
M00560100	Gold Probe Catheter	10	300	Standard Plug*	Each	
M00560220	Gold Probe Catheter	7	350	Standard Plug*	Each	
M00560071	Gold Probe Catheter	7	300	Standard Plug*	Box 5	
M00560101	Gold Probe Catheter	10	300	Standard Plug*	Box 5	

*Compatible with Endostat Bipolar-Monopolar Electrosurgical Generators.

For connections to other generators, please refer to the Electrosurgical Bipolar Cable Adapters, (see page 45).

Hemostasis

Needles




Interject™ Contrast Single-Use Injection Therapy Needle Catheters*



Order Number	Sheath Design	Needle Gauge	Maximum Needle Extension Length (mm)	Sheath O.D. (mm)	Working Length (cm)	Packaging
M005 1815 1	Contrast	23	4	1.8	200	Box 5
M005 1816 1	Contrast	25	4	1.8	200	Box 5
M005 1825 1	Contrast	23	6	1.8	200	Box 5
M005 1826 1	Contrast	25	6	1.8	200	Box 5
M005 1835 1	Contrast	23	4	2.3	240	Box 5
M005 1836 1	Contrast	25	4	2.3	240	Box 5

*With Star Catheter™ Technology.



Interject Clear Single-Use Injection Therapy Needle Catheters*




Order Number	Sheath Design	Needle Gauge	Maximum Needle Extension Length (mm)	Sheath O.D. (mm)	Working Length (cm)	Packaging
M005 1810 1	Clear	23	4	1.8	200	Box 5
M005 1811 1	Clear	25	4	1.8	200	Box 5
M005 1830 1	Clear	23	4	2.3	240	Box 5
M005 1831 1	Clear	25	4	2.3	240	Box 5

*With Star Catheter Technology.

Hemostasis

Needles (Continued)





Optiflo™ Single-Use Hemostasis Catheters



Order Number	Description	Needle Gauge	Sheath O.D. (mm)	Working Length (cm)	Working Channel (mm)	Packaging
M00511731	Optiflo Catheter	23	2.3	200	2.8	Box 5
M00511751	Optiflo Catheter	25	2.3	200	2.8	Box 5

Multiple Band Ligators




Speedband Superview Super 7™ Multiple Band Ligators



Order Number	Description	Required Working Channel (mm)	Scope O.D. (mm)	Bands / Ligating Unit	Packaging
M00542251	Speedband Superview Super 7 Band Ligator	2.8	8.6 – 11.5	7	Box 2
M00542253	Speedband Superview Super 7 Band Ligator	2.8	8.6 – 11.5	7	Box 4

Hemostasis

Generator



Endostat™ III Bipolar / Monopolar Electrosurgical Generator


Order Number	Description	Packaging
M00540200	Endostat III Bipolar / Monopolar Electrosurgical Generator <ul style="list-style-type: none"> - Power Cord w/ I.E.C. Molded Plug - Footswitch - Water Bottle - Pump Tubing (spare) - Fuse (spare) - Operating Manual 	Each

Endostat Generator Accessories

Order Number	Description	Packaging
M00540720	Water Bottle	Each
M00540240	Endostat III Dual Footswitch	Each
M00540760	Power Cord w/ I.E.C. Molded Plug	Each
M00540771	E-Z Load Pump Tubing – Endostat III Generator (for generators w/ E-Z Load Pump)	Box 5

Hemostasis

Electrosurgical Bipolar Cable Adapter



Electrosurgical Bipolar Cable Adapter*

Order Number	Description	Pin Spacing (mm)	Generator Compatibility**	Packaging
M00561440	Bipolar Cable Adapter	18	EMC (Coagulasem Dolley)	Each
M00561450	Bipolar Cable Adapter	22	ERBE, EMC Dolley, Olympus™ (ESG100 & PSD60)	Each
M00561460	Bipolar Cable Adapter	25	Lamidey (Optima 4)	Each
M00561470	Bipolar Cable Adapter	28	Valleylab™, Lamidey (Duo), Bovie (Aaron 3250™), Olympus (UES10, UES30 & Surgemaster), Bowa, Soring	Each
M00561480	Bipolar Cable Adapter	28 / Third Pin	ConMed™ (System 5000), Bovie (Aaron 2250™), Ethicon™	Each
M00561490	Bipolar Cable Adapter	32 / Third Pin	ConMed (Sabre 180, Excalibur, System 7550)	Each
M00561500	Bipolar Cable Adapter	39	Codman, Malis™	

*Bipolar Cable Adapters are used to connect Gold Probe™ and Injection Gold Probe™ Catheters with the above mentioned generators. Bipolar Cable Adapters are not required for the connection of Gold Probe and Injection Gold Probe to Endostat™ II or III Generators.

**For generators and their associated cable adaptors not listed above, please contact your Boston Scientific Representative.

Note: Monopolar active cords are in the Polypectomy Section, (see page 65).

Note: Some variations in generator faceplates may exist within a generator model. We recommend that you confirm pin spacing prior to ordering.

