

Pulmonary

Controlled Radial Expansion Pulmonary Balloon Dilators



CRE™ Single-Use Pulmonary Balloon Dilator*



Order Number	Balloon Inflated O.D.		Balloon Length (cm)	Inflation Pressures		Catheter Length (cm)	Packaging
	(mm)	(Fr)		(ATM)	(kPa)		
M005 50300	12-13.5-15	36-40.5-45	5.5	3-4.5-8	304-456-811	75	Each
M005 50310	15-16.5-18	45-49.5-54	5.5	3-4.5-7	304-456-709	75	Each
M005 50320	18-19-20	54-57-60	5.5	3-4.5-6	304-456-608	75	Each
M005 50330	8-9-10	24-27-30	3.0	3-5.5-9	304-557-912	75	Each
M005 50340	10-11-12	30-33-36	3.0	3-5-8	304-507-811	75	Each
M005 50350	12-13.5-15	36-40.5-45	3.0	3-4.5-8	304-456-811	75	Each

*Recommended Inflation Device: Alliance™ II Integrated Inflation System, Handle Order #5062, Sterile Syringe Order #5060 (see below).


Pulmonary Guidewires

Pulmonary Jagwire™ Pulmonary Guidewire



Order Number	Description	O.D. (in)	Length (cm)	Packaging
M005 15171	Pulmonary Jagwire Pulmonary Guidewire	.035	180	Box 2

Inflation System



Alliance II Inflation System

Order Number	Description	Packaging
M005 50620	Inflation Handle	Each
M005 50601	Single-Use Syringe / Gauge Assembly	Box 5