

Infection Prevention Portfolio

Compliance. Consistency. Care.

The Boston Scientific Infection Prevention portfolio of products is designed to provide the components needed to **comply with industry guidelines, mitigate cross-contamination risk, and reduce variability** during endoscope reprocessing.

OrcaPod™

Compliance EndoKit™

Blox™ Bite Blocks

CinchPad™

EndoArmor™ Gowns

EndoGlide+™
Lubricating Gel

HedgeHog™
Cleaning Brushes

Hydra™
Irrigation Products

Koala™ Sponge

Pure™ Enzymatic
Detergent



OrcaPod™

Seal™

Single Use Biopsy Valve

- Single Use to reduce the risk of cross-contamination between patients
- Designed to minimize leakage and maintains scope pressure while passing device

Orca™

Sterile Single Use Air/Water & Suction Valves

- Silicone free in accordance with manufacturer's recommendations
- Precision manufactured shaft and symmetrical center hole for a consistent fit, function and performance


Hydra™

Disposable Auxiliary Water Jet Connector


- Single Use – no lost components to replace
- One-way Valve designed to reduce the risk of cross-contamination

ORDERING INFORMATION

Order #	Description	Unit
SUV-617-50	Orca Sterile, Single Use Air/Water and Suction Valves	Box 50
SUV-628-50	OrcaPod Sterile, Single Use Air/Water and Suction Valves with Seal Biopsy Valve	Box 50
SUV-629-50	OrcaPod Sterile, Single Use Air/Water and Suction Valves, Seal Biopsy Valve and Hydra Water Jet Connector	Box 50
SBC-501-100	Seal Single Use Biopsy Valve	Box 100
SIT-361-100	Hydra Disposable Auxiliary Water Jet Connector	Box 100



Compliance EndoKit™




ORDERING INFORMATION

Name	SafeStart Kit	SafeStart CP	Compliance EndoKit	Compliance EndoKit	Compliance EndoKit	Compliance EndoKit	Compliance EndoKit	Compliance EndoKit
Item #	SIK-242	SIK-243	CEK-746	CEK-747	CEK-721	KS-3	KS-5	KS-16
Koala Sponge w/ Pure Imbibed	✓	✓	✓	✓	✓	✓	✓	✓
16 oz Bowl	✓	✓	✓	✓	✓	✓	✓	✓
Cinch Pad		✓	✓	✓	✓	✓	✓	✓
BioHazard Sticker		✓	✓	✓	✓	✓	✓	✓
Gauze (10)			✓	✓	✓	✓	✓	✓
EndoGlide+ 1.1 oz Lubricant			✓		✓			✓
EndoGlide+ 2.0 oz Lubricant						✓	✓	
HedgeHog Dual End Brush				✓	✓		✓	
HedgeHog Double End Brush						✓		✓



Work with your Boston Scientific Representative to find the right custom kit configuration for you.

Blox™ Bite Blocks



Designed for optimal access and patient comfort


- Available in three sizes to accommodate various patient and clinical needs
- Large lateral openings enhance digital guidance and provides easy access for suctioning
- Radiused dental rim fits comfortably to protect teeth, lips and gums
- Soft, adjustable strap, not made with natural rubber latex



ORDERING INFORMATION

Order #	Description	Unit
SBT-546-100	60Fr Blox Bite Block with Dental Rim	Box 100
SBT-114-100	54Fr Blox Bite Block with Dental Rim	Box 100
SBT-321-100	54Fr Blox Bite Block w/o Dental Rim	Box 100
SBP-236-40	38Fr Pediatric Blox Bite Block	Box 40

CinchPad™ Procedure and Endoscope Transport Pad



Designed to provide leak-proof containment to reduce bio-hazard risks and complies with society guidelines for endoscope transport

- Flexible drape size for full surface coverage for large and small cart surfaces
- Absorbent pad for fluid control and transportation protection
- Strong leak-proof liner to prevent leak through of fluids
- Cinching cords to convert the surface pad into your enclosed transport container




ORDERING INFORMATION

Order #	Description	Unit
SPX-017-25	CinchPad Endoscope Transport Pad	Box 25

EndoArmor™ Procedure Gowns

Designed by GI professionals to meet the comfort requirements for endoscopy procedures

- Breathable Microweave fabric
- Featherweight mesh back keeps staff cool
- MagicFit universal sizing




ORDERING INFORMATION

Order #	Description	Unit
SGE-368-100	EndoArmor GI Procedure Gown	Box 100
SGE-389-100	EndoArmor+ Procedure Gown	Box 100

Lubricants

EndoGlide+™ lubricating gel was specifically formulated by gastroenterologists for endoscopy


- 1.1 oz packet is designed to minimize waste and promote single use
- Viscosity facilitates scope passage in the anatomy



ORDERING INFORMATION


Order #	Description	Unit
SLT-394-100	1.1 oz EndoGlide+ Lubricating Gel	Box 100
SLT-395-50	2 oz EndoGlide+ Lubricating Gel	Box 50

HedgeHog™ Cleaning Brushes



Designed for maximum scrubbing and efficiency

- Available in various configurations to provide the right size and density for optimal scrubbing
- Stiffness designed to provide smooth scope passage
- Clinically evaluated to be as effective as OEM cleaning brushes¹



ORDERING INFORMATION


Order #	Description	Unit
SBS-227-50	HedgeHog Single-end Channel Brush	Box 50
SBD-228-50	HedgeHog Double-end Channel Brush	Box 50
SBD-289-50	HedgeHog Dual-end Channel/Valve Brush	Box 50
SBD-291-50	HedgeHog Sweeper Brush	Box 50
SBD-382-50	HedgeHog Double Bundle	Box 50

1. Data on file.

Hydra™ Irrigation System

Provides a comprehensive disposable water bottle and irrigation solution

- Carefully designed to address and reduce the risks associated with cross-contamination
- Designed to maintain optimal flow and performance
- Compatible with air or CO₂ insufflation



ORDERING INFORMATION

Hydra Irrigation Tubing (24-Hour Use)

Order #	Irrigation Pump Compatibility	Unit
SIT-355-15	Olympus™ AFU-100, ERBE™ EIP 2, Medivators™ EndoStratus™ & EndoGator™	Box 15
SIT-576-15	Olympus OFP2	Box 15

Hydra Disposable Auxiliary Water Jet Connectors (Single Use)

Order #	Scope Compatibility	Unit
SIT-361-100	Olympus 160, 180, & 190 Series	Box 100

Hydra Water Bottle Caps (24-Hour Use)

Order #	CO ₂	Air/Water Tubing Length	CO ₂ Tubing Length	Unit
SCT-350-10	No	16" (41 cm)	N/A	Box 10
SCT-351-10	Yes	16" (41 cm)	18" (46 cm)	Box 10

Koala™ Sponge

Designed with absorbent sponge cell technology to remove harmful bio-burden from endoscopes

- Ergonomically shaped design to help ensure full contact on insertion tubes
- Designed to protect the distal tip of an endoscope during transport




ORDERING INFORMATION

Order #	Description	Unit
SSB-210-50	Koala Cleaning Sponge	Box 50
SES-240-50	Koala Cleaning Sponge Imbibed w/ Pure	Box 50

Pure™ Enzymatic Detergent

Provides low foaming detergent in an environmentally friendly biodegradable solution

- Low foaming to be compliant with clinical guidelines
- 2 minute activation time for rapid cleaning power
- Environmentally friendly biodegradable solution
- 1/2 ounce per gallon concentration



ORDERING INFORMATION

Order #	Description	Unit
SEE-220-4	Pure Enzymatic Solution (spearmint scented)	1 Gallon Box 4
SEE-573-4	Pure Enzymatic Solution (fragrance free)	1 Gallon Box 4



**Boston
Scientific**
Advancing science for life™

Boston Scientific Corporation
300 Boston Scientific Way
Marlborough, MA 01752-1234
www.bostonscientific.com

©2018 Boston Scientific Corporation
or its affiliates. All rights reserved.

ENDO-515216-AC

All images taken by Boston Scientific.

Registered trademarks are the property of their respective owners.

Indications, contraindications, warnings and instructions for use can be found in the product labeling supplied with each device.

CAUTION: Federal (US) law restricts these devices to sale by or on the order of a physician.