



FARAPULSE™

Pulsed Field Ablation System






FARAWAVE™ Ordering Information

Reference Number	
M004PF41M401	31 mm (distal diameter)
M004PF41M402	35 mm (distal diameter)
Catheter Specifications	
Catheter type	Over-the-wire
Shaft diameter	12F
Usable length	115 cm
Number of electrodes	20
Compatibility	
Guidewire	0.035" diameter, 180 cm minimum length
Sheath	FARADRIVE™ (REF M004PF21M402)
Accessories	
FARASTAR catheter connection cable	REF M004PF41M434

FARAWAVE™ is sterile and single-use only





FARADrive™ Ordering Information

Reference Number	
M004PF21M402	Clear shaft
Sheath Specifications	
Usable length	74 cm
Sheath overall length	91 cm
Inner diameter	13F
Deflection	Unidirectional, 35 mm reach
Radio-opaque marker	Distal end
Dilator Specifications	
Total length	97 cm
Usable length	94 cm
Compatibility	
Guidewires	Maximum diameter 0.035"
Pulsed Field Ablation Catheter	FARAWAVE™ (REF M004PF41M401, M004PF41M402)
VersaCross™ integration	VersaCross™ Connect Access Solution for FARADrive (180 cm Pigtail, VXAK0041) (180 cm J-Tip, VXAK0045)

FARADrive™ Steerable Sheath and dilator are sterile and single-use only





FARASTAR™ Ordering Information

Reference Number	
M004PF61M401	FARASTAR Pulsed Field Ablation Generator
M004PF61M407	FARASTAR Recording System Module
Shipping Information – FARASTAR Pulsed Field Ablation Generator	
Weight	552 lbs (250 kg)
Dimensions	L 120 x H 200 x W 80 cm
Shipping Information – FARASTAR Recording System Module	
Weight	111 lbs (50 kg)
Dimensions	L 120 x H 70 x W 80 cm
Technical Specifications	
Input voltage	120 V / 60 Hz
Output voltage	Proprietary FARAPULSE biphasic, bipolar waveform
Weight	243 lbs (110 kg)
Dimensions – FARASTAR Pulsed Field Ablation Generator	L 79 x W 61 x H 117 cm (160 cm with monitor)
Dimensions – FARASTAR Recording System Module	L 41 x W 19 x H 46 cm
User Interface	Articulating touchscreen
Accessories	
M004PF61M404	FARASTAR Stimulation Module Cable
M004PF61M405	FARASTAR EGM Cable
M004PF61M409	FARASTAR Stimulation Module Male Cable
M004PF61M410	FARASTAR Stimulation Module Female Cable
M004PF61M411	FARASTAR Stimulation Module Y-Cable Long
M004PF61M412	FARASTAR Stimulation Module Y-Cable Short
M004PF61M413	FARASTAR Recording System Module Catheter Pin Cable
M004PF61M415	FARASTAR Recording System Module ECG Trunk Cable
M004PF61M416	FARASTAR Recording System Module ECG Output Module
M004PF61M417	FARASTAR Recording System Module EGM Input Module