

LithoVue™ System

Single-Use Digital
Flexible Ureteroscope

Design features

- 270° deflection in both directions
- 7.7F tip diameter
- 9.5F [$\leq 3.23\text{mm}$] outer diameter
- 3.6F ID working channel
- Digital CMOS imager
- Working distance of 2mm–50mm
- Light source built into the handle
- Flexible sheath
- Integrated camera head
– no secondary external attachments required
- Compatible with laser lithotripsy
– uses existing technologies and familiar surgical tools
- Mobile cart for portability
- All-in-one touchscreen PC
– includes monitor workstation, image processing software





Caution: Federal (U.S.) law restricts this device to sale by or on the order of a physician. Please refer to package insert provided with the product for complete indications for use, contraindications, warnings, precautions, adverse events, and instructions prior to using these products.

Bench test results on file with Boston Scientific. Bench test results may not necessarily be indicative of clinical performance.

All trademarks are the property of their respective owners.

Boston Scientific Corporation
300 Boston Scientific Way
Marlborough, MA 01752
www.bostonscientific.com

© 2015 by Boston Scientific Corporation
or its affiliates. All rights reserved.

URO-340302-AC DEC 2015