

Navigator™ HD

Ureteral Access Sheath

Boston
Scientific

Seamless transition
for ease of insertion



Navigator™ HD Ureteral Access Sheath

The Navigator HD Ureteral Access Sheath provides ureteral dilation and a working channel for the introduction of ureteroscopes and devices during ureteroscopic procedures. The new flexible tip, higher durometer dilator and sheath and additional sizes provide solutions for your access needs.

Benefits of a Ureteral Access Sheath

- ▶ Protect the ureter during repeated instrument exchanges, reducing the potential for trauma.
- ▶ The continuous channel also protects delicate instruments and the flexible ureteroscopes from damage, reducing costly repairs.
- ▶ Clinically proven to protect the urinary system, reduce O.R. time, and reduce the need for balloon dilators.
- ▶ Facilitates the navigation of tortuous anatomy.
- ▶ Aids in the washing out of stone fragments, facilitating complex ureteroscopic procedures and preventing over-distension of the collecting system.
- ▶ Reduces or eliminates need for balloon dilators.



Navigator™ HD Ureteral Access Sheath Features

Kink resistant sheath – reinforced with a stainless steel coil, designed to provide stability and strength through tortuous anatomy.

Hydrophilic coating – on sheath and dilator allows for ease of insertion. Improved coating is designed for durability in bilateral cases.

Radiopaque marker band – confirms placement of sheath tip under fluoroscopy.

Dual durometer dilator – stiffer body to facilitate placement in tortuous anatomy with flexible tip for reduced trauma during insertion and dilation.

Innovative hub design – secure single-latch mechanism that allows for one handed latch engagement without rotational alignment.



flexible

Flexible, single-tapered tip designed for reduced trauma during insertion and dilation

seamless

Seamless transition from dilator to sheath for ease of insertion

Navigator™ HD Ureteral Access Sheath

Order Number	Sheath Inner (Fr)	Outer Diameter (Fr)	Length (cm)
M0062502210	11	13	28
M0062502220	11	13	36
M0062502230	11	13	46
M0062502240	12	14	28
M0062502250	12	14	36
M0062502260	12	14	46
M0062502270	13	15	28
M0062502280	13	15	36
M0062502290	13	15	46



Boston Scientific
Advancing science for life™

Boston Scientific Corporation
One Boston Scientific Place
Natick, MA 01760-1537
www.bostonscientificstone.com

Ordering Information
1.888.272.1001

©2013 Boston Scientific Corporation
or its affiliates. All rights reserved.

URO-118610-AA XM 1/13

Navigator is a registered or unregistered trademark of Boston Scientific Corporation or its affiliates. All other trademarks are property of their respective owners.

Caution: Federal Law (USA) restricts this device to sale by or on the order of a physician. Refer to package insert provided with this product for complete Indications for Use, Contraindications, Warnings, Precautions, Adverse Events, and Instructions prior to using this product.