



**Boston
Scientific**
Advancing science for life™

UroMax Ultra™

Balloon Dilatation Catheter



Designed for radial dilatation of the urinary tract to relieve strictures and to facilitate the placement of larger surgical instruments.

The UroMax Ultra™ Balloon Dilatation Catheter features two lumen hubs for placement and infusion of contrast medium during balloon inflation.¹

- Low Profile — proprietary quadra-fold technology allows low profile during placement and withdrawal
- High Pressure Dilatation — controlled inflation pressure up to 20 ATM for dilating narrow strictures
- Non-Compliant Balloon — designed with specially treated polymers and processes to maintain its outer diameter during inflation
- Two lumen hubs for placement and infusion of contrast medium during balloon inflation



Why Radial Dilatation?



Single stage

Balloon dilatation may allow the procedure to be completed at a single stage and avoid additional medical costs and side effects of a staged, pre-stented approach.² Radial dilatation may help avoid additional procedures and subsequent financial burden to the hospital system and patient.



Focused Dilatation

Balloon dilatation of strictures under visualization may avoid the potential shearing forces associated with sequential rigid dilatation.³

Ordering Information

UPN	Size	Configuration
M0062251000	12F x 4cm	Single Balloon
M0062251010	15F x 4cm	Single Balloon
M0062251020	18F x 4cm	Single Balloon
M0062251030	21F x 4cm	Single Balloon
M0062251040	24F x 4cm	Single Balloon
M0062251050	30F x 4cm	Single Balloon
M0062251060	12F x 6cm	Single Balloon
M0062251070	15F x 6cm	Single Balloon
M0062251080	18F x 6cm	Single Balloon
M0062251090	21F x 6cm	Single Balloon
M0062251100	24F x 8cm	Single Balloon
M0062251110	30F x 8cm	Single Balloon
M0062251150	12F x 10cm	Single Balloon
M0062251160	15F x 10cm	Single Balloon
M0062251170	18F x 10cm	Single Balloon
M0062251180	21F x 10cm	Single Balloon

UPN	Size	Configuration
M0062251200	12F x 4cm	Kit
M0062251210	15F x 4cm	Kit
M0062251220	18F x 4cm	Kit
M0062251230	21F x 4cm	Kit
M0062251240	24F x 4cm	Kit
M0062251250	30F x 4cm	Kit
M0062251260	12F x 6cm	Kit
M0062251270	15F x 6cm	Kit
M0062251280	18F x 6cm	Kit
M0062251290	21F x 6cm	Kit
M0062251300	24F x 8cm	Kit
M0062251310	30F x 8cm	Kit
M0062251350	12F x 10cm	Kit
M0062251360	15F x 10cm	Kit
M0062251370	18F x 10cm	Kit
M0062251380	21F x 10cm	Kit

Packaging Configurations

- Available as a kit or a single balloon
- Kit includes one UroMax Ultra™ Balloon, Encore™ 26 Inflator, and two-way stopcock



Please consult your sales representative for more information and ordering details. To learn more visit www.bostonscientific.com/UroMax

1. The UroMax Ultra™ Balloon Dilatation Catheter device has been tested with passing results according to ASTM D6499-032.
2. Kuntz NJ, Neisius A, Tsivian M, et al. Balloon dilation of the ureter: A contemporary review of outcomes and complications. *J Urol.* 2015 Aug;194(2):413-7.
3. Stamatiou K, Christopoulos G, Moschouris H, et al. Treatment of urethral strictures with balloon dilation: A forgotten tale. *Arch Ital Urol Androl.* 2015 Sep 30;87(3):254-5.

Caution: U.S. Federal law restricts this device to sale by or on the order of a physician.

CAUTION: The law restricts these devices to sale by or on the order of a physician. Indications, contraindications, warnings and instructions for use can be found in the product labelling supplied with each device. Information for use only in countries with applicable health authority registrations. Material not intended for use in France.

Products shown for INFORMATION purposes only and may not be approved or for sale in certain countries. Please check availability with your local sales representative or customer service.

All images are the property of Boston Scientific. All trademarks are the property of their respective owners.

**Boston
Scientific**
Advancing science for life™

Boston Scientific Corporation
300 Boston Scientific Way
Marlborough, MA 01752
www.bostonscientific.com

© 2018 Boston Scientific Corporation
or its affiliates. All rights reserved.

URO-567313-AA OCT 2018