




Be a Leader in Bariatric and Gastrointestinal Procedures


OverStitch™
Endoscopic Suturing System



The OverStitch™ Endoscopic Suturing System enables advanced endoluminal surgery by allowing physicians to place **full-thickness sutures** through a flexible endoscope. With the capability to **customize suture patterns**, OverStitch™ facilitates a **broad range** of bariatric and gastrointestinal procedures.



Full thickness sutures incorporate muscle layer to help increase tissue approximation durability.



Customize suture patterns to address simple or complex gastrointestinal and bariatric procedures

Catalog Information

Product Family	Code	Description	Package Qty.
OverStitch™	ESS-G02-160	OverStitch™ Endoscopic Suturing System	3 units
	PLY-G02-020-APL	OverStitch™ 2-0 Polypropylene Suture	12 units
	CNH-G01-000	OverStitch™ Suture Cinch	6 units
Accessories	THX-165-028	Tissue Helix, 165 cm Length, 2.8mm OD	6 units
	OVT-027-160	OverTube™ Endoscopic Access System	3 units

To learn more, visit overstitch.apolloendo.com/learn-more

CAUTION: The law restricts these devices to sale by or on the order of a physician. Indications, contraindications, warnings and instructions for use can be found at www.apolloendo.com/dfus. Products shown for INFORMATION purposes only and may not be approved or for sale in certain countries. This material is not intended for use in France. 2024 Copyright © Boston Scientific Corporation or its affiliates. All rights reserved.

ENDO-1737407-AA

CE 2797