


ZERO EXCHANGE


to your therapy delivery
in a SINGLE SOLUTION



VersaCross Connect™ Access Solution for FARADRIE™

DESIGNED for use with the FARADRIE™ Steerable Sheath†

3-IN-1 RF WIRE
for use as a starter wire, puncture device, and exchange rail



†FARADRIE™ Steerable Sheath as shown in this brochure is sold separately.

*The therapy referred to here may be one of several left heart procedures that require transseptal access. Additional procedural steps and/or devices may be required to deliver the therapy. Before use, consult Instructions for Use for any devices accordingly.
†FARADRIE™ Steerable Sheath as shown in this brochure is sold separately.
‡The VersaCross Connect Transseptal Dilator is for use with a 13F (4.31 mm) ID FARADRIE™ Steerable Sheath which is 74cm in length, specifically, model: M004PF21M402 (USA), M004PFCE21M402 (Canada).

All trademarks are property of their respective owners. Patents Pending and/or issued. CAUTION: The law restricts this device to sale by or on the order of a physician. Indications, Contraindications, Warnings, and Instructions For Use can be found in the product labelling supplied with each device or at www.IFU-BSCI.com.

Products shown for INFORMATION purposes only and may not be approved or for sale in certain countries. This material not intended for use in France.

Boston Scientific is a Global Company. Please note that model numbers, indications, contraindications, warnings and specifications may differ depending on geographic region. Not all information displayed in this brochure may be licensed in accordance with Canadian law. Please contact your Boston Scientific representative for local labeling, product specifications and licensed model numbers.



For brief summaries, scan QR code or visit qrco.de/bs-vxafd

Boston Scientific
Advancing science for life™


Baylis Medical Company Inc.
5959 Trans-Canada Highway
Montreal, QC Canada H4T 1A1
www.bostonscientific.com

Customer Service
United States: 1-888-272-1001
Canada: 1-888-359-9691

© 2024 Boston Scientific Corporation or its affiliates. All rights reserved.
EP-1782905-AB

Boston Scientific
Advancing science for life™

VersaCross Connect™ Access Solution for FARADRIE™




Driving Peace of Mind

ZERO EXCHANGE
FARADRIE™ Steerable Sheath delivery
For use with the FARAPULSE™ Pulsed-Field Ablation System


INTEGRATED PREDICTABILITY

with the **VersaCross Connect™** Access Solution for **FARADrive™**


EXCHANGELESS 3-in-1 RF Wire



1 Start
by tracking up to the SVC



2 Access
with purpose-built RF




3 Deliver
with confidence


Predictable and Efficient Transseptal

for smooth **FARADrive™** Steerable Sheath delivery

DELIVER THERAPY ON TARGET



TRUform™
Shapeable Technology



Precise RF
Puncture Technology



OMNIVIZ™ Technology

Reliably locate the **VersaCross Connect™** solution on fluoroscopy and ultrasound. Positional markers indicate position of RF tip within dilator. Connect to DuoMode™ Cable to track and mark RF tip on mapping systems.

VersaCross Connect™

Access Solution for **FARADrive™**

Personalize your solution

1. Choose your VersaCross™ RF Wire



2. Connect with VersaCross Connect™ Transseptal Dilator



Ordering information

Part. No.	Wire	Wire length	Dilator curve	Dilator length
VXAK0041	Pigtail	180 cm	D1	93 cm
VXAK0045	J-Tip	180 cm	D1	93 cm